


Supplementary Figure S1

(a) **PC3**



(b) **HEK-293**




Figure S1. CCK-8 assays in PC3 and HEK-293 cells. (a,b) Relative absorbances in PC3 and HEK-293 cells. Cells were treated with auraptene (or vehicle) for 24 or 48 h at the indicated concentrations, and then absorbances at 450 nm were measured 3 h after adding CCK-8 solution. The absorbances (after background subtraction) relative to auraptene untreated cells were calculated ($n = 4$) (** $p < 0.01$, *** $p < 0.001$).