

Supplementary Table 1. Comparison between Spadaro's and manufacture's methods in wine matrix.

Clean-up method	Validation levels (µg/L)	Recovery (%) (n = 3)	RSD (%) (n = 3)	LOD (µg/L) (n = 3)	LOQ (µg/L) (n = 3)
Method of Spadaro <i>et al.</i>	2	84.03	19.22	0.14	0.48
	10	102.24	15.51		
	20	93.07	10.53		
Method of manufacturer	2	82.31	17.84	0.12	0.41
	10	99.15	12.26		
	20	93.96	13.95		