

Supplemental Materials

Berberine-Coated Biomimetic Composite Microspheres for Simultaneously Hemostatic and Antibacterial Performance

Xiaojian Zhang ¹, Kaili Dai ¹, Chenyu Liu ¹, Haofeng Hu ², Fulin Luo ², Qifan Qi ², Lei Wang ^{2,3},
Fei Ye ^{2,3}, Jia Jin ^{2,3,*}, Jie Tang ^{1,*} and Fan Yang ^{1,*}

¹ Shanghai Engineering Research Center of Molecular Therapeutics and New Drug Development, East China Normal University, Shanghai 200062, China; 51184300155@stu.ecnu.edu.cn (X.Z.); 15988810463@163.com (K.D.); 51194300147@stu.ecnu.edu.cn (C.L.)

² College of Life Sciences and Medicine, Zhejiang Sci-Tech University, Hangzhou 310018, China; 201920201014@mails.zstu.edu.cn (H.H.); 202020801049@mails.zstu.edu.cn (F.L.); 2018339901070@mails.zstu.edu.cn (Q.Q.); leiwang1986@hotmail.com (L.W.); yefei@zstu.edu.cn (F.Y.)

³ Zhejiang Provincial Key Laboratory of Silkworm Bioreactor and Biomedicine, Hangzhou 310018, China

* Correspondence: jjin@zstu.edu.cn (J.J.); jtang@chem.ecnu.edu.cn (J.T.); fyang@chem.ecnu.edu.cn (F.Y.)

Supplemental Materials

Table S1. The proportion of hemolytic experimental materials.

	NC (0.9%NaCl)	PC (ddH ₂ O)	TS (test sample)				
			BAG	BACG- 1B	BACG-3B	BACG-6B	CMPHP
2% RBCs/mL	0.5	0.5	0.5	0.5	0.5	0.5	0.5
0.9%NaCl/mL	0.5		0.44	0.44	0.44	0.44	0.44
Distilled water/mL		0.5					
Test group/mL			0.06	0.06	0.06	0.06	0.06

NC: negative control (0.9% NaCl solution); PC: positive control (distilled water).