

Figure S1. Exemplar agarose gel showing PCR-RFLP banding patterns for (A) the single enzyme (TaqI) assay, and (B) the double digest (BsaI + DdeI) assay, both of which can be used to distinguish between two *O. disjunctus* mtDNA clades. Lanes are labeled as follows: L, 100-bp ladder; U/C, uncut PCR product; 1–6: sample lanes that contain individuals of known clade membership (1–4 and 6 = clade A, 5 = clade B). Known fragment sizes (bp) derived from the ladder are labeled on the left, and expected fragment sizes are on the right.