

Protozoan parasite	Commercial kit	Molecular method	Positive results	References
<i>Cryptosporidium</i> spp	QIAamp DNA stool mini kit (Qiagen Hilden, Germany)	Nested PCR	93.8% (n=32)	[18]
<i>Giardia duodenalis</i>	Wizard magnetic DNA purification system (Promega Corporation, Madison, WI, USA)	Real-time PCR	100% (n=47)	
<i>Entamoeba</i> spp	PowerFecal DNA isolation kit (MoBio, Carlsbad, CA, USA)	Real-time PCR	93.8% (n=13)	
<i>Cryptosporidium hominis</i>	ZR Fecal DNA MiniPrep™ Kit (Zymo Research, Irvine, CA)	Real-time PCR	100% (n=12)	[19]
<i>Cryptosporidium parvum</i>	ZR Fecal DNA MiniPrep™ Kit (Zymo Research, Irvine, CA)	Real-time PCR	100% (n=11)	
<i>Cryptosporidium felis</i>	ZR Fecal DNA MiniPrep™ Kit (Zymo Research, Irvine, CA)	Real-time PCR	100% (n=3)	
<i>Blastocystis</i> sp.	ZR Fecal DNA MiniPrep™ Kit (Zymo Research, Irvine, CA)	PCR	94% (n=50)	[20]
<i>Enterocytozoon bienuesi</i>	QIAamp stool mini kit (Qiagen Hilden, Germany)	PCR	100% (n=290)	[21]
<i>Toxoplasma gondi</i>	RTA DNA Isolation Kit from stool (Istambul, Turkey)	LAMP	100% (n=41)	[22]