



Article

# Culture-Independent PCR Detection and Differentiation of *Mycobacteria* spp. in Antemortem Respiratory Samples from African Elephants (*Loxodonta Africana*) and Rhinoceros (*Ceratotherium Simum*, *Diceros Bicornis*) in South Africa

Wynand J. Goosen <sup>1,\*</sup>, Charlene Clarke <sup>1</sup>, Léanie Kleynhans <sup>1</sup>, Tanya J. Kerr <sup>1</sup>, Peter Buss <sup>2</sup> and Michele A. Miller <sup>1</sup>

- DSI-NRF Centre of Excellence for Biomedical Tuberculosis Research, South African Medical Research Council Centre for Tuberculosis Research, Division of Molecular Biology and Human Genetics, Faculty of Medicine and Health Sciences, Stellenbosch University, P.O. Box 241, Cape Town 8000, South Africa; cclarke@sun.ac.za (C.C.); leaniek@sun.ac.za (L.K.); tjkerr@sun.ac.za (T.J.K.); miller@sun.ac.za (M.A.M.)
- Veterinary Wildlife Services, Kruger National Park, South African National Parks, Skukuza 1350, South Africa; peter.buss@sanparks.org
- \* Correspondence: wjgoosen@sun.ac.za

**Abstract:** Since certain *Mycobacterium tuberculosis* complex (MTBC) members, such as *M. bovis*, are endemic in specific South African wildlife reserves and zoos, cases of clinically important nontuberculous mycobacteria (NTM) in wildlife may be neglected. Additionally, due to the inability of tests to differentiate between the host responses to MTBC and NTM, the diagnosis of MTBC may be confounded by the presence of NTMs. This may hinder control efforts. These constraints highlight the need for enhanced rapid detection and differentiation methods for MTBC and NTM, especially in high MTBC burden areas. We evaluated the use of the GeneXpert MTB/RIF Ultra, the Hain CM*direct* V1.0 line probe assay, and novel amplicon sequencing PCRs targeting the mycobacterial *rpoB* and *ku* gene targets, directly on antemortem African elephant (n = 26) bronchoalveolar lavage fluid (BALF) (n = 22) and trunk washes (n = 21) and rhinoceros (n = 23) BALF (n = 23), with known MTBC culture-positive and NTM culture-positive results. Our findings suggest that the Ultra is the most sensitive diagnostic test for MTBC DNA detection directly in raw antemortem respiratory specimens and that the *rpoB* PCR is ideal for *Mycobacterium* genus DNA detection and species identification through amplicon sequencing.

**Keywords:** African elephants; bronchioloalveolar lavage; GeneXpert MTB/RIF Ultra; Hain CM*direct* V1.0 LPA; *ku* PCR; *mycobacterium tuberculosis* complex; non-tuberculous mycobacteria; rhinoceros; *rpoB* PCR; trunk wash

# check for updates

Citation: Goosen, W.J.; Clarke, C.; Kleynhans, L.; Kerr, T.J.; Buss, P.; Miller, M.A. Culture-Independent PCR Detection and Differentiation of *Mycobacteria* spp. in Antemortem Respiratory Samples from African Elephants (*Loxodonta Africana*) and Rhinoceros (*Ceratotherium Simum*, *Diceros Bicornis*) in South Africa. *Pathogens* 2022, 11, 709. https://doi.org/10.3390/pathogens11060709

Academic Editors: Mitchell V. Palmer, Paola M. Boggiatto and Carly Kanipe

Received: 27 May 2022 Accepted: 18 June 2022 Published: 20 June 2022

**Publisher's Note:** MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2022 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https://creativecommons.org/licenses/by/4.0/).

#### 1. Introduction

Mycobacteria are a diverse group of microorganisms that can be found in practically every environmental niche [1]. They are divided into two categories: *Mycobacterium tuberculosis* complex (MTBC) and non-tuberculous mycobacteria (NTM) [2]. *Mycobacterium tuberculosis*, *M. bovis*, *M. orygis*, *M. bovis* bacillus Calmette–Guerin (BCG), *M. africanum*, *M. cannettii*, *M. pinnipedii*, *M. caprae*, and *M. microti* are all members of the MTBC. *Mycobacterium bovis* and *M. tuberculosis* are well known for causing chronic infectious disease in people and animals, including livestock and wildlife species [3–6]. Cases of *M. bovis* in different South African (SA) wildlife, including African buffaloes, African lions, African wild dogs, and rhinoceros, are sporadically reported. The disease control efforts for buffaloes include test-and-slaughter strategies, and quarantine for endangered species [3,5,7–9]. Many of the high TB burden countries are heavily dependent on animal-related industries, such as tourism and agriculture, leading to habitat encroachment and increased opportunities

Pathogens 2022, 11, 709 2 of 16

for disease transmission at the animal–human interface. This phenomenon has recently been highlighted by the unexpected discovery of a fatal *M. tuberculosis* infection in an African elephant (*Loxodonta africana*) in Kruger National Park (KNP), which is endemic for *M. bovis* [10,11].

Non-tuberculous mycobacteria, also known as environmental mycobacteria, are natural inhabitants of the environment and comprise more than 150 species listed on public bacterial databases (http://www.bacterio.net (accessed on 19 May 2022)). Interestingly, more than a third of these species have been implicated in diseases in livestock, wildlife, and humans [12]. Numerous animal and human NTM infections are reported globally, probably due to their ubiquitous presence and opportunistic nature [13,14]. The Mycobacterium avium complex (MAC) organisms, consisting of M. avium subsp. avium, M. avium subsp. paratuberculosis, M. avium subsp. hominissuis, M. intracellulare, M. sylvaticum, M. colombiense, M. bouchedurhonense, M. timonense, M. chimaera, M. arosiense, M. yongonense, and M. marseillense, are the most well-known opportunistic NTMs reported in animals and humans [15]. In addition, M. kansasii, M. marinum, and M. ulcerans commonly cause opportunistic NTM disease. Infections with M. kansasii can cause pathological signs resembling the M. bovis-associated disease in animals [12,16]. Another NTM, M. scrofulaceum has been found in the lymph nodes of cattle, buffaloes, farmed deer, swine, wild pigs, patas monkeys, fish, and mice [17]. These organisms can spread indirectly through contaminated environments, as is evident by the isolation of mycobacteria from cattle and mice faeces [17].

In Africa, most of the NTMs are reported as incidental findings in livestock with gross pathological lesions identified in slaughterhouses, when tissues are submitted for M. bovis surveillance [18,19]. There are substantially fewer reports of NTMs isolated from free-ranging wildlife. Therefore, the extent and distribution of NTM infections among the African wildlife species are largely unknown. In SA, the M. avium subspecies paratuberculosis, M. terrae, M. nonchromogenicum, M. vaccae/M. vanbaalenii, and unidentified species closely related to M. moriokense, have been isolated from African buffaloes, livestock, and their environments, indicating that NTMs may be exchanged at environment-animal interfaces [20]. Besides being opportunistic pathogens, some of the NTM species may colonize the host without development of disease, but instead priming the host's immune system, confounding the immuno-diagnosis of bovine tuberculosis (bTB) due to cross-reactivity to shared antigens [21,22]. The genetic investigations into NTM isolates from the samples collected from various wildlife species showed the presence of orthologous genes, such as ESAT-6 and CFP-10, situated in the ESX-1 to ESX-5 regions [2]. The occurrence of positive M. bovis immunological test results in these animals suggest that NTMs potentially interfere with bTB diagnostic assays. Since M. bovis is endemic in several of the SA wildlife reserves and zoos, the cases of NTM infection in wildlife may be neglected, due to the primary focus being on M. bovis' diagnosis. However, due to the inability of tests to differentiate the host responses towards MTBC and NTMs, due to the interference of NTMs, this may hinder control efforts. Unfortunately, correctly identifying a NTM requires mycobacterial culture with a minimum 6–8-week incubation period and further genetic speciation [23]. Previously, using elephants' and rhinoceros' respiratory samples, mycobacterial culture was shown to have a lower limit of detection of 1000 colony forming units (CFU) for culturing M. bovis, and 100 CFU for M. tuberculosis [23,24]. This constraint highlights the need for enhanced and rapid MTBC and NTM detection and differentiation assays, especially in the high *M. bovis* burden areas.

Assays using PCR have been used to rapidly detect and differentiate between MTBC and NTM. Recently, Cepheid's GeneXpert MTB/RIF Ultra assay (Ultra) has been shown to provide rapid detection of MTBC DNA in the tissue and respiratory samples collected from infected African buffaloes, African elephants, and rhinoceros [5,6]. The Hain GenoType CM*direct* VER 1.0 line probe assay (Hain LPA) can also lead to sensitive and specific detection of the DNA from the *Mycobacterium* genus organisms, MTBC (without differentiation), and the differentiation of more than 20 clinically relevant NTMs directly from human

Pathogens 2022, 11, 709 3 of 16

patient specimens [25]. Although 16S rRNA sequencing of Mycobacterium cultures has been used for speciation, there are increasing reports that this technique is sub-optimal [26]. More recent studies, using hundreds of strains of Mycobacterium spp. isolates, have demonstrated that the combined PCR amplification of highly conserved regions of the Mycobacterium ku and rpoB genes, using primers specifically designed for the Mycobacterium genus, have superior performance compared to the 16S rRNA amplicon sequencing for genus detection and the speciation of low-level mixed microbial populations [26-29]. Consequently, these findings have major implications for the field of targeted, next generation sequencing directly from clinical samples, especially for the improved surveillance of clinically important Mycobacteria spp. Therefore, the aims of this pilot study were to evaluate the use of (1) the Ultra and Hain LPA for rapid MTBC DNA detection; (2) the Hain LPA for NTM DNA detection and species differentiation without required sequencing; and (3) ku and rpoB amplicon sequencing for the *Mycobacterium* genus DNA detection and species differentiation, directly from antemortem respiratory samples. These samples were collected from known MTBC and NTM culture-confirmed positive African elephants and rhinoceros from KNP and a zoo in SA.


#### 2. Results

#### 2.1. Mycobacterial Culture Results

The respiratory samples from six animals (five elephants and one white rhinoceros) were defined as the culture-confirmed MTBC positive cohort for this study (Figure 1 and Table 1). Six of the eight elephant respiratory samples, three bronchoalveolar lavage fluid (BALF) samples (including one duplicate sample), and three trunk washes from five MTBC-infected elephants (free-ranging and zoo), were confirmed to contain viable MTBC through mycobacterial culture and speciation. The remaining two respiratory samples (TW and BALF) did not grow MTBC, but rather M. elephantis and M. stomatepiae were isolated (Table 1). The infected animals included: (1) two zoo elephants (18/85 and 18/177), from which M. tuberculosis was cultured from the first animal's (18/85) BALF sample, M. tuberculosis and M. africanum from the second elephant's (18/177) duplicate BALF samples and M. elephantis from its TW sample, respectively; and (2) three free-ranging KNP elephants (18/527, 18/533, and 18/538), from which M. bovis was cultured from each animal's TW sample (Table 1). Even though the KNP elephant 18/527 TW sample contained viable culturable M. bovis, its BALF sample did not (Table 1). Lastly, M. bovis was also cultured from a BALF sample collected from a white rhinoceros (19/46) from KNP (Table 1).

Respiratory samples from nine animals (eight elephants and one white rhinoceros) were defined as the culture-confirmed NTM positive cohort for this study (Figure 1 and Table 1). Ten respiratory samples (seven TW and three BALF) from eight elephants and a single (BALF) sample from one rhinoceros were confirmed to contain viable NTMs through mycobacterial culture (Table 1). The *Mycobacterium mantenii* was successfully isolated and characterized from a KNP elephant's (19/460) BALF sample and the *M. abscessus* strain from its TW sample. Similarly, *M. interjectum* was isolated and characterized from a KNP elephant's (18/530) BALF sample and the *M. avium* complex from its TW sample (Table 1). Non-tuberculous mycobacteria isolated from the TW samples collected from the remaining four KNP elephants and one zoo elephant included the *M. avium* complex strain (18/173), *M. mageritense* strain (18/532), *M. intracellulare* (18/534 and 18/539), and *M. fortuitum* strain (21/496). The elephant (18/176) and rhinoceros (18/31) BALF samples were culture positive for the *M. foliorum* and *M. scrofulaceum* strain, respectively (Table 1).

Pathogens 2022, 11, 709 4 of 16



**Figure 1.** Study method flow chart for African elephants (n = 26) and rhinoceros (n = 23) respiratory sample processing and PCR testing for mycobacterial identification. TW: Trunk wash; BALF: Bronchioalveolar lavage fluid; RD: Region of difference; MTBC: *Mycobacterium tuberculosis* complex; NTM: Non-tuberculous mycobacteria.

Pathogens **2022**, 11, 709 5 of 16

**Table 1.** Mycobacterial isolates identified by culture and virulence determination of n = 20 respiratory samples from 13 African elephants and 2 white rhinoceros, including the combinational use of four nucleic acid amplification tests directly on respiratory specimens for the prediction of culture outcome.

Mycobacterial spp. Present	Species (#Animals)	Sample Type (Animal ID)	Location (South Africa)	Mycobacterial Culture Result from Antemortem Respiratory Specimens		Nucleic Acid Amplification Test Results from Antemortem Respiratory Specimens			
				Combined MGIT and MGIT-TiKa Result	ESAT-6 and CFP-10 Virulence Determination	GeneXpert MTB/RIF Ultra Result	Hain CM <i>direct</i> V1.0 LPA Result	rpoB PCR Result	ku PCR Result
	African elephants n = 5	BALF (18/85)	Zoo	M. tuberculosis	Positive	MTB DETECTED Medium; RIF resistance indeterminate	MTBC	M. tuberculosis	M. tuberculosis
		BALF (18/177)	Zoo	M. tuberculosis	Positive	MTB TRACE DETECTED	Mixed NTM	M. tuberculosis	M. tuberculosis
		BALF (18/177)	Zoo	M. africanum	Positive	MTB TRACE DETECTED	M. avium and/or M. interjectum	M. africanum	M. africanum
		TW (18/177)	Zoo	M. elephantis strain	Positive	MTB NOT DETECTED	M. fortuitum group	M. elephantis	M. fortuitum group
Confirmed MTBC-positive		BALF (18/527)	KNP	M. stomatepiae	Negative	MTB NOT DETECTED	MTBC and/or M. fortuitum group	M. stomatepiae	mixed NTMs-M. smegmatis
		TW (18/527)	KNP	M. bovis	Positive	MTB TRACE DETECTED	Mixed NTM	M. bovis	M. bovis
		TW (18/533)	KNP	M. bovis	Positive	MTB TRACE DETECTED	MTBC and/or mixed NTM infection	M. bovis	M. bovis
		TW (18/538)	KNP	M. bovis	Positive	MTB TRACE DETECTED	M. fortuitum group	M. bovis	M. bovis
	White rhinoceros $n = 1$	BALF (19/46)	KNP	M. bovis	Positive	MTB TRACE DETECTED	MTBC and/or mixed NTM infection	M. bovis	M. bovis

Pathogens **2022**, 11, 709 6 of 16

Table 1. Cont.

Mycobacterial spp. Present	Species (#Animals)	Sample Type (Animal ID)	Location (South Africa)	Mycobacterial Culture Result from Antemortem Respiratory Specimens		Nucleic Acid Amplification Test Results from Antemortem Respiratory Specimens			
				Combined MGIT and MGIT-TiKa Result	ESAT-6 and CFP-10 Virulence Determination	GeneXpert MTB/RIF Ultra Result	Hain CM <i>direct</i> V1.0 LPA Result	rpoB PCR Result	ku PCR Result
	African elephants n = 8	TW (18/173)	KNP	M. avium complex strain	Positive	MTB TRACE DETECTED	MTBC and/or mixed NTM infection	M. avium complex	Mycobacteria spp.
		BALF (18/176)	Zoo	M. foliorum	Positive	MTB NOT DETECTED	M. fortuitum group	M. fortuitum group	M. fortuitum group
		BALF (19/460)	KNP	M. mantenii	Positive	MTB TRACE DETECTED	MTBC and/or mixed NTM infection	M. mantenii	M. mantenii
		TW (19/460)	KNP	M. abscessus strain	Positive	MTB NOT DETECTED	M. fortuitum group and/or M. abscessus	M. abscessus	M. abscessus
Confirmed		BALF (18/530)	KNP	M. interjectum strain	Positive	MTB NOT DETECTED	Mixed NTM	M. interjectum strain	mixed NTMs-M. avium complex
NTM-positive		TW (18/530)	KNP	M. avium strain	Negative	MTB NOT DETECTED	Mixed NTM	M. avium complex	M. elephantis strair
		TW (18/532)	KNP	M. mageritense strain	Negative	MTB NOT DETECTED	M. fortuitum group	Negative	Negative
		TW (18/534)	KNP	M. intracellulare	Positive	MTB TRACE DETECTED	MTBC and/or mixed NTM infection	M. avium complex	M. intracellulare
		TW (18/539)	KNP	M. intracellulare	Negative	MTB NOT DETECTED	M. fortuitum group	M. intracellulare	mixed NTMs-M. avium complex
		TW (21/496)	KNP	M. fortuitum strain	Positive	MTB NOT DETECTED	M. fortuitum group	M. fortuitum	M. elephantis strair
	White rhinoceros $n = 1$	BALF (18/31)	KNP	M. scrofulaceum strain	Positive	MTB NOT DETECTED	Mixed NTM	M. avium subsp. Paratuberculosis	mixed NTMs-M. avium subsp. Paratuberculosis

BALF: Bronchoalveolar lavage fluid; TW: Trunk wash; KNP: Kruger National Park; MGIT: Measurable growth in the *Mycobacterium* growth indicator tube (MGIT) detected by the BACTEC MGIT 960 mycobacterial detection system and strain typing (Warren et al., 2016); TiKa: modified MGIT system by alternative use of cationic D-enantiomer peptides for sample decontamination and the additional use of Supplement B during tube inoculation (Goosen et al., 2022); *ESAT-6/CFP-10*: PCR amplification of both virulence genetic targets from all culture; GeneXpert MTB/RIF Ultra: MTBC DNA detection through probe-based qPCR that simultaneously targets insertion elements IS6110 and IS1081 and, if positive, targets *rpoB* for drug resistance determination; Hain CM*direct* V1.0 line probe assay is a test system for the detection of *M. tuberculosis* complex and differentiation of more than 20 clinically relevant NTM directly from patient specimens; *rpoB* PCR uses *Mycobacterium* genus specific primers for the identification of genus DNA and subsequent speciation through amplicon sequencing; *Ku* PCR uses *Mycobacterium* genus specific primers for the identification of genus DNA and subsequent speciation through amplicon sequencing.

Pathogens 2022, 11, 709 7 of 16

Additionally, 25 respiratory samples (15 BALF and 10 TW) from 13 elephants and 21 BALF samples from 21 rhinoceros (20 white- and 1 black rhinoceros) were considered mycobacterial culture negative by the BACTEC 960 MGIT system and defined as the culture negative cohort in this study (Figure 1; File S1, Supplementary Materials).

# 2.2. Presence of ESAT-6/CFP-10 in All Mycobacterial Cultures

All of the MTBC culture isolates from seven respiratory samples (one rhinoceros and five elephants, including one duplicate BALF sample), were PCR-positive for *ESAT-6* and *CFP-10* (Table 1). Nine of the thirteen NTMs isolated by the respiratory sample culture were also PCR-positive for *ESAT-6* and *CFP-10* (one rhinoceros and ten African elephants) (Table 1). However, one of the PCR-positive NTM strains (*M. elephantis*) was cultured from a TW sample from an elephant (18/177), with *M. tuberculosis* and *M. africanum* isolated by culture in the BALF samples (Table 1).

Within the culture negative cohort, the presence of *ESAT-6* and *CFP-10* was detected in 7 out of the 46 respiratory sample cultures (13 elephants and 21 rhinoceros; File S1, Supplementary Materials). For these *ESAT-6/CFP-10* positive samples, the potential presence of the following *Mycobacterium* species was also identified by PCR amplicon sequencing directly from the respiratory samples: (1) *M. africanum*; (2) *M. bovis*; (3) *M. interjectum*; (4) *M. intracellulare*; (5) *M. avium* complex; and (6) *M. orygis* (File S1, Supplementary Materials). The target amplification was confirmed through Sanger sequencing and using the NCBI's Basic Local Alignment Search Tool for nucleotide (BLASTn) [30].

# 2.3. Nucleic Acid Amplification Test Results on Raw Respiratory Samples

#### 2.3.1. Ultra and Hain LPA for MTBC DNA Detection

The Ultra successfully identified all five of the infected elephants and the one infected rhinoceros as MTBC infected, which included four positive BALF and three elephant TW samples out of the nine respiratory samples (7/9) (Table 1). The Hain LPA also correctly identified three of the MTBC-infected elephants (18/85, 18/527, and 18/533), based on one positive TW sample and two positive BALF samples. The *M. bovis*-positive rhinoceros (19/46) was also identified as MTBC positive by the Hain LPA, using its BALF sample (Table 1). Three of the four positive MTBC results on the Hain LPA agreed with the Ultra results, with a discordant result (negative Ultra; MTBC and/or *M. fortuitum* group identified by the Hain LPA) for the BALF sample of the *M. bovis*-positive elephant (18/527). This BALF sample was *ESAT-6/CFP-10* negative and *M. stomatepiae* was isolated (Table 1). The TW sample from the *M. tuberculosis*-infected zoo elephant (18/177), based on the BALF culture, only contained culturable *M. elephantis*, which was negative on the Ultra, Hain LPA, and *ESAT-6/CFP-10* PCRs (Table 1).

When the Ultra and Hain LPA results were compared, using a two-tailed z-test, there was a significant difference (p < 0.00001) between these two tests for MTBC detection within the confirmed MTBC-infected cohort. However, agreement between the Ultra and Hain LPA for MTBC DNA detection from all of the specimens (regardless of culture outcome) was "substantial" ( $\kappa = 0.75$ , 95% CI 0.55–0.96: standard error (SE) = 0.10) [31]. Both the Ultra and Hain LPA identified the same three elephants (18/173, 19/460, and 18/534) as MTBC infected, based on one BALF and two TW samples, although the culture isolates categorized these individuals in the NTM positive cohort (Table 1). All three of these respiratory samples were also positive for ESAT-6/CFP-10 by PCR. Similarly, within the culture negative cohort, the Ultra and Hain LPA detected the same four elephants as MTBC infected, based on the results from the three BALF samples and one TW sample (all ESAT-6/CFP-10 positive) (File S1, Supplementary Materials).

A comparison between the Ultra and Hain LPA results with culture detection showed "substantial" agreement between the Ultra and culture for MTBC detection ( $\kappa=0.61$ , 95% CI 0.36–0.86, SE = 0.13), although lower than agreement between the Ultra and Hain LPA results ( $\kappa=0.75$ ). However, there was only "slight" agreement between the Hain LPA and culture results ( $\kappa=0.23$ , 95% CI-0.07–0.54: SE = 0.16). When the results of the Ultra and

Pathogens 2022, 11, 709 8 of 16

Hain LPA were combined and compared to culture for MTBC detection, there was only "moderate agreement" ( $\kappa = 0.58, 95\%$  CI 0.33–0.83: SE = 0.13).

# 2.3.2. Hain LPA NTM DNA Detection and Species Differentiation

The Hain LPA identified all eleven of the respiratory samples from the cultureconfirmed NTM positive cohort (eight elephants and one rhinoceros), as containing NTM DNA (11/11). Eight of the samples (8/11) were also positive for ESAT-6/CFP-10 amplification and three of those included samples from the elephants (18/173, 19/460, and 18/534) that were also identified as containing MTBC DNA by the Ultra and Hain LPA (Table 1). The Hain LPA detected mixed NTM DNA (>1 NTM strain) within the samples from these three KNP elephants, as well as an additional KNP elephant (18/530) and KNP rhinoceros (18/31) within the NTM positive cohort (Table 1). Similarly, within the MTBC-infected cohort, the Hain LPA also detected NTM DNA in eight out of nine respiratory samples, with five of the respiratory samples showing a mixture of NTM DNA (Table 1). Within the culture negative cohort, NTM DNA was detected by the Hain LPA in 19 of the 46 respiratory samples. These included four samples from the four elephants also identified as containing MTBC DNA by both the Ultra and Hain LPA. Three of these four elephant samples, and six additional elephant respiratory samples, were shown to have a mixture of NTM DNA based on the Hain LPA (File S1, Supplementary Materials). When the Hain LPA was compared to the culture results for the detection of NTM using all of the respiratory samples, only "fair" agreement was observed (Table 2). However, when the identified NTM species by the Hain LPA were compared with culture results, the Hain LPA only correctly identified one culture isolate (KNP elephant 19/460, TW, M. abscessus strain) (Table 1).

**Table 2.** The kappa ( $\kappa$ ), 95% confidence interval, standard error (SE) of agreement between four tests and culture for antemortem *Mycobacterium* genus DNA detection from respiratory samples from African elephants (n = 26) and rhinoceros (n = 23).

Test and Combinations	rpoB PCR	ku PCR	Hain LPA	Culture	rpoB/ku	<i>rpoB/</i> Hain LPA	<i>ku</i> /Hain LPA	rpoB/ku/ Hain LPA
rpoB PCR	1							
ku PCR	0.91 (0.81–1.00, 0.05)	1						
Hain LPA	0.94 (0.85–1.00, 0.04)	0.88 (0.76–0.99, 0.06)	1					
Culture	0.47 (0.29–0.66, 0.09)	0.50 (0.31–0.68, 0.09)	0.32 (0.16–0.48, 0.09)	1				
rpoB/ku PCRs	0.97 (0.91–1.00, 0.03)	0.94 (0.86–1.00, 0.04)	0.97 (0.91–1.00, 0.03)	0.45 (0.27–0.63, 0.09)	1			
rpoB/Hain LPA	0.94 (0.85–1.00, 0.04)	0.94 (0.85–1.00, 0.04)	1.00 (0.95–1.00, 0.03)	0.46 (0.29–0.64, 0.09)	0.94 (0.85–1.00, 0.04)	1		
ku/Hain LPA	0.94 (0.85–1.00, 0.04)	0.94 (0.85–1.00, 0.04)	1.00 (0.95–1.00, 0.03)	0.46 (0.29–0.64, 0.09)	0.94 (0.85–1.00, 0.04)	1.00 (0.95–1.00, 0.03)	1	
<i>rpoB/Ku/</i> Hain LPA	0.94 (0.85–1.00, 0.04)	0.94 (0.85–1.00, 0.04)	1.00 (0.95–1.00, 0.03)	0.46, (0.29–0.64, 0.09)	0.94 (0.85–1.00, 0.04)	1.00 (0.95–1.00, 0.03)	1.00 (0.95–1.00, 0.03)	1

*rpoB* PCR uses *Mycobacterium* genus specific primers for the identification of genus DNA and subsequent speciation through amplicon sequencing; *Ku* PCR uses *Mycobacterium* genus specific primers for the identification of genus DNA and subsequent speciation through amplicon sequencing; Hain CM*direct* V1.0 line probe assay is a test system for the detection of *M. tuberculosis* complex and differentiation of more than 20 clinically relevant NTM directly from patient specimens; Mycobacterial culture detection based on the combined isolation of *Mycobacteria* spp. through conventional BACTEC 960 MGIT system and MGIT-TiKa, followed by appropriate genetic speciation of all isolates (Warren et al., 2016). Landis Kappa scale: Kappa < 0: No agreement; Kappa between 0.00 and 0.20: Slight agreement; Kappa between 0.21 and 0.40: Fair agreement; Kappa between 0.41 and 0.60: Moderate agreement; Kappa between 0.61 and 0.80: Substantial agreement; Kappa between 0.81 and 1.00: Almost perfect agreement.

Pathogens 2022, 11, 709 9 of 16

2.3.3. Ku and rpoB Amplicon Sequencing for Mycobacterium Genus Detection and Speciation

Nine of the respiratory samples, confirmed to contain MTBC ( $M.\ tuberculosis$ ,  $M.\ africanum$ ,  $M.\ bovis$ ) by culture, showed an amplification in the ku and rpoB PCRs, which identified the presence of the Mycobacterium genus DNA (Table 1). Upon amplicon sequencing, both PCRs correctly identified the same MTBC species, directly from respiratory samples, as by culture (Table 1). Discordant culture/ku PCR results were observed for two of the elephant (18/177 and 18/527) TW and BALF samples. The elephant (18/177) TW was identified as containing  $M.\ elephantis$  (ESAT-6/CFP-10 positive) by culture and rpoB PCR, but as the " $M.\ fortuitum$  group" by the ku PCR (Table 1). Similarly, the elephant (18/527) BALF was identified as containing  $M.\ stomatepiae$  (ESAT-6/CFP-10 negative) by culture and rpoB PCR, but "mixed NTMs/ $M.\ smegmatis$ " by ku PCR (Table 1). When comparing the results from the ku and rpoB PCRs, there was no significant difference (two-tailed z-test, p=1) between these two tests for the Mycobacterium genus detection and the MTBC species identification within the culture-confirmed MTBC-infected cohort.

Both the ku and rpoB PCRs identified the presence of the Mycobacterium genus DNA in all but one of the eleven respiratory samples from the NTM-positive cohort. The TW sample from the KNP elephant (18/532) was culture positive for the M. mageritense strain (ESAT-6/CFP-10 negative) but identified as the M. fortuitum group by the Hain LPA (Table 1). The rpoB amplicon sequences predicted the same species of NTM as culture, but directly from the respiratory samples; the ku amplicon sequences identified five of the eleven NTM species assigned after culture (Table 1). Although there was no statistical difference (p = 1.0) between the two PCRs for Mycobacterium genus detection, there was a significant difference in the ability to identify NTMs (p < 0.02).

Within the culture negative cohort, the presence of the *Mycobacterium* genus was detected by both PCRs in the same 15 elephant respiratory samples (also all were Hain LPA positive). Two of the Hain LPA positive elephant (18/537 TW and 18/539 BALF) samples were *rpoB* PCR positive, but *ku* PCR negative; in addition, one elephant 18/157 TW sample that was also Hain LPA positive was *ku* PCR positive but *rpoB* PCR negative (File S1, Supplementary Materials). Three respiratory samples from different elephants (18/255, 18/534 and 18/536) were shown to have *M. africanum*, *M. bovis*, and *M. orygis* DNA, respectively, by both the *ku* and *rpoB* PCRs, although these samples were culture negative but all *ESAT-6/CFP-10* PCR positive. These three samples were also positive for MTBC DNA by both the Ultra and Hain LPA (File S1, Supplementary Materials).

To evaluate the best test or combination of tests, including the Hain LPA for *Mycobacterium* genus detection, agreement analysis was performed (Table 2). For the *Mycobacterium* genus detection, "fair" to "moderate" agreement were observed between the culture and individual tests (*rpoB-*, *ku* PCR, and Hain LPA) as well as for the combinations of these tests versus culture results (Table 2). "Almost perfect" agreement was reported between the individual tests (*rpoB*, *ku* PCRs and Hain LPA), as well as test combinations (Table 2).

For the mycobacterial species identification compared to culture, the Hain LPA only correctly identified a single NTM-positive sample (as previously mentioned), and it was incapable of differentiating between the MTBCs. Compared to the culture results, the rpoB and ku PCRs both correctly identified the same 6/7 MTBC culture positive samples. For all of the NTM positive samples, the rpoB PCR correctly identified 10/13 NTM culture positive samples, with 29/46 culture negative samples also negative on the rpoB PCR (Table 1; File S1, Supplementary Materials). The ku PCR correctly identified 4/13 NTM culture positive samples, and 30/46 of the culture negative samples were also negative on the ku PCR (Table 1; File S1, Supplementary Materials). The four NTMs identified by the ku PCR were also correctly identified by the rpoB PCR (Table 1).

# 3. Discussion

This pilot study described the successful culture isolation and genetic speciation of *M. tuberculosis*, *M. africanum* and *M. bovis* from antemortem respiratory samples collected from zoo and free-ranging African elephants and rhinoceros. Our recent culture isolation

Pathogens 2022, 11, 709 10 of 16

success was largely due to the combined use of conventional MGIT culture with a novel modified version called MGIT-TiKa [23,24], as well as improved isolate speciation through the simultaneous detection and sequencing of three different genetic markers (16S rRNA, rpoB, and hsp65). In the zoo elephant (18/85), M. tuberculosis had previously been isolated from the lung and lymph node tissue samples [5], and in this study, M. tuberculosis was also isolated from its antemortem BALF sample (Table 1). Notably, the elephant 18/85 BALF samples were collected before it was euthanized for tissue sample collection, removing the possibility of BALF contamination at necropsy. Moreover, in the same zoo, the simultaneous culture isolation of both M. tuberculosis and M. africanum from duplicate BALF samples was found in a single contact elephant 18/177. This was the only MTBC co-infection discovered among the elephants and rhinoceros in this study. Since this elephant was in a zoo environment in SA with known MTBC infections, the results were not that surprising [5]. The effects of MTBC co-infections on disease development are still largely unknown, especially in zoo settings, but certainly warrant further investigations. The successful culture isolation of M. bovis from three free-ranging African elephants and one rhinoceros in a M. bovis endemic wildlife reserve, such as KNP, were not unexpected, but it was surprising that it was cultured from paucibacillary antemortem respiratory samples. These findings suggest that possible MTBC shedding by elephants and rhinoceros may occur in zoos and wildlife reserves, however, it must also be noted that the isolation of M. bovis from elephant trunk wash samples could also have been due to contamination from other infected hosts, such as African buffaloes, shedding into the environment [32].

The isolation of various NTMs by culture from antemortem respiratory samples was also expected (Table 1), especially with the novel culturing approach used. It is noteworthy that some clinically important NTMs were isolated, along with supportive PCR evidence for the presence of ESAT-6 and CFP-10 virulence factors. These included M. abscessus, M. avium complex strains, M. interjectum, M. fortuitum, and M. scrofulaceum (Table 1). However, one limitation was that all of these animals were from MTBC endemic areas and the M. tuberculosis and M. bovis may have been present in small amounts below the detection threshold of 100 CFU for M. tuberculosis culture and 1000 CFU for M. bovis, but detectable by PCR [5,24]. This could result in positive ESAT-6/CFP-10 PCR, but culture negative results (Table 1). For many of these NTMs, opportunistic human and animal infections have been reported [2,33–35]. Two cases of atypical mycobacteriosis, caused by M. szulgai, have been reported in zoo elephants, as well as disease caused by M. kansasii infection in a bontebok herd, which were all positive for antibodies to ESAT-6/CFP-10 [16,36]. These cases highlight the diagnostic challenges around the detection of NTMs and differentiation from MTBC infections. The need to differentiate mycobacterial species is becoming increasingly recognized, despite the fact that many NTMs may not cause disease, but could prime the host's immune system, subsequently impeding the accurate diagnosis of MTBC infections, especially when using virulence factor proteins as test antigens [37].

An important finding, within the culture negative cohort, was the PCR detection and identification of MTBC DNA (*M. bovis*, *M. africanum*, *M. orygis*) in antemortem samples (File S1, Supplementary Materials). Without this result, any positive immunological assay results would have been classified as false-positive, rather than MTBC infected. These findings were supported by the *ESAT-6/CFP-10* positive PCR results, combined by the simultaneous species' identifications by three separate tests, the *rpoB* PCR, *ku* PCR, and the Hain LPA. However, a significant limitation is the inability of the PCRs to differentiate between live and dead bacteria, this may account for a positive PCR/culture negative result, which would be important for evaluating the transmission risk. Therefore, it is likely that shedding in MTBC-infected elephants and rhinoceros may occur more frequently than reported, due to the paucibacillary nature of antemortem respiratory samples and the overall sub-optimal sensitivity of culture [38,39].

The predictive ability of three candidate PCRs (Ultra, *rpoB* PCR, *ku* PCR) and one LPA, used directly on antemortem respiratory samples from culture-defined (MTBC, NTM, and negative) animal cohorts, was assessed to identify the most sensitive technique (Table 1).

Pathogens **2022**, 11, 709 11 of 16

For the MTBC DNA detection, the Ultra qPCR was significantly (p < 0.00001) more sensitive than the Hain LPA, with samples from the MTBC culture positive cohort. In addition, one advantage was that the Ultra used raw samples as input material whereas the Hain LPA required DNA extraction. Discordant results between the tests were reported for elephant 18/537 (M. bovis culture positive TW) where the Ultra detected MTBC DNA only in the TW sample, whereas the Hain LPA detected MTBC DNA in the BALF sample (from which M. stomatepiae was isolated by culture). These findings highlight the possibility that the MTBC may be present in samples from infected animals but may be unculturable due to their paucibacillary nature and sub-optimal sensitivity of culture [24]. A sample's bacterial load may also be reduced during sampling, sample handling, and storage. It is also possible that the TW sample reflected the overall respiratory load of bacilli and the possibility of not sampling the infected site using BAL. This highlights the importance of testing multiple samples.

Within the NTM culture positive and negative cohorts, MTBC DNA was also simultaneously detected by both the Ultra and Hain LPA, with positive amplification of the virulence factors for those samples. Again, this may have been dead bacteria or samples with small unculturable amounts of MTBC bacilli. Agreement between the Ultra, Hain LPA, and their combined use versus the culture results for all of the samples, showed that the Ultra was the most sensitive test for MTBC DNA detection. Combining the results from the Ultra and Hain LPA for MTBC DNA detection did not improve the overall detection, supporting the individual use of the Ultra for MTBC DNA detection directly from raw antemortem respiratory samples, which may identify infected animals that are culture negative.

The Hain LPA, for NTM DNA detection and species identification, identified all of the samples from the NTM positive cohort as containing NTM DNA (Table 1). Additionally, it also detected NTM DNA in the MTBC positive and culture negative cohort samples, which would be important for identifying the co-infections. The Hain LPA produced fair agreement compared to culture for the NTM detection, although there was no agreement with the culture for species' identification. This observation may be due to the use of hybridization technology by the Hain LPA producing unclear subjective results, especially when used directly on raw animal samples, or further evidence of the selective pressure introduced by the mycobacterial culture during sample processing [24]. Most of the samples were identified by Hain LPA as containing DNA from a mixture of NTM species, but the cultures were positive for only a single NTM species. This could indicate that the NTM that was successfully cultured most likely had outcompeted the rest during incubation, or that they were eradicated prior to inoculation during the decontamination process [24].

Through the PCR amplification and subsequent amplicon sequencing of the *rpoB* and ku gene targets, using Mycobacterium specific primers, both of the PCRs correctly detected and identified all of the MTBC species directly from the antemortem respiratory samples within the culture-confirmed MTBC positive cohort. No significant difference (p = 1) for genus detection or species identification was reported between these PCRs within the MTBC positive cohort. Similarly, both PCRs detected the Mycobacterium genus DNA from all of the samples in the NTM positive cohort, with no significant difference (p = 1), except for one TW sample from a KNP elephant. The negative result for both PCRs may have been due to damaged DNA during sample storage or handling, since the same DNA sample produced a positive M. fortuitum result in the Hain LPA. Otherwise, the amplicon sequencing of the rpoB target identified all of the cultured NTM species correctly with 100% accuracy, directly from respiratory samples, unlike the sequenced ku amplicons that only predicted 45% of the cultured species. A significant difference (p < 0.02) between the *rpoB* PCR and ku PCR was observed for the species identification. The success of the rpoB PCR for culture isolate prediction is based on a study performed by Adékambi et al. (2003), where the authors focused on a 723 bp variable region, exhibiting 83.9 to 97% interspecies similarity and 0 to 1.7% interspecies divergence, to design a primer pair for both PCR amplification and sequencing of this region for the identification of

Pathogens **2022**, 11, 709 12 of 16

rapidly growing mycobacteria [28]. Using these PCRs within the culture negative cohort, both the *rpoB* and *ku* PCRs identified the presence of *M. africanum*, *M. bovis*, and *M. orygis* DNA in three separate elephant samples. These samples were also all *ESAT-6/CFP-10*, Ultra, and Hain LPA positive, suggesting that these elephants were truly infected (File S1, Supplementary Materials). Agreement analysis for all of the samples between tests and test combinations versus culture for genus detection and species identification revealed that the *rpoB* PCR is the ideal individual test, compared to the rest. The Hain LPA detected slightly more samples containing *Mycobacterium* genus DNA, but it differentiated the species very poorly compared to culture. Therefore, combining the *rpoB* PCR with either the *ku* PCR and/or the Hain LPA slightly increased the genus detection, but not the species identification compared to culture. All of these findings have significant consequences for species management, with regards to screening animals prior to translocation, reducing risk of spillover into humans, improving veterinary disease control to avoid further transmission between species (including livestock), and protecting susceptible endangered wildlife species against infection [40–42].

The limitations of this study included the use of Sanger sequencing, as opposed to deep sequencing, resulting in shallow depth coverage, limited sample size, multiple freeze—thawing of samples prior to testing, the Hain LPA specialized equipment and training requirements, and possible underestimated sensitivities and overestimated specificities when using animal samples. Lastly, the cohort only included animals from *M. bovis* and *M. tuberculosis* endemic areas. Since the animals were naturally infected, the sporadic shedding of MTBC and the lack of corroborating evidence of NTM infection could have significantly affected the findings in this study. Therefore, additional studies are needed to assess the performance of culture-independent techniques for the rapid identification of MTBC and NTM infections in livestock and wildlife.

In conclusion, based on this pilot study, the Ultra appeared to be the optimal diagnostic test for MTBC DNA detection directly from raw antemortem respiratory specimens, and the *rpoB* PCR for *Mycobacterium* genus DNA detection and species identification through amplicon sequencing. Notably, the need for an NTM-specific test, such as the Hain LPA, could be circumvented by the combined use of the Ultra and *rpoB* PCR.

# 4. Materials and Methods

# 4.1. Animals

Throughout 2018, 2019, and 2021, the BALF samples (n = 23) were opportunistically collected from free-ranging white (n = 21) and black (n = 1) rhinoceros from M. bovisendemic KNP and zoo kept black (n = 1) rhinoceros. Trunk wash (n = 21) and BALF (n = 22) samples were also opportunistically collected from free-ranging M. bovis-exposed African elephants (n = 23) in KNP and zoo kept elephants (n = 3).

### 4.2. Antemortem Sample Collection and Mycobacterial Culture

The bronchoalveolar lavage fluid was collected endoscopically, and the TW samples as previously described [5,43,44]. Briefly, 150 mL respiratory samples were collected infield from each animal into a 500 mL sterile suction vacuum container, kept cool with ice bricks and transported to the laboratory within 4 h of collection. Thereafter, total volumes were aliquoted into three separate Corning<sup>TM</sup> Falcon<sup>TM</sup> 50 mL Conical Centrifuge Tubes (Waltham, MA, USA) per animal, concentrated by centrifugation at  $2000 \times g$  for 30 min, pooled into a 4 mL concentrate and frozen at -20 °C until further processing. For optimal *Mycobacteria* spp. isolation, all antemortem respiratory samples were processed for mycobacterial culture, using the conventional Mycobacteria Growth Indicator Tubes (MGIT, Becton Dickson, Franklin Lakes, NJ, USA), in parallel with a novel decontamination and mycobacterial culture technique (TiKa) for improved sensitivity, as previously reported [8,24]. Briefly, one aliquot of sample was decontaminated using BBL MycoPrep (Becton Dickinson, Franklin Lakes, NJ, USA) and inoculated into conventional MGIT tubes containing BBL MGIT PANTA-OADC enrichment media (Becton Dickinson, Franklin Lakes,

Pathogens 2022, 11, 709 13 of 16

NJ, USA), as described by the manufacturers. A second aliquot of the sample was decontaminated using TiKa-KiC (TiKa Diagnostics, London, UK) and inoculated into MGIT-TiKa tubes containing PANTA-OADC enrichment media and TiKa supplement B reagent (TiKa Diagnostics, London, UK). All of the MGITs (conventional and TiKa) were transferred to the BACTEC MGIT 960 mycobacterial detection system (Becton Dickinson). All of the liquid cultures with detected bacterial growth were sub-cultured onto blood agar plates to exclude contaminants and were subjected to Ziehl–Neelsen (ZN) acid fast staining [45]. All of the MGITs (conventional and TiKa) with positive bacterial growth underwent a further *Mycobacteria* spp. genetic speciation using the region of difference PCR, *16S* rRNA, *rpoB* and *hsp65* PCRs, and Sanger sequencing (Illumina Inc., San Diego, CA, USA), as previously described [28,29,46]. The animals with samples confirmed to contain either MTBC and/or NTM by culture and genetic speciation were defined as "MTBC positive" or "NTM positive". Furthermore, the presence or absence of immunogenic proteins in all of the bacterial growth-positive MGIT tubes were confirmed by PCR amplification and subsequent amplicon sequencing of *ESAT-6* and *CFP-10* genes, as previously described [47].

# 4.3. Nucleic Acid Amplification Tests (NAAT) for Mycobacteria spp. Detection and Differentiation

The GeneXpert MTB/RIF Ultra assay (Cepheid, Sunnyvale, CA, USA) was performed on the raw BALF and TW samples for the detection of MTBC DNA, as previously described [5]. Briefly, 700  $\mu$ L of samples were treated to chemical lysis, as prescribed by the manufacturer, and 2 mL of the solution aliquoted into the GeneXpert MTBC/RIF Ultra cartridge sample chamber. The result outputs were as follows: (1) "MTB not detected"; (2) "MTB trace detected"; and (3) "MTB detected high/medium/low/very low" [5].

Prior to performing the remaining NAATs, the total DNA was first extracted from all of the raw BALF and TW samples using the QIAGEN DNeasy® Blood and Tissue kit (Qiagen, Hilden, Germany), as previously described [48]. The DNA was used in the Genotype CMdirect version 1.0, line probe assay (Hain Lifescience, Germany) to determine the presence of the *Mycobacterium* genus, MTBC (without differentiation) and the subsequent detection and differentiation of more than 20 clinically relevant NTM species, according to manufacturer's instructions. Briefly, the multiplex PCRs were performed using biotinylated primers and proprietary amplification mixes in a MiniAmp $^{\text{TM}}$  Thermal Cycler (ThermoFisher Scientific, Waltham, MA, USA) with a ramp rate of  $\leq$ 2.2 °C/s. Thereafter, the amplicons were reverse hybridized in an automated GT-Blot 48 hybridization washer (Bruker, Billerica, MA, USA) onto precoated membrane strips and the hybridization results interpreted according to a template provided within the kit.

The extracted DNA was also subjected to conventional PCR using Mycobacteriumspecific primers for two highly conserved genetic regions, the ku (600 bp) and rpoB (740 bp) genes [26–29]. The amplicons were speciated by Sanger sequencing (Illumina, Inc.), as previously described [27–29]. Briefly, for both targets, a total reaction volume of 25 µL was prepared that consisted of 12.5 µL OneTaq Hot Start 2x master mix (New England Biolabs, Ipswich, MA, USA), 0.8 μM of each respective forward and reverse primer (Integrated DNA Technologies, Coralville, IA, USA), and 6.5 μL nuclease-free water and 5 μL DNA template. The positive controls which included 5  $\mu$ L 30 ng/ $\mu$ L extracted H37Rv M. tuberculosis and M. bovis DNA, as well as no template controls, were used in each PCR. Using a Veriti<sup>TM</sup> 96-well Thermal Cycler (Applied Biosystems, Waltham, MA, USA), the cycling conditions were 95 °C for 15 min, followed by 45 cycles at 95 °C for 1 min, 64 °C for 1 min, and 72 °C for 1 min, and a final elongation at 72 °C for 5 min. The amplicon presence, size, and intensity were confirmed by 1.5% agarose gel electrophoresis, followed by gel imaging using the ChemiDoc M.D. Universal Hood III Gel Documentation System (Bio-Rad, Hercules, CA, USA). All of the amplicons were sequenced by Sanger sequencing (Illumina, Inc., San Diego, CA, USA) through the Central Analytical Facility (CAF) at Stellenbosch University, SA. The sequence contigs were generated using Sequencher 5.1. software (Gene Codes, Ann Arbor, MI, USA) and blasted on NCBI's Basic Local Alignment Search Tool for nucleotide (BLASTn) [30].

Pathogens 2022, 11, 709 14 of 16

# 4.4. Data Analysis

All of the positive and negative NAAT and mycobacterial culture results are reported as proportions of the total number of animals and respiratory samples tested. For all of the known Mycobacteria spp. (MTBC and NTM) culture-positive specimens, each relevant NAAT and combination of NAATs with similar target species, the test-positive proportions were compared using a two-tailed z-test, where the z-scores and p-values were calculated (https://www.socscistatistics.com/tests/ztest/default2.aspx (accessed on 25 May 2022)). The p-values were considered statistically significant if p < 0.05. Thereafter, agreement analysis was performed between the relevant individual NAATs and test combinations for genus detection and differentiation, MTBC, and NTM detection and differentiation by calculating Cohen's kappa coefficient ( $\kappa$ ) using the online agreement calculator webtool (https://graphpad.com/quickcalcs/kappa1/ (accessed on 25 May 2022)). All of the NAATs and culture results for each sample are provided in File S1, Supplementary Materials.

**Supplementary Materials:** The following are available online at: https://www.mdpi.com/article/10.3390/pathogens11060709/s1, File S1: Supplementary material\_V1\_amended.xlsx.

**Author Contributions:** Conceptualization, W.J.G., M.A.M. and C.C.; methodology, W.J.G., C.C. and M.A.M.; formal analysis, W.J.G.; writing—original draft preparation, W.J.G.; writing—review and editing, W.J.G., T.J.K., P.B., L.K., C.C. and M.A.M.; supervision, M.A.M.; project administration, W.J.G. and T.J.K.; funding acquisition, W.J.G. and M.A.M. All authors have read and agreed to the published version of the manuscript.

**Funding:** This research was funded by the American Association of Zoo Veterinarians Wild Animal Health Fund [grant #S007355]. Sample collections were funded by the Elephant Sanctuary in Tennessee, U.S.A. [donation], the South African government through the South African Medical Research Council and the National Research Foundation South African Research Chair Initiative [grant #86949]. Personal funding for W.J.G. through Wellcome Foundation (grant #222941/Z/21/Z).

**Institutional Review Board Statement:** Ethical approval for this study was granted by the Stellenbosch University Animal Care and Use committee (ACU-2021–19019 and ACU-2022–21572). Permission to perform animal research in terms of section 20 of the Animal Diseases Act was granted by the South African Department of Agriculture, Land Reform and Rural Development (DALRRD), South Africa (12/11/1/7/2 and 12/11/1/7/6). All of the animals were handled by the Veterinary Wildlife Services veterinarians and game capture teams of the South African National Parks according to their standard operating procedures.

**Informed Consent Statement:** Not applicable.

Data Availability Statement: Data are available in the Supplementary Materials.

**Acknowledgments:** We thank Guy Hausler, Leana Freese, Tebogo Manamela, the Game Capture staff from Kruger National Park, and Drs. Jennie Hewlett and Robert Campbell from the National Zoological Gardens of South Africa for their assistance with the study.

Conflicts of Interest: The authors declare no conflict of interest.

# References

- 1. Johansen, M.D.; Herrmann, J.L.; Kremer, L. Non-tuberculous mycobacteria and the rise of *Mycobacterium abscessus*. *Nat. Rev. Microbiol.* **2020**, *18*, 392–407. [CrossRef] [PubMed]
- Gcebe, N.; Hlokwe, T.M. Non-tuberculous Mycobacteria in South African wildlife: Neglected pathogens and potential impediments for bovine tuberculosis diagnosis. Front. Cell. Infect. Microbiol. 2017, 7, 15. Available online: http://journal.frontiersin.org/article/10.3389/fcimb.2017.00015/full (accessed on 15 May 2022). [CrossRef] [PubMed]
- 3. Bernitz, N.; Kerr, T.J.; Goosen, W.J.; Chileshe, J.; Higgitt, R.L.; Roos, E.O.; Meiring, C.; Gumbo, R.; de Waal, C.; Clarke, C.; et al. Review of diagnostic tests for detection of *Mycobacterium bovis* infection in South African wildlife. *Front. Vet. Sci.* **2021**, *8*, 26. [CrossRef] [PubMed]
- 4. Goosen, W.J.; Kerr, T.; Kleynhans, L.; Buss, P.; Cooper, D.; Warren, R.M.; Van Helden, P.D.; Schröder, B.; Parsons, S.; Miller, M.A. The VetMAX<sup>TM</sup> *M. tuberculosis* complex PCR kit detects MTBC DNA in antemortem and postmortem samples from white rhinoceros (*Ceratotherium simum*), African elephants (*Loxodonta africana*) and African buffaloes (*Syncerus caffer*). *BMC Vet. Res.* **2020**, *16*, 220. [CrossRef]

Pathogens **2022**, 11, 709 15 of 16

 Goosen, W.J.; Kerr, T.; Kleynhans, L.; Warren, R.M.; Van Helden, P.D.; Persing, D.H.; Parsons, S.; Buss, P.; Miller, M.A. The Xpert MTB/RIF Ultra assay detects Mycobacterium tuberculosis complex DNA in white rhinoceros (Ceratotherium simum) and African elephants (Loxodonta africana). Sci. Rep. 2020, 10, 14482. [CrossRef]

- 6. Clarke, C.; Smith, K.; Goldswain, S.J.; Helm, C.; Cooper, D.V.; Kerr, T.J.; Kleynhans, L.; van Helden, P.D.; Warren, R.M.; Miller, M.A.; et al. Novel molecular transport medium used in combination with Xpert MTB/RIF ultra provides rapid detection of *Mycobacterium bovis* in African buffaloes. *Sci. Rep.* **2021**, *11*, 7061. [CrossRef]
- 7. Hlokwe, T.M.; van Helden, P.; Michel, A.L. Evidence of increasing intra and inter-species transmission of *Mycobacterium bovis* in South Africa: Are we losing the battle? *Prev. Vet. Med.* **2014**, *115*, 10–17. [CrossRef]
- 8. Goosen, W.; Miller, M.; Chegou, N.N.; Cooper, D.; Warren, R.; van Helden, P.D.; Parsons, S.D. Agreement between assays of cell-mediated immunity utilizing *Mycobacterium bovis*-specific antigens for the diagnosis of tuberculosis in African buffaloes (*Syncerus caffer*). Vet. Immunol. Immunopathol. 2014, 160, 133–138. [CrossRef]
- 9. Espie, I.W.; Hlokwe, T.M.; Van Pittius, N.C.G.; Lane, E.; Tordiffe, A.S.W.; Michel, A.L.; Müller, A.; Kotze, A.; Van Helden, P.D. Pulmonary infection due to *Mycobacterium bovis* in a black rhinoceros (*Diceros bicornis*) in South Africa. *J. Wildl. Dis.* **2009**, 45, 1187–1193. [CrossRef]
- 10. Miller, M.A.; Buss, P.; Roos, E.; Hausler, G.; Dippenaar, A.; Mitchell, E.; Van Schalkwyk, L.; Robbe-Austerman, S.; Waters, W.R.; Sikar-Gang, A.; et al. Fatal Tuberculosis in a free-ranging African elephant and One Health implications of human pathogens in wildlife. Front. Vet Sci. 2019, 6, 18. Available online: https://www.frontiersin.org/articles/10.3389/fvets.2019.00018/full (accessed on 8 February 2019). [CrossRef]
- 11. Miller, M.A.; Buss, P.; Parsons, S.D.C.; Roos, E.; Chileshe, J.; Goosen, W.J.; van Schalkwyk, L.; de Klerk-Lorist, L.-M.; Hofmeyr, M.; Hausler, G.; et al. Conservation of white rhinoceroses threatened by bovine tuberculosis, South Africa, 2016–2017. *Emerg. Infect. Dis.* 2018, 24, 2373–2375. [CrossRef] [PubMed]
- Okoi, C.; Anderson, S.T.B.; Antonio, M.; Mulwa, S.N.; Gehre, F.; Adetifa, I.M.O. Non-tuberculous Mycobacteria isolated from pulmonary samples in sub-Saharan Africa-systematic review and Meta Analyses. Sci. Rep. 2017, 7, 12002. [CrossRef] [PubMed]
- 13. Zulu, M.; Monde, N.; Nkhoma, P.; Malama, S.; Munyeme, M. Nontuberculous Mycobacteria in humans, animals, and water in Zambia: A systematic review. *Front. Trop Dis.* **2021**, 2. Available online: https://www.frontiersin.org/article/10.3389/fitd.2021.679501 (accessed on 17 June 2022). [CrossRef]
- 14. Falkinham, J.O. Ecology of Nontuberculous Mycobacteria. Microorganisms 2021, 9, 2262. [CrossRef] [PubMed]
- Kaczmarkowska, A.; Didkowska, A.; Kwiecień, E.; Stefańska, I.; Rzewuska, M.; Anusz, K. The Mycobacterium avium complex–an underestimated threat to humans and animals. Ann. Agric. Environ. Med. 2022, 29, 22–27. [CrossRef]
- 16. Miller, M.; Terrell, S.; Lyashchenko, K.; Greenwald, R.; Harris, B.; Thomsen, B.V.; Fontenot, D.; Stetter, M.; Neiffer, D.; Fleming, G. *Mycobacterium kansasii* infection in a bontebok (*Damaliscus pygaragus dorcas*) herd: Diagnostic challenges in differentiating from the *Mycobacterium tuberculosis* complex. *J. Zoo Wildl. Med. Off. Publ. Am. Assoc. Zoo Vet.* **2011**, 42, 468–472.
- 17. Bercovier, H.; Vincent, V. Mycobacterial infections in domestic and wild animals due to *Mycobacterium marinum*, M. fortuitum, M. chelonae, M. porcinum, M. farcinogenes, M. smegmatis, M. scrofulaceum, M. xenopi, M. kansasii, M. simiae and M. genavense. Rev. Sci. Tech. Int. Off. Epizoot. **2001**, 20, 265–290. [CrossRef]
- 18. Tingan, T.K.; Mensah, G.I.; Agyekum, E.B.; Amanor, I.B.; Addo, S.O.; Ayamdoo, Y.I.; Duah, M.S.; Mosi, L.; Addo, K.K. Non-tuberculous mycobacteria, not *Mycobacterium bovis*, are a significant cause of TB-like lesions observed in slaughtered cattle in Ghana. *IJID Reg.* **2022**, *3*, 8–14. [CrossRef]
- 19. Hannah, C.E.; Ford, B.A.; Chung, J.; Ince, D.; Wanat, K.A. Characteristics of Nontuberculous mycobacterial infections at a midwestern tertiary hospital: A retrospective study of 365 patients. *Open Forum. Infect. Dis.* **2020**, *7*, ofaa173. [CrossRef]
- 20. Gcebe, N.; Rutten, V.; Gey van Pittius, N.C.; Michel, A. Prevalence and distribution of non-tuberculous mycobacteria (NTM) in cattle, African buffaloes (*Syncerus caffer*) and their environments in South Africa. *Transbound. Emerg. Dis.* **2013**, *60* (Suppl. 1), 74–84. [CrossRef]
- 21. Vordermeier, H.M.; Brown, J.; Cockle, P.J.; Franken, W.P.J.; Drijfhout, J.W.; Arend, S.M.; Ottenhoff, T.H.M.; Jahans, K.; Hewinson, R.G. Assessment of cross-reactivity between *Mycobacterium bovis* and *M. kansasii ESAT-6* and *CFP-10* at the T-cell epitope level. *Clin. Vaccine Immunol. CVI* 2007, 14, 1203–1209. [CrossRef] [PubMed]
- 22. Schiller, I.; Oesch, B.; Vordermeier, H.M.; Palmer, M.V.; Harris, B.N.; Orloski, K.A.; Buddle, B.M.; Thacker, T.C.; Lyashchenko, K.P.; Waters, W.R. Bovine tuberculosis: A review of current and emerging diagnostic techniques in view of their relevance for disease control and eradication. *Transbound. Emerg. Dis.* 2010, *57*, 205–220. [CrossRef] [PubMed]
- Bull, T.J.; Munshi, T.; Mikkelsen, H.; Hartmann, S.B.; Sørensen, M.R.; Garcia, J.S.; Lopez-Perez, P.M.; Hofmann, S.; Hilpert, K.; Jungersen, G. Improved culture medium (TiKa) for *Mycobacterium avium* subspecies paratuberculosis (MAP) matches qPCR sensitivity and reveals significant proportions of non-viable MAP in lymphoid tissue of vaccinated MAP challenged animals. *Front Microbiol.* 2017, 7, 2112. Available online: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5209360/ (accessed on 26 March 2018). [CrossRef]
- Goosen, W.J.; Kleynhans, L.; Kerr, T.J.; van Helden, P.D.; Buss, P.; Warren, R.M.; Miller, M.A. Improved detection of *Mycobacterium tuberculosis* and *M. bovis* in African wildlife samples using cationic peptide decontamination and mycobacterial culture supplementation. *J. Vet. Diagn Investig. Off. Publ. Am. Assoc. Vet. Lab. Diagn. Inc.* 2022, 34, 61–67. [CrossRef] [PubMed]
- 25. World Health Organization. *Global Tuberculosis Report* 2016; World Health Organization: New York, NY, USA, 2016; 142p. Available online: https://apps.who.int/iris/handle/10665/250441 (accessed on 21 April 2022).

Pathogens 2022, 11, 709 16 of 16

26. Colman, R.E.; Schupp, J.M.; Hicks, N.D.; Smith, D.E.; Buchhagen, J.L.; Valafar, F.; Crudu, V.; Romancenco, E.; Noroc, E.; Jackson, L.; et al. Detection of low-level mixed-population drug resistance in *Mycobacterium tuberculosis* using high fidelity amplicon sequencing. *PLoS ONE* **2015**, *10*, e0126626. [CrossRef]

- 27. Zhou, L.; Ma, C.; Xiao, T.; Li, M.; Liu, H.; Zhao, X.; Wan, K.; Wang, R. A new single gene differential biomarker for *Mycobacterium tuberculosis* complex and non-tuberculosis Mycobacteria. *Front. Microbiol.* **2019**, *10*, 1887. [CrossRef]
- 28. Adékambi, T.; Colson, P.; Drancourt, M. *rpoB*-based identification of nonpigmented and late-pigmenting rapidly growing mycobacteria. *J. Clin. Microbiol.* **2003**, *41*, 5699–5708. [CrossRef]
- 29. de Zwaan, R.; van Ingen, J.; van Soolingen, D. Utility of *rpoB* gene sequencing for identification of nontuberculous Mycobacteria in the Netherlands. *J. Clin. Microbiol.* **2014**, *52*, 2544–2551. [CrossRef]
- 30. Altschul, S.F.; Gish, W.; Miller, W.; Myers, E.W.; Lipman, D.J. Basic local alignment search tool. J. Mol. Biol. 1990, 215, 403–410. [CrossRef]
- 31. Landis, J.R.; Koch, G.G. The measurement of observer agreement for categorical data. *Biometrics* 1977, 33, 159–174. [CrossRef]
- 32. Dwyer, R.A.; Witte, C.; Buss, P.; Goosen, W.J.; Miller, M. Epidemiology of tuberculosis in multi-host wildlife systems: Implications for black (*Diceros bicornis*) and white (*Ceratotherium simum*) rhinoceros. *Front. Vet. Sci.* 2020, 7, 580476. Available online: https://www.frontiersin.org/articles/10.3389/fvets.2020.580476/full?field=&journalName=Frontiers\_in\_Veterinary\_Science&id=580476 (accessed on 4 November 2020).
- 33. Gharbi, R.; Khanna, V.; Frigui, W.; Mhenni, B.; Brosch, R.; Mardassi, H. Phenotypic and genomic hallmarks of a novel, potentially pathogenic rapidly growing Mycobacterium species related to the *Mycobacterium fortuitum* complex. *Sci. Rep.* **2021**, *11*, 13011. [CrossRef] [PubMed]
- 34. Ferrell, K.C.; Johansen, M.D.; Triccas, J.A.; Counoupas, C. Virulence mechanisms of *Mycobacterium abscessus*: Current knowledge and implications for vaccine design. *Front Microbiol.* **2022**, *13*, 842017. Available online: https://www.frontiersin.org/article/10.3389/fmicb.2022.842017 (accessed on 16 May 2022). [CrossRef]
- 35. Ly, A.; Liu, J. Mycobacterial virulence factors: Surface-exposed lipids and secreted proteins. *Int. J. Mol. Sci.* **2020**, 21, 3985. [CrossRef] [PubMed]
- 36. Lacasse, C.; Terio, K.; Kinsel, M.J.; Farina, L.L.; Travis, M.A.; Greenwald, R.; Lyashchenko, K.P.; Miller, M.; Gamble, K.C. Two cases of atypical mycobacteriosis caused by *Mycobacterium szulgai* associated with mortality in captive African elephants (*Loxodonta africana*). *J. Zoo Wildl. Med. Off. Publ. Am. Assoc. Zoo Vet.* 2007, 38, 101–107. [CrossRef] [PubMed]
- 37. Gcebe, N.; Michel, A.; Gey van Pittius, N.C.; Rutten, V. Comparative genomics and proteomic analysis of four non-tuberculous Mycobacterium species and *Mycobacterium tuberculosis* complex: Occurrence of shared immunogenic proteins. *Front. Microbiol.* **2016**, 7, 795. Available online: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4894912/ (accessed on 27 November 2020). [CrossRef] [PubMed]
- 38. de la Rua-Domenech, R.; Goodchild, A.T.; Vordermeier, H.M.; Hewinson, R.G.; Christiansen, K.H.; Clifton-Hadley, R.S. Ante mortem diagnosis of tuberculosis in cattle: A review of the tuberculin tests, gamma-interferon assay and other ancillary diagnostic techniques. *Res. Vet. Sci.* **2006**, *81*, 190–210. [CrossRef]
- 39. Lorente-Leal, V.; Liandris, E.; Castellanos, E.; Bezos, J.; Domínguez, L.; De Juan, L.; Romero, B. Validation of a real-time PCR for the detection of *Mycobacterium tuberculosis* complex members in bovine tissue samples. *Front. Vet. Sci.* **2019**, *6*, 61. [CrossRef]
- 40. Clarke, C.; van Helden, P.; Miller, M.; Parsons, S. Animal-adapted members of the *Mycobacterium tuberculosis* complex endemic to the Southern African subregion. *J. S. Afr. Vet. Assoc.* **2016**, *87*, 1322. [CrossRef]
- 41. Love, D.M.; Garner, M.M.; Lyashchenko, K.P.; Sikar-Gang, A.; Bradway, D.S.; Robbe-Austerman, S.; Miller, M.; Ramer, J. Tuberculosis caused by *Mycobacterium orygis* in a greater one-horned rhinoceros (*rhinoceros unicornis*): First report in the western hemisphere. *J. Zoo Wildl. Med.* **2020**, *50*, 1000–1004. [CrossRef]
- 42. Thapa, J.; Paudel, S.; Sadaula, A.; Shah, Y.; Maharjan, B.; Kaufman, G.E.; McCauley, D.; Gairhe, K.P.; Tsubota, T.; Suzuki, Y.; et al. *Mycobacterium orygis*–associated tuberculosis in free-ranging rhinoceros, Nepal, 2015. *Emerg. Infect Dis.* **2016**, 22, 570–572. [CrossRef]
- 43. Michel, A.L.; Lane, E.P.; De Klerk-Lorist, L.-M.; Hofmeyr, M.; Van Der Heijden, E.M.D.L.; Botha, L.; Van Helden, P.; Miller, M.; Buss, P. Experimental *Mycobacterium bovis* infection in three white rhinoceroses (*Ceratotherium simum*): Susceptibility, clinical and anatomical pathology. *PLoS ONE* **2017**, 12, e0179943. [CrossRef] [PubMed]
- 44. Parsons, S.D.C.; Morar-Leather, D.; Buss, P.; Hofmeyr, J.; McFadyen, R.; Rutten, V.P.; Van Helden, P.D.; Miller, M.A.; Michel, A. The kinetics of the humoral and interferon-gamma immune responses to experimental *Mycobacterium bovis* Infection in the White Rhinoceros (*Ceratotherium simum*). *Front. Immunol.* **2017**, *8*, 1831. Available online: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5743672/ (accessed on 13 March 2019). [CrossRef] [PubMed]
- 45. Forbes, B.A.; Hall, G.S.; Miller, M.B.; Novak, S.M.; Rowlinson, M.-C.; Salfinger, M.; Somoskövi, A.; Warshauer, D.M.; Wilson, M.L. Practical guidance for clinical microbiology laboratories: Mycobacteria. *Clin. Microbiol. Rev.* **2018**, *31*, e00038-17. Available online: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5967691/ (accessed on 16 March 2021). [CrossRef] [PubMed]
- 46. Warren, R.M.; Van Pittius, N.C.G.; Barnard, M.; Hesseling, A.; Engelke, E.; De Kock, M.; Gutierrez, M.C.; Chege, G.K.; Victor, T.C.; Hoal, E.; et al. Differentiation of *Mycobacterium tuberculosis* complex by PCR amplification of genomic regions of difference. *Int. J. Tuberc. Lung Dis. Off. J. Int. Union Tuberc. Lung Dis.* 2006, 10, 818–822.
- 47. van Ingen, J.; de Zwaan, R.; Dekhuijzen, R.; Boeree, M.; van Soolingen, D. Region of Difference 1 in nontuberculous Mycobacterium species adds a phylogenetic and taxonomical character. *J. Bacteriol.* **2009**, *191*, 5865–5867. [CrossRef]
- 48. Stanton, J.J.; Nofs, S.A.; Peng, R.; Hayward, G.S.; Ling, P.D. Development and validation of quantitative real-time polymerase chain reaction assays to detect elephant endotheliotropic herpesviruses-2, 3, 4, 5, and 6. *J. Virol. Methods* **2012**, *186*, 73–77. [CrossRef]