



Remien

Rapid Antifungal Susceptibility Testing of Yeasts and Molds by MALDI-TOF MS: A Systematic Review and Meta-Analysis

Miriam Alisa Knoll 1,*, Hanno Ulmer 20 and Cornelia Lass-Flörl 10

- Institute of Hygiene and Medical Microbiology, Medical University of Innsbruck, 6020 Innsbruck, Austria; cornelia.lass-floerl@i-med.ac.at
- Department of Medical Statistics, Informatics and Health Economics, Medical University of Innsbruck, 6020 Innsbruck, Austria; hanno.ulmer@i-med.ac.at
- * Correspondence: miriam.knoll@i-med.ac.at; Tel.: +43-(0)512-9003-72573

Abstract: Due to the growing burden of fungal infections and a recent rise in antifungal resistance, antifungal susceptibility testing (AFST) is of increasing importance. The common methods of AFST have turnaround times of 24 to 48 h, and the available rapid methods are limited by applicability, costefficiency or accuracy. Given the urgency of adequate antifungal treatment in invasive mycoses, the need for the rapid and reliable detection of resistance is evident. In this systematic review and metaanalysis, we evaluated the diagnostic accuracy of AFST based on matrix-assisted laser desorption ionization time of flight mass spectrometry (MALDI-TOF MS). Twelve studies were reviewed, and data for the comparative analysis of their accuracy and methodology were systematically extracted. Compared to broth dilution as the gold standard, MALDI-TOF MS-based AFST reached a pooled sensitivity and specificity of 91% (95% Confidence Interval [CI], 84% to 96%) and 95% (95% CI, 90% to 98%), respectively. A comparative analysis showed that the sensitivity was higher for the semi-quantitative matrix-assisted laser desorption ionization Biotyper antibiotic susceptibility test rapid assay (MBT ASTRA) technique (96%) than for the correlate composite index (CCI) approach (85%), which is based on spectrum changes. Turnaround times below eight hours reached better diagnostic values than longer incubation periods, qualifying MALDI-TOF MS-based AFST as a rapid and accurate method for the detection of antifungal resistance.

Keywords: antifungal susceptibility testing; MALDI-TOF mass spectrometry; rapid AFST; detection of antifungal resistance

updates

check for

Citation: Knoll, M.A.; Ulmer, H.; Lass-Flörl, C. Rapid Antifungal Susceptibility Testing of Yeasts and Molds by MALDI-TOF MS: A Systematic Review and Meta-Analysis. J. Fungi 2021, 7, 63. https://doi.org/10.3390/jof7010063

Received: 12 November 2020 Accepted: 14 January 2021 Published: 18 January 2021

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2021 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https://creativecommons.org/licenses/by/4.0/).

1. Introduction

Fungal infections have been neglected in comparison to other infectious agents, despite their ubiquity and the rise of infections in the recent past. Worldwide, more than 1.5 million deaths per year are caused by fungal infections, which is comparable to the annual mortality rates of malaria and tuberculosis [1,2]. Advances in medicine have led to an increase of populations that are especially susceptible to fungal infections. In order to protect and treat those individuals, antifungal drugs are used in increasing amounts and for increasing durations, causing a selection of mutated fungal strains. Additionally, the extensive use of antifungals in agriculture enhances the development of antifungal resistance. Even with some new antifungal drugs in the pipeline, only five classes of antifungal drugs are available for patient treatment [3]. Due to intrinsic resistance of many fungal species and the limiting toxicities of some antifungals, the number of applicable drugs for fungal pathogens is even further narrowed. It is thus evident that any additional acquired resistance can have severe consequences for the affected individual, and the detection and surveillance of antifungal resistance are of increasing importance.

Antifungal susceptibility testing (AFST) is performed to investigate the response of a fungal pathogen to an antifungal drug. The gold standard for AFST is broth dilution, as specified by the current guidelines ((Clinical & Laboratory Standards Institute [CLSI]

J. Fungi **2021**, 7, 63 2 of 16

and European Committee on Antimicrobial Susceptibility Testing [EUCAST])). In order to perform those, 18 to 24 h cultures are required for EUCAST [4,5], while CLSI even requires double subculturing before testing a 24 h culture [6]. Before the results are available, broth dilution requires 24 to 48 h of incubation [4–6]. This long turnaround time represents a serious hindrance for the timely adjustment of an antifungal treatment and the initiation of disease control measures. The fact that broth dilution is still recommended as a reference method is due to the universal detection of resistance irrespective of the underlying mechanism. Fungal growth in the presence of an antifungal is the ultimate indicator for in vitro resistance, and visible growth quite simply takes time. Novel methods with shorter turnaround times have been developed, but all of those come with their own limitations of cost efficiency, applicability, or a limited number of detectable mechanisms [7–9].

Matrix-assisted laser desorption ionization time of flight mass spectrometry (MALDI-TOF MS) is a rapid, cost-effective and easy-to-handle method for microbial identification. This technology generates mass spectral profiles that are characteristic for each microorganism, which can be assigned to the respective species by comparison to a database. While it is well-established in microbial identification, attempts to utilize it in AFST have been relatively recent [10]. A variety of different approaches have been developed for the detection of resistance in bacteria [11], whereas AFST was investigated by merely two main approaches. The first method is based on the change of fungal spectrum profiles upon exposition to an antifungal substance. A Minimal Profile Change Concentration (MPCC) analogous to a minimum inhibitory concentration (MIC) value is inferred from a so-called composite correlation index (CCI) based on spectra from different concentrations of the antifungal substance. A second semi-quantitative technique named MBT ASTRA was originally developed for bacteria [12]. The growth of the infective agent is determined for different concentrations of antifungals and compared to a drug-free control, which provides a ratio of relative growth (RG). By defining a cut-off value, the RG at the corresponding antifungal concentration can be used to classify the quantity of the pathogen's growth as susceptible or resistant. The turnaround times for those two approaches vary from three to thirty hours, with a median of six hours for the reviewed studies.

As the time to adequate antifungal therapy is a major factor for patient outcomes in invasive mycoses [13], there is a demand for novel rapid technologies in AFST. In this systematic review and meta-analysis, we evaluate the existing data on MALDI-TOF MS applications in AFST, and assess the diagnostic accuracy of this method. We comparatively analyze the different technical approaches and evaluate the influence of test constellations as the tested species, antifungal substances, and incubation time on the performance of the assay.

To our knowledge, this is the first analysis evaluating the diagnostic accuracy of MALDI-TOF MS applications in AFST, and the first study both to statistically compare the different approaches and to investigate their suitability for rapid diagnostic appliance.

2. Materials and Methods

2.1. Study Design

We performed a systematic review and meta-analysis of MALDI-TOF MS-based assays for the detection of antifungal resistance in yeasts and molds regarding azoles and echinocandins. The Preferred Reporting Items for Systematic Reviews and Meta-Analyses (PRISMA 2009) guidelines were used as the methodological support. The statistical analysis was conducted using MedCalc® Statistical Software, version 19.5.3 (MedCalc Software Ltd., Ostend, Belgium; https://www.medcalc.org; 2020) and IBM SPSS Statistics for Windows, Version 25.0 (IBM Corporation, Armonk, NY, USA; 2019).

2.2. Search Strategy

Original articles published in English up to end of September 2020 were searched with PubMed and Google Scholar. 'MALDI-TOF MS-based susceptibility testing' and 'mass spectrometry rapid detection resistance' were combined with 'fungi', 'Candida' and

J. Fungi **2021**, 7, 63 3 of 16

'Aspergillus'. New links displayed beside the abstracts and articles mentioned in former reviews were followed. Finally, the bibliographies of each article were carefully reviewed, and relevant articles were also retrieved.

2.3. Inclusion and Exclusion Criteria

Only studies that evaluated antifungal susceptibility testing or the rapid detection of resistance towards antifungals through the use of MALDI-TOF MS technology were included. Additionally, the study reports must have had extractable data to fill the 4 cells of a 2×2 table for diagnostic tests (true resistant, false resistant, false susceptible and true susceptible). Additionally, broth dilution according to CLSI or EUCAST had to be used as reference standard or, if sequencing was originally mentioned as a gold standard, the data obtained by broth dilution had to be extractable.

2.4. Data Extraction

The full texts of the included articles were retrieved from the university library. The data from study reports was extracted twice. In order to ensure the reproducibility and completeness of the data extraction, an Excel spreadsheet (Microsoft Corp.) compiling all of the variables to be extracted was used. The data items included the authors, year of publication, tested antifungal drugs, technical approach, incubation time, sample size, reference method for broth dilution, sources of breakpoints or epidemiological cut-off (ECOFF) values, and values of true resistant, false resistant, false susceptible and true susceptible. When more than one antifungal or species was tested, these data were extracted separately. If possible, the raw data were used for the verification of the results. Otherwise, data for a 2×2 table for diagnostic tests and comparisons of agreement were extracted from the main text or the results section of the respective report. Strains classified as dose-dependently or intermediately susceptible by CLSI or EUCAST broth dilution were included in the calculations for essential and categorical agreement, but excluded from the preparation of the 2×2 tables. When more than one antifungal or species was tested, these data were extracted separately.

2.5. Breakpoints

The classification of the strains into susceptible (S), susceptible dose-dependent (SDD), intermediate (I) or resistant (R) was performed by applying the breakpoints (NCCLS, CLSI, EUCAST) or ECOFFs mentioned in the original report.

2.6. Data Analysis

Essential agreement was defined as \pm 2 dilution steps around the MIC value obtained by broth microdilution. Categorical agreement was defined as concordant categorization (susceptible, susceptible dose-dependent, intermediate or resistant) by both broth dilution and MALDI-TOF MS. All of the results, irrespective of categorization according to EUCAST or CLSI broth dilution (S, SDD, I or R), were included in the analysis for essential and categorical agreement.


For the comparative analysis, the studies were divided into groups based on their technical approach, antifungal substance, tested species, and incubation time. The sensitivity and specificity proportions of the selected studies were summarized using the meta-analysis for a proportion algorithm of the MedCalc® Statistical Software Version 19.5.3 (MedCalc Software Ltd, Ostend, Belgium). This algorithm uses a Freeman–Tukey transformation (arcsine square root transformation; Freeman & Tukey, 1950 [14]) to calculate the weighted summary proportion under the fixed and random effects model (DerSimonian & Laird, 1986 [15]). Heterogeneity was significant in 69% of the sub-analyses. Therefore, the results were decided using the random effects approach. The results of the meta-analysis were visualized with forest plots. Sensitivity was defined as the proportion of drug resistant fungal strains correctly identified as resistant by MALDI-TOF MS. Specificity was defined as the proportion of susceptible strains correctly identified as susceptible by MALDI-TOF

J. Fungi **2021**, 7, 63 4 of 16


MS. Strains classified as susceptible dose-dependent or intermediate by CLSI or EUCAST broth dilution were excluded from the sensitivity, specificity and forest plots. An exception to this was Saracli et al., 2015 [16], in which the SDD for *C. glabrata* was defined as S, and thereby made comparable.

3. Results

In total, 424 articles were screened, and 17 full-text articles were assessed for eligibility. Twelve reports fulfilled the inclusion criteria for the meta-analysis. The study characteristics of the twelve included studies are summarized in Table 1. The collective sensitivity and 100-specificity for each study compared to the gold standard broth dilution (CLSI or EUCAST) are shown in Figure 1.



Overall sensitivity of MALDI-TOF MS-based AFST.



Overall specificity of MALDI-TOF MS-based AFST.

Figure 1. Overall sensitivity and specificity of MALDI-TOF MS-based AFST compared to the gold standard of broth dilution (CLSI or EUCAST).

J. Fungi **2021**, 7, 63 5 of 16

Table 1. Characteristics of the included studies (n = 12).

	Antifungal	Species	Incubation Time	Essential Agreement	Categorical Agreement	Sensitivity	95% CI	Specificity	95% CI	True Resistant	False Resistant	True Susceptible	False Susceptible	Reference Method/Breakpoints
CCI Approach														
Azoles														
Marinach 2009 [17]	Fluconazole	C. albicans	15 h	100%	88%	80%	28% to 99%	100%	63% to 100%	4	0	8	1	CLSI M27-A2 [18]
Saracli 2015 [16]	Fluconazole	C. albicans	16 h	n.a.	80%	86% *	68% to 97%	88% *	47% to 100%	21	1	7	3	CLSI M27-A3 [19]
		C. glabrata	16 h	n.a.	86%	0%	0% to 98%	88% *	73% to 97%	0	4	30	1	CLSI M27-A3 [19]
		C. tropicalis	16 h	n.a.	65%	86%	42% to 100%	57% *	37% to 76%	6	12	16	1	CLSI M27-A3 [19]
	Voriconazole	C. albicans	16 h	n.a.	91%	93% *	68% to 100%	87% *	60% to 98%	14	2	13	1	CLSI M27-A3 [19]
		C. glabrata	16 h	n.a.	66%	100%	40% to 100%	61% *	42% to 78%	4	12	19	0	CLSI M27-A3 [19]
		C. tropicalis	16 h	n.a.	84%	80% *	28% to 99%	82% *	63% to 94%	4	5	23	1	CLSI M27-A3 [19]
	Posaconazole	C. albicans	16 h	n.a.	69%	55% *	32% to 76%	92% *	64% to 100%	12	1	12	10	CLSI M27-A3 [19]
		C. glabrata	16 h	n.a.	100%	n.a.	n.a.	100% *	90% to 100%	0	0	35	0	CLSI M27-A3 [19]
		C. tropicalis	16 h	n.a.	84%	75% *	19% to 99%	85% *	68% to 95%	3	5	28	1	CLSI M27-A3 [19]
Paul 2017 [20]	Fluconazole	C. tropicalis	4 h	100%	94%	100%	75% to 100%	100%	82% to 100%	13	0	19	0	CLSI M27-A3 [19]
Gitman 2017 [21]	Voriconazole	Aspergillus spp.	30 h	100%	100%	100%	59% to 100%	100%	75% to 100%	7	0	13	0	Pfaller 2009 [22]
Vella 2017 [23]	Fluconazole	C. glabrata	3 h	n.a.	79%	95%	83% to 99%	96%	81% to 100%	37	1	26	2	CLSI M27-S4 [24]
Delavy 2020 [25]	Fluconazole **	C. albicans	3 h	86%	n.a.	83%	52% to 98%	89%	52% to 100%	10	1	8	2	EUCAST v.9 [26]
Echinocandins														
De Carolis 2012 [27]	Caspofungin	C. albicans	15 h	100%	100%	100%	69% to 100%	100%	40% to 100%	10	0	4	0	Pfaller 2011 [28]
		C. glabrata	15 h	100%	100%	100%	63% to 100%	100%	40% to 100%	8	0	4	0	Pfaller 2011 [28]
		C. parapsilosis	15 h	100%	100%	n.a.	n.a.	100%	40% to 100%	0	0	4	0	Pfaller 2011 [28]
		C. krusei	15 h	100%	50%	100%	16% to 100%	0%	0% to 84%	2	2	0	0	Pfaller 2011 [28]
		A. fumigatus	15 h	100%	100%	100%	3% to 100%	100%	48% to 100%	1	0	5	0	Espinel-Ingroff 2011 [29

J. Fungi **2021**, 7, 63 6 of 16

 Table 1. Cont.

	Antifungal	Species	Incubation Time	Essential Agreement	Categorical Agreement	Sensitivity	95% CI	Specificity	95% CI	True Resistant	False Resistant	True Susceptible	False Susceptible	Reference Method/Breakpoints
		A. flavus	15 h	100%	100%	n.a.	n.a.	100%	40% to 100%	0	0	4	0	Espinel-Ingroff 2011 [29]
Vella 2013 [30]	Caspofungin	C. albicans	3 h	n.a.	98%	91%	59% to 100%	100%	93% to 100%	10	0	50	1	Pfaller 2011 [28]
Vella 2017 [23]	Anidulafungin	C. glabrata	3 h	n.a.	88%	67%	41% to 87%	100%	94% to 100%	12	0	58	6	CLSI M27-S4 [24]
Roberto 2020 [31]	Anidulafungin	C. parapsilosis complex	3 h	90%	95%	0%	0% to 98%	100%	82% to 100%	0	0	19	1	CLSI M27-A3 & -S4 [19,24]
	Caspofungin	C. parapsilosis complex	3 h	95%	95%	n.a.	n.a.	95%	75% to 100%	0	1	19	0	CLSI M27-A3 & -S4 [19,24]
	Micafungin	C. parapsilosis complex	3 h	95%	100%	n.a.	n.a.	100%	83% to 100%	0	0	20	0	CLSI M27-A3 & -S4 [19,24]
MBT ASTRA														
Echinocandins														
Vatanshenassan 2018 [32]	Caspofungin	C. albicans	6 h	91%	100%	100%	85% to 100%	100%	87% to 100%	22	0	25	0	CLSI M60 [33]
		C. glabrata	6 h	93%	61%	94%	80% to 99%	80%	28% to 99%	31	1	4	2	CLSI M60 [33]
Vatanshenassan 2019 [34]	Anidulafungin	C. glabrata	6 h	90%	89%	96%	80% to 100%	98%	91% to 100%	25	1	62	1	CLSI M60 [33]
Vatanshenassan 2019 [35]	Anidulafungin	C. auris	6 h	80%	98%	100%	54% to 100%	98%	88% to 100%	6	1	43	0	Tentative BPs [36–38]
	Micafungin	C. auris	6 h	86%	96%	100%	54% to 100%	96%	85% to 99%	6	2	42	0	Tentative BPs [36–38]
	Caspofungin	C. auris	6 h	68%	82%	100%	80% to 100%	73%	54% to 87%	17	9	24	0	Tentative BPs [36–38]


^{*} After the application of a 5% tolerance rate; ** \pm Cyclosporine A; Abbreviations: CI: Confidence Interval; n.a.: not applicable; BPs: breakpoints.

J. Fungi **2021**, 7, 63 7 of 16

Overall, the MALDI-TOF MS-based detection of resistance achieved 91.1% sensitivity (95% Confidence Interval [CI], 84% to 96%) and 95.1% specificity (95% CI, 90% to 98%) compared to the gold standard broth dilution (CLSI or EUCAST).


3.1. Sub-Analysis by Approach

The studies were divided into two groups based on their technical approach. Nine studies were based on the detection of altered spectrum profiles (CCI approach), and three studies were assigned to the MBT ASTRA technique. When more than one antifungal was tested, these study arms were analyzed separately. The overall sensitivity was higher for the MBT ASTRA approach, with 96% vs. 85.3% for the CCI approach, while the specificity was similar for both approaches (93.2% vs. 94.2%, respectively; Figure 2).



Overall·sensitivity·for·the·CCI·approach.

$Overall \cdot sensitivity \cdot for \cdot the \cdot MBT \cdot ASTRA \cdot approach.$



 $Overall \cdot sensitivity \cdot for \cdot the \cdot CCI \cdot approach.$

 $Overall \cdot sensitivity \cdot for \cdot the \cdot MBT \cdot ASTRA \cdot approach.$

Figure 2. Sensitivity and specificity by technical approach.

3.2. Sub-Analysis by Antifungal Substance

The studies were divided into two groups based on the tested antifungals. Six studies tested with azoles, and seven studies tested with echinocandins. When more than one antifungal was tested, these study arms were analyzed separately.


The largest number of isolates was tested for caspofungin, with 104 resistant and 157 susceptible *Candida* and *Aspergillus* strains of different species. The highest sensitivity was reached for micafungin; however, only one study—with merely six resistant isolates—was

suitable for analysis (Table 2). Micafungin also showed the best combined sensitivity and specificity. The highest specificity was achieved for anidulafungin. The lowest sensitivity was shown for posaconazole, and the lowest specificity was shown for voriconazole.


Antifungal	Number of Resistant Isolates	Sensitivity (%)	Number of Susceptible Isolates	Specificity (%)
Fluconazole	101	89.3	133	91.4
Posaconazole	26	57.7 *	81	92.6 *
Voriconazole	31	92.0	87	88.0
Caspofungin	104	96.4	157	92.1
Anidulafungin	51	79.6	184	98.4
Micafungin	6	100 *	64	96.2


Table 2. Sensitivity and specificity by antifungal agent.

Both the sensitivity and specificity were higher for echinocandins than for azoles, despite an outlier with a very low sensitivity for anidulafungin (Figure 3). The sensitivities for echinocandins and azoles were 92% and 86.4%, respectively. The specificities for echinocandins and azoles were 95.9% and 90.1%, respectively.



Overall sensitivity for azoles.





Overall sensitivity for azoles.

Overall sensitivity for echinocandins.

Figure 3. Sensitivity and specificity by the group of the antifungal substance.

A subclassification into the two technical approaches was applicable for echinocandins. Again, the MBT ASTRA approach achieved better sensitivity than the CCI approach (96% and 93.2%, respectively), while better specificity was obtained by CCI (98.1%, compared to

^{*} Only one study was included in the analysis.

J. Fungi 2021, 7, 63 9 of 16


> 79% for MBT ASTRA). Inversely, analyzing only studies based on the CCI approach, no significant difference in sensitivity between the two antifungal groups could be observed (84.4% for azoles vs. 83.3% for echinocandins). The specificity was still higher in the echinocandin group (97.3% vs. 92.3% for azoles).

3.3. Sub-Analysis by Species


The studies were divided into groups based on the tested species. When more than one antifungal was tested, these studies were analyzed separately. The largest number of samples per individual species was tested for C. glabrata, followed by C. albicans. Altogether, six studies analyzed C. albicans, five studies analyzed C. glabrata, five studies included other Candida species, and two studies analyzed different Aspergillus species.

A subclassification into the two approaches was applicable for C. glabrata. Again, the MBT ASTRA approach achieved a better sensitivity than the CCI approach (93.6% and 92.2%) and a lower specificity (83.9% and 92.6%, respectively).

The individually tested best sensitivity and specificity for yeasts was reached for C. auris, with 97.7% and 90%, respectively (Table 3). Altogether, 30 Aspergillus strains were included in the analysis, and both sensitivity and specificity were high (95.1% and 97.9%, respectively; Figures 4 and 5).



Overall sensitivity for C. albicans.



De Carolis 2012 - Caspofungin Saracli 2015 - Fluconazole

Vatanshenassan 2019 - Anidulafungin Vatanshenassan 2019 - Micafungin

Vatanshenassan 2019 - Caspofungin Roberto 2020 - Anidulafungin

Saracli 2015 - Voriconazole Saracli 2015 - Posaconazole

Paul 2017 - Fluconazole

Total (fixed effects) Total (random effects)

Overall sensitivity for other Candida species.

0.0 0.2 0.4 0.6


Proportion

Overall sensitivity for Aspergillus species.

Figure 4. Sensitivity by the group of the antifungal substance.

Overall sensitivity for C. glabrata.

J. Fungi 2021, 7, 63 10 of 16



Overall specificity for other Candida species.

3

species

0.0 0.2 0.40.6 0.8

Overall specificity for Aspergillus species.

0.0 0.2 0.4 0.6

4

100 *

Figure 5. Specificity by the group of the antifungal substance.

Species Number of Resistant Isolates Sensitivity (%) **Number of Susceptible Isolates** Specificity (%) 94.7 C. albicans 121 133 C. glabrata 129 88.8 257 92.5 29 108 C. tropicalis 86.7 82.7 0 * C. parapsilosis 1 63 97.2 C. krusei 100 * 0 * 2 2 97.7 * C. auris 29 90 * 121 Aspergillus 5 100 18 100 fumigatus Other Aspergillus

Table 3. Sensitivity and specificity by species.

100 *


3.4. Sub-Analysis by Time to Result

The studies were divided into groups based on the time to result. A cut-off of eight hours was defined to determine if the results were available within the same working day, and could thus be considered 'rapid results'.

The pooled sensitivity and specificity for the turnaround times below eight hours were 91.4% and 96.1%, respectively, while the times above eight ours were both noticeably lower (86.3% and 87.5%; Figure 6). Analyzing only studies based on the CCI approach, still better specificity was achieved below eight hours (98.2% vs. 87.5% above eight hours),


^{*} Only one study was included in the analysis.

while this difference became less marked regarding sensitivity (87.5% below vs. 86.3% above eight hours).



Overall·sensitivity·for·turnaround·times·below·8·h.

 $Overall \cdot sensitivity \cdot for \cdot turn around \cdot times \cdot above \cdot 8 \cdot h.$



Overall-specificity-for-turnaround-times-below-8-h.

Overall·specificity·for·turnaround·times·above·8·h.

Figure 6. Sensitivity and specificity by turnaround time.

4. Discussion

In invasive mycoses, early antifungal therapy is crucial for patient survival [13]. Inadequate antifungal therapy—and even short delays to the initiation of an adequate therapy—increase mortality rates in patients with candidemia [39–41]. Given the recent rise in the drug resistance of fungal pathogens, antifungal susceptibility testing (AFST) is becoming increasingly relevant [42–48]. As the common methods of AFST have turnaround times of 24 to 48 h, new rapid methods to detect resistances are necessary in order to ensure the fast and adequate adaption of antifungal therapy when needed.

The current methods for the rapid detection of resistance include molecular-based analyses [49,50], flow cytometry [51] or porous aluminum oxide (PAO)-based cultures [52,53], which deliver results within a few hours. Furthermore, the improvement of the conventional methods by the direct inoculation of standard tests with positive blood cultures [54] or shortening the time to result by automated approaches such as as the Vitek 2 yeast panel [55–57] have been investigated. However, all of these methods have their own

limitations of applicability, diagnostic accuracy or time to result. MALDI-TOF MS has been investigated for application in AFST, as the technology is well established for species identification in routine microbiology. Many diagnostic laboratories are already equipped with a MALDI-TOF MS system, which would facilitate the utilization upon the development of an appropriate assay for AFST. In this meta-analysis, we provide an overview of the recent studies utilizing MALDI-TOF MS for AFST. The overall sensitivity and specificity of MALDI-TOF MS-based AFST has been evaluated at 91% and 95%, respectively. Given the further-expandable amount of data, these results suggest that MALDI-TOF MS is a promising new tool for rapid AFST.

Two main approaches for the utilization of MALDI-TOF MS in AFST have been established. The first method was developed in 2009 [17], and is based on a change of fungal spectrum profiles upon exposure to an antifungal substance. Spectra from different concentrations of the antifungal substance are used to create a composite correlation index (CCI) and thereby infer a Minimal Profile Change Concentration (MPCC) analogous to an MIC value. Some studies also evaluated a simplified three-point assay, in which fungal spectra at the breakpoint concentration were compared to a maximum concentration as well as a drug-free control, and were then classified into resistant or susceptible depending upon their similarity to one of the extremes [16,21,23,25,30]. However, while the hands-on time is somewhat shortened, the time to result mainly depends on the incubation time, and is thereby not abbreviated.

Originally developed for the rapid antibiotic susceptibility testing of bacteria, a second semi-quantitative technique named MBT ASTRA was recently adapted for AFST [33]. Instead of focusing on changes in the spectra of fungi, this method is based on the determination of the growth of cells in the absence or presence of an antifungal substance. The relative growth compared to a drug-free control is calculated for each concentration, and a cut-off is defined in order to classify the relative growth values as resistant or susceptible.

Comparing these two approaches, MBT ASTRA showed, overall, a better sensitivity in our meta-analysis. While some studies also obtained high diagnostic values with the CCI approach, these results were rather heterogeneous. In contrast, consistently high sensitivity was achieved by MBT ASTRA. In regard to specificity, no significant differences could be observed between the two approaches. Considering the other results of this meta-analysis, one could deduce that this difference in sensitivity may be caused by the choice of more suitable antifungal drugs and species for analysis. While Vatanshenassan et al. tested with echinocandins only, many of the CCI studies were performed using azoles, for which generally-worse diagnostic results were obtained. In order to evaluate whether the choice of the tested antifungal had an influence on these results, a sub-analysis divided by the approach was performed. Only studies using echinocandins were compared, and still a markedly-higher pooled sensitivity was obtained by MBT ASTRA than by CCI. The same trend could be shown upon direct comparison on a species level, in which -among studies testing C. glabrata—those using the MBT ASTRA technique obtained higher sensitivities than the CCI approach. In conclusion, the better diagnostic results for the semi-quantitative MBT ASTRA seem to be partly due to a better selection of the tested bug-drug combinations, but upon close examination of the data, the measure of relative growth seems to obtain better results than the comparison of spectra, independently of the test conditions.

Upon division into groups based on the antifungal substance, better diagnostic values for echinocandins become apparent. However, only one study tested members of both antifungal classes, thereby enabling a direct comparison by the same setup of experiments [23]. In this study, no clear difference between the antifungals is observable. In order to exclude a bias by the implemented technical approach, only studies working with the CCI technique were compared. In doing so, no significant difference in sensitivity between the two antifungal groups could be observed, whilst the specificity was still better for echinocandins. Not enough data were available for the performance of this sub-analysis on

MBT ASTRA-based studies. Given that the CCI approach is based on spectrum changes, those differences might be due to a different influence on the proteome by the respective antifungal. However, since MALDI-TOF MS-based AFST is fairly new, this is—to our knowledge—the first report on different diagnostic values between antifungal classes. In conclusion, the overall analysis favors diagnostic results for echinocandins over azoles, but the reason might be partly due to a bias in the technical approach. More data with direct comparisons of two different classes of antifungal substances would be necessary. Nevertheless, the goal of developing new methods for rapid AFST should be to obtain valid and reliable results regardless of species or antifungal substance.

Similarly, upon the separation of the studies on a species level, differences in sensitivity and specificity between different yeast and mold species can be recognized. However, the only species with an overall case number above 100 tested isolates were C. albicans and glabrata. These species have been tested in several studies, and an overall sensitivity of 86.6% for *C. albicans* and 88.8% for *C. glabrata* was obtained. The specificity was above 90% for both species, which was lower than the overall values (91% sensitivity, 95% specificity). Good results were obtained for Aspergilli, which were included in two studies [21,27]. While AFST for molds commonly requires at least 24-48 h, De Carolis et al. were able to obtain a sensitivity and specificity of 100% within merely 15 h for different Aspergilli. In a second study, the incubation time was prolonged to 30 h for valid results [21], which still represents a reduction of the time to result compared to common methods. Only a small number of isolates has been tested, and not enough data is available for a conclusion, but nevertheless, it appears to be a promising tool for rapid AFST of Aspergillus sp. For other species, so far, not enough data is available to reliably compare the diagnostic accuracy. Hence, the very high or low diagnostic accuracies in this meta-analysis for some species have to be attributed to the small number of tested isolates, and are therefore not comparable. We have to assume that those values will approximate to the general values.

The eventual objective of the development of new methods for AFST is to obtain reliable results within the shortest possible time. However, shortening the time to result should not be traded against accuracy, but rather reach the same reliability as the reference method. Since this is a balancing act, one would suppose that shorter incubation times naturally come with lower diagnostic accuracy. Remarkably, in this meta-analysis, turnaround times below 8 h achieved better pooled sensitivity and specificity compared to studies with incubation times above 8 h. One confounding factor could be the distribution of the two technical approaches, since all of the studies based on MBT ASTRA with better sensitivity than CCI based reports were in the group with shorter incubation times. Analyzing only the CCI based studies, the difference becomes smaller, and only the specificity remains markedly better in the group of shorter times to result. This could be due to the fact that newer studies tended to apply shorter incubation times, and that expertise and technical know-how have increased over the years. Nonetheless, it is substantial to note that MALDI-TOF MS-based AFST can achieve satisfactory diagnostic results in under 8 h, which qualifies the method for future rapid applications in routine diagnostics.

In summary, the target of MALDI-TOF MS-based AFST is the rapid and reliable detection of resistance, in order to ensure a prompt modification of therapy if it is proven necessary. Vatanshenassan et al. [35] recently evaluated AFST via MALDI-TOF MS directly from positive blood cultures, which will presumably become one of the future adaptions to investigate. The reviewed articles cover a highly-divergent setup of experiments with different fungal species, antimycotics, microdilution methods and breakpoints. Despite this limitation, some conclusions can be drawn from the currently-available literature. Given the results of this analysis, reaching short incubation times without decreasing diagnostic accuracy is feasible. Better sensitivity has been achieved by the application of the semi-quantitative MBT ASTRA technique compared to the CCI approach based on spectrum changes. With a sensitivity of 91% and specificity of 95% overall, MALDI-TOF MS-based AFST has proven to be a promising new method for the rapid detection of resistance in

fungi. Further investigations using a more homogenous methodology with a larger number of isolates from a range of species and antifungal drugs will be necessary.

Author Contributions: Conceptualization, M.A.K. and C.L.-F.; methodology, M.A.K. and H.U.; software, H.U.; validation, M.A.K., H.U. and C.L.-F.; formal analysis, M.A.K. and H.U.; investigation, M.A.K.; resources, M.A.K.; data curation, M.A.K.; writing—original draft preparation, M.A.K.; writing—review and editing, C.L.-F. and H.U.; visualization, M.A.K.; supervision, C.L.-F.; project administration, C.L.-F. All authors have read and agreed to the published version of the manuscript.

Funding: This research received no external funding.

Conflicts of Interest: The authors declare no conflict of interest.

References

1. Brown, G.D.; Denning, D.W.; Gow, N.A.; Levitz, S.M.; Moorlag, S.J.; White, T.C. Hidden Killers: Human Fungal Infections. *Sci. Transl. Med.* **2012**, *4*, 165rv13. [CrossRef]

- 2. Bongomin, F.; Gago, S.; Oladele, R.; Denning, D.W. Global and Multi-National Prevalence of Fungal Diseases—Estimate Precision. *J. Fungi* 2017, 3, 57. [CrossRef]
- 3. Gintjee, T.J.; Donnelley, M.A.; Iii, G.R.T. Aspiring Antifungals: Review of Current Antifungal Pipeline Developments. *J. Fungi* **2020**, *6*, 28. [CrossRef]
- 4. EUCAST. EUCAST Antifungal MIC Method for Yeasts. Method for the Determination of Broth Dilution Minimum Inhibitory Concentrations of Antifungal Agents for Yeasts; (E.DEF 7.3.2); EUCAST Development Laboratory for Fungi: Copenhagen, Denmark, 2020; Available online: https://www.eucast.org/astoffungi/methodsinantifungalsusceptibilitytesting/susceptibility_testing_of_yeasts/(accessed on 27 August 2020).
- 5. EUCAST. EUCAST Antifungal MIC Method for Conidia Forming Moulds. Method for the Determination of broth Dilution Minimum Inhibitory Concentrations of Antifungal Agents for Conidia Forming Moulds; (E.DEF9.3.2); EUCAST Development Laboratory for fungi: Copenhagen, Denmark, 2020; Available online: https://www.eucast.org/fileadmin/src/media/PDFs/EUCAST_files/AFST/Files/EUCAST_E_Def_9.3.2_Mould_testing_definitive_revised_2020.pdf (accessed on 27 August 2020).
- 6. Clinical and Laboratory Standards Institute. *Performance Standards for Antifungal Susceptibility Testing of Yeasts*, 2nd ed.; CLSI Supplement M60; Clinical and Laboratory Standards Institute: Wayne, PA, USA, 2020; ISBN 978-1-68440-082-9.
- 7. Sanguinetti, M.; Posteraro, B. Susceptibility Testing of Fungi to Antifungal Drugs. J. Fungi 2018, 4, 110. [CrossRef] [PubMed]
- 8. Berkow, E.L.; Lockhart, S.R.; Ostrosky-Zeichner, L. Antifungal Susceptibility Testing: Current Approaches. *Clin. Microbiol. Rev.* **2020**, 33, e00069-19. [CrossRef] [PubMed]
- 9. Knabl, L.; Lass-Flörl, C. Antifungal susceptibility testing in Candida species: Current methods and promising new tools for shortening the turnaround time. *Expert Rev. Anti Infect. Ther.* **2020**, *18*, 779–787. [CrossRef] [PubMed]
- 10. Florio, W.; Tavanti, A.; Ghelardi, E.; Lupetti, A. MALDI-TOF MS Applications to the Detection of Antifungal Resistance: State of the Art and Future Perspectives. *Front. Microbiol.* **2018**, *9*, 2577. [CrossRef] [PubMed]
- 11. Oviaño, M.; Bou, G. Matrix-Assisted Laser Desorption Ionization–Time of Flight Mass Spectrometry for the Rapid Detection of Antimicrobial Resistance Mechanisms and Beyond. *Clin. Microbiol. Rev.* **2018**, 32. [CrossRef] [PubMed]
- 12. Lange, C.; Schubert, S.; Jung, J.; Kostrzewa, M.; Sparbier, K. Quantitative Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry for Rapid Resistance Detection. *J. Clin. Microbiol.* **2014**, 52, 4155–4162. [CrossRef]
- 13. Garey, K.W.; Rege, M.; Pai, M.P.; Mingo, D.E.; Suda, K.J.; Turpin, R.S.; Bearden, D.T. Time to Initiation of Fluconazole Therapy Impacts Mortality in Patients with Candidemia: A Multi-Institutional Study. Clin. Infect. Dis. 2006, 43, 25–31. [CrossRef]
- 14. Freeman, M.F.; Tukey, J.W. Transformations Related to the Angular and the Square Root. *Ann. Math. Stat.* **1950**, 21, 607–611. [CrossRef]
- 15. DerSimonian, R.; Laird, N. Meta-analysis in clinical trials. Control. Clin. Trials 1986, 7, 177–188. [CrossRef]
- 16. Saracli, M.A.; Fothergill, A.W.; Sutton, D.A.; Wiederhold, N.P. Detection of triazole resistance amongCandidaspecies by matrix-assisted laser desorption/ionization-time of flight mass spectrometry (MALDI-TOF MS). *Med. Mycol.* **2015**, *53*, 736–742. [CrossRef] [PubMed]
- 17. Marinach, C.; Alanio, A.; Palous, M.; Fekkar, A.; Brossas, J.-Y.; Brun, S.; Snounou, G.; Hennequin, C.; Sanglard, D.; Datry, A.; et al. MALDI-TOF MS-based drug susceptibility testing of pathogens: The example of Candida albicans and fluconazole. *Proteomics* **2009**, *9*, 4627–4631. [CrossRef] [PubMed]
- 18. NCCLS. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Approved Standard, 2nd ed.; NCCLS: Wayne, PA, USA, 2002; ISBN 1-56238-469-4.
- Clinical and Laboratory Standards Institute. Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts, 3rd ed.; Approved Standard; CLSI Document M27-A3; Clinical and Laboratory Standards Institute: Wayne, PA, USA, 2008; Volume 28, pp. 1–25.
- 20. Paul, S.; Singh, P.; AS, S.; Rudramurthy, S.M.; Chakrabarti, A.; Ghosh, A. Rapid detection of fluconazole resistance in Candida tropicalis by MALDI-TOF MS. *Med. Mycol.* **2017**, *56*, 234–241. [CrossRef] [PubMed]

 Gitman, M.R.; McTaggart, L.; Spinato, J.; Poopalarajah, R.; Lister, E.; Husain, S.; Kus, J.V. Antifungal Susceptibility Testing of Aspergillus spp. by Using a Composite Correlation Index (CCI)-Based Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry Method Appears to Not Offer Benefit over Traditional Broth Microdilution Testing. J. Clin. Microbiol. 2017, 55, 2030–2034. [CrossRef] [PubMed]

- Pfaller, M.; Diekema, D.J.; Ghannoum, M.A.; Rex, J.H.; Alexander, B.D.; Andes, D.; Brown, S.D.; Chaturvedi, V.; Espinel-Ingroff, A.; Fowler, C.L.; et al. Wild-Type MIC Distribution and Epidemiological Cutoff Values for Aspergillus fumigatus and Three Triazoles as Determined by the Clinical and Laboratory Standards Institute Broth Microdilution Methods. J. Clin. Microbiol. 2009, 47, 3142–3146. [CrossRef]
- 23. Vella, A.; De Carolis, E.; Mello, E.; Perlin, D.S.; Sanglard, D.; Sanguinetti, M.; Posteraro, B. Potential Use of MALDI-ToF Mass Spectrometry for Rapid Detection of Antifungal Resistance in the Human Pathogen Candida glabrata. *Sci. Rep.* **2017**, 7, 1–9. [CrossRef]
- 24. Clinical and Laboratory Standards Institute. *Reference Method for Broth Dilution Antifungal Susceptibility Testing of Yeasts. Fourth Informational Supplement*; CLSI Document M27-S4; Clinical and Laboratory Standards Institute: Wayne, PA, USA, 2012; ISBN 1-56238-864-9.
- Delavy, M.; Cerutti, L.; Croxatto, A.; Prod'Hom, G.; Sanglard, D.; Greub, G.; Coste, A.T. Machine Learning Approach for Candida albicans Fluconazole Resistance Detection Using Matrix-Assisted Laser Desorption/Ionization Time-of-Flight Mass Spectrometry. Front. Microbiol. 2020, 10, 3000. [CrossRef]
- 26. EUCAST. European Committee on Antimicrobial Susceptibility Testing Antifungal Agents Breakpoint Tables for Interpretation of MICs. Version 9.0, Valid from 2018-02-12. 2018. Available online: http://www.eucast.org/fileadmin/src/media/PDFs/EUCAST_files/AFST/Clinical_breakpoints/Antifungal_breakpoints_v_9.0_180212.pdf (accessed on 1 September 2020).
- De Carolis, E.; Vella, A.; Florio, A.R.; Posteraro, P.; Perlin, D.S.; Sanguinetti, M. Use of Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry for Caspofungin Susceptibility Testing of Candida and Aspergillus Species. *J. Clin. Microbiol.* 2012, 50, 2479–2483. [CrossRef]
- 28. Pfaller, M.; Diekema, D.; Andes, D.; Arendrup, M.; Brown, S.; Lockhart, S.; Motyl, M.; Perlin, D. Clinical breakpoints for the echinocandins and Candida revisited: Integration of molecular, clinical, and microbiological data to arrive at species-specific interpretive criteria. *Drug Resist. Updat.* **2011**, *14*, 164–176. [CrossRef] [PubMed]
- 29. Espinel-Ingroff, A.; Fothergill, A.; Fuller, J.; Johnson, E.; Pelaez, T.; Turnidge, J. Wild-Type MIC Distributions and Epidemiological Cutoff Values for Caspofungin and Aspergillus spp. for the CLSI Broth Microdilution Method (M38-A2 Document). *Antimicrob. Agents Chemother.* **2011**, *55*, 2855–2859. [CrossRef] [PubMed]
- 30. Vella, A.; De Carolis, E.; Vaccaro, L.; Posteraro, P.; Perlin, D.S.; Kostrzewa, M.; Posteraro, B.; Sanguinetti, M. Rapid Antifungal Susceptibility Testing by Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry Analysis. *J. Clin. Microbiol.* **2013**, *51*, 2964–2969. [CrossRef] [PubMed]
- 31. Roberto, A.E.M.; Xavier, D.E.; Vidal, E.E.; Vidal, C.F.D.L.; Neves, R.P.; De Lima-Neto, R.G. Rapid Detection of Echinocandins Resistance by MALDI-TOF MS in Candida parapsilosis Complex. *Microorganisms* **2020**, *8*, 109. [CrossRef] [PubMed]
- 32. Vatanshenassan, M.; Boekhout, T.; Lass-Flörl, C.; Lackner, M.; Schubert, S.; Kostrzewa, M.; Sparbier, K. Proof of Concept for MBT ASTRA, a Rapid Matrix-Assisted Laser Desorption Ionization–Time of Flight Mass Spectrometry (MALDI-TOF MS)-Based Method to Detect Caspofungin Resistance in Candida albicans and Candida glabrata. *J. Clin. Microbiol.* 2018, 56, e00420-18. [CrossRef] [PubMed]
- 33. Clinical and Laboratory Standards Institute. *Performance Standards for Antifungal Susceptibility Testing of Yeast*, 1st ed.; CLSI Supplement M60; Clinical and Laboratory Standards Institute: Wayne, PA, USA, 2017; ISBN 1-56238-829-0.
- 34. Vatanshenassan, M.; Arastehfar, A.; Boekhout, T.; Berman, J.; Lass-Flörl, C.; Sparbier, K.; Kostrzewa, M. Anidulafungin Susceptibility Testing of Candida glabrata Isolates from Blood Cultures by the MALDI Biotyper Antibiotic (Antifungal) Susceptibility Test Rapid Assay. *Antimicrob. Agents Chemother.* **2019**, *63*, e00554-19. [CrossRef]
- 35. Vatanshenassan, M.; Boekhout, T.; Meis, J.F.; Berman, J.; Chowdhary, A.; Ben-Ami, R.; Sparbier, K.; Kostrzewa, M. Candida auris Identification and Rapid Antifungal Susceptibility Testing Against Echinocandins by MALDI-TOF MS. Front. Cell. Infect. Microbiol. 2019, 9, 20. [CrossRef]
- 36. Ben-Ami, R.; Berman, J.; Novikov, A.; Bash, E.; Shachor-Meyouhas, Y.; Zakin, S.; Maor, Y.; Tarabia, J.; Schechner, V.; Adler, A.; et al. Multidrug-ResistantCandida haemuloniiandC. auris, Tel Aviv, Israel. *Emerg. Infect. Dis.* **2017**, 23, 195–203. [CrossRef]
- 37. Spivak, E.S.; Hanson, K.E. Candida auris: An Emerging Fungal Pathogen. J. Clin. Microbiol. 2017, 56, e01588-17. [CrossRef]
- 38. Tsay, S.; Kallen, A.; Jackson, B.R.; Chiller, T.M.; Vallabhaneni, S. Approach to the Investigation and Management of Patients with Candida auris, an Emerging Multidrug-Resistant Yeast. *Clin. Infect. Dis.* **2018**, *66*, 306–311. [CrossRef]
- 39. Morrell, M.; Fraser, V.J.; Kollef, M.H. Delaying the Empiric Treatment of Candida Bloodstream Infection until Positive Blood Culture Results Are Obtained: A Potential Risk Factor for Hospital Mortality. *Antimicrob. Agents Chemother.* **2005**, *49*, 3640–3645. [CrossRef] [PubMed]
- 40. Tumbarello, M.; Posteraro, B.; Trecarichi, E.M.; Fiori, B.; Rossi, M.; Porta, R.; Donati, K.D.G.; La Sorda, M.; Spanu, T.; Fadda, G.; et al. Biofilm Production by Candida Species and Inadequate Antifungal Therapy as Predictors of Mortality for Patients with Candidemia. *J. Clin. Microbiol.* 2007, 45, 1843–1850. [CrossRef] [PubMed]

41. Savage, R.D.; Fowler, R.A.; Rishu, A.H.; Bagshaw, S.M.; Cook, D.; Dodek, P.; Hall, R.; Kumar, A.; Lamontagne, F.; Lauzier, F.; et al. The Effect of Inadequate Initial Empiric Antimicrobial Treatment on Mortality in Critically Ill Patients with Bloodstream Infections: A Multi-Centre Retrospective Cohort Study. *PLoS ONE* **2016**, *11*, e0154944. [CrossRef]

- 42. Castanheira, M.; Messer, S.A.; Rhomberg, P.R.; Pfaller, M.A. Antifungal susceptibility patterns of a global collection of fungal isolates: Results of the SENTRY Antifungal Surveillance Program (2013). *Diagn. Microbiol. Infect. Dis.* **2016**, *85*, 200–204. [CrossRef] [PubMed]
- 43. Alexander, B.D.; Johnson, M.D.; Pfeiffer, C.D.; Jiménez-Ortigosa, C.; Catania, J.; Booker, R.; Castanheira, M.; Messer, S.A.; Perlin, D.S.; Pfaller, M.A. Increasing Echinocandin Resistance in Candida glabrata: Clinical Failure Correlates with Presence of FKS Mutations and Elevated Minimum Inhibitory Concentrations. *Clin. Infect. Dis.* 2013, 56, 1724–1732. [CrossRef]
- 44. Farmakiotis, D.; Tarrand, J.J.; Kontoyiannis, D.P. Drug-ResistantCandida glabrataInfection in Cancer Patients. *Emerg. Infect. Dis.* **2014**, 20, 1833–1840. [CrossRef]
- Guinea, J.; Zaragoza, Ó.; Escribano, P.; Martín-Mazuelos, E.; Pemán, J.; Sánchez-Reus, F.; Cuenca-Estrella, M. Molecular Identification and Antifungal Susceptibility of Yeast Isolates Causing Fungemia Collected in a Population-Based Study in Spain in 2010 and 2011. Antimicrob. Agents Chemother. 2013, 58, 1529–1537. [CrossRef]
- 46. Van Der Linden, J.; Arendrup, M.; Warris, A.; Lagrou, K.; Pelloux, H.; Hauser, P.; Chryssanthou, E.; Mellado, E.; Kidd, S.; Tortorano, A.; et al. Prospective Multicenter International Surveillance of Azole Resistance in Aspergillus fumigatus. *Emerg. Infect. Dis.* **2015**, 21, 1041–1044. [CrossRef]
- 47. Vallabhaneni, S.; Cleveland, A.A.; Farley, M.M.; Harrison, L.H.; Schaffner, W.; Beldavs, Z.G.; Derado, G.; Pham, C.D.; Lockhart, S.R.; Smith, R.M. Epidemiology and Risk Factors for Echinocandin Nonsusceptible Candida glabrata Bloodstream Infections: Data from a Large Multisite Population-Based Candidemia Surveillance Program, 2008–2014. *Open Forum Infect. Dis.* 2015, 2, ofv163. [CrossRef]
- 48. Alanio, A.; Denis, B.; Hamane, S.; Raffoux, E.; De Latour, R.P.; Menotti, J.; Amorim, S.; Touratier, S.; Bergeron, A.; Bretagne, S. Azole Resistance of Aspergillus fumigatusin Immunocompromised Patients with Invasive Aspergillosis. *Emerg. Infect. Dis.* **2016**, 22, 157–158. [CrossRef]
- 49. Dudiuk, C.; Gamarra, S.; Jiménez-Ortigosa, C.; Leonardelli, F.; Macedo, D.; Perlin, D.S.; Garcia-Effron, G. Quick Detection of FKS1 Mutations Responsible for Clinical Echinocandin Resistance in Candida albicans. *J. Clin. Microbiol.* **2015**, *53*, 2037–2041. [CrossRef] [PubMed]
- 50. Zhao, Y.; Nagasaki, Y.; Kordalewska, M.; Press, E.G.; Shields, R.K.; Nguyen, M.H.; Clancy, C.J.; Perlin, D.S. Rapid Detection of FKS-Associated Echinocandin Resistance in Candida glabrata. *Antimicrob. Agents Chemother.* **2016**, *60*, 6573–6577. [CrossRef]
- 51. Chaturvedi, V.; Ramani, R.; Pfaller, M.A. Collaborative Study of the NCCLS and Flow Cytometry Methods for Antifungal Susceptibility Testing of Candida albicans. *J. Clin. Microbiol.* **2004**, *42*, 2249–2251. [CrossRef] [PubMed]
- 52. Ingham, C.; Boonstra, S.; Levels, S.; De Lange, M.; Meis, J.F.; Schneeberger, P.M. Rapid Susceptibility Testing and Microcolony Analysis of Candida spp. Cultured and Imaged on Porous Aluminum Oxide. *PLoS ONE* **2012**, 7, e33818. [CrossRef] [PubMed]
- 53. Ingham, C.; Schneeberger, P.M. Microcolony Imaging of Aspergillus fumigatus Treated with Echinocandins Reveals Both Fungistatic and Fungicidal Activities. *PLoS ONE* **2012**, *7*, e35478. [CrossRef]
- 54. Oz, Y.; Gokbolat, E. Evaluation of direct antifungal susceptibility testing methods of Candida spp. from positive blood culture bottles. *J. Clin. Lab. Anal.* **2017**, 32, e22297. [CrossRef]
- 55. Cuenca-Estrella, M.; Gomez-Lopez, A.; Alastruey-Izquierdo, A.; Bernal-Martinez, L.; Cuesta, I.; Buitrago, M.J.; Rodriguez-Tudela, J.L. Comparison of the Vitek 2 Antifungal Susceptibility System with the Clinical and Laboratory Standards Institute (CLSI) and European Committee on Antimicrobial Susceptibility Testing (EUCAST) Broth Microdilution Reference Methods and with the Sensititre YeastOne and Etest Techniques for In Vitro Detection of Antifungal Resistance in Yeast Isolates. *J. Clin. Microbiol.* **2010**, 48, 1782–1786. [CrossRef]
- 56. Astvad, K.M.; Perlin, D.S.; Johansen, H.K.; Jensen, R.H.; Arendrup, M.C. Evaluation of Caspofungin Susceptibility Testing by the New Vitek 2 AST-YS06 Yeast Card Using a Unique Collection of FKSWild-Type and Hot Spot Mutant Isolates, Including the Five Most Common Candida Species. Antimicrob. Agents Chemother. 2012, 57, 177–182. [CrossRef]
- 57. Pfaller, M.A.; Diekema, D.J.; Procop, G.W.; Wiederhold, N.P. Multicenter Evaluation of the New Vitek 2 Yeast Susceptibility Test Using New CLSI Clinical Breakpoints for Fluconazole. *J. Clin. Microbiol.* **2014**, *52*, 2126–2130. [CrossRef]