


A.



B.




Figure S1. Immunoblot analysis of total RNA samples using anti-dsRNA antibody

Antibody against dsRNA was used to estimate the amount of dsRNA in total RNA isolated from Vero cell samples or mouse tissue samples. Immunoblotting was performed as described in [22, 25] (a) Blot of native acrylamide gel containing 2 μ g of total RNA from Vero cell samples or dilution series of known amount of in vitro transcribed 0.9kb dsRNA [22]. (b) 1 μ g of total RNA from each mouse total RNA sample was loaded in a slot blot and compared to a dilution series of in vitro transcribed 0.9kb dsRNA.