

Table S1. TP, TN, FP and FN cases of the Bcl-2 test (qualitative and quantitative), of the cIg LC test and of the two tests together in groups with inconclusive sIg LC.

Groups with inconclusive sIg LC	Test	TP	TN	FP	FN
Negative (N=159)	Bcl-2	50	93	1	15
	cIg LC	8	4	1	146
	Bcl-2 and/or cIg LC	51	93	1	14
Dual positive (N=47)	Bcl-2	37	4	0	6
	cIg LC	30	1	0	16
	Bcl-2 and/or cIg LC	41	4	0	2
Difficult interpretation (N=14)	Bcl-2	9	3	0	2
	cIg LC	4	0	0	10
	Bcl-2 and/or cIg LC	10	3	0	1

sIg LC surface immunoglobulin light chains, cIg LC cytoplasmic immunoglobulin light chains, TP true positives, TN true negatives, FP false positives, FN false negatives, N number of cases