

Supplementary Information

FRCaMP, Red Fluorescent Genetically Encoded Calcium Indicator Based on Calmodulin from *Schizosaccharomyces Pombe* Fungus

Oksana M. Subach, Natalia V. Barykina, Elizaveta S. Chefanova, Anna V. Vlaskina, Vladimir P. Sotskov, Olga I. Ivashkina, Konstantin V. Anokhin, and Fedor V. Subach

Supplementary Tables and Figures

Table S1	2
Table S2	3
Figure S1	4
Figure S2	5
Figure S3	6
Figure S4	7
Figure S5	8
Figure S6	9
Figure S7	10

Table S1. List of primers.

Primer	Primer sequence (5'-3')
RSP-BglII	CTGAGATCTATGAGGAAACCGTTCCGTG
CaMSP-HindIII-r	CCAAGCTTCTACTTCCAAGAAATGACACGAGAGAATTG
MSP-3-5	ggTAGATCTatgAGAAAAnnsnnnsnnsACTGCCTATAACGCTGTAC
MSP-6-8	ggTAGATCTatgAGAAAAACGTTCCGTnnsnnsnnsAACGCTGTACGTGCTTTC
MSP-9-11	ggTAGATCTatgAGAAAAACGTTCCGTACTGCCTATnnsnnsnnsCGTGCTTTC AACACTTGG
dRSP-BglII2	cgagatctGTGCCGCCGGTTCCGAG
BamHI-bJun	TGGggatccgcaccatggTGAAGGCGGAGAGGAAGCGCATG
bJun-NheI-r	CAGgctagcGGTGGCGATGGATCTTCTAG
NheI-RSPC	ACCgctagcCTGGTAAGCAAGGGCGAGG
BamHI-bFos	TGGggatccgcaccatggTGGGTCTGCGCAGTCCATCGG
bFos-NheI-r	CATgctagcGTGGTTCATGACTTCTG
NheI-RSPN	CACgctagcATGAGGAAACCGTTCCGTG
RSPN-HindIII-r	GATAagttCTACTTGTACAGCGCGTCCGTG

Table S2. In vitro $\Delta F/F$ response of truncated versions (with deleted M13-like peptide) of the purified FRCaM and GCaM6s indicators to the saturating calcium ion concentrations.

Indicator	$\Delta F/F$		
	0-39 μM ^a	0-820 μM ^a	0-2000 μM ^a
FRCaM	0.17 \pm 0.03	0.00 \pm 0.04	0.04 \pm 0.01
GCaM6s	0.16 \pm 0.03	0.18 \pm 0.02	0.05 \pm 0.03

^a 39, 820, and 2000 μM free calcium concentration corresponds to the 30 mM MOPS, 100 mM KCl, pH 7.20 buffer supplemented with either 10 mM CaEGTA, or 10 mM CaNTA or 2 mM CaCl₂, respectively. 0 free calcium concentration corresponds to the 30 mM MOPS, 100 mM KCl, pH 7.20 buffer supplemented with either 10 mM EGTA, or 10 mM NTA or 0 mM CaCl₂, respectively. Data were averaged across 8 repeats. SD is shown.

NES-FRCaMP-stop gene:

ATGCTTCAACTCCTCTTGAACGTCTTACTCTTCGAGATCT**ATGAGGAAACCGTTCCGTGGC**
GCGGGCAACGCTGTGCGTCTTCAGCACTGGAAAAAGCTAGTGCCGCGGTTCCGAGTGGAT
GTACCCCGAGGACGGGCCCTGAAGAGCGAGATCAGGAAGGGCTGAGGCTGAAGGACGGCG
CCACTACGCCGCCAGGTCAAGACCACCTACAAGGCCAAGAACGCCGTGCAGCTGCCGGCG
CTACATCGTCGACATCGAGTTGGACATCTTGTCCCACAACGAGGACTACACCATCGTGGAACAGT
GCGAACGCCCGTGGGCCACCCCCACCGGTGGCACGGACGGCTGTACAAGGGAGGTACAG
GCTCCGGCGGGAGTCTGTAAGCAAGGGCGAGGAGGATAACATGCCATCATCAAGGAGTTCAT
GCGCTCAAGGTGCACATGGAGGGCTCCGTGAACGGCCACGAGTCGAGATCGAGGGCGAGGG
CGAGGCCGCCCTACGAGGCCTCCAGACCGCTAGGCTGAAGGTGACCAAGGGTGGCCCCCTG
CCCTCGCCTGGACATCCTGTCCCCTCAGTTCATGTACGGCTCCAAGGCCTACATTAAGCACCC
AGCCGACATCCCCGACTACTTCAAGCTGTCCCTCCCCGAGGGCTTCAGGTGGAGCGCGTGTATGA
GCTTCGAGGACGGCGGCATTATTACGTAAATCAGGACTCCTCCCTGCAGGACGGCGTATTATC
TACAAGGTGAAGCTGCTCGGCACCAACTCCCCCGACGGCCCCGTAATGCAGAAAGAGACCA
TGGGCTGGGAGGGCTCCTACGAGATGACTACCGTAACCTACAGATGAGCAGATTGCGGAGTT
CGTGAGGCCTTCGCTGCTGATCGTACAGGATGAAATATCAGTCCAAATGAATTGGGTGT
GGTTATGAGGTGTTAGGTCTATGCCACTGCGCCGAATTACAAGATATGATTAATGAGGTG
ATGCCGATGGAATGGCACAATTGATTACGAATTGACAATGATGGCCCGAAAATGAAG
GATACCGACGACGAAGAGGAAGTCGCGAAGCCTTAAAGTCTCGATAAAGATGGAAGTGGAT
ACATTACAGTCGAGGAGCTGACTCATGTTCTACAAGTCTCGGTGAACGTTGTCTCGAGAAGAA
GTAGCCGATGTGATACGTGAAGCCGACTCGATGGCGATGGTGTAACTACGAAGAATTCT
CTCGTGTCAATTCTCCAAGTAA

Figure S1. Nucleotide sequence of NES-FRCaMP protein.




Figure S2. Calcium titration curves for R-GECO1.




Figure S3. Photochromism in the FRCaMP_{sat} indicator. 5 cycles of continuous metal halide lamp illumination were composed of photobleaching with yellow light (550/25BP) for 4 seconds followed by 30 seconds of darkness between cycles (marked as grey lines). Data were averaged across 6 drops in oil from 3 fields of view.




Figure S4. Calcium imaging of non-specific (spontaneous) activity of neuronal cultures co-expressing the FRCaMP and NCaMP7 calcium indicators. **(a, c)** Confocal images of two fields of view for neuronal cultures co-expressing the FRCaMP and NCaMP7 indicators. Scale bar, 50 μ m. **(b, c)** Examples of $\Delta F/F$ traces for the 8 cells are shown for each of two fields of view. Neuronal cultures co-expressing the NES-FRCaMP and NES-NCaMP7 indicators were imaged on DIV 15th. Neuronal cultures were transduced on DIV 4th with the mixture of rAAVs carrying NES-FRCaMP and NES-NCaMP7.




Figure S5. Puncta-like localization of the red FRCaMP indicator co-expressing with green NCaMP7 in cultured neurons. Confocal images of neurons co-expressing green NCaMP7 indicator with even distribution and red indicator FRCaMP with uneven puncta-like distribution. Neuronal cultures were imaged on DIV 22th. Scale bar, 20 μ m.




Figure S6. Localization and time-lapse response of the split-version of the FRCaMP indicator to Ca^{2+} variations in HeLa cells depending on the presence of the heterodimerizing (bJun-FRCaMPC/bFos-FRCaMPN) or non-heterodimerizing (bJun-FRCaMPC/bFOS Δ Zip-FRCaMPN) pair. **(a)** Zoomed area of the cell selected on Figure 7a. **(b)** The graph illustrates $\Delta F/F$ changes over time in red fluorescence of the split-version of the FRCaMP indicator in response to the addition of 2.5 μM of ionomycin depending on the presence of the heterodimerizing (bJun-FRCaMPC/bFos-FRCaMPN) or non-heterodimerizing (bJun-FRCaMPC/bFOS Δ Zip-FRCaMPN) pair. The changes on the graph correspond to the area indicated, as a numbered circle on the panel a, Figure 7 or panel a, Figure S5. Time of ionomycin addition is shown by arrow.




Figure S7. Brightness of the split-version of the FRCaMP indicator at low Ca^{2+} concentrations in HeLa cells depending on the presence of the heterodimerizing (bJun-FRCaMPC / bFos-FRCaMPN) or non-heterodimerizing (bJun-FRCaMPC / bFOS Δ Zip-FRCaMPN) pair. Comparison of the averaged brightness (for the heterodimerizing bJun-FRCaMPC and bFos-FRCaMPN split calcium indicator and its control non-heterodimerizing bJun-FRCaMPC and bFos-FRCaMPN split calcium indicator in HeLa cells at physiological calcium concentration before ionomycin addition. Error bars are standard deviations across twenty one cells (three cultures). p values show statistical difference between the respective values. ****, p - value is lower than 0.0001.