

Table S1. Primers used in ISSR-PCR analysis of *Dracocephalum forrestii* shoot treatments.

Primer Code	Primer Sequence 5' → 3'	Annealing Temperature	Number of bands for treatment (Growth regulator combination)							
			C	C-BI	TRIA 2.5	TRIA 5	TRIA 10	TRIA 2.5-BI	TRIA 5-BI	TRIA 10-BI
ISSR-X2	CTCCTCCTCCTCRC*	45°C	3	3	3	3	3	3	3	3
UBC808	AGAGAGAGAGAGAGAGC	48°C	6	6	6	6	6	6	6	6
UBC809	AGAGAGAGAGAGAGAGG	48°C	5	5	5	5	5	5	5	5
UBC840	GAGAGAGAGAGAGAGAYT*	48°C	7	7	7	7	7	7	7	7
UBC 818	CACACACACACACACAG	48°C	5	5	5	5	5	5	5	5
UBC 834	AGAGAGAGAGAGAGAGYT*	48°C	7	7	7	7	7	7	7	7

* Y = C or T; R = A or G.