

Table S1. The primer used for Selfie-Digital PCR

Gene	Amplicon size (bp)	Forward primer sequence 5'-3'	Reverse primer sequence 5'-3'
<i>GAPDH</i>	74	GATGGCCCCTCCGGGAAACT	GCGCCAGTAGAGGCAGGGAT
<i>PPIA</i>	70	ACCGCCGAGGAAAACCGTGT	CGTCGACGGCAATGTCGAAGAA
<i>RPL30</i>	92	GTGGGAAGTACGTCCTGGGGT	AGCTGGGCAGTTGTTAGCGA
<i>DYDC2</i>	61	GCTTTGGAAATTGCCTGGCCC	GGGTCACTGGGCCGAACCTT
<i>TCEA3</i>	96	AGGAAGAAGAAAGAGGAAGCGCTGA	ACCATCCCACTCAGACAGAGTTCA
<i>GRID1</i>	91	CTCTCGGCCCTGGAGATGGG	AGGTGCTGACCGAGAGCTGG
<i>HKDC1</i>	100	TGGGGGAGATTGTGCGGCAG	TCGAAGATGCCCCTGGTCCG
<i>SEL1L3</i>	87	GCTGTTCTGGGGGCAGCAAG	GGATCCTCCGTCTCCAGGGC
<i>DHRS2</i>	85	GGGCTGGTTTCATCCCTGTGC	TACCCGGTTAGCCAGGACGC
<i>FABP3</i>	74	CCTTGACGCCTGCTCTCTT	GTGCCCAGGAAAGCGTCCAC
<i>RYS3</i>	65	CGACCAGGTGGACCCCTTCC	AAGAAGAGGCCCCAGCCCACC
<i>GRM2</i>	82	GGCCATGGGATCGCTGCTTG	CAGCACCTTCTTGGCTGGGC
<i>DPYSL4</i>	74	TGGTCCTTCCTGGTGGCGTT	GTCGTCAGCCGGTGTCTATGC
<i>ABCG1</i>	80	ATGACGGAGCCCAAGTCGGT	AGCAGGTCCGTCTCAGTGGC