

**Supplemental Material 1.** Marginal gap evaluation using optical microscopy, of the width interfaces  $w$  of each surface and segment of each sample of the two considered groups: Group 1, conventional veneers (CO) and Group 2, crenelated veneers (CR).

Sample	Surface	Mean value of the interface width $w$ ( $\mu\text{m}$ ) of each segment	
		Group 1, classical veneers (CO)	Group 2, crenelated veneers (CR)
1	Incisal	mesial 160; middle 90; distal 170	mesial 50; middle 60; distal 40
	Cervical	mesial 180; middle 210; distal 190	mesial 60; middle 40; distal 40
	Mesial	incisal 160; middle 220; cervical 210	incisal 20; middle 40; cervical 90
	Distal	incisal 190; middle 220; cervical 230	incisal 80; middle 90; cervical 100
2	Incisal	mesial 140; middle 180; distal 140	mesial 100; middle 80; distal 60
	Cervical	mesial 40; middle 110; distal 90	mesial 80; middle 30; distal 80
	Mesial	incisal 120; middle 140; cervical 150	incisal 30; middle 30; cervical 40
	Distal	incisal 230; middle 190; cervical 340	incisal 80; middle 30; cervical 50
3	Incisal	mesial 140; middle 80; distal 80	mesial 180; middle 140; distal 140
	Cervical	mesial 120; middle 130; distal 170	mesial 40; middle 60; distal 40
	Mesial	incisal 180; middle 200; cervical 80	incisal 80; middle 80; cervical 80
	Distal	incisal 140; middle 90; cervical 240	incisal 100; middle 40; cervical 80
4	Incisal	mesial 160; middle 130; distal 140	mesial 80; middle 60; distal 60
	Cervical	mesial 100; middle 120; distal 80	mesial 60; middle 70; distal 80
	Mesial	incisal 150; middle 170; cervical 100	incisal 40; middle 30; cervical 50
	Distal	incisal 180; middle 110; cervical 280	incisal 50; middle 40; cervical 80
5	Incisal	mesial 170; middle 200; distal 230	mesial 80; middle 60; distal 50
	Cervical	mesial 320; middle 290; distal 330	mesial 70; middle 70; distal 90
	Mesial	incisal 180; middle 210; cervical 270	incisal 100; middle 90; cervical 90
	Distal	incisal 200; middle 230; cervical 300	incisal 80; middle 70; cervical 60
6	Incisal	mesial 160; middle 140; distal 190	mesial 80; middle 70; distal 90
	Cervical	mesial 210; middle 240; distal 260	mesial 100; middle 70; distal 90
	Mesial	incisal 150; middle 180; cervical 210	incisal 50; middle 30; cervical 40
	Distal	incisal 180; middle 200; cervical 220	incisal 40; middle 40; cervical 30
7	Incisal	mesial 250; middle 230; distal 290	mesial 70; middle 90; distal 90
	Cervical	mesial 310; middle 270; distal 240	mesial 70; middle 80; distal 80
	Mesial	incisal 190; middle 230; cervical 210	incisal 40; middle 30; cervical 40
	Distal	incisal 210; middle 250; cervical 230	incisal 30; middle 30; cervical 40
8	Incisal	mesial 190; middle 170; distal 220	mesial 50; middle 70; distal 70
	Cervical	mesial 290; middle 250; distal 280	mesial 60; middle 50; distal 40
	Mesial	incisal 330; middle 310; cervical 290	incisal 40; middle 30; cervical 30
	Distal	incisal 380; middle 310; cervical 350	incisal 50; middle 30; cervical 50
9	Incisal	mesial 260; middle 200; distal 240	mesial 100; middle 80; distal 90
	Cervical	mesial 290; middle 260; distal 230	mesial 110; middle 90; distal 80
	Mesial	incisal 280; middle 240; cervical 300	incisal 70; middle 60; cervical 60
	Distal	incisal 320; middle 290; cervical 310	incisal 70; middle 50; cervical 60

<b>10</b>	Incisal	mesial 420; middle 380; distal 430	mesial 60; middle 80; distal 70
	Cervical	mesial 440; middle 390; distal 410	mesial 80; middle 70; distal 90
	Mesial	incisal 390; middle 350; cervical 410	incisal 50; middle 30; cervical 30
	Distal	incisal 400; middle 370; cervical 410	incisal 40; middle 40; cervical 40
<b>11</b>	Incisal	mesial 330; middle 290; distal 290	mesial 90 middle 70; distal 60
	Cervical	mesial 350; middle 300; distal 310	mesial 80; middle 60; distal 40
	Mesial	incisal 370; middle 350; cervical 390	incisal 30; middle 40; cervical 50
	Distal	incisal 380; middle 340; cervical 400	incisal 40; middle 50; cervical 50
<b>12</b>	Incisal	mesial 210; middle 190; distal 220	mesial 100; middle 120; distal 110
	Cervical	mesial 220; middle 180; distal 240	mesial 90; middle 110; distal 100
	Mesial	incisal 280; middle 210; cervical 230	incisal 50; middle 70; cervical 50
	Distal	incisal 290; middle 250; cervical 310	incisal 60; middle 40; cervical 50