

Table S1 : sera titrated using the BioPro ELISA kit versus sera titrated using the FAVN test

		BioProRabies ELISA Ab kit – sera					
Assay	Result	Dogs- Experimental conditions (France)			Dogs- Field conditions (Tunisia)		
		Positive (≥70%)	Negative (<70%)	<i>Overall agreement</i>	Positive (≥70%)	Negative (<70%)	<i>Overall agreement</i>
FAVN test sera	Positive (≥0.5 IU/mL)	63	1		39	8	
	– Negative (<0.5 IU/mL)	16	53	87.2%	7	29	81.9%
	Total	79	54		46	37	