

Figure S1. Gating strategy for analysis of CD4+IL-17+ T cells and CD4+CD25+Foxp3+ T cells.

Cell were gated based on size using FSC-A and SSC-A to eliminate dead cells and debris. Single cells were gated using FSC-A and FSC-H. $CD4^+$ T cells were selected for further analysis. Subsequently, $CD4^+$ T cell was followed by examination of expression of IL-17 PerCP-CyTM 5.5 and double expression of CD25 PE and Foxp3 Alexa Flour 488.