

Supplementary Materials: An On-Site Simultaneous Semi-quantification of Aflatoxin B1, Zearalenone, and T-2 Toxin in Maize- and Cereal-based Feed via Multicolor Immunochromatographic Assay

Lin Xu, Zhaowei Zhang, Qi Zhang, Wen Zhang, Li Yu, Du Wang, Hui Li and Peiwu Li

Table S1. Mass spectrum conditions of LC-MS/MS analysis.

Compound	Mode	Precursor ion	Product ion(m/z)	Collision energy (eV)
AFB1	ESI ⁺	313[M+H] ⁺	285.0	35
			241.1	35
ZEN	ESI ⁻	317[M-H] ⁻	130.8	30
			175.9	25
T-2	ESI ⁺	489[M+H] ⁺	245.2	26
			387.3	20