

Table S2. Results of dilution scheme evaluation for ATX-fortified pooled human urine samples showing average (avg), standard deviation (Stdev), percent accuracy, and relative standard deviation (RSD). Fortified samples were prepared in triplicate.

Dilutions for ATX-fortified Human Urine Samples				
	2000 ng/mL ATX Spike 1:100 Dilution	2000 ng/mL ATX Spike 1:50 Dilution	2000 ng/mL ATX Spike 1:20 Dilution	2000 ng/mL Spike ATX 1:5 Dilution
Avg (ng/mL)	17.2	44.7	98.0	358
Calculated Conc (ng/mL)	1723	2235	1959	1790
Stdev (ng/mL)	1.21	2.96	6.26	8.19
RSD (%)	7.02	6.62	6.39	2.29
Accuracy (%)	86.2	112	98.0	89.5