

**Table S6.** The clones of Dako antibodies used for immunohistochemistry analyses.

<b>MoAb</b>	<b>CD3</b>	<b>CD5</b>	<b>CD10</b>	<b>CD20</b>	<b>BCL2</b>	<b>BCL6</b>	<b>MUM1</b>	<b>Ki-67</b>
Clone	F7.2.38	SP19	56C6	L26	124	PG-B6p	MUM1p	MIB-1
Dilution	1:50	F	F	F	F	F	F	F

MoAb, monoclonal antibody; Dako, (Glostrup, Denmark); F - FLEX, EnVision FLEX. High-sensitivity visualization system, intended for use in immunohistochemistry together with Dako Omnis autostainer was used.