

Table S1

Antibody	Company	Clone	Species	Reference number	Dilution	IHC incubation time (min)	Detection method
CD044	Ventana	SP37	Rb	06364985001	RTU	32	Optiview
MIB 1	Dako	MIB1	Mo	M7240	1/50	32	Optiview
OTP	Sigma-Aldrich	Polyclonal	Rb	HPA059342	1/50	64	Optiview
RB	BD Biosciences	G3-245	Mo	554136	1/1000*	32	Optiview
p16	Ventana	E6H4	Mo	06695248001	1/2	32	Optiview

Table S1. Characteristics of antibodies used for immunohistochemistry; IHC: immunohistochemistry; RTU: ready-to-use; \* in Dako Antibody Diluent.