

Supplementary Table S1. Primers used in PCR amplifications from FFPE tumor sections for Sanger sequencing.				
Primer	Gene	Coding Exon	Sequence	Amplicon size (bp)
KEAP1Ex1F	KEAP1	1	cctcatccagccctgtcttc	303
KEAP1Ex1R			ctcgatcacgtagaagacctt	
KEAP1Ex2-2F	KEAP1	2	caaggactacctggcgaagtc	302
KEAP1Ex2-2R			gggtgtaacagtcaggggc	
KEAP1Ex3F	KEAP1	3	tccacgaaggcagctataatgg	254
KEAP1Ex3R			gcttcggatgggttcattgc	
KEAP1Ex4F	KEAP1	4	gagtcaccttctctgcatgg	248
KEAP1Ex4R			gatgggctagtcaggactctt	
KEAP1Ex5F	KEAP1	5	gctgcatctctctttctgtc	239
KEAP1Ex5R			gtacagttctgctggcgaac	
NRF2Ex2F	NRF2	2	accatcaacagtggcataatgtg	397
NRF2Ex2R			ctgccataacttccaagaac	