

Table S1. The list of primers used for Real-time PCR.

Name	Forward sequence: (5' to 3')	Reverse sequence: (5' to 3')
<i>INTS2</i>	CCTCTTGCCCTGCCTGGTGA	GCTCTCTCCACTGCCTCCTCC
<i>PIGN</i>	TGATCCAAGCCTGTCCCTGGA	TGAGAAACGGCCAGCCTGCA
<i>MROH6</i>	GCCACTCGCGCTCTTGGAGA	AGGCCTCCACTGTGCAGCTG
<i>NEB</i>	TGGGATCCAGCACGCCAAGG	CCCTCCTTGGCGGCGTTGAT
<i>GMFB</i>	CCTGGTGCTGGATGAGGA	TGGTTAGCTCGGCAGTCTGGA C
<i>HEG1</i>	CACCGTCACAACCCAGCCA	TCACACGTGGCTCCCTTGGC
<i>INHBA</i>	CAAGGCGGCGCTTCTCAACG	GCCACACTTCTGCACGCTCC