

Table S2. Primers used in qPCRs to amplified specific targets.

Name	Forward primer	Reverse primer
uc.84	CTGCTTGATGTAACTTGATTAGCTC	AAATCTACATCTATCATTTTGCCAAT
uc.96	ACCCATCATGCATGTAAACG	AATGGAGTTTGCCGTGAGTT
uc.110	CCGTCCTCCCTGCACTAA	CAGCCAAAGGGGAAGTGTAT
uc.183	CTGCCCTATCAACTTTCGATGGTAG	CCGTTTCTCAGGCTCCCTCTC
pre-miR-221	GAGAACATGTTTCCAGGTAG	CAGGTCTGGGGCATGAAC
CCNB1	AGGCGAAGATCAACATGGCA	GGGCTTGGAGAGGCAGTATC
E2F1	CTTGGAGGGACCAGGGTTTC	GTTACCTTCATTCCCCGGT
TP53	TCACTTCCTTTTTCCGCAGT	TTCTGGCACTTCATCACAGC
CDKN1A	GCTTCATGCCAGCTACTTCC	CTGTGCTCACTTCAGGGTCA
CDKN1B	GTGCGAGTGTCTAACGGGAG	AGCCGGAGCCCCAATTAAAG
FBXW11	CAGCTAGAGTGGGTGCGAT	GGCCTCTTCTGGAGACGAT
B-Act	CATTGCCGACAGGATGCA	GCTGATCCACATCTGCTGGA