

Current progress on marine microplastics pollution: A review on pollution occurrence, detection and environmental effects

Fei-fei Liu ^{1,*}, Su-chun Wang ¹, Zhi-lin Zhu ¹, Guang-zhou Liu ¹

¹ Institute of Marine Science and Technology, Shandong University, Qingdao 266237, PR China; liufeifei@sdu.edu.cn (F.L.); suchunw@163.com (S.W.); zhuzhilinsdu@163.com (Z.Z.); liuguangzhou@sdu.edu.cn (G.L.)

* Correspondence: liufeifei@sdu.edu.cn




Figure S1. Sampling technologies and methods of processing.

Table S1. Advantages and disadvantages of the identification and quantification methods for microplastics.

Method	Advantages	Disadvantages
Microscopy	<ul style="list-style-type: none"> ▪ Simple, fast, and easy 	<ul style="list-style-type: none"> ▪ High possibility of false positive ▪ High possibility of missing small and transparent plastic particles
FTIR	<ul style="list-style-type: none"> ▪ Suitable for aliphatic compounds and polyesters ▪ Short measurement time 	<ul style="list-style-type: none"> ▪ Very small particles (<10 µm) are not detectable ▪ Spectra quality depends on sample purity, particle size and thickness
Raman	<ul style="list-style-type: none"> ▪ Suitable for microplastics made of aliphatic, aromatic or C=C compounds ▪ Detection of down to 1 µm microplastics in size 	<ul style="list-style-type: none"> ▪ Long measurement time ▪ Interference by sample purity and measurement parameters
Py-GC-MS	<ul style="list-style-type: none"> ▪ Simultaneously identify polymer types of microplastics particles and associated organic plastic additives 	<ul style="list-style-type: none"> ▪ Only a single particle can be analyzed in each run and not suitable for complex sample mixtures ▪ A few polymer identification
TGA-TDS-GC-MS	<ul style="list-style-type: none"> ▪ Allow the identification of high amount of environmental samples without any presorting or sample preparation 	<ul style="list-style-type: none"> ▪ Interference by TGA experiment conditions
alkaline hydrolysis LC-MS/MS	<ul style="list-style-type: none"> ▪ High recovery and analyzing efficiency ▪ Separation of microplastics from the samples is unnecessary ▪ Low loss of microplastics caused by picking 	<ul style="list-style-type: none"> ▪ Not suitable to quantify the different microplastics which contain the same building block compound