



**Figure S1. Arbovirus diagnostic algorithm.** Diagnostic algorithm followed by the Department of Arbovirology and Hemorrhagic Fevers of Instituto Evandro Chagas, currently the national reference for arbovirus diagnostics in Brazil.

**Table S1.** Contingency table with the RVNT assays diagnostics predictive values. All RVNT assays were evaluated in comparison to the reference assay PRNT.

Query	Tests	PPV	NPV	Sensitivity (%)	Specificity (%)	Accuracy (%)	McNemar value
PRNT90	DENV RVNT <sub>90</sub> Rluc	91.2	100	100	78.3	93.3	0.218 (p>0,05)
	RVNT <sub>90</sub> Rluc	32.0	100.0	100.0	74.6	77.3	0.000 (p<0,05)
	RVNT <sub>90</sub> mCherry	33.3	98.1	88,00	79.10	80.0	0.001 (p<0,05)
ZIKV Rluc	RVNT <sub>90</sub> mCherry	100.00	92.59	84.0	100.0	94.7	0.125 (p>0,05)