

Table S1. Micro-tensile bond strength test groups.

Group	Adhesive System	Application Technique	Teeth (<i>n</i>)	Slabs/24h	Slabs/6 months
1	CSE	A	5	30	30
2	CSE	APP	5	30	30
3	CSE	AAP	5	30	30
4	CSE	AAA	5	30	30
5	OBFL	A	5	30	30
6	OBFL	APP	5	30	30
7	OBFL	AAP	5	30	30
8	OBFL	AAA	5	30	30
9	OBU	A	5	30	30
10	OBU	APP	5	30	30
11	OBU	AAP	5	30	30
12	OBU	AAA	5	30	30
13	PBU	A	5	30	30
14	PBU	APP	5	30	30
15	PBU	AAP	5	30	30
16	PBU	AAA	5	30	30
Total:			80	480	480

Clearfil SE Bond (CSE); OptiBond FL (OBFL); OptiBond Universal (OBU); Prime&Bond Universal (PBU); Single active application (A); triple application, Active–Passive–Passive (APP); triple application, Active–Active–Passive (AAP); triple application, Active–Active–Active (AAA).