

Supplementary Materials

Article

Response Surface Methodology for the Optimisation of Electrochemical Biosensors for Heavy Metals Detection

Giuseppe Egidio De Benedetto ¹, Sabrina Di Masi ^{2,*}, Antonio Pennetta ¹ and Cosimino Malitesta ²

¹ Dipartimento di Beni Culturali, Università del Salento, Via D. Birago 64, 73100 Lecce, Italy; giuseppe.debenedetto@unisalento.it (G.E.D.B.); antonio.pennetta@unisalento.it (A.P.)

² Dipartimento di Scienze e Tecnologie Biologiche ed Ambientali, Via per Monteroni 1, 73100 Lecce, Italy; cosimino.malitesta@unisalento.it

* Correspondence: sabrina.dimasi@unisalento.it; Tel.: +39-0832-29-7145

Received: 11 January 2019; Accepted: 09 February 2019; Published: 13 February 2019

Contents

1. Table S1
2. Analysis of Variance of the model for Al³⁺ and Bi³⁺
3. Optimisation of responses
4. Selectivity of the biosensor

1. Table S1

Table S1. Variables and levels considered for the design of experiment (DOE).

Name (Factor)	Units	Low	High
Enzyme concentration (X1)	U/mL	50	800
Number of CV cycles (X2)	-	10	30
Flow rate (X3)	mL/min	0.3	1.0

2. Analysis of Variance of the Model for Al³⁺ and Bi³⁺

The analysis of variance with the calculated F and p values were obtained for all the tested metal ions.

Response Surface Regression: Sensitivity towards Al³⁺ ions versus [Enzyme] (U/mL); Flow rate (mL/min); Number of cycles

Analysis of Variance.

Source	DF	Adj SS	Adj MS	F-Value	P-Value
Model	9	0.128467	0.014274	1.79	0.188
Linear	3	0.071472	0.023824	3.00	0.082
[Enzyme] (U/mL)	1	0.049279	0.049279	6.20	0.032
Flow rate (mL/min)	1	0.000127	0.000127	0.02	0.902
Number of cycles	1	0.022066	0.022066	2.77	0.127
Square	3	0.004253	0.001418	0.18	0.909
[Enzyme] (U/mL)*[Enzyme] (U/mL)	1	0.003239	0.003239	0.41	0.538

Flow rate (mL/min)*Flow rate (mL/min)	1	0.000005	0.000005	0.00	0.981
Number of cycles*Number of cycles	1	0.000657	0.000657	0.08	0.780
2-Way Interaction	3	0.052743	0.017581	2.21	0.150
[Enzyme] (U/mL)*Flow rate (mL/min)	1	0.006555	0.006555	0.82	0.385
[Enzyme] (U/mL)*Number of cycles	1	0.045753	0.045753	5.75	0.037
Flow rate (mL/min)*Number of cycles	1	0.000435	0.000435	0.05	0.820
Error	10	0.079532	0.007953		
Lack-of-Fit	5	0.056098	0.011220	2.39	0.180
Pure Error	5	0.023433	0.004687		
Total	19	0.207999			


Model Summary.

S	R-sq	R-sq(adj)	R-sq(pred)
0.0891806	61.76%	27.35%	0.00%

Fits and Diagnostics for Unusual Observations.

Obs	S_Al3+	Fit	Resid	Std	Resid
3	0.0260	0.1535	-0.1275	-2.28	R
7	0.4170	0.3128	0.1042	2.03	R

R: Large residual.

Residual Plots for S_Al3+**Response Surface Regression: Sensitivity towards Bi3+ versus [Enzyme] (U/mL); Flow rate (mL/min); Number of cycles.****Analysis of Variance.**

Source	DF	Adj SS	Adj MS	F-Value	P-Value
Model	9	2.93908	0.32656	11.04	0.000
Linear	3	1.61191	0.53730	18.16	0.000
[Enzyme] (U/mL)	1	1.58134	1.58134	53.44	0.000
Flow rate (mL/min)	1	0.01054	0.01054	0.36	0.564
Number of cycles	1	0.02003	0.02003	0.68	0.430
Square	3	1.14477	0.38159	12.90	0.001
[Enzyme] (U/mL)*[Enzyme] (U/mL)	1	1.12542	1.12542	38.03	0.000
Flow rate (mL/min)*Flow rate (mL/min)	1	0.05734	0.05734	1.94	0.194
Number of cycles*Number of cycles	1	0.02296	0.02296	0.78	0.399
2-Way Interaction	3	0.18240	0.06080	2.05	0.170
[Enzyme] (U/mL)*Flow rate (mL/min)	1	0.00108	0.00108	0.04	0.852

[Enzyme] (U/mL)*Number of cycles	1	0.17731	0.17731	5.99	0.034
Flow rate (mL/min)*Number of cycles	1	0.00401	0.00401	0.14	0.721
Error	10	0.29591	0.02959		
Lack-of-Fit	5	0.22835	0.04567	3.38	0.104
Pure Error	5	0.06757	0.01351		
Total	19	3.23499			


Model Summary.

S	R-sq	R-sq(adj)	R-sq(pred)
0.172022	90.85%	82.62%	43.41%

Fits and Diagnostics for Unusual Observations.

Obs	S_Bi3+	Fit	Resid	Std	Resid
7	1.675	1.455	0.220	-2.23	R
16	1.642	1.890	-0.248	2.30	R

R: Large residual.

Residual Plots for S_Bi3+




Figure S1. Response Surface for Al^{3+} (A), and Bi^{3+} (B): the sensitivities were improved when low concentrations of enzyme and a high number of cycles were employed during the synthesis of biosensor.

3. Optimisation of the responses

Response Optimization: S_Bi³⁺; S_Al³⁺.

Response	Goal	Lower	Target	Upper	Weight	Importance
S_Bi ³⁺	Maximum	0.438	1.675		1	1
S_Al ³⁺	Maximum	0.001	0.417		1	1

Optimization Plot




Figure S2. Numerical optimisation performed by the software.

Solution.

Solution	[Enzyme](U/mL)	Flow rate (mL/min)	Number of cycles	S_Bi3+	S_Al3+	Composite Desirability
1	50	0.3	30	2.68352	0.574429	1

4. Selectivity of biosensor




Figure S3. Electrochemical responses of PPD/GOx biosensor after the injection of 0.8 μM of atrazine and 36.8 μM of 2,4-DB. Baseline: Glucose 1 mM.