



# Article Using a Quantitative High-Throughput Screening Platform to Identify Molecular Targets and Compounds as Repurposing Candidates for Endometriosis

Molly L. Churchill <sup>1</sup>, Sarah J. Holdsworth-Carson <sup>1,2</sup>, Karla J. Cowley <sup>3</sup>, Jennii Luu <sup>3</sup>, Kaylene J. Simpson <sup>3,4</sup>, Martin Healey <sup>1,5</sup>, Peter A. W. Rogers <sup>1</sup> and J. F. Donoghue <sup>1,\*</sup>

- <sup>1</sup> Gynaecology Research Centre, Department of Obstetrics and Gynaecology, University of Melbourne and The Royal Women's Hospital, Parkville, VIC 3052, Australia; molly.churchill@outlook.com (M.L.C.); sarah.carson@epworth.org.au (S.J.H.-C.); parogers@unimelb.edu.au (P.A.W.R.)
- <sup>2</sup> Julia Argyrou Endometriosis Centre, Epworth HealthCare, Richmond, VIC 3121, Australia
- <sup>3</sup> Victorian Centre for Functional Genomics, Peter MacCallum Cancer Centre, Parkville, VIC 3010, Australia; kaylene.simpson@petermac.org (K.J.S.)
- <sup>4</sup> Sir Peter MacCallum Department of Oncology, University of Melbourne, Parkville, VIC 3010, Australia
- <sup>5</sup> Gynaecology Endometriosis and Pelvic Pain Unit, Royal Women's Hospital, Parkville, VIC 3052, Australia
- Correspondence: jacqueline.donoghue@unimelb.edu.au

Abstract: Endometriosis, defined as the growth of hormonally responsive endometrial-like tissue outside of the uterine cavity, is an estrogen-dependent, chronic, pro-inflammatory disease that affects up to 11.4% of women of reproductive age and gender-diverse people with a uterus. At present, there is no long-term cure, and the identification of new therapies that provide a high level of efficacy and favourable long-term safety profiles with rapid clinical access are a priority. In this study, quantitative high-throughput compound screens of 3517 clinically approved compounds were performed on patient-derived immortalized human endometrial stromal cell lines. Following assay optimization and compound criteria selection, a high-throughput screening protocol was developed to enable the identification of compounds that interfered with estrogen-stimulated cell growth. From these screens, 23 novel compounds were identified, in addition to their molecular targets and in silico cellsignalling pathways, which included the neuroactive ligand-receptor interaction pathway, metabolic pathways, and cancer-associated pathways. This study demonstrates for the first time the feasibility of performing large compound screens for the identification of new translatable therapeutics and the improved characterization of endometriosis molecular pathophysiology. Further investigation of the molecular targets identified herein will help uncover new mechanisms involved in the establishment, symptomology, and progression of endometriosis.

**Keywords:** endometriosis; endometrial stromal cells; estrogen-signalling pathways; high-throughput screening; high-content imaging

# 1. Introduction

Endometriosis is an estrogen-dependent gynaecological disease that affects up to 11.4% women of reproductive age [1,2] as well as trans and gender-diverse people (for whom statistics are currently unavailable). In extremely rare conditions such as Rokitansky–Küster–Hauser (MRKH) syndrome (congenital aplasia of the vagina, cervix, and uterus), endometriomas have been identified, most likely due to congenital remnants of endometrial tissue [3]. Further to this, endometriosis of the bladder, peritoneum, and genitals has been described in men [4]. While more substantial evidence of endometriotic lesions in males is needed, these rare cases may be demonstrations of remnant embryological cell growth [4] or mesothelium metaplasia [5] in response to elevated estrogen. Endometriosis is characterized by lesions comprising endometrial-like cells (stromal and epithelial), which grow most commonly within the peritoneal cavity, stimulated by estrogen [6,7]. Symptoms



Citation: Churchill, M.L.; Holdsworth-Carson, S.J.; Cowley, K.J.; Luu, J.; Simpson, K.J.; Healey, M.; Rogers, P.A.W.; Donoghue, J.F. Using a Quantitative High-Throughput Screening Platform to Identify Molecular Targets and Compounds as Repurposing Candidates for Endometriosis. *Biomolecules* 2023, *13*, 965. https:// doi.org/10.3390/biom13060965

Academic Editor: Kazuhiro Tamura

Received: 27 March 2023 Revised: 30 May 2023 Accepted: 2 June 2023 Published: 8 June 2023



**Copyright:** © 2023 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). of endometriosis commonly include dysmenorrhea, chronic pelvic pain, subfertility, and reduced quality of life. At present, there is no long-term cure, with an average 5-year recurrence rate following surgery of 20.5 to 43.5% [8–10]. Current treatments involve surgery to remove lesions or medical approaches including non-steroidal anti-inflammatories, analgesics, hormonal contraceptives, gonadotropin-releasing hormone (GnRH) agonists, and others, such as synthetic androgens to reduce systemic estrogen levels [11]. Patient responsiveness to these treatments is variable, with serious and undesirable consequences such as hypestrogenism, early menopause, and infertility [11]. It is therefore a priority that new therapeutic options are identified that do not directly inhibit estrogen or the estrogen receptor (ER).

Endometriosis is dependent upon estrogen signalling, which controls several physiological processes by directly and indirectly regulating gene transcription via the genomic and non-genomic estrogen signalling pathways. The genomic pathway is mediated by two estrogen receptors (ERs), ER $\alpha$  and Er $\beta$ , which either directly translocate to the nucleus and bind estrogen response elements (EREs) in or near the promoters of target genes or indirectly via transcription factors that interact with gene promoters [12,13]. A third member of the ER family is membrane-bound G-protein-coupled ER1 (GPER1), formally known as GPR30. This receptor, in addition to membrane-bound ER $\alpha$  and ER $\beta$ , induces non-genomic estrogen signalling pathways that rapidly respond to estrogen to generate signal transduction cascades that activate PI3K/AKT, MAPK, cAMP, and others [14–16]. Further, estrogen also signals through GPR30 independent of ER and transactivates tyrosine kinase receptors such as EGFR to indirectly regulate gene expression through the activation of cell-signalling factors and transcription factors [17]. Intracellular signalling cascades also mobilise Ca<sup>2+</sup>, cAMP production, NO production, and cytoskeletal rearrangements [13,18]. Taken together, these intracellular signalling cascade elements downstream of estrogen and its receptors are therefore ideal candidates for the indirect targeting of estrogen signalling pathways via therapeutic compound screens.

The high-throughput screening (HTS) of compounds is an increasingly popular approach to in vitro drug discovery as it reduces the time and cost of drug development compared to conventional target-based discovery [19,20]. HTS involves the rapid investigation of thousands of compounds with the use of 384- or 1536-well plates and the automation of cell and compound dispensing [19]. HTS also provides an opportunity to repurpose pre-existing compounds due to previously established clinical safety parameters, which substantially fast-tracks clinical trials, drug approval, and patient access [21–24]. Therefore, applying HTS to thousands of compounds in eutopic endometrial stromal cell lines may result in the discovery of compounds that can be repurposed for a condition lacking an efficacious drug choice. The HTS may also provide new mechanistic insights into the pathogenesis of endometriosis and support the development of new diagnostics and therapeutics. As endometriosis is an estrogen-dependent disease, the inhibition of cell proliferation in an HTS compound screen must be performed under estrogenic conditions to best mimic the endometriosis environment and to ensure the compounds identified are targeting the downstream effectors of estrogen stimulation and not estrogen or its receptor directly. Therefore, the aim of this study was to develop an HTS protocol using endometriosis-patient-derived immortalized endometrial stromal cell lines for the identification of non-estrogen-receptor-targeting compounds that have translational potential as therapeutics for endometriosis.

#### 2. Materials and Methods

#### 2.1. Reagents and Kits

Phenol-red-free DMEMF/12 (PRF-DMEM/F12), charcoal-stripped FBS (CS-FBS), penicillin–streptomycin, GlutaMAX-I, insulin-transferrin-selenium (ITS), Dulbecco's phosphate-buffered saline (DPBS), TrypLE express, paraformaldehyde (PFA), 4',6-Diamidino-2-Phenylindole, Dihydrochloride (DAPI), and Cell Tracker<sup>TM</sup> Green 5-chloromethylfluorescein diacetate (CMFDA) were purchased from Thermo Fisher Scientific, Waltham, MA, USA.

Dimethyl sulphoxide (DMSO), Triton X-100, and estradiol-17 $\beta$  (estrogen) were purchased from Sigma-Aldrich, Castle Hill NSW 2154, Australia. The Real Time-Glo<sup>TM</sup> MT Cell Viability Assay was purchased from Promega, Alexandria NSW 2015, Australia. Salinomycin, panobinostat, fulvestrant, tamoxifen, and Y26732 were purchased from Selleck Chem, Houston, TX 77014 USA.

# 2.2. Cell Lines and Culture

Two hTERT immortalized endometriosis-patient-derived hESC cell lines were screened [25]. These cell lines were generated from endometrial tissue collected via curet-tage from women who had provided informed consent prior to undergoing laparoscopy at the Royal Women's Hospital, Melbourne. The ethical approval was provided by the Human Research Ethics Committee (HREC 10–43 and HREC 11–24) at the Royal Women's Hospital, Melbourne, Australia. All experiments were performed in accordance with relevant guide-lines and regulations. The 1455 hESC line was derived from a stage II endometriosis patient with no pain symptomology and no reported recurrent disease, while the 1458 hESC line was derived from a stage IV endometriosis patient with severe pain and recurrent disease. The cell lines were confirmed to be free of mycoplasma prior to their use and were cultured in PRF-DMEM/F12/10% CS-FBS/1 × penicillin–streptomycin/2 mM GlutaMAX-I/1× ITS at 37 °C in a 5% CO<sub>2</sub> humidified incubator. For high HTS, large stocks of cells were grown (P10-P20) and frozen in batches of 20 vials of  $1 \times 10^6$  cells per vial. This was carried out to ensure that each compound plate screened was performed on the same passage of cells over several weeks.

#### 2.3. Compound Libraries

Assay-ready, pre-stamped 384-well compound library plates were sourced from an open access compound library from Compounds Australia (Griffith University, Nathan QLD 4111, Australia). These plates were from the CA-FDA Collection, the CA-Epigenetics Collection, and the CA-Kinase Collection. Each plate left columns 23 and 24 reserved empty for the addition of in-house quality control (QC) compounds. The compound library contained 3517 compounds selected from the FDA clinically approved drug and compound library and the Selleck Kinase and Epigenetics drug libraries. Each pre-stamped compound plate contained 50  $\mu$ M (5× final concentration) of compound in a 650 nL aliquot. On the day of their use, pre-stamped compound plates were reconstituted with 65  $\mu$ L/well of PRF-DMEM/F12 + 5% CS-FBS using a BioTek EL406 dispenser (5  $\mu$ L cassette, high flow rate) (Biotek<sup>TM</sup> Winooski, VT, 05404-1356, USA) (Supplementary Figure S1A).

## 2.4. High-Throughput Screening (HTS)

Cells (800 cells/well) were seeded in 384-well black-walled clear-base plates (40  $\mu$ L/well) (Corning Incorporated, Corning, NY 14831, USA) in PRF-DMEM/F12 + 10% CS-FBS using a Biotek<sup>TM</sup> EL406 dispenser. The seeded plates were pulse-centrifuged briefly and incubated for 24 h in a LiCONiC high-throughput incubator (LiCONiC Instruments, Liechtenstein). Following cell adhesion, spent media were removed using a BioTek<sup>TM</sup> wash-aspirator manifold and replaced with fresh media (40  $\mu$ L) supplemented with 5% CS-FBS with 0.3% DMSO (vehicle) or 10<sup>-5</sup> M estradiol-17 $\beta$  (estrogen) (5  $\mu$ L cassette, high flow rate; aspirator coordinates: Z = 36). The compounds were then added as described below.

#### 2.5. Quality Control (QC) Compound Screen

QC positive (QC<sup>pos</sup>) control compounds *inhibit* cell survival (the outcome desired for the HTS), while the QC negative (QC<sup>neg</sup>) control compounds *promote* cell survival (the outcome not desired for the HTS). In plates pre-seeded with the 1455 or 1458 cells, the initial QC compounds, panobinostat (HDAC inhibitor), salinomycin, (antibiotic), fulvestrant, tamoxifen (selective estrogen receptor modulators), and Y26732 (ROCK inhibitor) were serially diluted (1:2 for 10 dose points from 50  $\mu$ M to 0.05  $\mu$ M) in triplicate. The cells were cultured for 3 days under both vehicle and estrogen conditions using the HTS protocol as

described above (panobinostat, salinomycin, and Y26732), while the Real Time-Glo<sup>TM</sup> MT Cell Viability Assay was used for fulvestrant and tamoxifen according to the manufacturer's instructions (Promega). Acceptable QC<sup>neg</sup> control compounds exhibited no inhibition of cell growth during estrogen treatment, while the QC<sup>pos</sup> controls exhibited >70% inhibition of cell growth at 10 µM during estrogen treatment [26,27].

### 2.6. Compound Library Screen

Pre-stamped compound plates were reconstituted (described above), and the QC compounds (50  $\mu$ M) were added manually (65  $\mu$ L/well) into columns 23 and 24 in a randomized layout including 8 wells of DMSO (1.5%), 4 wells of media supplemented with 5% CS-FBS, and 4 wells of each QC compound. All compounds were then dispensed (10  $\mu$ L/well) using the ALH 3000 Advanced Liquid Handling system (Caliper Life Sciences, Hopkinton, MA 01748, USA). The final compound concentrations were 10  $\mu$ M in DMSO (0.3%) in a total volume of 50  $\mu$ L.

Cell plates and compounds were incubated for 3 additional days, after which the media were removed, and the cells were washed with DPBS. Using the BioTek EL406 plate wash/dispense automation, the cells were then live-cell stained with 4  $\mu$ M of green CMFDA in serum-free media (25  $\mu$ L/well, 30 min, 37 °C), fixed with 2% PFA in DPBS (25  $\mu$ L/well, 15 min, room temperature [RT]), washed with DPBS (50  $\mu$ L/well), and then stained with DAPI nuclear stain (0.005 mg/mL in 50 mM Tris/PBS + 10% Triton X-100; 25  $\mu$ L/well, 20 min, RT). The cells were then washed twice and mounted in DPBS. Each plate was thermo-sealed with foil (Plateloc, Agilent, CA, USA) and was thus ready for imaging.

## 2.7. High-Content Imaging (HCI)

The wells were imaged with cell membranes and nuclei detected via segmentation; the cell nuclei were quantified using a fully automated Cellomics ArrayScan<sup>TM</sup> VTi High Content Imager (Thermo Fisher Scientific). An image acquisition software (Cellomics Scan software, v7.6.2.1, Thermo Fisher Scientific) was used to capture 16 fields per well using a  $10 \times$  objective to ensure >90% of cells were counted (Supplementary Figure S1B). The "Target Activation V4" bioapplication (Cellomics, Thermo Fisher Scientific) was used with thresholding to identify valid nuclei, which included all objects with a total area between 60 and 1400 pixels and an average intensity of less than 4095. The images were manually checked for acquisition errors. This included assessing the image sharpness, nuclei tracing, low DMSO control cell counts, and artefacts. If any of these errors were detected, the wells were removed from the analysis.

#### 2.8. Selection of Compound Hits

Cell counts were normalized to the mean cell counts of the vehicle control wells on a per-plate basis. Outlier control wells were removed from analysis if they were 3 standard deviations (SD) above or below the mean of all vehicle control wells per plate. If the cell viability in the vehicle test wells fell below 50% of the vehicle-alone wells (an indication of compound toxicity in the absence of estradiol-17 $\beta$ ), the wells were removed from all further analysis. The percentage of estrogen-driven cell growth inhibition (%E-DI) was determined by the following equation: %E-DI = 100 ( $\beta/\alpha$ ), where  $\beta$  is the normalized cell count of estrogentreated cells, and  $\alpha$  is the normalized cell count of vehicle-treated cells (Supplementary Figure S1C). The %E-DI cut-off for compound hit determination was calculated based on the average %E-DI + 2.5 standard deviations (SDs) of all samples tested [28].

The compound "hits" that were identified in the initial HTS for confirmation were individually (cherry) picked and stamped onto 384-well compound plates by Compounds Australia, leaving columns 23 and 24 blank for QC compound additions as previously described. Each compound was serially diluted (1:5 for 5 dose points from 10  $\mu$ M to 0.02  $\mu$ M) to determine the IC<sub>50</sub> (Supplemental Figure S1D). The plates were reconstituted in 65  $\mu$ L of media and dispensed onto 1455 and 1458 pre-seeded cell plates and cultured for 3 days, followed by staining, fixation, imaging, and quantification as previously described.

To determine the IC<sub>50</sub> of each compound, the % survival was plotted against the log compound concentration (log  $\mu$ M) using a nonlinear regression in GraphPad Prism 7. Compounds with an IC<sub>50</sub>  $\leq$  1  $\mu$ M were designated lead compounds.

#### 2.9. In Silico Analysis of Compound Hits

Lists of molecular targets for each compound hit were assembled using public bioinformatics databases: Drug Bank [29] and SwissTargetPrediction [30]. Protein–protein interaction (PPI), Gene Ontology (GO), and Kyoto Encyclopedia of Genes and Genomes (KEGG) pathway analyses were performed using the online String Pathway Database [31].

#### 2.10. Data and Statistical Analysis

After normalization, the assay quality and robustness were determined using the Z' factor and %CV. To determine the Z' factor for the antagonist/inhibition assay, the following equation was used: Z' factor =  $1 - 3(\delta n + \delta p)/(\mu n - \mu p)$ , where  $\delta$  = standard deviation,  $\mu$  = mean, n = negative control (MAX value), and p = positive control (MIN value), respectively. A Z' factor = 1 indicates an ideal assay in which the SD dynamic range is extremely large. A Z' score between 0.5 and 1 indicates an excellent assay as the SD dynamic range was large enough to consider test compounds as real hits. A Z' factor between 0 and 0.5 indicates a marginal but acceptable assay, and a Z' < 0 indicates an unacceptable assay [26]. This metric was determined for all control wells during the analysis. The %CV is a measurement of sample variability that measures the ratio between the standard deviation and the mean of a sample. The lower the value, the less variation occurs, demonstrating reliable reproducibility. To determine the %CV (coefficient of variance) [26,27], the following equation was used: %CV = 100 ( $\delta/\mu$ ), where  $\delta$  is the standard deviation, and  $\mu$  is the mean of the sample. An arbitrary threshold was set at %CV < 25 for an HTS to be considered acceptable, while a %CV > 25 was considered unacceptable.


#### 3. Results

#### 3.1. Establishing the HTS Quality Control Parameters for Estrogen-Stimulated Cell Survival

Two immortalized hESC cell lines were used for the HTS compound screen. The 1455 and 1458 hESC lines were previously described and are known to be hormonally responsive [25]. For this screen, the hESC cells were plated out into 384-well clear-bottom plates (800 cells/well) and after 24 h, DMSO, estradiol-17 $\beta$ , and a series of quality control (QC) compounds were added in replicate. The QC control compounds were selected based upon a strong reliable inhibition or promotion of cell survival over several wells and plates under estrogen and vehicle conditions. QC control compounds that caused estrogen-specific responses were excluded. This requirement was essential to ensure that the HTS could be reliably performed over several weeks due to the size of the library used (3517 compounds) and the different conditions required for the screen. At a single dose of 10  $\mu$ M, panobinostat and salinomycin were identified as positive QC (QCpos) control compounds as they inhibited cell survival—the outcome required for HTS compound hits (Supplementary Figure S2A–D). Panobinostat is a pan-HDAC inhibitor that inhibits epigenetic activity, while salinomycin is an antibiotic that interferes with membrane ion concentrations, resulting in osmotic pressure imbalances and cell death [32]. Both compounds reliably inhibited cell survival in a dose-dependent manner across multiple wells and plates, including the HTS dose of 10 µM.

At 10  $\mu$ M, fulvestrant was identified as the only reliable negative QC (QC<sup>neg</sup>) control compound compared to tamoxifen and Y26732, as it reliably promoted cell survival—the unwanted outcome for HTS compound hits (Supplementary Figure S2E–J). While fulvestrant is an anti-estrogenic compound that reduces cell survival by degrading the estrogen receptor, at 10  $\mu$ M, cell survival was not compromised for these two cell lines. Tamoxifen, another anti-estrogenic compound selectively inhibited cell survival during estrogen treatment compared to the vehicle and was therefore excluded. While Y26732 is a ROCK inhibitor that promotes cell survival across different concentrations, it was not

reliably mitogenic across several wells and plates for these cell lines. For these reasons, only panobinostat, salinomycin, and fulvestrant were used as reliable QC compounds for the initial HTS (Figure 1).




**Figure 1.** High-content imaging of QC wells. Images taken of 1458 (**A**) and 1455 (**B**) stromal cells exposed to DMSO, media, fulvestrant, panobinostat, salinomycin, tamoxifen, and Y26732. Cells are labelled with CFMDA (green) and DAPI (blue).

#### 3.2. Demonstrating HTS Assay Robustness

Following the normalization of the QC wells, the Z-prime (Z') factor and percentage coefficient of variance (%CV) were determined. The Z' factors for all plates, conditions, and cell lines were  $\geq$ 0.5 (Supplemental Table S1), demonstrating good separation between QC<sup>pos</sup> and QC<sup>neg</sup> [26]. The only exceptions to this were the 1455 tamoxifen wells, which fell below 0.3, indicating that tamoxifen was not suitable as a QC control compound for HTS, confirming previous findings. The %CVs of the QC<sup>neg</sup> control compounds (DMSO, fulvestrant, and Y26732) for all plates, conditions, and cell lines were <12% (Supplemental Table S1), well below our set threshold of <25% and the commonly published thresholds of 10–20% [33,34]. The %CVs of the QC<sup>pos</sup> control compounds (Panobinostat and salinomycin) for all plates, conditions, and cell lines were <22%, also below our set threshold of <25%. Again, tamoxifen proved to be an unsuitable QC control compound as the %CV ranged from 19.29 to 27.74%. These data established the initial HTS QC<sup>neg</sup> and QC<sup>pos</sup> compounds and a reliably robust HTS protocol.

#### 3.3. Identification of Compounds That Specifically Inhibited Estrogen-Stimulated Cell Survival

Samples with more than 50% reduced cell viability in vehicle wells were excluded from further analysis. This initial triage was used to ensure that only compounds that were inhibitory during estrogen treatment were considered for analysis. Following this, 283 compounds (8.1%) were removed from analysis, leaving 3234 compounds (91.9%) for analysis. To further triage this large list of compounds, a threshold of 70% inhibition of cell survival was set. This threshold was derived from the averaged percentage of estrogen-dependent inhibition (%E-DI) + 2.5 standard deviations. At this threshold, 19 compounds were identified for the 1455 cell line (Figure 2A), and 36 compounds were identified for the



1458 cell line (Figure 2B). Together, the HTS identified 55 compounds (1.7%) of interest for further validation.

**Figure 2.** Initial qHTS, new QC compounds, and hits. The hESC lines 1455 (**A**) and 1458 (**B**) were screened against 3517 compounds in an HTS. Hits were identified as compounds that induced >70% cell growth inhibition (above red dashed line). New negative QC compounds SNO1004380 and SNO1006318 were selected from cell counts imaged from 1455 (**C**) and 1458 (**D**) to improve assay precision through the reliable and strong separation between positive and negative QC compound cell counts (**E**). Data are represented as means  $\pm$  standard deviations.

#### 3.4. Identification of Lead Compounds Using Quantitative HTS (qHTS)

A quantitative HTS (qHTS) was used to serially dilute the 55 compounds identified in a 5-point dose response to generate a half-maximal inhibitory concentration (IC<sub>50</sub>). The first screen was without replicate samples, while the second screen was performed with duplicates. An average of the two independent screens was then used for all calculations. To ensure assay confidence, two extra QC<sup>neg</sup> control compounds were included in this screen. These two compounds were identified from the initial HTS as they maintained cell viability during the vehicle and estradiol-17 $\beta$  treatments more consistently than Y26732. These compounds were SNO1004380 (beta-adrenergic agonist) and SNO1006318 (steroid sapogenin) (Figure 2C,D). Therefore, the qHTS included the following control wells: 0.03% DMSO (compound solute), SNO1004380, SNO1006318, and fulvestrant as QC<sup>neg</sup> control compounds and salinomycin as the QC<sup>pos</sup> control compound. As salinomycin was consistently inhibitory to cell growth, panobinostat was no longer required. The plating efficiency and cell counts for each replicate were good with small SEM, and good separation between positive and negative QC control compounds was observed (Figure 2E).

The Z' factors for all plates, conditions, and cell lines were  $\geq 0.5$  except for the 1455 fulvestrant vehicle wells (0.43) (Supplemental Table S2). The %CV values of the QC<sup>neg</sup> control compounds (DMSO, fulvestrant, SNO1004380, and SNO1006318) for all plates, conditions, and cell lines were <13% (threshold set to <25%) (Supplemental Table S2). The two new QC<sup>neg</sup> control compounds had Z' factors and %CVs consistent with the DMSO data from the initial HTS, making these two new QC compounds suitable for future HTS. The %CV values of the QC<sup>pos</sup> control compound (salinomycin) for all plates, conditions, and cell lines were <20%, also below our set threshold of <25%. Together, these data indicate the qHTS was performing within the previously established assay parameters, enabling the identification of lead compounds from this qHTS.

From the first qHTS (samples without replicates), lead compounds were triaged using the calculated IC<sub>50</sub>, which was derived from the 5-point dilution curves. An IC<sub>50</sub> (<1  $\mu$ M) was employed as an arbitrary cut-off to reduce the identification of non-specific toxic compounds that would have no clinical or translational relevance. Twenty-three (23) compounds were identified from the initial 5-point dilution curves and used in the second qHTS (samples in duplicate) to confirm IC<sub>50</sub> data. From this screen, four compounds were identified as hits for the 1455 hESC cell line only (Aminothiazole, Hydroxyzine Pamoate, Ketotifen Fumarate, and Spectinomycin HCl), six compounds were identified as hits for the 1458 hESC cell line only (Benazepril HCl, Ceftibuten, Chloroquinalol, Editol, Repaglinide, and Silenafil Citrate), and seven compounds were identified as hits for both cell lines (7,8-Dimethoxyflavone, Chlordiazepoxide, Cytidine triphosphate disodium, Indoprofen, Pantoprazole, Pregabalin, and Promazine HCl) (Table 1). Six compounds were excluded due to a non-determination (ND) of the IC<sub>50</sub> (or IC<sub>50</sub> > 1  $\mu$ M) (Ethynodiol diacetate, Baicalin, Penicillamine, Fomepizole, Ergotamine, and Iodoquinol).

Compound Name	Generic Name	Compound Group	IC <sub>50</sub> [μM] 1455	IC <sub>50</sub> [μM] 1458
SN01006330	7,8- Dimethoxyflavone	anti-inflammatory	0.08	0.19
SN01005561	Aminothiazole	anti-microbial	0.12	2.08
SN01005320	Benazepril HCl	anti-hypertensive	9.37	0.06
SN01005071	Ceftibuten	anti-microbial	2.29	0.07
SN00852779	Chlordiazepoxide	anti-depressant	0.02	0.08
SN01005451	Chloroquinalol	anti-microbial	>20	0.08
SN01006117	Cytidine triphosphate disodium	anti-inflammatory	0.02	0.08
SN01005419	Editol	anti-inflammatory	ND	0.87
SN01004366	Hydroxyzine Pamoate	anti-histamine	0.02	ND
SN01004587	Indoprofen	anti-inflammatory	0.02	0.001

Table 1. Lead compounds identified with IC<sub>50</sub> from qHTS.

Compound Name	Generic Name	Compound Group	IC <sub>50</sub> [μM] 1455	IC <sub>50</sub> [μM] 1458
SN01004583	Ketotifen Fumarate	anti-histamine	0.22	ND
SN01005061	Pantoprazole	proton pump inhibitor	0.19	0.37
SN01005391	Pregabalin	anti-inflammatory	0.08	0.08
SN01004486	Promazine HCl	anti-psychotic	0.02	0.15
SN01005316	Repaglinide	anti-diabetic	>20	0.07
SN01005445	Sildenafil Citrate	anti-inflammatory	1.88	0.96
SN01004511	Spectinomycin HCl	anti-microbial	0.02	ND

Table 1. Cont.


## 3.5. In Silico Analysis of Compound Molecular Targets and Pathways

The final 17 lead compounds identified were entered into the online drug and compound databases (see methods section), to determine which molecular targets (proteins/receptors, etc.) were known for each compound. Lists of the molecular targets were uploaded into the STRING protein pathway analysis database to determine protein interactions and the biological pathways utilized the by cell lines.

From the four compounds identified as hits for the 1455 hESC line, nine molecular targets were identified. To perform the string analysis, extra intuitive interactive targets were generated by the STRING program (Figure 3; Supplemental Table S3). The most significant functional enrichment pathways were the GO:0050999 regulation of nitric oxide synthase activity (FDR  $4.44 \times 10^{-10}$ ), the GO:0009051 pentose-phosphate shunt, oxidative branch (FDR  $1.42 \times 10^{-07}$ ), the GO:0055114 oxidation reduction process (FDR  $7.75 \times 10^{-06}$ ), and HSA1430728 metabolism (FDR  $2.40 \times 10^{-05}$ ). The most significant KEGG signalling pathways were prolactin signalling (FDR  $8.54 \times 10^{-06}$ ), estrogen signalling (FDR 0.00069), and VEGF signalling (FDR0.0413).


From the six compounds identified as hits for the 1458 line only, 92 molecular targets were identified (Figure 4; Supplemental Table S3). The most significant functional enrichment pathways were HSA143078 metabolism (FDR  $5.39 \times 10^{-09}$ ), HSA162582 signal transduction (FDR  $1.42 \times 10^{-06}$ ), HSA168249 innate immune system (FDR 0.0061), and apoptosis (FDR 0.0427). The most significant KEGG signalling pathways were the cAMP signalling (FDR  $3.73 \times 10^{-05}$ ), PI3K-Akt signalling (FDR 0.00083), calcium signalling (FDR 0.00081), p53 signalling (FDR 0.00067), cGMP-PKG signalling (FDR 0.00039), HIF-1 signalling (FDR 0.0039), insulin signalling (FDR 0.0046), estrogen signalling (FDR 0.0046), TNF signalling (FDR 0.0131), FoxO signalling (FDR 0.0186), PPAR signalling (FDR 0.0278), and IL-17 signalling (0.0435) pathways.

From the seven compounds identified as hits for both stromal cell lines, 40 molecular targets were identified (Figure 5; Supplemental Table S3). The most significant functional enrichment pathways were the hsa04080 neuroactive ligand–receptor interactions (FDR  $1.19 \times 10^{-16}$ ), hsa00591 linoleic acid metabolism (FDR  $2.84 \times 10^{-05}$ ), hsa01100 metabolic pathways (FDR 0.0220), and hsa04010 MAPK signalling (FDR 0.0342) pathways. The most significant KEGG signalling pathways were the calcium signalling pathway (FDR  $1.34 \times 10^{-11}$ ), cAMP signalling pathway (FDR 0.0013), MAPK signalling pathway (FDR 0.0342), and cGMP-PKG signalling pathway (FDR 0.0496).




	<b>Biological Pathways</b>	FDR
	GO:0050999 Regulation of nitric-oxide synthase activity	4.44 x 10 <sup>-10</sup>
$\bigcirc$	GO:0009051 Pentose-phosphate shunt, oxidative branch	1.42 x 10 <sup>-07</sup>
$\bigcirc$	GO:0055114 Oxidation-reduction process	7.75 x 10 <sup>-06</sup>
$\bigcirc$	HSA-1430728 Metabolism	2.40 x 10 <sup>-05</sup>
Ο	Others (undefined)	

**Figure 3.** Molecular pathways identified for the hESC 1455 line. The 1455 cell line exhibited significant compound sensitivity to the regulation of nitric oxide synthase activity pathways, with an FDR  $\leq 4.44 \times 10^{-10}$ .



	<b>Biological Pathways</b>	FDR
	HSA-143078Metabolism	5.39 x 10 <sup>-09</sup>
$\bigcirc$	HSA-162582 Signal Transduction	1.42 x 10 <sup>-06</sup>
$\bigcirc$	HSA-168249 Innate Immune System	0.0061
$\bigcirc$	HSA-109581 Apoptosis	0.0427
0	Others (undefined)	

**Figure 4.** Molecular pathways identified for hESC 1458. The 1458 cell line exhibited significant compound sensitivity to metabolism pathways, with an FDR  $\leq 0.5.39 \times 10^{-09}$ .



	<b>Biological Pathways</b>	FDR
	HSA-04080 Neuroactive ligand receptor interaction	1.19 x 10 <sup>-16</sup>
$\bigcirc$	HSA-00591 Linoleic acid metabolism	2.84 x 10 <sup>-05</sup>
$\bigcirc$	HSA-01100 Metabolism	0.0220
$\bigcirc$	HSA-04010 MAPK Signalling	0.0342
Ο	Others (undefined)	

**Figure 5.** Molecular pathways common to both cell lines. The 1455 and 1458 hESC cells exhibited significant compound sensitivity to neuroactive ligand–receptor pathways, with a false discovery rate (FDR)  $\leq 1.19 \times 10^{-16}$ .

# 4. Discussion

This study presents the first high-throughput screen of large compound libraries using patient-derived endometrial stromal cells to identify new clinically translatable therapeutics for endometriosis. While the compounds identified from this HTS show promise, further evaluation will be required before clinical translation to demonstrate their efficacy and specificity as anti-endometriosis agents. HTS of large compound libraries is proving to be a valuable tool for drug discovery and drug repurposing in a variety of diseases [19]. For endometriosis, the development of HTS technologies has been limited to small-molecule drug/compound screens using the in silico analysis of predicted endometriosis targets [35–37] and high-throughput RNA sequencing [38,39] to identify new endometriosis "druggable targets". The use of patient-derived endometrial stromal cell lines in this HTS was made possible through protocol miniaturization and the identification of robust QC compounds. This protocol also benefited from the adaptable optimization of each screen [40]. This process allowed for the identification of new negative QC compounds for inclusion in the final confirmation HTS to improve assay robustness. Using adaptable

optimization, the Z' factor and the %CV were improved across all assay plates and reduced the risk of false discovery.

The qHTS was performed under estrogen stimulation to mimic the diseased state of endometriosis. In doing so, changes in gene expression induced by estrogen [41] result in compound targets that may be selectively expressed in ectopic lesions where estrogen signalling is aberrantly upregulated [42]. From this study, 17 compounds were identified that targeted a diverse range of cell surface receptors, cell-signalling molecules, and metabolic pathways, and the compound classifications ranged from antibiotic, antidepressant, anti-psychotic, anti-histamine, anti-diabetic, and anti-convulsive, to other compounds. Several of the compounds identified targeted molecules involved in endometriosisassociated processes, including downstream estrogen-signalling molecules (PI3K/AKT, NO, cAMP) (7–8 Dimethoxyflavone, and Sildenafil citrate), inflammation (Aminothiazole, Cytidine triphosphate disodium, 7–8 Dimethoxyflavone, Ketotifen Fumarate, and Sildenafil Citrate), and neurotransmitter pathways (Hydroxyzine Pamoate, Pregabalin, Promazine, and Chlordiazepoxide). Together, these findings suggest that the endometrial stromal cells may contribute to the inflammation, neuropathic pain, and nociception associated with endometriosis through the expression of pro-inflammatory cytokines [43] and neuroreceptors [44]. At present, there are no reports of these compounds being used to treat endometriosis.

From the in-silico data, the most significant KEGG pathway identified was the neuroactive ligand–receptor interaction that was identified from compounds that suppressed the growth of both cell lines during estrogen treatment. These neuroactive receptors included histamine receptors (HRH1), serotonin receptors (HTR<sub>2A</sub>) and dopamine receptors (DRD1–3). The expression of neuroactive receptors on non-neuronal cells is not new [44,45]. Histamine receptors have previously been reported on myometrial smooth muscle cells [46], endothelial cells [47], and reproductive tissues [48]. While histamine (HRH ligand) has been shown to stimulate estrogen production in granulosa cells [49], estrogen also increases histamine secretion and the release of proinflammatory cytokines from mast cells [50–52]. Further, the placental mast cell release of histamines promotes myometrial contraction during labour via the HRH1 receptors [53]. Together, the histamine receptor (and histamines) may have roles in proinflammatory responses and vascular tone in endometriosis.

Serotonin and its receptor (HTR) are well-described neurotransmitters. However, over 90% of serotonin is found outside the central nervous system (CNS) [54], and 5-HT<sub>2A</sub> is expressed in several non-neuronal tissues, including vascular smooth muscle cells [55], the myometrial smooth muscle cells of pregnant human myometrium [45], and the uterine artery [56]. In addition, 5-HT<sub>2A</sub> increases cell adhesion and cytoskeletal remodelling [57] and plays a significant role in platelet function [58]. Serotonin is also a key modulator of spinal nociceptive transmission [59] and may therefore have a role in endometriosis aetiology, pathophysiology, and pain perception, although the expression of these receptors in ectopic lesions has not been described.

Dopamine and its receptor (*DRD1–3*) are also neurotransmitters that have been associated with increased endometrial cancer severity and reduced progression-free survival [60], with the *DRD2* antagonist ONC201 shown to reduce tumour growth in vivo [61], including the highly aggressive *H3K27M* mutant gliomas [62]. *DRD2* is also expressed in human eutopic and ectopic endometrium [44], with *DRD2* polymorphisms identified as endometriosis candidate genes in women with peritoneal moderate/severe endometriosis [63]. Further to this, cabergoline, (a dopamine receptor agonist used for the treatment of hyperprolactinemia), was shown to reduce the size of endometriosis lesions via the inhibition of *VEGFR2* in a mouse model of endometriosis [64,65]. Cabergoline also induced a "lax stroma characteristic to atrophic or degenerative tissue". In a more recent study, *DRD2* was upregulated in secretory phase tissue, with cabergoline inducing rapid stromal cell decidualization [66]. Taken together, the *DRD2* agonist cabergoline may terminally differentiate endometrial stromal cells via decidualization, leading to stromal cell death [67]. While a drug like cabergoline may have potential as an endometriosis therapeutic, it targets more than just *DRD2*, making a more selective *DRD2* inhibitor such as ONC201 more suitable for targeting *DRD2* in endometriosis. ONC201 was not part of this current screening library as it is not currently approved by the FDA and is listed as investigational.

This study was designed to demonstrate the suitability of using high-throughput compound screens to identify new therapeutics and targets for endometriosis. A limitation to this protocol was the use of eutopic endometrial stromal cells. As this protocol has generated significant robust data, future screens using ectopic lesion stromal cells can now be performed. Another limitation to this study is the use of stromal cells derived from the eutopic endometria of women diagnosed with endometriosis. To further improve and validate our findings and confirm that "normal" endometrial cells are not impacted by these compounds, the inclusion of eutopic endometrial stromal cells from healthy women with no pain, no endometriosis, and no other pathology should be included in future studies and screens.

#### 5. Conclusions

This study has established the first quantitative high-throughput compound screening protocol for endometriosis and has identified several new compounds for translational consideration. While many of these compound hits and their molecular targets are well documented, there is some risk that the compounds will have off-target effects. The results of this study therefore establish a robust HTS process for evaluating novel molecules, peptides, and compounds as-endometriosis specific therapeutics. This study also demonstrates novel pathways and new molecular targets that may be significant contributors to endometriosis pathogenesis, inflammation, and nociception pain pathways and provides several new avenues of investigation for future endometriosis pathophysiological studies.

**Supplementary Materials:** The following supporting information can be downloaded at: https://www.mdpi.com/article/10.3390/biom13060965/s1, Figure S1: Quantitative High Throughput Screening (qHTS) experimental set up; Figure S2: Dose response curves for QC compounds; Table S1: Z' factors and %CV for confirmation HTS; Table S2: Z' factors and %CV for confirmation HTS; Table S3: Target genes identified from qHTS used for protein-protein-interaction.

**Author Contributions:** Study design: J.F.D., M.L.C., K.J.S. and K.J.C. Experimentation and data interpretation: J.F.D., M.L.C., S.J.H.-C., K.J.S., K.J.C., J.L., M.H. and P.A.W.R. Manuscript preparation: M.L.C. and J.F.D. with the help of all authors. J.F.D. supervised the study. All authors have read and agreed to the published version of the manuscript.

**Funding:** J.F.D. was supported by the Norman Beischer Medical Research Foundation, and the University of Melbourne, Department of Obstetrics and Gynaecology Innovation Grant. The Victorian Centre for Functional Genomics (K.J.S.) is funded by the Australian Cancer Research Foundation (ACRF), Phenomics Australia through funding from the Australian Government's National Collaborative Research Infrastructure Strategy (NCRIS) program, the Peter MacCallum Cancer Centre Foundation, and the University of Melbourne Research Collaborative Infrastructure Program.

**Institutional Review Board Statement:** The study was conducted in accordance with the Declaration of Helsinki, and approved by the Human Ethics Committee of the Royal Women's Hospital, Melbourne, Australia (HREC 10–43 and HREC 11–24; 18 September 2020).

Informed Consent Statement: Informed consent was obtained from all subjects involved in the study.

**Acknowledgments:** We also thank members from the Victorian Centre for Functional Genomics at the Peter MacCallum Cancer Centre for their technical assistance.

Conflicts of Interest: There are no conflict of interest.

#### References

- 1. Burney, R.O.; Giudice, L.C. Pathogenesis and pathophysiology of endometriosis. *Fertil. Steril.* **2012**, *98*, 511–519. [CrossRef]
- 2. AIHW. Endometriosis in Australia: Prevalence and Hospitalisations. Available online: https://www.aihw.gov.au/reports/ chronic-disease/endometriosis-prevalence-and-hospitalisations/contents/table-of-contents. (accessed on 29 August 2019).
- Konrad, L.; Dietze, R.; Kudipudi, P.K.; Horne, F.; Meinhold-Heerlein, I. Endometriosis in MRKH cases as a proof for the coelomic metaplasia hypothesis? *Reproduction* 2019, 158, R41–R47. [CrossRef]

- 4. Rei, C.; Williams, T.; Feloney, M. Endometriosis in a Man as a Rare Source of Abdominal Pain: A Case Report and Review of the Literature. *Case Rep. Obstet. Gynecol.* 2018, 2018, 2083121. [CrossRef]
- Zamecnik, M.; Hostakova, D. Endometriosis in a mesothelial cyst of tunica vaginalis of the testis. Report of a case. *Ceskoslov.* Patol. 2013, 49, 134–136.
- Bulun, S.E.; Zeitoun, K.; Takayama, K.; Noble, L.; Michael, D.; Simpson, E.; Johns, A.; Putman, M.; Sasano, H. Estrogen production in endometriosis and use of aromatase inhibitors to treat endometriosis. *Endocr. Relat. Cancer* 1999, *6*, 293–301. [CrossRef]
- Kitawaki, J.; Noguchi, T.; Amatsu, T.; Maeda, K.; Tsukamoto, K.; Yamamoto, T.; Fushiki, S.; Osawa, Y.; Honjo, H. Expression of aromatase cytochrome P450 protein and messenger ribonucleic acid in human endometriotic and adenomyotic tissues but not in normal endometrium. *Biol. Reprod.* 1997, 57, 514–519. [CrossRef]
- 8. Guo, S.W. Recurrence of endometriosis and its control. Hum. Reprod. Update 2009, 15, 441-461. [CrossRef]
- 9. Schliep, K.C.; Mumford, S.L.; Peterson, C.M.; Chen, Z.; Johnstone, E.B.; Sharp, H.T.; Stanford, J.B.; Hammoud, A.O.; Sun, L.; Buck Louis, G.M. Pain typology and incident endometriosis. *Hum. Reprod.* **2015**, *30*, 2427–2438. [CrossRef]
- Zondervan, K.T.; Becker, C.M.; Koga, K.; Missmer, S.A.; Taylor, R.N.; Vigano, P. Endometriosis. *Nat. Rev. Dis. Prim.* 2018, *4*, 9. [CrossRef]
  Giudice, L.C. Clinical practice. Endometriosis. *N. Engl. J. Med.* 2010, *362*, 2389–2398. [CrossRef]
- 12. Enmark, E.; Gustafsson, J.A. Estrogen receptors—An overview. J. Intern. Med. 1999, 246, 133–138. [CrossRef]
- 13. Bjornstrom, L.; Sjoberg, M. Mechanisms of estrogen receptor signaling: Convergence of genomic and nongenomic actions on target genes. *Mol. Endocrinol.* 2005, *19*, 833–842. [CrossRef]
- Titolo, D.; Mayer, C.M.; Dhillon, S.S.; Cai, F.; Belsham, D.D. Estrogen facilitates both phosphatidylinositol 3-kinase/Akt and ERK1/2 mitogen-activated protein kinase membrane signaling required for long-term neuropeptide Y transcriptional regulation in clonal, immortalized neurons. J. Neurosci. 2008, 28, 6473–6482. [CrossRef]
- 15. Kelly, M.J.; Levin, E.R. Rapid actions of plasma membrane estrogen receptors. Trends Endocrinol. Metab. 2001, 12, 152–156. [CrossRef]
- Vasudevan, N.; Pfaff, D.W. Membrane-initiated actions of estrogens in neuroendocrinology: Emerging principles. *Endocr. Rev.* 2007, 28, 1–19. [CrossRef]
- Maggiolini, M.; Vivacqua, A.; Fasanella, G.; Recchia, A.G.; Sisci, D.; Pezzi, V.; Montanaro, D.; Musti, A.M.; Picard, D.; Ando, S. The G protein-coupled receptor GPR30 mediates c-fos up-regulation by 17beta-estradiol and phytoestrogens in breast cancer cells. J. Biol. Chem. 2004, 279, 27008–27016. [CrossRef]
- 18. Edwards, D.P. Regulation of signal transduction pathways by estrogen and progesterone. *Annu. Rev. Physiol.* **2005**, *67*, 335–376. [CrossRef]
- Hughes, J.P.; Rees, S.; Kalindjian, S.B.; Philpott, K.L. Principles of early drug discovery. *Br. J. Pharmacol.* 2011, 162, 1239–1249. [CrossRef]
- Ashburn, T.T.; Thor, K.B. Drug repositioning: Identifying and developing new uses for existing drugs. *Nat. Rev. Drug Discov.* 2004, 3, 673–683. [CrossRef]
- Hernandez, J.J.; Pryszlak, M.; Smith, L.; Yanchus, C.; Kurji, N.; Shahani, V.M.; Molinski, S.V. Giving Drugs a Second Chance: Overcoming Regulatory and Financial Hurdles in Repurposing Approved Drugs as Cancer Therapeutics. *Front. Oncol.* 2017, 7, 273. [CrossRef]
- 22. Fox, S.; Farr-Jones, S.; Sopchak, L.; Boggs, A.; Nicely, H.W.; Khoury, R.; Biros, M. High-throughput screening: Update on practices and success. J. Biomol. Screen. 2006, 11, 864–869. [CrossRef] [PubMed]
- Badr, C.E.; Wurdinger, T.; Tannous, B.A. Functional drug screening assay reveals potential glioma therapeutics. Assay Drug Dev. Technol. 2011, 9, 281–289. [CrossRef] [PubMed]
- 24. Roberts, J.P.; Tesdorpf, U.J. Phenotypic Drug Discovery with High Content Screening; Perkin Elmer: Waltham, MA, USA, 2015.
- Holdsworth-Carson, S.J.; Colgrave, E.M.; Donoghue, J.F.; Fung, J.N.; Churchill, M.L.; Mortlock, S.; Paiva, P.; Healey, M.; Montgomery, G.W.; Girling, J.E.; et al. Generation of immortalized human endometrial stromal cell lines with different endometriosis risk genotypes. *Mol. Hum. Reprod.* 2019, 25, 194–205. [CrossRef] [PubMed]
- Zhang, J.H.; Chung, T.D.; Oldenburg, K.R. A Simple Statistical Parameter for Use in Evaluation and Validation of High Throughput Screening Assays. J. Biomol. Screen. 1999, 4, 67–73. [CrossRef] [PubMed]
- Saunders, D.N.; Falkenberg, K.J.; Simpson, K.J. High-throughput approaches to measuring cell death. *Cold Spring Harb. Protoc.* 2014, 591–601. [CrossRef]
- 28. Dragiev, P.; Nadon, R.; Makarenkov, V. Two effective methods for correcting experimental high-throughput screening data. *Bioinformatics* **2012**, *28*, 1775–1782. [CrossRef]
- 29. Wishart, D.S.; Feunang, Y.D.; Guo, A.C.; Lo, E.J.; Marcu, A.; Grant, J.R.; Sajed, T.; Johnson, D.; Li, C.; Sayeeda, Z.; et al. DrugBank 5.0: A major update to the DrugBank database for 2018. *Nucleic Acids Res.* **2018**, *46*, D1074–D1082. [CrossRef]
- Daina, A.; Michielin, O.; Zoete, V. SwissTargetPrediction: Updated data and new features for efficient prediction of protein targets of small molecules. *Nucleic Acids Res.* 2019, 47, W357–W364. [CrossRef]
- Szklarczyk, D.; Gable, A.L.; Lyon, D.; Junge, A.; Wyder, S.; Huerta-Cepas, J.; Simonovic, M.; Doncheva, N.T.; Morris, J.H.; Bork, P.; et al. STRING v11: Protein-protein association networks with increased coverage, supporting functional discovery in genome-wide experimental datasets. *Nucleic Acids Res.* 2019, 47, D607–D613. [CrossRef]
- 32. Mitani, M.; Yamanishi, T.; Miyazaki, Y.; Otake, N. Salinomycin effects on mitochondrial ion translocation and respiration. *Antimicrob. Agents Chemother.* **1976**, *9*, 655–660. [CrossRef]

- 33. Iversen, P.W.; Beck, B.; Chen, Y.F.; Dere, W.; Devanarayan, V.; Eastwood, B.J.; Farmen, M.W.; Iturria, S.J.; Montrose, C.; Moore, R.A.; et al. HTS Assay Validation. In Assay Guidance Manual; Sittampalam, G.S., Grossman, A., Brimacombe, K., Arkin, M., Auld, D., Austin, C., Baell, J., Bejcek, B., Caaveiro, J.M.M., Chung, T.D.Y., et al., Eds.; Eli Lilly & Company and the National Center for Advancing Translational Sciences: Bethesda, MD, USA, 2004.
- Sittampalam, G.S.; Kahl, S.D.; Janzen, W.P. High-throughput screening: Advances in assay technologies. *Curr. Opin. Chem. Biol.* 1997, 1, 384–391. [CrossRef] [PubMed]
- Kuohung, W.; Burnett, M.; Mukhtyar, D.; Schuman, E.; Ni, J.; Crowley, W.F.; Glicksman, M.A.; Kaiser, U.B. A high-throughput small-molecule ligand screen targeted to agonists and antagonists of the G-protein-coupled receptor GPR54. *J. Biomol. Screen.* 2010, 15, 508–517. [CrossRef] [PubMed]
- 36. Baurle, S.; Nagel, J.; Peters, O.; Brauer, N.; Ter Laak, A.; Preusse, C.; Rottmann, A.; Heldmann, D.; Bothe, U.; Blume, T.; et al. Identification of a Benzimidazolecarboxylic Acid Derivative (BAY 1316957) as a Potent and Selective Human Prostaglandin E2 Receptor Subtype 4 (hEP4-R) Antagonist for the Treatment of Endometriosis. *J. Med. Chem.* 2019, 62, 2541–2563. [CrossRef] [PubMed]
- 37. Liu, S.; Yang, Y.; Li, W.; Tian, X.; Cui, H.; Zhang, Q. Identification of small-molecule ligands that bind to MiR-21 as potential therapeutics for endometriosis by screening ZINC database and in-vitro assays. *Gene* **2018**, *662*, 46–53. [CrossRef]
- Rekker, K.; Saare, M.; Eriste, E.; Tasa, T.; Kukuskina, V.; Roost, A.M.; Anderson, K.; Samuel, K.; Karro, H.; Salumets, A.; et al. High-throughput mRNA sequencing of stromal cells from endometriomas and endometrium. *Reproduction* 2017, 154, 93–100. [CrossRef] [PubMed]
- Zhao, L.; Gu, C.; Ye, M.; Zhang, Z.; Li, L.; Fan, W.; Meng, Y. Integration analysis of microRNA and mRNA paired expression profiling identifies deregulated microRNA-transcription factor-gene regulatory networks in ovarian endometriosis. *Reprod. Biol. Endocrinol.* 2018, 16, 4. [CrossRef]
- 40. An, W.F.; Tolliday, N. Cell-based assays for high-throughput screening. Mol. Biotechnol. 2010, 45, 180–186. [CrossRef]
- Ren, C.E.; Zhu, X.; Li, J.; Lyle, C.; Dowdy, S.; Podratz, K.C.; Byck, D.; Chen, H.B.; Jiang, S.W. Microarray analysis on gene regulation by estrogen, progesterone and tamoxifen in human endometrial stromal cells. *Int. J. Mol. Sci.* 2015, 16, 5864–5885. [CrossRef]
- Huhtinen, K.; Desai, R.; Stahle, M.; Salminen, A.; Handelsman, D.J.; Perheentupa, A.; Poutanen, M. Endometrial and endometriotic concentrations of estrone and estradiol are determined by local metabolism rather than circulating levels. *J. Clin. Endocrinol. Metab.* 2012, 97, 4228–4235. [CrossRef]
- 43. Malutan, A.M.; Drugan, T.; Costin, N.; Ciortea, R.; Bucuri, C.; Rada, M.P.; Mihu, D. Pro-inflammatory cytokines for evaluation of inflammatory status in endometriosis. *Cent. Eur. J. Immunol.* **2015**, *40*, 96–102. [CrossRef]
- Novella-Maestre, E.; Carda, C.; Ruiz-Sauri, A.; Garcia-Velasco, J.A.; Simon, C.; Pellicer, A. Identification and quantification of dopamine receptor 2 in human eutopic and ectopic endometrium: A novel molecular target for endometriosis therapy. *Biol. Reprod.* 2010, *83*, 866–873. [CrossRef] [PubMed]
- 45. Cordeaux, Y.; Pasupathy, D.; Bacon, J.; Charnock-Jones, D.S.; Smith, G.C. Characterization of serotonin receptors in pregnant human myometrium. *J. Pharmacol. Exp. Ther.* **2009**, *328*, 682–691. [CrossRef] [PubMed]
- 46. Willets, J.M.; Taylor, A.H.; Shaw, H.; Konje, J.C.; Challiss, R.A. Selective regulation of H1 histamine receptor signalling by G protein-coupled receptor kinase 2 in uterine smooth muscle cells. *Mol. Endocrinol.* **2008**, *22*, 1893–1907. [CrossRef] [PubMed]
- 47. Heltianu, C.; Simionescu, M.; Simionescu, N. Histamine receptors of the microvascular endothelium revealed in situ with a histamine-ferritin conjugate: Characteristic high-affinity binding sites in venules. *J. Cell Biol.* **1982**, *93*, 357–364. [CrossRef]
- 48. Dey, S.K.; Villanueva, C.; Abdou, N.I. Histamine receptors on rabbit blastocyst and endometrial cell membranes. *Nature* **1979**, 278, 648–649. [CrossRef] [PubMed]
- 49. Bodis, J.; Tinneberg, H.R.; Schwarz, H.; Papenfuss, F.; Torok, A.; Hanf, V. The effect of histamine on progesterone and estradiol secretion of human granulosa cells in serum-free culture. *Gynecol. Endocrinol.* **1993**, *7*, 235–239. [CrossRef]
- Cocchiara, R.; Albeggiani, G.; Di Trapani, G.; Azzolina, A.; Lampiasi, N.; Rizzo, F.; Diotallevi, L.; Gianaroli, L.; Geraci, D. Estradiol enhances in vitro the histamine release induced by embryonic histamine-releasing factor (EHRF) from uterine mast cells. *Hum. Reprod.* 1992, 7, 1036–1041. [CrossRef]
- Zaitsu, M.; Narita, S.; Lambert, K.C.; Grady, J.J.; Estes, D.M.; Curran, E.M.; Brooks, E.G.; Watson, C.S.; Goldblum, R.M.; Midoro-Horiuti, T. Estradiol activates mast cells via a non-genomic estrogen receptor-alpha and calcium influx. *Mol. Immunol.* 2007, 44, 1977–1985. [CrossRef]
- 52. Zierau, O.; Zenclussen, A.C.; Jensen, F. Role of female sex hormones, estradiol and progesterone, in mast cell behavior. *Front. Immunol.* **2012**, *3*, 169. [CrossRef]
- 53. Needham, K.; Fadia, M.; Dahlstrom, J.E.; Harrington, K.; Shadbolt, B.; Robson, S.J. Significance of mast cell distribution in placental tissue and membranes in spontaneous preterm birth. *J. Inflamm. Res.* **2016**, *9*, 141–145. [CrossRef]
- 54. Cote, F.; Fligny, C.; Bayard, E.; Launay, J.M.; Gershon, M.D.; Mallet, J.; Vodjdani, G. Maternal serotonin is crucial for murine embryonic development. *Proc. Natl. Acad. Sci. USA* **2007**, *104*, 329–334. [CrossRef]
- Hayashi, T.; Sumi, D.; Matsui-Hirai, H.; Fukatsu, A.; Arockia Rani, P.J.; Kano, H.; Tsunekawa, T.; Iguchi, A. Sarpogrelate HCl, a selective 5-HT2A antagonist, retards the progression of atherosclerosis through a novel mechanism. *Atherosclerosis* 2003, *168*, 23–31. [CrossRef] [PubMed]

- 56. Karlsson, C.; Bodelsson, G.; Bodelsson, M.; Stjernquist, M. 5-Hydroxytryptamine contracts human uterine artery smooth muscle predominantly via 5-HT2 receptors. *Hum. Reprod.* **1997**, *12*, 361–367. [CrossRef] [PubMed]
- 57. John Jayakumar, J.A.K.; Panicker, M.M.; Basu, B. Serotonin 2A (5-HT(2A)) receptor affects cell-matrix adhesion and the formation and maintenance of stress fibers in HEK293 cells. *Sci. Rep.* 2020, *10*, 21675. [CrossRef] [PubMed]
- Duerschmied, D.; Suidan, G.L.; Demers, M.; Herr, N.; Carbo, C.; Brill, A.; Cifuni, S.M.; Mauler, M.; Cicko, S.; Bader, M.; et al. Platelet serotonin promotes the recruitment of neutrophils to sites of acute inflammation in mice. *Blood* 2013, 121, 1008–1015. [CrossRef] [PubMed]
- Courteix, C.; Dupuis, A.; Martin, P.Y.; Sion, B. 5-HT2A Receptors and Pain. In 5-HT2A Receptors in the Central Nervous System; The Receptors; Guiard, B., Di Giovanni, G., Eds.; Humana Press: Cham, Switzerland, 2018; Volume 32, pp. 339–352. [CrossRef]
- 60. Pierce, S.R.; Fang, Z.; Yin, Y.; West, L.; Asher, M.; Hao, T.; Zhang, X.; Tucker, K.; Staley, A.; Fan, Y.; et al. Targeting dopamine receptor D2 as a novel therapeutic strategy in endometrial cancer. *J. Exp. Clin. Cancer Res.* **2021**, *40*, 61. [CrossRef]
- 61. Ishizawa, J.; Kojima, K.; Chachad, D.; Ruvolo, P.; Ruvolo, V.; Jacamo, R.O.; Borthakur, G.; Mu, H.; Zeng, Z.; Tabe, Y.; et al. ATF4 induction through an atypical integrated stress response to ONC201 triggers p53-independent apoptosis in hematological malignancies. *Sci. Signal.* **2016**, *9*, ra17. [CrossRef]
- Chi, A.S.; Tarapore, R.S.; Hall, M.D.; Shonka, N.; Gardner, S.; Umemura, Y.; Sumrall, A.; Khatib, Z.; Mueller, S.; Kline, C.; et al. Pediatric and adult H3 K27M-mutant diffuse midline glioma treated with the selective DRD2 antagonist ONC201. *J. Neurooncol.* 2019, 145, 97–105. [CrossRef]
- 63. Bilibio, J.P.; Matte, U.; de Conto, E.; Genro, V.K.; Souza, C.A.; Cunha-Filho, J.S. Dopamine receptor D2 genotype (3438) is associated with moderate/severe endometriosis in infertile women in Brazil. *Fertil. Steril.* **2013**, *99*, 1340–1345. [CrossRef]
- 64. Delgado-Rosas, F.; Gomez, R.; Ferrero, H.; Gaytan, F.; Garcia-Velasco, J.; Simon, C.; Pellicer, A. The effects of ergot and non-ergotderived dopamine agonists in an experimental mouse model of endometriosis. *Reproduction* **2011**, *142*, 745–755. [CrossRef]
- Novella-Maestre, E.; Carda, C.; Noguera, I.; Ruiz-Sauri, A.; Garcia-Velasco, J.A.; Simon, C.; Pellicer, A. Dopamine agonist administration causes a reduction in endometrial implants through modulation of angiogenesis in experimentally induced endometriosis. *Hum. Reprod.* 2009, 24, 1025–1035. [CrossRef] [PubMed]
- Yu, J.; Berga, S.L.; Meng, Q.; Xia, M.; Kohout, T.A.; van Duin, M.; Taylor, R.N. Cabergoline Stimulates Human Endometrial Stromal Cell Decidualization and Reverses Effects of Interleukin-1beta In Vitro. *J. Clin. Endocrinol. Metab.* 2021, 106, 3591–3604. [CrossRef] [PubMed]
- Das, S.K. Cell cycle regulatory control for uterine stromal cell decidualization in implantation. *Reproduction* 2009, 137, 889–899.
  [CrossRef] [PubMed]

**Disclaimer/Publisher's Note:** The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.