

Supplementary Material

# Metabolites of Prickly Rose: Chemodiversity and Digestive-Enzyme-Inhibiting Potential of *Rosa acicularis* and the Main Ellagitannin Rugosin D

Daniil N. Olennikov <sup>1,\*</sup>, Vladimir V. Chempsov <sup>2</sup> and Nadezhda K. Chirikova <sup>2</sup>

- <sup>1</sup> Laboratory of Medical and Biological Research, Institute of General and Experimental Biology, Siberian Division, Russian Academy of Science, 670047 Ulan-Ude, Russia  
<sup>2</sup> Department of Biology, Institute of Natural Sciences, North-Eastern Federal University, 677027 Yakutsk, Russia; soulkage94@gmail.com (V.V.C.), hofnung@mail.ru (N.K.C.)  
\* Correspondence: olennikovdn@mail.ru; Tel.: +8-902-160-06-27

---

## Content

**Table S1.** Reference standards and internal standards used for the qualitative and quantitative analysis by HPLC-MS and HPLC-DAD.






**Table S2.** Optimized MRM transitions of compounds **1–123** in HPLC-MS/MS analysis.


**Table S3.** Regression equations, correlation coefficients, standard deviation, limits of detection, limits of quantification and linear ranges for 34 reference standards.

**Figure S1.** HPLC-DAD chromatogramms of gastro-intestinal digestion markers after digestion of artificial substrate mixture.

**Table S1.** Reference standards and internal standards used for the qualitative and quantitative analysis by HPLC-MS and HPLC-DAD.

Compound (No in Table 2), Formula	Purity ( $\geq$ ), %	Manufacturer (cat. No) or Isolation Reference *	Used for Qualitative (QL) and Quantitative (QT) Analysis of Compounds (No in Table 2)
1-O-Galloyl- $\beta$ -D-glucopyranose ( <b>1</b> ) 	90	Sigma (69288)	QL: <b>1</b> . QT: <b>1</b>
Gallic acid ( <b>2</b> ) 	97	Sigma (G7384)	QL: <b>2</b> . QT: <b>2</b>
1-O-Caffeoylquinic acid ( <b>3</b> ) 	98	ChemFaces (CFN99121)	QL: <b>3</b> . QT: <b>3, 7</b>
2-Pyrone-4,6-dicarboxylic acid ( <b>8</b> ) 	95	Lab collection/isolated from <i>Comarum palustre</i> [87]	QL: <b>8</b> . QT: <b>8</b>
Epigallocatechin ( <b>9</b> ) 	98	Extrasynthese (0979S)	QL: <b>9</b> . QT: <b>9</b>
5-O-Caffeoylquinic acid ( <b>10</b> ) 	95	Sigma (C3878)	QL: <b>10</b> . QT: <b>10</b>
Catechin ( <b>11</b> ) 	99	Extrasynthese (0976S)	QL: <b>11</b> . QT: <b>11</b>
1,6-Di-O-galloyl glucose ( <b>13</b> ) 	92	Toronto (D293195)	QL: <b>13</b> . QT: <b>12, 13</b>
Tellimagrandin I <sub>1</sub> ( <b>14</b> ) 	No data	Lab collection/isolated from <i>Filipendula ulmaria</i> [88]	QL: <b>14</b>
Epicatechin ( <b>15</b> ) 	99	Extrasynthese (0977S)	QL: <b>15</b> . QT: <b>15</b>


<b>Rugosin B<sub>1</sub> (16)</b>	No data	Lab collection/isolated from <i>Filipendula ulmaria</i> [88]	QL: 16
			
<b>1,3,6-Tri-O-galloyl glucose (18)</b>	90	Sigma (78864)	QL: 18. QT: 17, 18, 20
			
<b>Tellimagrandin I<sub>2</sub> (19)</b>	No data	Lab collection/isolated from <i>Filipendula ulmaria</i> [88]	QL: 19
			
<b>Tellimagrandin II (21)</b>	No data	Lab collection/isolated from <i>Filipendula ulmaria</i> [88]	QL: 21
			
<b>Rugosin B<sub>2</sub> (22)</b>	No data	Lab collection/isolated from <i>Filipendula ulmaria</i> [88]	QL: 22
			



No data      Lab collection/isolated  
from *Filipendula ulmaria*

QL: 24


[88]



No data      Lab collection/isolated  
from *Filipendula ulmaria*

QL: 25


[88]



No data      Lab collection/isolated  
from *Filipendula ulmaria*

QL: 26



[88]



No data      Lab collection/isolated  
from *Filipendula ulmaria*

QL: 27










[88]











No data      Lab collection/isolated  
from *Punica granatum*







QL: 28

[40]

<b>Ellagic acid (29)</b>		95	Sigma (E2250)	QL: 29. QT: 6, 29, 30, 31, 36, 39, 43, 44, 90, 91, 94, 110, 112, 113
<b>1,2,3,6-Tetra-O-galloyl glucose (33)</b>		98	BioBioPha (BBPP05354)	QL: 33. QT: 32, 33
<b>1,2,3,4,6-Penta-O-galloyl glucose (47)</b>		96	Sigma (67548)	QL: 47. QT: 46, 47, 48, 49, 50, 92
<b>Rutin = quercetin-3-O-rutinoside (65)</b>		95	Sigma (R2303)	QL: 65. QT: 52, 54, 57, 61, 66, 67, 65
<b>Miquelianin = quercetin-3-O-β-D-glucuronopyranoside (68)</b>		95	Sigma (00310590)	QL: 68. QT: 68
<b>Kaempferol-3-O-β-D-glucuronopyranoside (71)</b>		97	Sigma (79273)	QL: 71. QT: 58, 60, 62, 69, 70, 71
<b>Quercitrin = quercetin-3-O-α-L-rhamnopyranoside (72)</b>		98	Sigma (740580)	QL: 72. QT: 72
<b>Helichrysinoside = quercetin-3-O-(6''-O-p-coumaroyl)-β-D- glucopyranoside (73)</b>		No data	Lab collection/isolated from <i>Spiraea salicifolia</i> [89]	QL: 73, 78, 79, 80, 81
<b>Juglanin = kaempferol-3-O-α-L-arabinofuranoside (74)</b>		98	MCE (HY-N3442)	QL: 74. QT: 74

Tiliroside =			
kaempferol-3-O-(6''-O-p-coumaroyl)-β-D-glucopyranoside ( <b>75</b> )	98	Sigma (79257)	QL: <b>75</b> . QT: <b>75, 82, 83, 84, 85, 86</b>
			
Kaempferol-3-O-(6''-galloyl)-β-D-glucopyranoside ( <b>76</b> )	98	ChemFaces (CFN92383)	QL: <b>76</b> . QT: <b>59, 76, 77</b>
			
Peonidin 3,5-di-O-β-D-glucopyranoside ( <b>93</b> )	90	Extrasynthese (0918)	QL: <b>93</b> . QT: <b>93</b>
			
Procyanidin B <sub>3</sub> ( <b>97</b> )	95	Carbosynth (FP65542)	QL: <b>97</b> . QT: <b>97</b>
			
Procyanidin B <sub>1</sub> ( <b>98</b> )	99	MCE (HY-N0795)	QL: <b>98</b> . QT: <b>34, 35, 37, 38, 40, 41, 42, 45, 95, 96, 98, 99, 101, 102, 103, 104, 114, 115, 116, 117, 118, 119, 120, 121, 122, 123</b>
			
Procyanidin B <sub>4</sub> ( <b>100</b> )	98	ChemFaces (CFN91171)	QL: <b>100</b> . QT: <b>100</b>
			
Procyanidin B <sub>2</sub> ( <b>105</b> )	90	Sigma (42157)	QL: <b>105</b> . QT: <b>105</b>
			
Procyanidin B <sub>2</sub> 3-O-gallate ( <b>106</b> )	95	TransMIT (C6201-P044)	QL: <b>106</b> . QT: <b>106, 108</b>
			

Procyanidin B <sub>2</sub> 3''-O-gallate ( <b>107</b> )		97	ALB (ALB-RS-9598)	QL: <b>107</b> . QT: <b>107</b>
Procyanidin B <sub>2</sub> 3,3''-di-O-gallate ( <b>109</b> )		97	BOC (79907-44-1)	QL: <b>109</b> . QT: <b>109</b>
Amritoside = 1-O-ellagoyl-gentiobiose ( <b>110</b> )		No data	Lab collection/isolated from <i>Punica granatum</i> [40]	QL: <b>110</b>
Epicatechin gallate ( <b>111</b> )		97.5	Extrasynthese (0978S)	QL: <b>111</b> . QT: <b>111</b>
Eschweilenol C = ellagic acid 4-O- $\alpha$ -L-rhamnopyranoside ( <b>113</b> )		98	ChemFaces (CFN95363)	QL: <b>113</b> . QT: <b>113</b>
Corilagin		98	Sigma (G0424)	QT: <b>4, 5, 14, 16, 19, 21, 22, 23, 24, 25, 26, 27, 87, 88, 89</b>
Taxifolin = dihydroquercetin		95	Sigma (PHL89284)	QT: <b>51, 53, 55, 63, 64</b>
Quercetin-3-O-(6''-galloyl)- $\beta$ -D-galactopyranoside		98	MCE (HY-N7024)	QT: <b>56</b>
3',4'-Di-O-acetyl-cis-khellactone		No data	Lab collection/isolated from <i>Phlojodicarpus sibiricus</i> [90]	Internal standard 1

Isovitexin-7,2"-di-O-glucoside 	No data	Lab collection/isolated from <i>Silene dioica</i> [91]	Internal standard 2
Isojaceoside 	No data	Lab collection/isolated from <i>Artemisia frigida</i> [92]	Internal standard 3
6-O-Sinapoyl catalpol 	No data	Lab collection/isolated from <i>Veronica longifolia</i> [93]	Internal standard 4
Neomangiferin 	No data	Lab collection/isolated from <i>Gentiana daurica</i> [94]	Internal standard 5
<i>p</i> -Nitrophenol 	99	Sigma (241326)	Marker of ethylene- <i>p</i> -nitrophenyl- $\alpha$ -D-maltoheptose digestion
7-Amino-4-methylcoumarin 	99	Sigma (257370)	Marker of Na-benzoyl-L-arginine-7-amino-4-methylcoumarin hydrochloride digestion
4-Methylumbelliferon 	98	Sigma (M1381)	Marker of 4-methylumbelliferyl heptanoate digestion

\* Manufacturers list: ALB—ALB Materials Inc (Henderson, NV, USA); BioBioPha—BioBioPha (Kunming, Yunnan, PRC); BOC—BOC Sciences (Shirley, NY, USA); Carbosynth—Carbosynth Ltd (Compton, Great Britain); ChemFaces—ChemFaces (Wuhan, Hubei, PRC); Extrasynthese—Extrasynthese (Lyon, France); MCE—MCE Med Chem Express (Monmouth, NJ, USA); Sigma—Sigma-Aldrich (St. Louis, MO, USA); Toronto—Toronto Research Chemicals (North York, ON, Canada); TransMIT—TransMIT GmbH (Gießen, Germany).

**Table S2.** Optimized MRM transitions of compounds **1–123** in HPLC-MS/MS analysis.

No	Compound	Ionization	Precursor	Quan-tifier
1	1-O-Galloyl glucose (glucogallin)	Negative	331	169
2	Gallic acid	Negative	169	125
3	1-O-Caffeoylquinic acid	Negative	353	135
4	Galloyl-hexahydroxydiphenoyl-di-O-hexoside	Negative	795	301
5	Galloyl-hexahydroxydiphenoyl-di-O-hexoside	Negative	795	301
6	Ellagic acid tri-O-hexoside	Negative	787	301
7	Caffeic acid O-hexoside	Negative	341	179
8	2-Pyrone-4,6-dicarboxylic acid	Negative	183	111
9	Epigallocatechin	Negative	305	137
10	5-O-Caffeoylquinic acid	Negative	353	165
11	Catechin	Negative	289	137
12	Di-O-galloyl hexose	Negative	483	169
13	1,6-Di-O-galloyl glucose	Negative	483	125
14	Tellimagrandin I <sub>1</sub>	Negative	392	301
15	Epicatechin	Negative	289	137
16	Rugosin B <sub>1</sub>	Negative	476	301
17	Tri-O-galloyl hexose	Negative	635	169
18	1,3,6-Tri-O-galloyl glucose	Negative	635	169
19	Tellimagrandin I <sub>2</sub>	Negative	392	301
20	Tri-O-galloyl hexose	Negative	635	169
21	Tellimagrandin II <sub>1</sub>	Negative	468	301
22	Rugosin B <sub>2</sub>	Negative	476	301
23	Tellimagrandin II <sub>2</sub>	Negative	468	301
24	Rugosin E <sub>1</sub>	Negative	860	301
25	Rugosin E <sub>2</sub>	Negative	860	301
26	Rugosin A	Negative	552	301
27	Rugosin D	Negative	623	301
28	1-O-Ellagoyl glucose	Negative	463	301
29	Ellagic acid	Negative	301	185
30	Ellagic acid methyl ester O-hexoside	Negative	477	301
31	Ellagic acid methyl ester O-hexoside	Negative	477	301
32	Tetra-O-galloyl hexose	Negative	787	169
33	1,2,3,6-Tetra-O-galloyl glucose	Negative	787	169
34	Epicatechin / catechin trimer	Negative	865	451
35	Epicatechin / catechin trimer	Negative	865	451
36	Ellagic acid dimethyl ester O-hexoside	Negative	491	329
37	Epicatechin / catechin trimer	Negative	865	577
38	Epicatechin / catechin trimer	Negative	865	577
39	Di-ellagoyl methyl ester O-hexoside	Negative	761	477
40	Epicatechin / catechin tetramer	Negative	1153	577
41	Epicatechin / catechin tetramer	Negative	1153	577


42	Epicatechin / catechin tetramer	Negative	1153	577
43	Di-ellagoyl dimethyl ester O-hexoside	Negative	775	477
44	Ellagic acid trimethyl ester O-hexoside	Negative	505	343
45	Epicatechin / catechin tetramer	Negative	1153	863
46	Penta-O-galloyl hexose	Negative	939	787
47	1,2,3,4,6-Penta-O-galloyl glucose	Negative	939	787
48	Hexa-O-galloyl hexose	Negative	1091	787
49	Hexa-O-galloyl hexose	Negative	1091	787
50	Hexa-O-galloyl hexose	Negative	1091	787
51	Dihydroquercetin di-O-hexuronoside-di-O-hexoside	Negative	979	817
52	Quercetin di-O-hexuronoside-tri-O-hexoside	Negative	1139	815
53	Dihydroquercetin O-hexuronoside-di-O-hexoside	Negative	803	641
54	Quercetin O-hexuronoside-tri-O-hexoside	Negative	963	801
55	Dihydroquercetin O-hexuronoside-O-hexoside	Negative	641	479
56	Quercetin O-hexuronoside-tri-O-hexoside-O-gallate	Negative	1115	639
57	Quercetin di-O-hexuronoside-di-O-hexoside	Negative	977	653
58	Kaempferol O-hexuronoside-tri-O-hexoside	Negative	947	623
59	Kaempferol O-hexuronoside-tri-O-hexoside-O-gallate	Negative	1099	623
60	Kaempferol di-O-hexuronoside-di-O-hexoside	Negative	961	461
61	Quercetin O-hexuronoside-di-O-hexoside	Negative	801	477
62	Kaempferol O-hexuronoside-di-O-hexoside	Negative	785	461
63	Dihydroquercetin O-hexoside	Negative	465	303
64	Dihydroquercetin O-hexuronoside	Negative	479	303
65	Quercetin-O-rutinoside (rutin)	Negative	609	463
66	Quercetin O-hexuronoside-O-hexoside	Negative	639	477
67	Quercetin O-hexuronoside-O-pentoside	Negative	609	477
68	Quercetin-3-O-glucuronide (miquelianin)	Negative	477	301
69	Kaempferol O-hexuronoside-O-hexoside	Negative	623	461
70	Kaempferol O-hexuronoside-O-pentoside	Negative	593	461
71	Kaempferol-3-O-glucuronide	Negative	461	285
72	Quercetin-3-O-rhamnoside (quercitrin)	Negative	447	301
73	Quercetin-3-O-(6''-O-p-coumaroyl)-glucoside (helichryssoside)	Negative	609	463
74	Kaempferol-3-O-arabinoside (juglanin)	Negative	417	285
75	Kaempferol-3-O-(6''-O-p-coumaroyl)-glucoside (tiliroside)	Negative	593	447
76	Kaempferol-3-O-(6''-O-galloyl)-glucoside	Negative	599	447
77	Kaempferol O-hexuronoside-O-gallate	Negative	613	461
78	Quercetin O-hexoside-O-p-coumarate-O-gallate	Negative	761	609
79	Quercetin O-hexuronoside-O-p-coumarate-O-gallate	Negative	775	623
80	Quercetin O-hexoside-di-O-p-coumarate	Negative	755	463
81	Quercetin O-hexuronoside-di-O-p-coumarate	Negative	769	477
82	Kaempferol O-hexuronoside-O-p-coumarate-O-gallate	Negative	759	607
83	Kaempferol O-hexuronoside-O-hexoside-di-O-p-coumarate	Negative	915	461
84	Kaempferol O-hexoside-di-O-p-coumarate	Negative	739	447

85	Kaempferol <i>O</i> -hexuronoside-di- <i>O</i> - <i>p</i> -coumarate	Negative	753	461
86	Kaempferol <i>O</i> -hexuronoside-tri- <i>O</i> - <i>p</i> -coumarate	Negative	899	461
87	Hexahydroxydiphenoyl-tri- <i>O</i> -galloyl-hexose (tellimagrandin II isomer)	Negative	468	301
88	Valoneoyl-tri- <i>O</i> -galloyl-hexose (rugosin A isomer)	Negative	552	301
89	Hexahydroxydiphenoyl-valoneoyl-tetra- <i>O</i> -galloyl-di- <i>O</i> -hexose (rugosin D isomer)	Negative	623	765
90	Ellagic acid <i>O</i> -hexoside	Negative	463	301
91	Epicatechin / catechin tetramer	Negative	1153	577
92	Hexa- <i>O</i> -galloyl hexoside	Negative	1091	635
93	Peonidin 3,5-di- <i>O</i> -glucoside	Positive	625	301
94	Ellagic acid tri- <i>O</i> -hexoside	Negative	787	463
95	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
96	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
97	Procyanidin B <sub>3</sub>	Negative	577	407
98	Procyanidin B <sub>1</sub>	Negative	577	407
99	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
100	Procyanidin B <sub>4</sub>	Negative	577	407
101	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
102	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
103	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
104	Epicatechin / catechin dimer <i>O</i> -hexoside	Negative	739	451
105	Procyanidin B <sub>2</sub>	Negative	577	407
106	Procyanidin B <sub>2</sub> 3- <i>O</i> -gallate	Negative	729	441
107	Procyanidin B <sub>2</sub> 3''- <i>O</i> -gallate	Negative	729	441
108	Epicatechin / catechin dimer <i>O</i> -gallate	Negative	729	441
109	Procyanidin B <sub>2</sub> 3,3''-di- <i>O</i> -gallate	Negative	881	559
110	1-O-Ellagoyl-gentiobiose (amritoside)	Negative	625	301
111	Epicatechin gallate	Negative	441	289
112	Ellagic acid di- <i>O</i> -desoxyhexoside	Negative	593	301
113	Ellagic acid 4- <i>O</i> -rhamnoside (eschweilenol C)	Negative	447	301
114	Epicatechin / catechin tetramer	Negative	1153	577
115	Epicatechin / catechin tetramer	Negative	1153	577
116	Epicatechin / catechin tetramer	Negative	1153	577
117	Epicatechin / catechin tetramer	Negative	1153	577
118	Epicatechin / catechin pentamer	Negative	1441	863
119	Epicatechin / catechin pentamer	Negative	1441	863
120	Epicatechin / catechin pentamer	Negative	1441	863
121	Epicatechin / catechin pentamer	Negative	1441	863
122	Epicatechin / catechin pentamer	Negative	1441	863
123	Epicatechin / catechin pentamer	Negative	1441	863

**Table S3.** Regression equations, correlation coefficients ( $r^2$ ), standard deviation ( $S_{YX}$ ), limits of detection (LOD), limits of quantification (LOQ) and linear ranges for 34 reference standards.

Compound	Ionization <sup>a</sup>	CE <sup>b</sup> (eV)	Regression Equation <sup>c</sup>		$r^2$	$S_{YX}$	LOD/ LOQ ( $\mu\text{g/mL}$ )	Linear Range ( $\mu\text{g/mL}$ )
			$a$	$b \cdot 10^6$				
1-O-Caffeoylquinic acid	N	-15	2.5394	-1.2360	0.9994	0.45·10 <sup>-2</sup>	0.006/0.02	0.02–300.0
1-O-Galloyl-β-D-glucopyranose	N	-20	1.3586	-0.0663	0.9987	9.69·10 <sup>-2</sup>	0.24/0.71	0.8–100.0
1,6-Di-O-galloyl glucose	N	-20	1.7552	-0.0569	0.9982	8.89·10 <sup>-2</sup>	0.18/0.51	0.6–100.0
1,3,6-Tri-O-galloyl glucose	N	-25	2.1064	-0.0499	0.9985	8.74·10 <sup>-2</sup>	0.14/0.42	0.5–100.0
1,2,3,6-Tetra-O-galloyl glucose	N	-25	2.0168	-0.0629	0.9973	9.00·10 <sup>-2</sup>	0.15/0.45	0.5–100.0
1,2,3,4,6-Penta-O-galloyl glucose	N	-25	2.4561	-0.0171	0.9979	12.33·10 <sup>-2</sup>	0.17/0.50	0.6–100.0
2-Pyrone-4,6-dicarboxylic acid	N	-20	1.8535	0.0742	0.9970	4.55·10 <sup>-2</sup>	0.08/0.25	0.3–100.0
5-O-Caffeoylquinic acid	N	-15	2.9021	-1.4184	0.9998	0.39·10 <sup>-2</sup>	0.004/0.01	0.02–300.0
Catechin	N	-20	0.9562	-0.0521	0.9971	7.79·10 <sup>-2</sup>	0.27/0.82	0.90–100.0
Corilagin	N	-35	0.9634	-0.3518	0.9804	7.34·10 <sup>-2</sup>	0.25/0.76	0.8–100.0
Gallic acid	N	-10	2.6538	-0.1376	0.9990	1.17·10 <sup>-2</sup>	0.01/0.04	0.10–100.0
Epicatechin	N	-20	1.0828	-0.0456	0.9973	6.85·10 <sup>-2</sup>	0.21/0.63	0.70–100.0
Epicatechin gallate	N	-20	1.5152	-0.0523	0.9979	12.67·10 <sup>-2</sup>	0.28/0.84	0.90–100.0
Epigallocatechin	N	-35	1.2824	-0.1440	0.9987	14.44·10 <sup>-2</sup>	0.37/1.12	1.2–100.0
Ellagic acid	N	-30	0.9114	-0.6312	0.9887	6.37·10 <sup>-2</sup>	0.23/0.70	0.7–100.0
Eschweilenol C	N	-30	1.2820	-0.9634	0.9697	11.64·10 <sup>-2</sup>	0.30/0.91	1.0–100.0
Helichryside	N	-30	2.0319	-0.3615	0.9811	10.09·10 <sup>-2</sup>	0.17/0.52	0.6–100.0
Juglanin	N	-20	2.0384	-0.3640	0.9975	2.02·10 <sup>-2</sup>	0.03/0.10	0.10–350.0
Kaempferol-3-O-β-D-glucuronopyranoside	N	-30	2.2126	-0.5160	0.9987	8.11·10 <sup>-2</sup>	0.12/0.37	0.4–100.0
Kaempferol-3-O-(6''-galloyl)-β-D-glucopyranoside	N	-10	1.2416	-0.3615	0.9901	3.02·10 <sup>-2</sup>	0.08/0.24	0.3–100.0
Miquelianin	N	-20	1.6705	-0.4374	0.9988	12.79·10 <sup>-2</sup>	0.25/0.77	0.8–100.0
Peonidin	P	+20	0.9634	-0.8634	0.9832	10.37·10 <sup>-2</sup>	0.36/1.07	1.1–100.0
3,5-di-O-β-D-glucopyranoside								
Procyanidin B <sub>1</sub>	N	-30	1.3722	-0.0829	0.9973	9.93·10 <sup>-2</sup>	0.24/0.72	0.80–100.0
Procyanidin B <sub>2</sub>	N	-30	1.3620	-0.0820	0.9961	9.91·10 <sup>-2</sup>	0.21/0.72	0.80–100.0
Procyanidin B <sub>2</sub> 3-O-gallate	N	-35	0.9361	-0.4518	0.9870	9.35·10 <sup>-2</sup>	0.32/1.00	1.0–100.0
Procyanidin B <sub>2</sub> 3''-O-gallate	N	-35	1.5963	-0.3518	0.9910	14.63·10 <sup>-2</sup>	0.30/0.92	1.0–100.0
Procyanidin B <sub>2</sub> 3,3''-di-O-gallate	N	-25	1.4632	-0.0524	0.9953	8.12·10 <sup>-2</sup>	0.18/0.56	0.6–100.0
Procyanidin B <sub>3</sub>	N	-20	1.6341	-0.4283	0.9900	15.02·10 <sup>-2</sup>	0.30/0.92	1.0–100.0
Procyanidin B <sub>4</sub>	N	-25	1.0634	-0.0933	0.9902	10.01·10 <sup>-2</sup>	0.31/0.94	1.0–100.0
Quercitrin	N	-20	1.9871	-0.6871	0.9984	5.63·10 <sup>-2</sup>	0.09/0.28	0.40–400.0
Quercetin-3-O-(6''-galloyl)-β-D-galactopyranoside	N	-25	1.1492	-0.6010	0.9980	4.68·10 <sup>-2</sup>	0.14/0.41	0.50–400.0
Rutin	N	-20	1.2716	-0.7389	0.9897	9.14·10 <sup>-2</sup>	0.23/0.72	0.80–400.0
Tiliroside	N	-30	2.3312	-0.4563	0.9803	14.92·10 <sup>-2</sup>	0.21/0.64	0.7–100.0
Taxifolin	N	-15	1.1105	-0.3211	0.9937	4.18·10 <sup>-2</sup>	0.12/0.38	0.4–100.0

<sup>a</sup> Ionization mode : N—negative ; P—positive. <sup>b</sup> CE—collision energy. <sup>c</sup> Regression equation:  $y = a \cdot x + b$ .



**Figure S1.** HPLC-DAD chromatograms of gastro-intestinal digestion markers (*p*-nitrophenol, 7-amino-4-methylcoumarin, 4-methylumbellifrone) after digestion of artificial substrate mixture. Experimental groups: reference standards (*p*-nitrophenol, 20 µg/mL; 7-amino-4-methylcoumarin, 10 µg/mL; 4-methylumbellifrone, 15 µg/mL); acarbose group; trypsin-chymotrypsin inhibitor group; orlistat group; *R. acicularis* leaf extract, 1000 µg/mL, group; rugosin D, 100 µg/mL, group; pure gastro-intestinal final juice without addition of plant extracts and inhibitors.