



Article

# Development of a Colloidal Gold Immunochromatographic Assay Strip Using a Monoclonal Antibody for the Rapid Detection of Ofloxacin

Xiaolan Li <sup>1</sup>, Jin Huang <sup>1</sup>, Na Li <sup>1</sup>, Mahmoud Salah <sup>2</sup>, Shuoning Guan <sup>1</sup>, Wenwen Pan <sup>1</sup>, Ziyi Wang <sup>1</sup>, Xinghua Zhou <sup>1</sup>,\* <sup>1</sup> and Yun Wang <sup>1</sup>,\*

- School of Food and Biological Engineering, Jiangsu University, Zhenjiang 212013, China; 2222318040@stmail.ujs.edu.cn (X.L.); 3220910011@stmail.ujs.edu.cn (J.H.)
- Department of Environmental Agricultural Science, Faculty of Graduate Studies and Environmental Research, Ain Shams University, Cairo 11566, Egypt
- \* Correspondence: zxh2012@ujs.edu.cn (X.Z.); wangy1974@ujs.edu.cn (Y.W.); Tel.: +86-150-0528-2057 (X.Z.); +86-139-5128-0586 (Y.W.)

**Abstract:** The livestock industry uses ofloxacin, an antibiotic, to prevent several animal diseases; however, the overdose of ofloxacin used in animal farming treatments may appear in food products and cause some adverse human health effects. Hence, there is an immediate need to develop a method suitable for on site large-scale detection of ofloxacin residues in animal-derived foods. This study aimed to prepare a monoclonal antibody with high sensitivity and affinity for ofloxacin by re-synthesizing the ofloxacin hapten and synthesizing the corresponding complete antigen. The  $IC_{50}$  of the enzyme-linked immunosorbent assay (ic-ELISA) was 0.13 ng/mL, and the detection limit was 0.033 ng/mL. The visual detection limit of the established colloidal gold immunochromatographic test strip, for the visual detection of actual samples, was 1 ng/g. In summary, this work establishes a rapid detection method of ofloxacin residues on the basis of colloidal gold immunochromatography that is suitable for actual detection.

Keywords: ofloxacin; monoclonal antibody; immunochromatographic strip; rapid detection



Citation: Li, X.; Huang, J.; Li, N.; Salah, M.; Guan, S.; Pan, W.; Wang, Z.; Zhou, X.; Wang, Y. Development of a Colloidal Gold Immunochromatographic Assay Strip Using a Monoclonal Antibody for the Rapid Detection of Ofloxacin. *Foods* 2024, 13, 4137. https://doi.org/ 10.3390/foods13244137

Academic Editor: Ronaldo Vagner Thomatieli-Santos

Received: 30 October 2024 Revised: 15 December 2024 Accepted: 19 December 2024 Published: 20 December 2024



Copyright: © 2024 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https://creativecommons.org/licenses/by/4.0/).

#### 1. Introduction

Ofloxacin (OFL, see Figure 1) is a quinolone antibacterial drug with broad-spectrum antimicrobial effects and is, therefore, often used to treat respiratory tract infections, genitourinary infections, and infections caused by sensitive bacteria [1]. It is an important antibiotic for domestic animals [2]. However, with its widespread use, the problems caused by its residues have attracted more and more attention. The long-term low intake of animal-derived food and water containing residues of OFL by humans will lead to the development of drug resistance in humans in addition to potential health hazards, such as toxic effects (central nervous system toxicity), peripheral neuropathy, allergic reactions, etc. [3]. In 2015, China issued the Ministry of Agriculture Announcement No. 2292, stating that food should not contain OFL due to the increasing safety risk of such residues. Despite this, cases of ofloxacin residues exceeding the sampling limit still occur. Hence, there is an immediate need to develop a method suitable for on site large-scale detection of ofloxacin residues [4].

The common methods for detecting ofloxacin residues include instrument detection and immunoassay. The instrument detection method includes HPLC [5–8], HPLC-MS/MS [9–12], and surface-enhanced Raman spectroscopy (SERS) [13]. Although the instrumental detection method is characterized by high sensitivity and accuracy, its high cost and complicated operational procedures render it unsuitable for on site detection. The high sensitivity and simple operation of immunoassay lead to its widespread use in detecting antibiotic residues. Huet et al. [14] established a dc-ELISA universal detection

method for 15 FQs antibiotics in which the estimated detection limit of OFL was 25  $\mu$ g/kg. Zhang et al. [15] developed an ic-ELISA method to detect 13 FQs drugs synchronously and the detection limit for OFL was set at 1.15  $\mu$ g/L. This study primarily employs colloidal gold immunochromatography, a method known for its simple operation, low detection cost, and suitability for rapid screening in the field [16]. Yuping Wang [17] developed a colloidal gold immunochromatographic assay to detect several fluoroquinolones. The detection limits (LODs) for OFX and MBF in milk were 3.5–8.9 ng/mL. Additionally, Byzova et al. [1] developed a competitive colloidal gold immunochromatographic assay for detecting OFL residues with a detection limit of 30 ng/mL.

$$\begin{array}{c|c}
 & O \\
 & N \\
 & N \\
 & O \\$$

**Figure 1.** Chemical structure of OFL.

In this work, monoclonal antibodies with high sensitivity and high affinity for ofloxacin were prepared. On this basis, the ic-ELISA method and the colloidal gold immune test strip method were established. The IC $_{50}$  of ic-ELISA was 0.13 ng/mL, and the detection limit was 0.033 ng/mL. The visual detection line of the colloidal gold test strip was 1 ng/g, which met the national requirements and was lower than that reported in the previous literature.

# 2. Materials and Methods

#### 2.1. Reagents and Equipment

Ovalbumin (OVA), bovine serum albumin (BSA), horseradish peroxidase-labeled goat anti-mouse IgG, 4,4′-Bi-2,6-xylidine, Carbodiimides (EDC), Freund's incomplete adjuvant, Freund's complete adjuvant, Citric acid trisodium salt, Chloroaquinic acid (HAuCl<sub>3</sub>·4H<sub>2</sub>O) were purchased from Sigma (St. Louis, MO, USA). Fetal bovine serum (FBS), trisodium citrate and thymidine medium, hypoxanthine–thymidine medium, and RPMI-1640 medium were purchased from GIBCO (Carlsbad, CA, USA). All other reagents were acquired from the National Pharmaceutical Group Chemical Reagent Co., Ltd. (Shanghai, China). The microplate reader (HF4500) was purchased from Huaan Maike Biotechnology Co., Ltd. (Beijing, China). The CT300 CNC strip cutting machine and HM3030 XYZ three-dimensional dispensing platform were obtained from Shanghai Kinbio Tech Co., Ltd. (Shanghai, China). Ethics and statement: animals (6- to 8-week-old female BALB/c mice) for the experiments were obtained from Kavans Laboratory Animal Co., Ltd., Changzhou, China (Certificate No. 202122406).

# 2.2. Modification of Ofloxacin Semi-Antigen

To improve the sensitivity of the assay, a new semi-antigen (Figure 2) was formed by resynthesizing the ofloxacin—which lacks a methyl group compared with the original ofloxacin—and using this to synthesize the corresponding envelope, thus improving the sensitivity of the assay with the advantage of the heterologous envelope [18]. The method of modifying the OFL hapten structure is as follows [19]:

A total of 172 mg of 1,4-diazacyclohexane was weighed and added to a solution of DMSO (3 mL) containing 281 mg of 9,10-difluoro-2,3-dihydro-3-methyl-7-oxo-7H-pyrido[1,2,3-DE]-1,4-benzoxazine-6-carboxylic acid. The suspension was stirred at 95 °C for 12 h. Then, 10 mL

of acetone was added to the mixture. The brown residue obtained was ground in acetone (10 mL) to give the expected compound as an off-white powder.

Figure 2. Modification of the hapten of OFL.

# 2.3. Synthesis of Complete Antigens

The specific steps for the preparation of complete antigen (as shown in Figure 3) by the activated ester method [20] are as follows: the immunogen was prepared by coupling the unmodified OFL with BSA, and the coating antigen was prepared by coupling the modified OFL with OVA. Amounts of 6.4 mg of semi-antigen, 4.1 mg of N-hydroxysuccinimide (NHS), and 6.8 mg of 1-ethylcarbodiimide hydrochloride (EDC) were weighed and dissolved in 0.6 mL of dimethylformamide (DMF) and stirred under room-temperature conditions overnight to achieve activation. The activated semi-antigen solution was dropped into a carbonate buffer (10 mL) containing 25.5 mg BSA or 20 mg OVA, stirred at 25 °C, and reacted overnight. Afterward, dialysis was performed in phosphate buffer at 4 °C for 72 h, and the dialysate was replaced every 8 h. After dialysis, the samples were stored at -20 °C for later use.

Figure 3. Preparation of immunogen and coating antigen ((a) OFL immunogen; (b) OFL coating antigen).

#### 2.4. Preparation of Monoclonal Antibodies

The initial immunization was performed with an immunogen dose of  $100~\mu g$  diluted in 0.9% NaCl solution and fully mixed with an equal volume of Freund's complete adjuvant. Then, one booster immunization was performed every two weeks at a dose of 50~ug, again diluted with 0.9% NaCl solution, and fully mixed with the equal volume of Freund's incomplete adjuvant. Starting from the third immunization, ic-ELISA analyzed the sera of each group of mice.

The titer and inhibition rate of the sera are generally the two main considerations in the screening process of mice, in which the titer represents the affinity of the prepared Foods **2024**, 13, 4137 4 of 16

antibody; the inhibition rate is the main indicator of the specific recognition of the antibody, and its calculation formula is as follows:

Inhibition rate = 
$$\left(1 - \frac{B}{B_0}\right) \times 100\%$$

where  $B_0$  is the absorbance value of the negative control group; and B is the absorbance value of the detection group.

For better cell fusion, the mouse spleen cells with the highest titer and inhibition rate were mixed with Sp 2/0 myeloma cells at a ratio of 1:7 and fused with 1 mL of PEG 1500 [21,22]. The fused cells were cultured with the HAT medium for 8 days and then replaced with the HT medium. Ic-ELISA screened positive wells.

The hybridoma cells with the highest titer and the best sensitivity were injected intraperitoneally into mice using a limited dilution method with 3 subcloning, and then the mouse ascites were collected after one to two weeks. The purification method of the ascites was ammonium sulfate precipitation, and the purified monoclonal antibody was stored at  $-20\,^{\circ}\text{C}$ .

# 2.5. Determination of Antibody Affinity, Sensitivity, and Specificity

The type of monoclonal antibody was detected by the kit and the affinity constant (Ka) was determined by ic-ELISA [23,24]. The formula for calculating Ka is as follows:

$$Ka = \frac{n-1}{2(n[Ab]t - [Ab']t)}$$

where n = [Ag]t/[Ag']t, wherein "[Ag]t" and "[Ag']t" are two different coating antigens concentrations, and "[Ab']t" and "[Ab']t" are the concentration of OFL-mAb at 50% OD<sub>450</sub> nm at these two coating antigens concentrations.

The sensitivity of the antibody is usually expressed by half inhibitory concentration (IC $_{50}$ ). The specificity was expressed as the cross-reactivity (CR), and the cross-reactivity was calculated according to the following formula:

$$CR(\%) = \left(\frac{IC_{50} \text{ of OFL}}{IC_{50} \text{ of analogs}}\right) \times 100\%$$

# 2.6. Selection of the ic-ELISA Working Conditions

In the process of ic-ELISA, there are many factors affecting it. Therefore, to obtain the best detection results, the choice of containment conditions (37 °C containment for 1 h, 37 °C containment for 2 h, 4 °C overnight containment, 4 °C containment for 1 h, 4 °C containment for 2 h), competition time (0.125 h, 0.25 h, 0.5 h, 1 h), standard dilution pH (5, 6, 7.4, 8.6, 9.6), and ionic strength (0%, 0.5%, 1%, 1.5%, 2%) were optimized. Then, standard curves were established under these different conditions. IC $_{50}$  and A $_{max}$  (absorbance at zero concentration at 450 nm) were used to evaluate the performance of mAb in ic-ELISA [25–30].

#### 2.7. Preparation of Colloidal Gold-Labeled Antibody

According to the previous research literature, a sodium citrate reduction method was used for the preparation of colloidal gold [31,32], described as follows: Take 100 mL of 1% aqueous chloroauric acid in a flask and heat it until boiling, then add 4 mL of 1% trisodium citrate while stirring. In this process, when the reaction color of the solution changes to burgundy, the heating is stopped, allowing the solution to cool down to 20–25  $^{\circ}$ C and maintaining it at 4  $^{\circ}$ C.


Afterward, the pH of the colloidal gold solution was adjusted to 8.5 with 0.1 M  $K_2CO_3$  solution; then, 6 uL of antibody was added to each milliliter of colloidal gold and stirred in the solution for 45 min. Then, 100 uL of 10% bovine serum albumin solution was added to each milliliter and stirring was continued for 2 h at 20–25 °C, followed by centrifugation at

Foods **2024**, 13, 4137 5 of 16

8000 rpm for 15 min. Then, the precipitate was washed the obtained. Soluble impurities were removed by centrifugation, and the washing steps were repeated twice. Finally, the precipitate was resuspended and stored at  $4\,^{\circ}\text{C}$  in a dark place.

# 2.8. Preparation of Colloidal Gold Immunochromatographic Assay Strips

The immunochromatographic strip mainly consists of a sample pad, a backing plate (PVC), a nitrocellulose membrane (NC membrane), and an absorbent pad, as shown in Figure 4A. The coating antigens formed the detection line (T line), originally immobilized on the NC membrane, and the quality control line (C line) was formed by the sheep anti-mouse antibody (secondary antibody), immobilized on the NC membrane.



**Figure 4.** The basic principle of the colloidal gold immunochromatographic assay: **(A)** structure and **(B)** principle.

# 2.9. Principle of the Colloidal Gold Immunochromatographic Assay Strips

Figure 4B demonstrates the fundamental mechanism of a colloidal gold immunochromatographic test strip. The test strip is inserted into a microwell and the sample, mixed with gold-labeled antibody, migrates from the sample pad toward the absorbent end. Initially, the gold-labeled antibody binds to the T line, while the unbound portion continues to attach to the secondary antibody at the C line. In the case of a negative sample, the test strip exhibits two red bands of comparable intensity, or the T line may appear slightly darker than the C line. Conversely, if the sample is positive, the gold-labeled antibody preferentially binds to the target substance in the sample, resulting in reduced binding to the T line, and the intensity of the T line is weakened. The point where the T line starts to brighten signifies the visual limit of detection (vLOD) of the test strip.

Foods **2024**, 13, 4137 6 of 16

## 2.10. Optimization of Working Conditions of Immunochromatographic Strip

The amount of antibody and the pH of the working environment during the preparation of gold-labeled antibodies are two key factors that will affect the levels of detection.

# 2.10.1. Optimization of the Gold-Labeled Antibody pH

A total of 1 mL colloidal gold solution was added to the centrifuge tube, and then, 0, 2, 4, 6, 8, 10, 12, 14, and 16  $\mu$ L 0.1 M  $K_2CO_3$  solution was added, respectively. After mixing, 10 uL 0.5 mg/mL monoclonal antibody was added. Then, after 30 min the color of the mixed solution was observed. The group with unchanged color and the least amount of  $K_2CO_3$  solution was selected as the best dosage.

# 2.10.2. Optimization of Gold-Labeled Antibody Concentration

Based on the optimal  $K_2CO_3$  dosage, antibodies of 4, 6, 8, 10, and 12 µg were added to the centrifuge tubes with 1 mL of colloidal gold solution. After standing at 25 °C for 5 min, the group that had an unchanged color and the smallest number of antibodies displayed the optimal antibody concentration.

# 2.10.3. Optimization of the Concentration of Coating Antigen in T Line

The coating antigens with concentrations of 0.2, 0.5, and 1 ng/mL were coated on the T line using an instrument, and the color reaction was carried out after drying. The group with the most stable color rendering and the lowest visual detection limit was used as the optimal coating concentration.

### 2.11. The Actual Sample Detection Experiment

The actual samples included pork, fish, and chicken. Amounts of 0, 1, 2, and 4 ng/g were used as the added concentration of the OFL standard for standard addition detection, and the visual detection line of the colloidal gold immunochromatographic test strip was judged through megascopic examination.

# 3. Results and Discussion

## 3.1. Identification of Ofloxacin Semi-Antigen

Figures 5 and 6 show the mass spectra and spectroscopy of the modified OFL semi-antigen. There is a molecular ion peak of 349.22 m/z, which corresponds to the molecular weight of the modified OFL (347.33); therefore, the semi-antigen was successfully modified.





Figure 5. Mass spectrometry of the modified OFL.

Foods **2024**, 13, 4137 7 of 16



**Figure 6.** The <sup>1</sup>H NMR spectrum for the modified OFL.

#### 3.2. Identification of Complete Antigen

#### 3.2.1. UV Characterization of OFL Complete Antigen

As shown in Figure 7a,b, OFL has a characteristic absorption peak between 275 and 285 nm, and the maximum absorption peaks of OFL-BSA and OFL-OVA are at 285 nm and 280 nm, respectively, which are shifted compared with those of the carrier proteins. This indicates that the coupling of the OFL immunogen and encapsulant was successful [20].




Figure 7. UV-vis spectrum of OFL immunogen (a) and coating antigen (b).

# 3.2.2. Determination of Complete Antigen Concentration

The measurements of immunization, antibody potency, and inhibition rate are based on the concentrations of the immunogen and encapsulant. Figure 8 shows the established protein standard curve; the concentration of the immunogen (OFL-BSA) is 5.23 mg/mL according to the standard curve, and the concentration of the coating antigen (OFL (engineered)-OVA) is 5.25 mg/mL.




Figure 8. Protein standard curve.

#### 3.3. Selection of Hybridoma Cells

Three BALB/c mice of the same batch with similar body weights were selected and immunized with the prepared immunogen five times, and the anti-serum was detected by ic-ELISA. The serum test results obtained by the last immunization are shown in Table 1. The coating antigen concentration of 0.3  $\mu g$  / mL and the anti-serum dilution of 3000 times were selected for serum detection. After detection, the effect of the heterologous coating antigen was better than that of the homologous coating antigen; therefore, the heterologous coating antigen was used for the following series of experiments. In three mice, the titer and inhibition rate of No. 3 mice were not as good as those of No. 1 and No. 2; therefore, No. 3 mice were not considered. In contrast, the inhibition of No. 1 and No. 2 mice was not much different but the titer of No. 1 was slightly higher than that of No. 2. Therefore, the No. 1 mice were selected for the subsequent experiments.

**Table 1.** Serum test results after the fifth immunization.

Immunogen Mouse Number		OFL-BSA	
	1	2	3
OFL concentration (ppb)	50	50	50
$Titer_1 (OD_{450} nm)$	$2.113 \pm 0.014$	$2.102 \pm 0.11$	$1.479 \pm 0.037$
Titer <sub>2</sub> ( $OD_{450}$ nm)	$2.201 \pm 0.020$	$1.985 \pm 0.087$	$1.525 \pm 0.033$
Inhibition ratio <sub>1</sub> (%)	$24 \pm 0.94$	$19 \pm 2.16$	$15 \pm 1.21$
Inhibition ratio <sub>2</sub> (%)	$58 \pm 3.77$	$50\pm1.06$	$42\pm1.95$

Note. Titer<sub>1</sub> ( $OD_{450}$  nm), Inhibition ratio<sub>1</sub> (%): using a homologous coating antigen. Titer<sub>2</sub> ( $OD_{450}$  nm), Inhibition ratio<sub>2</sub> (%): using a heterologous coating antigen.

After the fusion of Sp 2/0 cells and B cells, the detected positive cells need to be sub-cloned three times using a limited dilution method to screen out the hybridoma cell lines with high purity and meet the requirements. After subcloning screening, five OFL hybridoma cell lines were successfully obtained. The results are shown in Figure 9 (The detailed data are presented in Supplementary Materials). According to the comprehensive data, the titers of antibodies produced by cell line 4A5 and cell line 8B9 were not much different but the IC50 and  $A_{max}/IC50$  values of 4A5 were better than those of 8B9. Therefore, selecting cell line 4A5 prepared the OFL monoclonal antibody and the establishment of related detection methods.

Foods 2024, 13, 4137 9 of 16





Figure 9. The screening results of hybridoma cells of OFL.

# 3.4. Identification of mAb

The subtype of monoclonal antibody was measured by the kit C060101-L (from Luoyang Baiaotong Experimental Material Center, Luoyang, China). As shown in Figure 10, the heavy chain subtype of OFL antibody was IgG2a, and the light chain subtype was kappa. The affinity constants of mAb were determined by selecting three coating antigen concentrations of 0.1, 0.03, and 0.01 µg/mL. Figure 11 displays the result of the affinity constant of the OFL antibody. According to the curve regression equation, the concentration of monoclonal antibody corresponding to  $OD_{max/2}$  was obtained. According to the affinity constant formula,  $Ka_1 = 7.73 \times 10^9 \text{ L/mol}$ ,  $Ka_2 = 1.32 \times 10^{10} \text{ L/mol}$ , and  $Ka_3 = 8.63 \times 10^9 \text{ L/mol}$  were calculated. The affinity constant of the OFL antibody is  $Ka = (Ka_1 + Ka_2 + Ka_3) / 3 = 9.85 \times 10^9 L/mol$ . The results showed that the prepared OFL antibody had a high affinity.



**Figure 10.** Subtype determination of mAb.





Figure 11. Affinity constant result of OFL antibody.

# 3.5. Optimization Results of ic-ELISA Working Conditions

The optimal working point of the ic-ELISA method was obtained by the chessboard method, that is, the concentration of coating antigen was 0.03  $\mu g$  / mL, and the concentration of antibody was 0.03  $\mu g$ /mL.

Based on the optimal working point, the optimization of the OFL antibody ic-ELISA method is shown in Figure 12 and the specific data are shown in Table 2 as well as Table 3. The optimal containment conditions for OFL's ic-ELISA method were 2 h of containment at 37  $^{\circ}$ C, a best competitive reaction time of 0.5 h, an optimal NaCl content of 1.5% for the standard dilutions, as well as an optimal pH of 7.4.



**Figure 12.** The optimization results of ELISA working conditions: blocking conditions ( $\mathbf{a}$ ), competition time ( $\mathbf{b}$ ), NaCl content ( $\mathbf{c}$ ), and pH ( $\mathbf{d}$ ).

	Blocking Condition					Competition Time/h			
•	37 °C 1 h	37 °C 2 h	4 °C Overnight	4 °C 1 h	4 °C 2 h	0.125	0.25	0.5	1
A <sub>max</sub>	1.635	1.49	1.521	1.771	1.678	1.36	1.42	1.53	1.56
$IC_{50}$	0.34	0.22	0.34	0.38	0.36	0.39	0.56	0.25	0.27
$A_{\text{max}}/IC_{50}$	4.81	6.77	4.27	4.66	4.66	3.49	2.54	6.12	5.78

Table 2. Screening results of optimal blocking conditions and competition time of ic-ELISA.

**Table 3.** The screening results of the optimal NaCl content and pH in ic-ELISA.

	NaCl Content					рН				
	0%	0.5%	1%	1.5%	2%	5	6	7.4	8.6	9.6
A <sub>max</sub>	1.69	1.63	1.56	1.56	1.51	1.62	1.63	1.52	1.5	1.51
$IC_{50}$	0.23	0.41	0.21	0.14	0.27	0.42	0.25	0.13	0.46	0.52
$A_{\text{max}}/IC_{50}$	7.35	3.98	7.43	11.14	5.59	3.86	6.52	11.69	3.26	2.9

# 3.6. Establishment of Standard Curve

In this work, IC $_{20}$ ~IC $_{80}$  is defined as the linear range of detection, and IC $_{20}$  is defined as the minimum detection limit of the ic-ELISA. The ic-ELISA standard inhibition curve of OFL is shown in Figure 13. The corresponding equation is  $y = -0.0086 + 1.303 / \left[1 + \left(\frac{x}{0.13}\right)^{1.01}\right]$ ,  $R^2 = 0.994$ . The IC $_{50}$  is 0.13 ng/mL, the detection limit is 0.033 ng/mL, and the linear detection range is 0.033~0.51 ng/mL.




Figure 13. The standard inhibition curves of OFL.

# 3.7. Determination of Antibody Sensitivity and Specificity


Table 4 shows the cross-reactivities of the OFL monoclonal antibody with five other common structural analogs, among which the cross-reactivities with four structural analogs—lomefloxacin (LOM), pefloxacin (PEF), ciprofloxacin (CIP), and norfloxacin (NOR)—was low at less than 0.1% but the cross-reactivity with marbofloxacin (MBF) was higher, reaching 86.67%; therefore, the immunological assay established with this monoclonal antibody can make a preliminary judgment of OFL.

Chemical Compound	IC <sub>50</sub> (ng/mL)	CR (%)
OFL	0.13	100
MBF	0.15	86.67
LOM	>500	<0.1
PEF	>500	<0.1
CIP	>500	<0.1
NOR	>500	<0.1

**Table 4.** The cross-reactivity of OFL monoclonal antibody.

# 3.8. Identification of Colloidal Gold-Labeled mAb

As shown in Figure 14, the color of the prepared colloidal gold solution was burgundy, clarified, and translucent without impurities, while the UV scanning graph showed that the prepared colloidal gold nanoparticles had the maximum absorbance value at 522 nm and the peak width was narrow.



**Figure 14.** The characteristic images of colloidal gold particles: **(a)** UV–visible spectra, and **(b)** transmission electron microscopy images.

# 3.9. Optimization of Parameters Related to Test Strips

The color of the gold-labeled antibody will change due to different pH values, and the naked eye will observe the color change, the results are summarized in Table 5, which show that the optimal amount of  $K_2CO_3$  solution for synthesizing OFL gold-labeled antibodies is 12  $\mu$ L.

**Table 5.** The color of gold-labeled antibody solution at different pH values.

		Volume of 0.1 M $K_2CO_3$ Used for 1 mL of Colloidal Gold ( $\mu$ L)							
	0	4	6	8	12	14	16	20	
Solution color	+	+	-++	-++	+++	+++	+++	+++	

Note: "--+": blue-violet; "-++": purple; "+++": red.

The antibody concentration also affects the color of the gold-labeled antibody solution, and the naked eye can observe the color change; the results are summarized in Table 6, which show that the optimal concentration of antibody for preparing gold-labeled antibodies is  $8 \, \mu g/mL$ .

<b>Table 6.</b> The color	of gold-labeled antibod	y solution at different antibod	y concentrations.

	The Amount of Antibody Used in 1 mL Colloidal Gold (µg)							
	4	6	8	10	12	14	16	18
Solution color	+	-++	+++	+++	+++	+++	+++	+++

Note: "--+": blue-violet; "-++": purple; "+++": red.

Table 7 displays the optimization results of the concentrations of coating antigen and antibody of the OFL colloidal gold test strip. When the concentration of coating antigen is 0.5 mg/mL, the color of the T line gradually deepens with the increase in the amount of antibody. When 12  $\mu$ L of the gold-labeled antibody is added, the gray value (by ImageJ software (Version: 1.52)) of the T line reaches 1004, and the vLOD has reached stability at 1 ng/mL. Therefore, the best detection performance of the colloidal gold test strip was established with 0.5 mg/mL of coating concentration and  $12 \mu$ L of the gold-labeled antibody.

**Table 7.** The optimization for the concentrations of antigen and antibody on the colloidal gold test strip.

Antigen Concentration (mg/mL)		0.2			0.5			1	
Gold-labeled antibody (μL)	4	8	12	4	8	12	4	8	12
Negative T line gray value	243	478	668	507	899	1004	789	1328	1811
vLOD (ng/mL)	-	-	0.5	0.5	1	1	1	3	4

Note: "-" indicates that the T-line signal cannot be judged by the naked eye.


# 3.10. The Actual Samples Detected by the ic-ELISA Method and Colloidal Gold Immunochromatographic Test Strip

This study investigated three types of meat: pork, chicken, and fish. Table 8 provides the determination results of the recovery rate of OFL by ic-ELISA; these show that the recovery rate of OFL in pork samples was 91.0%~103.2%, the recovery rate of OFL in fish was 82.0%~99.6%, and the recovery rate of OFL in chicken was 86.0%~104.4%. The recovery rates for all three samples fell within the range of 80%~120% and the average coefficient of variation was below 5%. These results indicate that the colloidal gold test strip method is feasible for detecting OFL residues applied to the actual detection of these three types of meat products.

**Table 8.** Recovery rates of OFL (n = 3).

Samples	Added Standard Concentration (ng/g)	Recovery (%)	CV (%)
	0.00	-	-
pork	0.25	$91.00 \pm 1.00$	2.56
-	0.50	$103.20 \pm 6.40$	6.12
	0.00	-	-
fish	0.25	$82.00 \pm 5.00$	4.11
	0.50	$99.60 \pm 2.40$	3.12
	0.00	-	-
chicken	0.25	$86.00 \pm 3.00$	2.37
	0.50	$104.40 \pm 5.60$	3.34

The vLOD of the colloidal gold immunochromatographic test strip was determined. In samples with no OFL present, the T-line color would be similar to the C-line color. As the OFL concentration in the sample increased, the T-line color became lighter. The results of testing three meat samples using the colloidal gold immunochromatographic test strip method are shown in Figure 15. From the diagram, it can be observed that the vLOD of this method in the three actual samples is 1 ng/g.



**Figure 15.** The images of colloidal gold immunochromatographic strip tests for OFL in (a) pork, (b) fish, and (c) chicken samples. (Note: the additional amounts of OFL standard in these samples are 1: 0 ng/g, 2: 1 ng/g, 3: 2 ng/g, 4: 4 ng/g).

The detection limit and detection time of the method used in this study were compared with other methods. The results are shown in Table 9. The comparison results show that the method used in this study is sensitive, rapid and simple.

Reference	Method	LOD (ng/mL)	Analysis Time (h)	Detection Object
This work	Colloidal gold immunochromatographic assay	1 (ng/g)	0.25	pork, fish, and chicken
	ELISA	0.033	5	1 , ,
[5]	HPLC	20	1–2	influent, effluent, and
[5]	TH EC	20	1-2	surface waters
[6]	HPLC	3.6	1–2	human urine
[7]	HPLC-FLD	0.7	1–2	chicken meat
[8]	HPLC-MEPS-FLD	0.05  (ng/g)	1–2	milk
[0]	HPLC-MS/MS	0.005-0.123	1–2	surface water and
[9]	TIFLC-IVIS/ IVIS	0.003-0.123	1-2	hospital wastewater
[10]	HPLC-MS/MS	4.30  (ng/g)	1–2	milk
[11]	SPE- HPLC-MS/MS	0.0023-0.0107	1–2	environmental water
[12]	HPLC-MS/MS	0.1-3  (ng/g)	1–2	bean sprouts
[13]	SERS	42.6-49.1	1–2	aquatic samples
	ELISA	2F ( / )	5	kidney, marine products,
[14]	ELISA	25 (ng/g)	3	eggs, and muscle
[15]	ELISA	1.15	5	rana catesbeianus
[17]	immunochromatographic assay	3.5-8.9	0.25	milk

30

**Table 9.** Comparison of this method with other measurement methods.

## 4. Conclusions

immunochromatographic assay

[1]

In this study, a colloidal gold immunochromatographic test strip based on mAb was used to detect ofloxacin residues in meat products. The affinity of the prepared monoclonal antibody was  $9.85 \times 10^9$  L/mol, and the cross-reactivity with the other four analogs except MBF was less than 0.1%, with high affinity and good specificity. The IC<sub>50</sub> of the established ic-ELISA method was 0.13 ng/mL, the detection limit was 0.033 ng/mL, and the recovery rate was 80%~120%. The visual detection line of the colloidal gold

0.25

aquatic samples

immunochromatographic test strip in three meat products was 1 ng/g. The above results indicate that the developed immunological detection method can effectively detect the residues of ofloxacin.

**Supplementary Materials:** The following supporting information can be downloaded at https://www.mdpi.com/article/10.3390/foods13244137/s1, Table S1: Positive cell screening results of titer of No. 1 microplate; Table S2: Positive cell screening results of titer of No. 2 microplate; Table S3: Positive cell screening results of titer of No. 3 microplate; Table S4: Positive cell screening results of titer of No. 4 microplate; Table S5: Positive cell screening results of titer of No. 5 microplate; Table S6: Positive cell screening results of titer of No. 6 microplate; Table S7: Positive cell screening results of titer of No. 7 microplate; Table S8: Positive cell screening results of titer of No. 8 microplate; Table S9: Preliminary screening results of positive cells.

**Author Contributions:** Conceptualization, S.G.; Methodology, X.L.; Software, Z.W.; Validation, J.H.; Investigation, X.L. and J.H.; Resources, W.P.; Data curation, X.L. and N.L.; Writing—original draft, X.L.; Writing—review and editing, M.S.; Project administration, X.Z. and Y.W.; Funding acquisition, X.Z. and Y.W. All authors have read and agreed to the published version of the manuscript.

**Funding:** This research was supported by Xinghua Zhou grant number 8421360001.

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

**Data Availability Statement:** The original contributions presented in this study are included in the article/Supplementary Materials. Further inquiries can be directed to the corresponding author.

Conflicts of Interest: The authors declare no conflicts of interest.

#### References

 Byzova, N.A.; Smirnova, N.I.; Zherdev, A.V.; Eremin, S.A.; Shanin, I.A.; Lei, H.-T.; Sun, Y.; Dzantiev, B.B. Rapid immunochromatographic assay for ofloxacin in animal original foodstuffs using native antisera labeled by colloidal gold. *Talanta* 2014, 119, 125–132. [CrossRef]

- Chen, X.; Li, Z.; Guo, J.; Li, D.; Gao, H.; Wang, Y.; Xu, C. Simultaneous screening for marbofloxacin and ofloxacin residues in animal-derived foods using an indirect competitive immunoassay. *Food Agric. Immunol.* 2017, 28, 489–499. [CrossRef]
- 3. Zheng, R.Q.; Yang, J.; Zhou, L.; Liao, S.D.; Wang, J.; Lyu, S.M.; Tan, A.J. Ecological toxicity of antibiotics on the germination of Medicago sativa. *J. South. Agric.* **2021**, *52*, 2457–2464.
- 4. Jiang, C.; Wu, T.; He, X.; Wang, Y.; Lian, H.Z. Preparation of Thermo-Sensitive Molecular Imprinted SERS Substrate with Robust Recyclability for Detection of Ofloxacin. *Chemosensors* **2022**, *10*, 437. [CrossRef]
- 5. Tong, C.; Zhuo, X.; Guo, Y. Occurrence and risk assessment of four typical fluoroquinolone antibiotics in raw and treated sewage and in receiving waters in Hangzhou, China. *J. Agric. Food Chem.* **2011**, *59*, 7303–7309. [CrossRef]
- 6. Vakh, C.; Pochivalov, A.; Andruch, V.; Moskvin, L.; Bulatov, A. A fully automated effervescence-assisted switchable solvent-based liquid phase microextraction procedure: Liquid chromatographic determination of ofloxacin in human urine samples. *Anal. Chim. Acta* 2016, 907, 54–59. [CrossRef] [PubMed]
- 7. Timofeeva, I.; Timofeev, S.; Moskvin, L.; Bulatov, A. A dispersive liquid-liquid microextraction using a switchable polarity dispersive solvent. Automated HPLC-FLD determination of ofloxacin in chicken meat. *Anal. Chim. Acta* 2017, 949, 35–42. [CrossRef] [PubMed]
- 8. Zhang, X.; Wang, C.; Yang, L.; Zhang, W.; Li, C. Determination of eight quinolones in milk using immunoaffinity microextraction in a packed syringe and liquid chromatography with fluorescence detection. *J. Chromatogr. B Anal. Technol. Biomed. Life Sci.* 2017, 1064, 68–74. [CrossRef] [PubMed]
- 9. Van Hoi, B.; Vu, C.-T.; Phung-Thi, L.-A.; Nguyen, P.T.; Mai, H.; Le, P.-T.; Duong, D.T.; Thi, H.N.; Le-Van, D.; Chu, D.B. Determination of pharmaceutical residues by UPLC-MS/MS method: Validation and application on surface water and hospital wastewater. *J. Anal. Methods Chem.* **2021**, 2021, 6628285. [CrossRef]
- 10. Chang, L.; Du, S.; Wu, X.; Zhang, J.; Gan, Z. Analysis, occurrence and exposure evaluation of antibiotic and anthelmintic residues in whole cow milk from China. *Antibiotics* **2023**, *12*, 1125. [CrossRef] [PubMed]
- 11. Wang, J.; Ye, K.; Tian, Y.; Liu, K.; Liang, L.; Li, Q.; Huang, N.; Wang, X. Simultaneous determination of 22 antibiotics in environmental water samples by solid phase extraction-high performance liquid chromatography-tandem mass spectrometry. *Chin. J. Chromatogr.* 2023, 41, 241–249. [CrossRef] [PubMed]
- 12. Feng, J.; Jiang, H.; Wang, J.; Jing, Z.; Zhang, F.; Tan, T.; He, F.; Jiang, L.; Li, H.; Chang, S.; et al. Simultaneous determination of 40 plant growth regulators, fungicides, insecticides, and antibiotics in bean sprouts by QuEChERS-high performance liquid chromatography-tandem mass spectrometry. *Chin. J. Chromatogr.* 2022, 40, 843–853. [CrossRef] [PubMed]

13. Yu, Z.; Huang, L.; Zhang, Z.; Li, G. Simultaneous and accurate quantification of multiple antibiotics in aquatic samples by surface-enhanced Raman scattering using a Ti<sub>3</sub>C<sub>2</sub>T<sub>x</sub>/DNA/Ag membrane substrate. *Anal. Chem.* **2021**, *93*, 13072–13079. [CrossRef] [PubMed]

- 14. Huet, A.-C.; Charlier, C.; Tittlemier, S.A.; Singh, G.; Benrejeb, S.; Delahaut, P. Simultaneous determination of (fluoro)quinolone antibiotics in kidney, marine products, eggs, and muscle by enzyme-linked immunosorbent assay (ELISA). *J. Agric. Food Chem.* **2006**, *54*, 2822–2827. [CrossRef]
- 15. Zhang, B.; Lang, Y.; Guo, B.; Cao, Z.; Cheng, J.; Cai, D.; Shentu, X.; Yu, X. Indirect competitive enzyme-linked immunosorbent assay based on broad-spectrum antibody for simultaneous determination of thirteen fluoroquinolone antibiotics in Rana catesbeianus. *Foods* **2023**, *12*, 2530. [CrossRef] [PubMed]
- 16. Biao, Z.; Danfeng, C.; Tianze, Y.; Yihan, L.; Yanyong, F.; Xuping, S.; Xiaoping, Y. Construction and application of chloramphenicol enzyme-linked immunosorbent assay in aquatic products. *Chin. J. Food Sci.* **2023**, *23*, 352–360.
- 17. Wang, Y. Application of Multi-Target Fluorescence Immunochromatographic Assay for the Detection of Three Types of Antibiotics. Ph.D. Theis, China Jiliang University, Hangzhou, China, 2022.
- 18. Li, R. Synthesis of pyrrole-2,5-diacetamide by N-hydroxysuccinimide activated ester method at room temperature. *J. Jiangxi Norm. Univ. (Nat. Sci. Ed.)* **2005**, *6*, 509–511.
- 19. Li, L.; Wan, Z.; Wen, X. Synthesis and identification of artificial antigens of eight fluoroquinolones. *Agric. Technol.* **2016**, 36, 17–18+48.
- 20. Guo, J.; Liu, L.; Xue, F.; Xing, C.; Song, S.; Kuang, H.; Xu, C. Development of a monoclonal antibody-based immunochromatographic strip for cephalexin. *Food Agric. Immunol.* **2015**, *26*, 282–292. [CrossRef]
- 21. Li, J.; Zhi, A.; Jia, G.; Nie, H.; Ai, L. Development of an ic-ELISA and immunochromatographic strip for detection of sparfloxacin in honey. *Food Agric. Immunol.* **2018**, 29, 147–158. [CrossRef]
- 22. Liu, B.-H.; Hung, C.-T.; Lu, C.-C.; Chou, H.-N.; Yu, F.-Y. Production of Monoclonal Antibody for Okadaic Acid and Its Utilization in an Ultrasensitive Enzyme-Linked Immunosorbent Assay and One-Step Immunochromatographic Strip. *Agric. Food Chem.* **2014**, *62*, 1254–1260. [CrossRef] [PubMed]
- 23. Wang, B.; Xiao, K.; Meng, L.; Shou, C. Generation and characterization of THBS2-specific monoclonal antibodies. *Chin. Immunol.* **2020**, *36*, 588–593.
- 24. Zhou, X.; Li, N.; Sun, C.; Zhang, X.; Zhang, C.; Zhou, J.; Guan, S.; Xiao, X.; Wang, Y. Development of a colloidal gold immunochromatographic strip for the rapid detection of pefloxacin in grass carp with a novel pretreatment method. *J. Environ. Sci. Health* 2022, 57, 517–525. [CrossRef] [PubMed]
- 25. Wang, L.; Zhang, Y.; Gao, X.; Duan, Z.; Wang, S. Determination of Chloramphenicol Residues in Milk by Enzyme-Linked Immunosorbent Assay: Improvement by Biotin-Streptavidin-Amplified System. *J. Agric. Food Chem.* **2010**, *58*, 3265–3270. [CrossRef] [PubMed]
- 26. Jiang, X.X.; Shi, H.Y.; Wu, N.; Wang, M.H. Development of an enzyme-linked immunosorbent assay for diniconazole in agricultural samples. *Food Chem.* **2011**, *125*, 1385–1389. [CrossRef]
- 27. Chen, X.; Xu, L.; Ma, W.; Liu, L.; Kuang, H.; Peng, C.; Wang, L.; Xu, C. Development of an Enzyme-Linked Immunosorbent Assay for Cyhalothrin. *Immunol. Investig.* **2013**, 42, 493–503. [CrossRef] [PubMed]
- 28. Nemzek, J.A.; Siddiqui, J.; Remick, D.G. Development and optimization of cytokine ELISAs using commercial antibody pairs. *J. Immunol. Methods* **2001**, 255, 149–157. [CrossRef] [PubMed]
- 29. Székács, A.; Le, H.T.; Szurdoki, F.; Hammock, B.D. Optimization and validation of an enzyme immunoassay for the insect growth regulator fenoxycarb. *Anal. Chim. Acta* **2003**, 487, 15–29. [CrossRef]
- 30. Rezapkin, G.; Dragunsky, E.; Chumakov, K. Improved ELISA test for determination of potency of Inactivated Poliovirus Vaccine (IPV). *Biologicals* **2005**, *33*, 17–27. [CrossRef]
- 31. Yao, J.; Sun, Y.; Li, Q.; Wang, F.; Teng, M.; Yang, Y.; Deng, R.; Hu, X. Colloidal gold-McAb probe-based rapid immunoassay strip for simultaneous detection of fumonisins in maize. *J. Sci. Food Agric.* **2017**, *97*, 2223–2229. [CrossRef]
- 32. Wu, W.-D.; Li, M.; Chen, M.; Li, L.-P.; Wang, R.; Chen, H.-L.; Chen, F.-Y.; Mi, Q.; Liang, W.-W.; Chen, H.-Z. Development of a colloidal gold immunochromatographic strip for rapid detection of Streptococcus agalactiae in tilapia. *Biosens. Bioelectron.* **2017**, 91, 66–69.

**Disclaimer/Publisher's Note:** The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.