

Figure S1: Gating strategies for T cell population frequencies in blood samples from inoculated and control animals




Figure S1: Gating strategies for estimating frequencies of T cells, here shown for one sample as an example. CD3⁺ cells were stained with FITC, CD4⁺ cells with PE, and CD8⁺ cells with SPRD detected with the PerCP filter. R4 corresponds to the population of CD3⁺ cells (plot 7). In the plot to the right (plot 13), the lower right (LR) quadrant indicates the frequency of CD4⁺CD8[−] cells (T helper cells, 15.7%), the upper left (UL) quadrant indicates the frequency of CD4[−]CD8⁺ cells (cytotoxic T cells, 4.3%), and the upper right (UR) quadrant indicates the frequency of CD4⁺CD8⁺ (double-positive T cells, 1.5%).