



# A Detailed Protocol for Constructing a Human Single-Chain Variable Fragment (scFv) Library and Downstream Screening via Phage Display

Ziyi Liu<sup>1,2,†</sup>, Dokyun Kim<sup>1,2,†</sup>, Seokmin Kang<sup>1,2</sup> and Jae U. Jung<sup>1,2,\*</sup>

- <sup>1</sup> Cancer Biology Department, Infection Biology Program, Lerner Research Institute, Cleveland Clinic, Cleveland, OH 44195, USA; liuz5@ccf.org (Z.L.); kimd9@ccf.org (D.K.); kangs5@ccf.org (S.K.)
- <sup>2</sup> Global Center for Pathogen Research and Human Health, Lerner Research Institute, Cleveland Clinic, Cleveland, OH 44195, USA
- \* Correspondence: jungj@ccf.org; Tel.: 1-216-444-4694
- <sup>+</sup> These authors contributed equally to this work.

**Abstract:** The development of monoclonal antibodies (mAbs) represents a significant milestone in both basic research and clinical applications due to their target specificity and versatility in therapeutic and diagnostic applications. The innovative strategy of mAb screening, utilizing phage display, facilitates the in vitro screening of antibodies with high affinity to target antigens. The single-chain variable fragment (scFv) is a subset of mAb derivatives, known for its high binding affinity and smaller size—just one-third of that of human IgG. This report outlines a detailed and comprehensive procedure for constructing a scFv phagemid library derived from human patients, followed by screening via phage display affinity selection. The protocol utilizes 348 primer combinations spanning the entire human antibody repertoire to minimize sequence bias and maintain library diversity during polymerase chain reaction (PCR) for scFv generation, resulting in a library size greater than  $1 \times 10^8$ . Furthermore, we describe a high-throughput phage display screening protocol using enzyme-linked immunosorbent assay (ELISA) to evaluate more than 1200 scFv candidates. The generation of a highly diverse scFv library, coupled with the implementation of a phage display screening methodology, is expected to provide a valuable resource for researchers in pursuit of scFvs with high affinity for target antigens, thus advancing both research and clinical endeavors.

Keywords: single-chain variable fragment (scFv); phage display library; antibody screening

# 1. Introduction

Monoclonal antibodies (mAbs) are pivotal in immunotherapeutic strategies for a broad spectrum of diseases, including autoimmune disorders, chronic inflammatory diseases, cancer, and infectious diseases [1,2]. Their clinical efficacy is evidenced by the success of Trastuzumab (Herceptin) in treating HER2-positive breast cancer [3], Adalimumab (Humira) for chronic rheumatoid arthritis and ulcerative colitis [4], and the Tixagevimab/Cilgavimab cocktail (Evusheld) targeting the SARS-CoV-2 Spike protein to combat COVID-19 [5]. Currently, there are more than 140 FDA-approved or regulatory-reviewed mAb-based drugs in the United States alone [2,6].

Motivated by this clinical success, derivatives of mAbs have been developed to expand the utilization of mAb-based diagnostics and therapeutics for various diseases, while preserving their target specificity and high binding affinity [7]. The conventional antibody, or immunoglobulin (Ig), is a "Y"-shaped tetrameric protein, formed by two heavy-chain–lightchain heterodimers connected by disulfide bonds (Figure 1). Comprising three functional domains, it includes two identical antigen-binding fragments (Fabs) and one crystallizable fragment (Fc). The variable fragment (Fv) located at the tip of the variable regions of the heavy chain (V<sub>H</sub>) and light chain (V<sub>L</sub>) is responsible for mediating target specificity. The



Citation: Liu, Z.; Kim, D.; Kang, S.; Jung, J.U. A Detailed Protocol for Constructing a Human Single-Chain Variable Fragment (scFv) Library and Downstream Screening via Phage Display. *Methods Protoc.* **2024**, *7*, 13. https://doi.org/10.3390/mps7010013

Academic Editor: Philip Hublitz

Received: 4 December 2023 Revised: 23 January 2024 Accepted: 30 January 2024 Published: 1 February 2024



**Copyright:** © 2024 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). combination of various kappa ( $\kappa$ ) or lambda ( $\lambda$ ) light chains with the heavy chain confers specificity against a wide range of antigens with high binding affinity. Consequently, V<sub>H</sub> and V<sub>L</sub> have been rearranged into diverse formats, including Fab, single-chain Fv (scFv), and variable heavy domains of heavy chain (V<sub>H</sub>H) antibodies [8].



**Figure 1.** mAb and its derivatives. IgG is composed of two identical heavy chains (V<sub>H</sub>) and two light chains (V<sub>L</sub>), where the light chains are either kappa ( $\kappa$ ) or lambda ( $\lambda$ ) light chains. Antibody derivatives are notably smaller than IgG. The scFv includes only the variable regions, comprising the antigen-binding domain.

scFv, first developed by Bird et al. in 1988 [9], is a recombinant polypeptide where  $V_H$  is joined to  $V_L$  via a flexible linker composed of glycine and serine amino acids [10]. With a molecular weight of approximately 27 kDa, scFvs are much smaller than conventional whole Ig-based monoclonal antibodies (~150 kDa), enabling more effective tissue penetration and access to epitope crevices, while retaining high antigen-binding affinity [11]. Their rapid systemic clearance has led to widespread use in therapeutic and imaging applications in clinical settings. Furthermore, the lack of an Fc region diminishes the immunogenicity of scFvs, broadening their therapeutic application scope. Additionally, their cost-effective bacterial production enhances the accessibility of scFvs as antibody-based therapeutics [10]. Due to these benefits, scFvs are increasingly utilized in clinical scenarios where conventional Ig and Fab are suboptimal, such as in Chimeric Antigen Receptor (CAR) T-cell therapy [12].

The discovery of antibodies involves genetic manipulation of recombinant antibody fragments. Historically, antibody development was revolutionized by hybridoma technology, established by Köhler and Milstein in 1975, enabling large-scale production of specific antibodies [13]. Nevertheless, this technique has limitations. For instance, most mAbs produced via hybridoma technology are of murine origin and may induce human antimouse antibody responses in clinical settings [14]. To mitigate immunogenicity, antibodies from non-human species undergo a humanization process, either through chimerization or by grafting complementarity-determining regions (CDRs) on scaffolds from human Ig. However, these modifications can compromise the antibodies' binding affinity or stability, potentially reducing their therapeutic efficacy [15,16]. Consequently, screening antibodies from a human library is deemed the most effective approach for developing clinically applicable antibodies [17].

Among various antibody libraries, such as synthetic, semi-synthetic, naïve, immune libraries, and predefined CDR libraries tailored to specific requirements, the immune library offers a distinct advantage due to its repertoire specifically elicited against the target antigen [17–21]. Affinity maturation through somatic hypermutation and repeated exposure via immunization or natural infection enhances the specificity and efficacy of the antibodies. Antibodies sourced from immune libraries exhibit superior binding affinities compared to those from synthetic or naïve libraries [17,22]. The human immune library, boasting its inherent diversity and eliminating the need for humanization, is an ideal source for library construction. Nonetheless, constructing a human immune library presents

challenges, such as ethical considerations and the impracticality of deliberate immunization with the target antigen. Hence, the creation of a human immune scFv library typically relies on donated samples from individuals naturally exposed to the antigen or infected with the target pathogen [17,18,20,23].

The innovative development of phage display by George P. Smith has revolutionized the field, offering an efficient technique for generating antibodies against specific antigens [24]. Phage display screening, or affinity selection, is characterized by its rapidity, the ability to screen vast libraries, superior efficiency compared to conventional hybridoma methods, and cost-effectiveness [25–27]. This method allows for the display of diverse scFv repertoires on the surface of filamentous phages, facilitating their binding to the desired antigen. Through repeated selection rounds, an in vitro evolution of the displayed antibodies occurs, enriching for scFvs with high affinity for the antigen [28].

This protocol delineates a comprehensive workflow for generating a highly diverse human scFv library from multiple donors. It includes an extensive set of primers for amplifying  $V_H$  and  $V_L$  (both  $\kappa$  and  $\lambda$  light chains), ensuring maximal library diversity. Moreover, the protocol describes subsequent phage display affinity selection procedures aimed at identifying scFvs that specifically recognize the target antigen with high binding affinity, utilizing high-throughput monoclonal ELISA (Figure 2). To date, this protocol has effectively created libraries specific to pathogens and identified numerous binders specific to various antigens. The reliability of the results has been established across a wide array of protein types. In summary, this protocol provides a thorough and well-devised procedure for generating scFv libraries, as well as an extensive account of the associated phage display screening process. The technique specified has demonstrated reliability and repeatability, holding promise for diverse applications within molecular biology and therapeutic development.



**Figure 2.** Workflow of scFv library generation and phage display affinity selection. (**A**) Repertoire amplification: variable heavy ( $V_H$ ) and variable light ( $V_L$ ) regions (including both  $\kappa$  and  $\lambda$  chains) are

PCR-amplified from the cDNA of human peripheral blood mononuclear cells (PBMCs) using primers containing *Sfil* restriction sites (indicated by dots). During amplification, short or long overlapping linker sequences (highlighted in green) are introduced upstream of  $V_H$  and downstream of  $V_L$ . Overlap extension: a subsequent PCR promotes the annealing of  $V_H$  to  $V_L$  with either short or long linkers, resulting in the formation of scFv constructs. These constructs are digested with *Sfil* and ligated into the *Sfil*-linearized pComb3XSS vector, creating the scFv phagemid library for subsequent phage display affinity selection toward specific antigens. (**B**,**C**) Structural predictions of scFvs: AlphaFold-predicted structures of scFvs with a short linker ( $V_H$ –S– $V_L$ ) (**B**) and a long linker ( $V_H$ –L– $V_L$ ) (**C**). The  $V_H$ -S- $V_L$  is predisposed to form bispecific diabodies, whereas the  $V_H$ –L– $V_L$  likely maintains a structure akin to the natural antibody paratope. (**D**) Phage display affinity selection procedure: the scFv phagemid library is introduced into F-pilus-bearing bacteria and subsequently infected with helper phages, resulting in the production of phages that display scFvs. These phages are then subjected to successive rounds of screening to bind the target antigen, as illustrated in steps 1–6. Typically, three to five rounds of affinity selection are necessary to isolate scFvs that bind with high affinity to the target antigen.

### 2. Experimental Design

2.1. Materials

- 1. SuperScript<sup>TM</sup> III First-Strand Synthesis System (Invitrogen, Waltham, MA, USA; Cat. no.: 18080051).
- 2. Oligo(dT)<sub>12-18</sub> Primer (Invitrogen, Waltham, MA, USA; Cat. no.: 18418012).
- Platinum<sup>TM</sup> Hot Start PCR Master Mix (2×) (Invitrogen, Waltham, MA, USA; Cat. no.: 13000012).
- UltraPure<sup>TM</sup> DNase/RNase-Free Distilled Water (Invitrogen, Waltham, MA, USA; Cat. no.: 10977015).
- 5. Gel loading dye, purple,  $6 \times$  (New England Biolabs, Ipswich, MA, USA; Cat. no.: B7024S).
- 6. SeaKem LE agarose (Lonza, Basel, Switzerland; Cat. no.: 50004).
- 7. Ethyl alcohol, pure (Sigma-Aldrich, St. Louis, MO, USA; Cat. no.: E7023).
- 8. Sodium acetate, 3 M, pH 5.2, molecular biology grade (Millipore Sigma, Burlington, MA, USA; Cat. no.: 567422).
- 9. Zymoclean Gel DNA Recovery kit (capped columns) (Zymo Research, Irvine, CA, USA; Cat. no.: D4007/D4008).
- 10. pComb3XSS plasmid vector (Addgene, Watertown, MA, USA; Cat. no.: 63890).
- 11. Sfil restriction enzyme (New England Biolabs, Ipswich, MA, USA; Cat. no.: R0123S).
- 12. T4 DNA ligase (New England Biolabs, Ipswich, MA, USA; Cat. no.: M0202S).
- 13. T4 DNA ligase reaction buffer (New England Biolabs, Ipswich, MA, USA; Cat. no.: B0202S).
- 14. XL1-Blue Electroporation-competent cells (pUC18 plasmid included) (Agilent, Santa Clara, CA, USA; Cat. no.: 200228).
- 15. GlycoBlue Co-precipitant (15 mg/mL) (Invitrogen, Waltham, MA, USA; Cat. no.: AM9516).
- 16. Electroporation cuvettes, 0.2 cm (Invitrogen, Waltham, MA, USA; Cat. no.: P45050).
- 17. 15 mL Falcon conical centrifuge tube (Fisher scientific, Hampton, NH, USA; Cat. no.: 14-959-53A).
- 18. 50 mL Falcon high-clarity conical centrifuge tube (Fisher scientific, Hampton, NH, USA; Cat. no.: 14-432-22).
- 19. SOC medium (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 15544034).
- 20. GeneJET Plasmid miniprep kit (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: K0502).
- 21. Falcon 14 mL round-bottom test tubes with cap (Fisher scientific, Hampton, NH, USA; Cat. no.: 14-959-11B).

- 22. CM13 interference-resistant helper phage (Antibody design labs, San Diego, CA, USA; Cat. no.: PH020L).
- 23. Nunc MaxiSorp ELISA plate (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 442404).
- 24. Immulon 4HBX ELISA plate (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 3855).
- 25. Glucose (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 15023021).
- 26. Ampicillin sodium salt (Millipore Sigma, Burlington, MA, USA; Cat. no.: A9518).
- 27. Carbenicillin disodium salt (Millipore Sigma, Burlington, MA, USA; Cat. no.: C3416).
- 28. Tetracycline (Thermo Fisher scientific, Waltham, MA, USA; Cat. No.: J61714-14).
- 29. Kanamycin disulfate salt (Millipore Sigma, Burlington, MA, USA; Cat. no.: K1876).
- 30. LB broth base powder (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 12780052).
- 31. LB agar powder (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 2270025).
- 32. Sulfuric acid (H<sub>2</sub>SO<sub>4</sub>) (Millipore Sigma, Burlington, MA, USA; Cat. no.: 339741).
- 33. Hydrochloric acid (HCl) (Millipore Sigma, Burlington, MA, USA; Cat. no.: 320331).
- 34. Sodium chloride (NaCl) (Millipore Sigma, Burlington, MA, USA; Cat. no. 13423).
- 35. Yeast extract (Millipore Sigma, Burlington, MA, USA; Cat. no.: Y1625).
- 36. Bacto-tryptone (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 211705).
- A 0.2 μm vacuum filter system (Millipore Sigma, Burlington, MA, USA; Cat. no.: Z358193).
- 38. Polyethylene glycol (PEG)-6000 (Millipore Sigma, Burlington, MA, USA; Cat. no.: 807491).
- 39. Tween-20 (Millipore Sigma, Burlington, MA, USA; Cat. no.: P1379).
- 40. Phosphate-buffered saline (PBS) (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: 10010023).
- 41. AEBSF, hydrochloride (Millipore Sigma, Burlington, MA, USA; Cat. no.: 101500)
- 42. Trypsin (Millipore Sigma, Burlington, MA, USA; Cat. no.: T4799).
- 43. M13 phage coat protein monoclonal antibody (clone A5B3), HRP (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: MA5-36125).
- 44. OptEIA<sup>TM</sup> TMB substrate reagent set (BD Biosciences, Franklin Lakes, NJ, USA, Cat. no.: 555214).
- 45. Polycarbonate Erlenmeyer Flask (battled base), 250 mL (Thomas scientific, Chadds Ford Township, PA, USA; Cat. no.: FBC0250S).
- 46. TG1 phage competent cells (Antibody design labs, San Diego, CA, USA; Cat. no.: PC001).
- 47. Falcon 96-well, non-treated round-bottom plate (Fisher scientific, Hampton, NH, USA; Cat. no.: 08-772-54).
- 48. PBMCs isolated from patients' blood were stored in RNAlater solution at -20 °C.

#### 2.2. Equipment

- C1000 Touch Thermal Cycler with Dual 48/48 Fast Reaction Module (Bio-Rad, Hercules, CA, USA; Cat. no.: 1851148).
- Gene Pulser II Electroporation system (Bio-Rad, Hercules, CA, USA; Cat. no.: 165–2110).
- Vac-Man Laboratory vacuum manifold (Promega, Madison, WI, USA; Cat. no.: A7231).
- FilterMax F5 Multi-Mode Microplate reader (Molecular devices, San Jose, CA, USA; Cat. no.: F5).
- Heidolph Titramax 101 Vibrating platform shaker (Marshall scientific, Hampton, NH, USA; Cat. no.: H036130080).
- Thermo MaxQ 6000 Incubated / Refrigerated shaker (Thermo Fisher scientific, Waltham, MA, USA, Cat. no.: SHKE6000-8CE).
- Thermo MaxQ 8000 Incubated / Refrigerated shaker (Thermo Fisher scientific, Waltham, MA, USA; Cat. no.: SHKE8000-8CE).

# 3. Procedure

#### 3.1. First Strand cDNA Synthesis (Reverse Transcription, 1 Day)

This protocol begins with the reverse transcription of variable heavy (V<sub>H</sub>) and variable light (V<sub>L</sub>) chain genes from whole blood or peripheral blood mononuclear cells (PBMCs) from human donors. Researchers may choose any RNA extraction kit or protocol that is optimal for the sample type to achieve the desired diversity of the final scFv library. While the quantity of human donor samples used greatly depends on the availability of resources and the specific aims of the research, we extracted total RNA from PBMCs from 19 donors. The extracted total RNA quality needs to be assessed via gel electrophoresis before proceeding to reverse transcription to generate V<sub>H</sub> and V<sub>L</sub> cDNA. High-quality, intact RNA should exhibit distinct bands at 4.2 kb and 1.8 kb, corresponding to the 28S and 18S rRNA subunits, in a 2:1 intensity ratio, respectively.

- 1. Quantify the RNA concentration in the purified RNA. If the RNA concentration is low, the concentration of the RNA is recommended prior to the next steps. Aim for a final concentration of 40  $\mu$ g of total RNA in 128  $\mu$ L of ultrapure water. If the concentration exceeds this, dilute the RNA accordingly.
- 2. To the RNA, add 16  $\mu$ L of 50  $\mu$ M oligo(dT)<sub>12-18</sub> and 16  $\mu$ L of 10 mM dNTP mix to achieve a final volume of 160  $\mu$ L. Incubate at 65 °C for 5 min, then place on ice for at least 1 min. Briefly centrifuge to collect any liquid from the cap or side of the tube.
- 3. In a separate RNase-free microcentrifuge tube, prepare a reverse transcription mix by combining 32  $\mu$ L of 10× RT buffer, 64  $\mu$ L of 25 mM MgCl<sub>2</sub>, 32  $\mu$ L of 100 mM DTT, 16  $\mu$ L of 40 U/ $\mu$ L RNaseOUT, and 16  $\mu$ L of 200 U/ $\mu$ L SuperScript III reverse transcriptase for a total of 160  $\mu$ L.
- 4. Combine the reverse transcription mix with the RNA/Oligo/dNTP mixture to reach a total volume of 320 μL.
- 5. Incubate the mixture at 50 °C for 50 min, then at 85 °C for 5 min. Cool on ice for at least 1 min.
- Add 16 μL of 2 U/μL E. coli RNase H and incubate for 20 min at 37 °C. Centrifuge briefly and pool the resultant first-strand cDNAs from different tubes.

**D PAUSE STEP** The first-strand cDNA may be stored at  $-20 \degree$ C for several weeks. For extended storage (several months to years), add 0.1 volume (33.6 µL) of 3 M sodium acetate (pH 5.2) and 2.2 volume (739.2 µL) of EtOH, vortex to mix well, and store at  $-80 \degree$ C.

#### 3.2. scFv Construction (6–10 Days)

#### 3.2.1. V Gene Amplification (3–7 Days)

The initial stage in scFv generation involves the amplification of the variable region (Figure 2A). The human genome contains three loci associated with antibodies: one for the heavy chain and two for the light chains, kappa ( $\kappa$ ) and lambda ( $\lambda$ ). The heavy chain consists of the variability (V), diversity (D), joining (J), and constant (C) segments. In contrast, both  $\kappa$  and  $\lambda$  light chains comprise only V and J segments (Table 1). The V, D, and J segments are further divided into various families. To encompass the entire immune repertoire, we designed 348 primer combinations that target all V families (forward primers) and J families (reverse primers). These primers, based on Rader et al.'s work and subsequently modified in accordance with Barbas et al. and the V base database [29,30], have been tailored to include sequences for both short and long linkers. The sequences of these primers are detailed in Appendix A.

	Heavy Chain	Light Chain (V <sub>L</sub> )						
	(V <sub>H</sub> )	к Light	λ Light					
V segment	40	40	30					
D segment	25	-	-					
J segment	6	5	4					

	Table 1. Composition of V	, D, and	I segments in the human ger	nome [31]
--	---------------------------	----------	-----------------------------	-----------

The scFv molecules can be assembled in two configurations:  $V_H$ -linker- $V_L$  or  $V_L$ -linker- $V_H$ . This protocol employs the  $V_H$ -linker- $V_L$  orientation, which has been demonstrated to exhibit higher binding affinity for the target antigen due to a more exposed  $V_H$  CDR3 [8,32].  $V_H$  and  $V_L$  regions are connected by either a short glycine and serine-rich flexible linker (GGSSRSS) or a long linker (SSGGGGGGGGGGGGGGSSRSS) [30]. ScFvs with the short linker are predisposed to forming bispecific diabodies (Figure 2B), while those with the long linker are more likely to mimic the natural elbow structure of the paratope, as predicted by AlphaFold (Figure 2C) [32–34]. Empirically, we have determined that long linker constructs are amplified more efficiently when positioned downstream of the  $V_H$  rather than upstream of  $V_L$ . Consequently, we produce four distinct pools of PCR products:  $V_H$  with a short linker ( $V_H$ -S),  $V_H$  with a long linker ( $V_H$ -L), kappa light chain ( $V_{\kappa}$ ), and lambda light chain ( $V_{\lambda}$ ), each derived from the corresponding antibody fragment isolated from the reactions. To mitigate nonspecific amplification and primer dimer formation during PCR, we utilize a Touch Down PCR protocol in lieu of a conventional PCR approach. The resulting amplified V region genes should be approximately 400 bp in size.

1. Prepare 348 unique PCR master mixes, with each PCR primer combination, and conduct Touch Down PCR as follows (Table 2):

Reaction Mixture		Touch Down PCR Condition						
Platinum Hot Start Taq $2 \times$ Master Mix	12.5 μL	94 °C	2 min					
Autoclaved H <sub>2</sub> O	7 μL	94 °C	30 s <sup>β</sup>					
First-strand cDNA	0.5 μL	80 °C to 68 °C	30 s (–1 °C/cycle) $^{\beta}$					
V region forward primer (10 $\mu$ M) (V segment) $^{\alpha}$	2.5 μL	72 °C	1.5 min <sup>β</sup>					
V region reverse primer (10 μM) (I segment) <sup>α</sup>	2.5 μL	94 °C	30 s <sup>x</sup>					
0 - 0 0		50 °C 72 °C 72 °C	30 s <sup>x</sup> 1.5 min <sup>x</sup> 10 min					
		105 °C	Lid temp					

Table 2. Conditions for V gene amplification PCR.

<sup> $\alpha$ </sup> Total primer combinations: **V**<sub>H</sub>-**S**: 19 of V<sub>H</sub> Forward (F') primers × 6 of V<sub>H</sub> Reverse (R') primers carrying short linker = 114; **V**<sub>H</sub>-**L**: 19 of V<sub>H</sub> F' primers × 6 of V<sub>H</sub> R' primers carrying long linker = 114; **V**<sub> $\kappa$ </sub> (light chain): 12 of V<sub>K</sub> F' primers × 5 of V<sub> $\kappa$ </sub> R' primers with partial overlap with short linker = 60; **V**<sub> $\lambda$ </sub> (light chain): 20 of V<sub> $\lambda$ </sub> F' primers × 3 of V<sub> $\lambda$ </sub> R' primers with partial overlap with long linker = 60. <sup> $\beta$ </sup> Perform 12 cycles of PCR with conditions of 30 s at 94 °C, 1 °C reduction/cycle from 80 °C to 68 °C, followed by 1.5 min at 72 °C. <sup>×</sup> Perform 22 cycles of PCR with conditions of 30 s at 94 °C, 30 s at 50 °C, 1.5 min at 72 °C.

- 2. Verify the PCR amplification success on a 1.5% agarose gel. Combine 3  $\mu$ L of the PCR product with 2  $\mu$ L of 6× DNA gel loading dye and perform gel electrophoresis at 135 V for 20 min.
- 3. Proceed with ethanol precipitation for PCR products that show successful amplification.

**PAUSE STEP** The PCR product can be stored at 4 °C for up to 7 days if not proceeding immediately.

3.2.2. Ethanol (EtOH) Precipitation of Pooled PCR Reactions (1 Day)

EtOH precipitation is employed to purify PCR products by removing enzymes and salts. Note that dNTPs and primers might co-precipitate, which could interfere with the Abs<sub>260nm</sub> and subsequent DNA quantification by UV spectrophotometry. To minimize this, replace sodium acetate with ammonium acetate in the protocol below.

1. Combine the PCR products for each V region into a separate 15 mL Falcon tube, ensuring you prepare four tubes for the four V regions. Mix thoroughly, then aliquot 330  $\mu$ L into individual 1.5 mL microcentrifuge tubes.

**CRITICAL STEP** It is essential to ensure a consistent mix in the 15 mL Falcon tubes before aliquoting, to maintain diversity and homogeneity, as some DNA might be lost during the EtOH precipitation process.

- 2. Add 0.1 volume (33  $\mu$ L) of 3 M sodium acetate to each tube, followed by 3 volumes (1089  $\mu$ L) of 100% EtOH. Incubate at -20 °C overnight to precipitate the DNA.
- 3. Centrifuge at  $16,000 \times g$  for 10 min at room temperature (RT) and aspirate most of the EtOH without disturbing the pellet.
- 4. Wash the DNA pellet with 1 mL of ice-cold 70% EtOH (in ultrapure water) and centrifuge at  $16,000 \times$  g for 2 min at RT. Carefully aspirate the ethanol without disturbing the pellet. Repeat this wash.

**CRITICAL STEP** Ice-cold 70% EtOH is recommended to prevent DNA loss, as DNA may dissolve in the aqueous portion of RT 70% EtOH.

- 5. After the second wash, centrifuge at  $16,000 \times g$  for 2 min at RT and remove all the EtOH.
- 6. Air-dry the pellet for 15 min to an hour with the tube lid open to ensure it is completely dry.
- 7. Resuspend the DNA pellet from each tube in 30 μL ultrapure water (preheated to 55 °C), and pool the resuspensions into four 4 tubes (one for each V fragment).
- Run 25 μL of each resuspended V fragment DNA on a 1.5% agarose gel at 100 V for 40 min, using multiple lanes as needed to fit all samples.

**A CRITICAL STEP** Employ low voltage for gel electrophoresis to prevent heatinduced DNA damage.

9. Extract DNA from the gel using a Zymoclean Gel DNA Recovery kit and elute in water heated to 55 °C to maximize yield. Determine DNA using a spectrophotometer.

**CRITICAL STEP** This step is crucial to diminish primer interference from the first PCR in subsequent amplifications.

**D** PAUSE STEP Store the purified DNA at -20 °C for long-term storage.

3.2.3. Secondary Overlapping Extension (SOE) PCR (2 Days)

Following the amplification of V regions,  $V_H$ -S and  $V_H$ -L amplicons are designed to contain sequences that overlap with  $V_{\kappa}$  and  $V_{\lambda}$  for fusion. After fusion, these constructs are further amplified via SOE PCR, resulting in four different scFv constructs:  $V_H$ -S-V $\kappa$ ,  $V_H$ -S-V $_{\lambda}$ ,  $V_H$ -L-V $_{\kappa}$ , and  $V_H$ -L-V $_{\lambda}$ . A minimum of 15 µg of each construct is necessary for subsequent experiments. Typically, conducting 12 SOE PCR reactions for each scFv combination (totaling 48 reactions) yields a sufficient DNA concentration for the next step.

1. Use 50 ng of each  $V_H$  and  $V_L$ , along with 5  $\mu$ L of primer mix, for the reactions. We recommend performing 12 SOE reactions for each construct to ensure an adequate quantity of DNA is produced. Execute the PCR according to the protocol outlined in Table 3.

Reaction M	ixture	Touch Down PCR Condition				
Platinum Hot Start Taq 2× Master Mix	12.5 μL	94 °C	2 min			
V <sub>H</sub> -short or V <sub>H</sub> -long	50 ng	94 °C	30 s <sup>β</sup>			
$V_{\kappa}$ or $V_{\lambda}$	50 ng	56 °C	30 s <sup>β</sup>			
SOE primer—F $^{\alpha}$	2.5 μL	72 °C	30 s <sup>β</sup> , <b>20 cycles</b>			
		12 °C X	~			
SOE primer—R $^{\alpha}$	2.5 μL	94 °C	30 s <sup>δ</sup>			
Autoclaved H <sub>2</sub> O	To 25 μL	68 °C	30 s <sup>δ</sup>			
		72 °C	1 min $\delta$			
		72 °C	10 min			
		105 °C	Lid temp			

Table 3. Conditions for SOE PCR.

 $^{\alpha}$  We recommend adding the primers after 20 PCR cycles to ensure adequate annealing of V<sub>H</sub>-V<sub>L</sub> annealing.  $^{\beta}$  Perform the initial 20 cycles at 94 °C for 30 s, and at 56 °C for 30 s to allow the overlapping sequences to anneal.  $^{\times}$  At this step, remove the tubes from the thermocycler and add 2.5 µL of forward primer (F') and 2.5 µL of reverse primer (R'). Subsequently, return the tubes to the thermocycler and proceed with the 'next step.'  $^{\delta}$  Perform an additional 15 cycles at 94 °C for 30 s, at 68 °C for 30 s, and at 72 °C for 1 min.

**CRITICAL STEP** Performing 20 cycles prior to an infinite hold at 12 °C is crucial. Fewer than 20 cycles may result in indistinct scFv bands, while more than 20 cycles can impair polymerase activity.

- 2. Assess the SOE PCR success by electrophoresing 3  $\mu$ L from each of the 48 tubes on a 1.5% agarose gel.
- 3. Pool the remaining SOE PCR products into four separate microcentrifuge tubes.
- Repeat DNA EtOH precipitation as described in Section 3.2.2, "EtOH precipitation of pooled PCR reactions." Combine approximately 360 μL of the reaction mixture from 12 PCR reactions per construct, adding 0.1 volume (30 μL) of 3 M sodium acetate, followed by 3 volumes (1080 μL) of 100% EtOH.
- 5. Precipitate DNA overnight at -20 °C.
- 6. The following day, wash the precipitated DNA as outlined in Section 3.2.2, "EtOH precipitation of pooled PCR reactions."
- 7. Utilize the required amount for digestion as detailed in Section 3.3, "Cloning scFv constructs to pComb3XSS by *Sfil* digestion and T4 DNA ligase reaction," or store the remainder of the DNA.

**PAUSE STEP** The precipitated DNA can be stored at -20 °C for several weeks.

#### 3.3. SfiI Digestion of pComb3XSS and scFv Constructs (2 Days)

The amplified scFv constructs include *Sfil* digestion sites for sticky end formation compatible with the pComb3XSS vector for cloning. 500 ng of each digested scFv construct is needed for library generation.

- Digest 50 μg of the pComb3XSS vector with 10 μL of *SfiI* enzyme (4 units/μg of DNA), bringing the total volume to 250 μL with water (50 μL for every 10 μg of DNA).
- 2. Digest 15 μg of the scFv construct precipitated from the prior step with 3 μL of *Sfil* enzyme, adjusting the final volume to 150 μL with water (50 μL for every 5 μg of DNA).
- 3. Incubate the digestion reactions at 37 °C overnight.
- 4. **OPTIONAL** Despite the manufacturer's recommendations for a 60 °C reaction temperature, 37 °C may be used for stable and prolonged digestion.
- 5. Confirm complete digestion by comparing the digested fragments with the undigested vector on a 1.5% agarose gel.
- 6. Use a Zymoclean Gel DNA Recovery kit to gel-extract the fully digested pComb3Xss (3379 bp), the stuffer fragment (1673 bp—removed from the pComb3XSS), and scFv constructs (700–800 bp). Elute with water preheated to 55 °C to maximize yield.

# **O PAUSE STEP** Store the eluted DNA at -20 °C for several weeks.

### 3.4. Test Library Generation (2 Days)

To ensure the generation of a library with a size greater than  $1 \times 10^7$ , a test library should be constructed to determine the necessary number of ligation reactions. T4 DNA ligase is used to ligate the digested scFv into the pComb3XSS vector. While alternative methods like Golden Gate assembly may serve the same purpose, our protocol specifically details the use of the traditional T4 DNA ligase reaction for its broad applicability. The test library is produced by direct electroporation of the ligation product into XL1-Blue electrocompetent cells. While this protocol specifically utilizes XL1-Blue, any bacterial strain that harbors the F-factor, such as TG1 or XL1-Blue, is suitable for generating libraries and subsequent affinity selection due to their ability to support filamentous phage infection.

3.4.1. Single T4 DNA Ligase Reaction of scFv into pComb3xSS (Adapted from [14], 1 Day)

- 1. Prepare four ligation reactions as follows, which include three controls and the scFv ligation reactions for each construct:
  - Control #1: without ligase;
  - Control #2: without insert;
  - Control #3: with stuffer insert;
  - Four separate scFv ligations for each scFv construct
  - The composition and volumes for these reactions are detailed in the table provided (Table 4).

Reaction Mixture	Control 1 (No Ligase)	Control 2 (No Insert)	Control 3 (with Stuffer Insert)	scFv	T4 DNA Lig Con	gase Reaction dition
T4 10× ligation buffer	2 µL	2 μL	2 µL	2 μL	16 °C	16 h
Digested scFv	-	-	-	70 ng	65 °C	10 min
Stuffer fragment	-	-	70 ng	-	4 °C	Hold
Digested pComb3XSS	140 ng	140 ng	140 ng	140 ng	105 °C	Lid temp
T4 ligase	-	1 μL (400 U)	1 μL (400 U)	1 μL (400 U)		
Autoclaved H <sub>2</sub> O	Το 20 μL	Το 20 μL	To 20 μL	Το 20 μL		

#### Table 4. Conditions for the T4 DNA ligase reaction.

- 2. Carry out the T4 DNA ligase reaction according to the protocol specified in the subsequent table (Table 4).
- 3. On the following day, combine 1  $\mu$ L from each of the four scFv construct ligation reactions (totaling 20  $\mu$ L) into a new microcentrifuge tube, resulting in a total volume of 4  $\mu$ L. This mixture will be electroporated into 50  $\mu$ L of XL1-Blue competent cells.
- 4. **OPTIONAL STEP** For quality control of individual scFv constructs, heat-shock TOP10 competent cells with 1 μL of each ligation product (four scFvs separately) and perform a dilution series inoculation on LB agar plates supplemented with ampicillin to assess diversity size.

3.4.2. Test Library Electroporation (1 Day)

- 1. In a new 1.5 mL microcentrifuge tube, combine 1 µL from each of the four scFv ligation products for electroporation into 50 µL of XL1-Blue competent cells.
- 2. Pre-warm SOC media to 37  $^\circ C$  and prepare 2 mL of pre-warmed SOC media in 15 mL Falcon tubes.
- 3. Thaw 250 μL of XL1-Blue competent cells, allocating 50 μL for each of the five electroporation samples. These include three negative controls (no ligase, no insert control,

and stuffer insert) and the combined scFv sample (4  $\mu L$  total), along with a pUC18 electroporation control.

- 4. Aliquot 50  $\mu$ L of competent cells into each pre-cooled microcentrifuge tube: mix 50  $\mu$ L of competent cells with 4  $\mu$ L of the ligation product for each control and scFv, along with 1  $\mu$ L of 0.1 ng/ $\mu$ L pUC18 plasmid as a positive control. Mix gently by pipetting and transfer into 0.2 cm gap electroporation cuvettes, avoiding bubble formation.
- 5. Set the Gene Pulser II apparatus to 25 µF capacitance, 2.5 kV, and 200 Ohm.
- 6. Tap the cuvette gently to remove bubbles and ensure all cells settle at the bottom, free from bubbles.

**A CRITICAL STEP** The pulse should result in 12.5 kV/cm with a time constant of 4–5 milliseconds.

7. Immediately add 950 μL of pre-warmed SOC to the electroporated cuvette, and transfer the mixture to the prepared 15 mL Falcon tubes containing 2 mL of SOC, making a total volume of 3 mL. Incubate at 37 °C with shaking at 225 rpm for 1 h. Meanwhile, pre-warm four LB agar plates supplemented with ampicillin to 37 °C.

**OPTIONAL** Given carbenicillin's enhanced stability compared to ampicillin, it can be used as a preferable alternative at working concentrations ranging from 50 to 100  $\mu$ g/mL.

- 8. After 1 h of recovery, take 10  $\mu$ L from the transformed cells and from the pUC18transformed cells to estimate the library size. Perform tenfold serial dilutions up to  $10^{-8}$  in a 96-well plate. Drop 5  $\mu$ L from each dilution onto a pre-warmed LB agar plate with ampicillin using a multichannel pipette or a 96-pin pronger. Plate the three ligation controls (no ligase, no insert, and stuffer insert) at a 1:100 dilution on separate LB agar plates.
- 9. Incubate all plates at 37 °C overnight (approximately 16 h).
- 10. The following day, calculate the library size from the test library using the formula:

# *Estimated library size from single reaction =*

Number of colonies on plate × dilution factor ×  $\frac{volume \text{ of ligation product } (=20 \text{ } \mu\text{L})}{volume \text{ of electroporation product } (=4 \text{ } \mu\text{L})}$ 

3.5. Direct (Large-Scale) Library Generation (3 Days)

3.5.1. Large-Scale Ligation (1 Day)

Calculate the number of ligation reactions needed to achieve a library size of  $10^7$  or  $10^8$ , based on test library results. For direct ligation, it is necessary to precipitate the ligated DNA to minimize salt interference, which can decrease electroporation efficiency and cell viability. For example, if a single ligation yields a library size of  $5 \times 10^6$ , then four ligation reactions are needed to surpass  $1 \times 10^7$ . The following steps apply to a single reaction; adjust quantities accordingly for the desired library size. After EtOH precipitation, combine all scFv constructs separately into four tubes. We recommend performing an equal number of ligation reactions for each scFv construct to ensure a balanced representation of each in the pooled ligation product.

- 1. Conduct multiple ligation reactions to reach the target library  $> 10^7$  size.
- Follow the T4 DNA ligase reaction protocol as outlined in Section 3.4.1, "Single T4 DNA ligase reaction of scFv into pComb3XSS."

3.5.2. EtOH Precipitation of the Ligated scFv Phagemid Library (1 Day)

1. Pool the ligation reactions up to a maximum of 330  $\mu$ L per tube. Add 0.1 volume of 3 M sodium acetate and 3 volumes of EtOH, as described in Section 3.2.2, "EtOH precipitation of pooled PCR reactions."

**OPTIONAL** For smaller numbers of ligation reactions or low DNA concentrations, add a co-precipitant like GlycoBlue or Glycogen at a 1:300 dilution during EtOH precipitation to enhance precipitation efficiency.

2. The next day, centrifuge the precipitated DNA and wash the pellet following the steps in Section 3.2.2, "EtOH precipitation of pooled PCR reactions." Resuspend the pellet in a total of 200  $\mu$ L of ultrapure water warmed to 55 °C.

#### 3.5.3. Electroporation into XL1-Blue Competent Cells (1 Day)

Electroporate 250 ng of the precipitated ligation product into 50  $\mu$ L of XL1-Blue competent cells. Using an excessive amount of DNA relative to the number of cells might result in the introduction of multiple DNA molecules into a single cell. Note that spectrophotometer readings may be inaccurate due to salt. The theoretical yield of precipitated DNA is calculated as follows: if 20 ligation reactions are precipitated, with each digested vector (3379 bp) yielding 171 ng of ligated scFv-carrying phagemid, the total yield will be 3420 ng in 200  $\mu$ L water (17.1 ng/ $\mu$ L). Therefore, 14.6  $\mu$ L of the precipitated DNA (equivalent to 250 ng) is needed for each 50  $\mu$ L aliquot of competent cells. In the actual library generation, the three negative ligation controls are unnecessary if using the same batch of digested vectors as in the test library.

- 1. Pre-cool the electroporation adaptor, a 15 mL Falcon tube, a microcentrifuge tube, and electroporation cuvettes (0.2 cm) on ice.
- Pre-warm 50 mL Falcon tubes containing 3 mL of SOC medium and 5 mL of SOC medium to 37 °C.
- 3. Prepare the required volume of 50 μL competent cell mixed with 14.6 μL DNA (based on calculations to reach the desired library size), mixing slowly for thorough integration.
- 4. Prepare a control mixture with pUC18 and competent cells in a separate microcentrifuge tube.
- 5. Aliquot the scFv master mix into the electroporation cuvettes, ensuring the absence of air bubbles. Wipe the cuvettes' metal surfaces to remove moisture.
- 6. Repeat the electroporation process as in Section 3.4.2, "Test library electroporation."
- Immediately add 935 μL of warmed SOC media to each cuvette and transfer the bacteria to the pre-warmed 50 mL Falcon tubes. Ensure each tube does not exceed 10 mL of total volume for adequate aeration.

**A CRITICAL STEP** The volume of bacteria in the Falcon tube should not exceed 20% of the tube's volume to aerate bacteria for growth. Thus, add only up to 7 electroporated bacteria (each 1 mL) to each 50 mL Falcon tube containing 3 mL of SOC media.

- 8. Shake the electroporated bacteria at 37 °C and 225 rpm for 1 h for recovery.
- 9. For library size calculation, dilute 100  $\mu$ L from the Falcon tube for serial dilutions to  $10^{-8}$ , then plate 5  $\mu$ L from each dilution. Store the remaining dilutions at 4 °C.
- 10. Spin down the rest of the bacteria and resuspend in a total of 4 mL. Incubate 1 mL onto each of several 245 mm square LB agar plates (with 2% w/v glucose and  $50 \mu g/mL$  ampicillin).
- 11. Incubate the plates overnight at 37  $^{\circ}$ C.
- 12. Scrape the bacteria from the plates using 4 mL of LB media, rinsing with an additional 2 mL of LB. Be gentle to avoid scraping the agar.
- 13. Add glycerol to a final concentration of 20%. Dilute 100  $\mu$ L of the scraped bacteria in 900  $\mu$ L LB and measure the OD at 600 nm (OD<sub>600</sub>). An OD<sub>600</sub> corresponds to 1  $\times$  10<sup>8</sup> bacteria per mL.

**CRITICAL STEP** Dilute the bacteria before measuring OD<sub>600</sub> to avoid exceeding the spectrophotometer's detection limit.

14. Aliquot the bacteria to achieve  $5 \times 10^9$  bacteria per cryovial for storage at -80 °C.

**PAUSE STEP** The scFv bacteria library can be stored at -80 °C for years.

15. **OPTIONAL STEP** Adjust the number of bacteria per vial based on library diversity. A higher number of bacteria based on library diversity can be used so that each library

stock contains 48 or more copies per scFv repertoire [35]. For example, if the diversity size is  $8 \times 10^8$ ,  $4 \times 10^{10}$  bacteria can be aliquoted to each cryovial so that each stock contains 50 copies of each scFv in the library. However, avoid exceeding  $1 \times 10^{12}$  bacteria per vial to prevent spontaneous phage precipitation.

#### 3.6. Library Quality Assessment (2–3 Days)

Quality control (QC) of the final library can include methods such as colony PCR, Sanger sequencing, *SfiI* digestion, or *BstNI* scFv fingerprinting. We recommend Sanger sequencing for identifying scFv sequences and assessing library diversity.

- 1. Thaw a vial of the library stock and prepare a serial dilution series to  $10^{-8}$ .
- 2. Inoculate 100  $\mu$ L of each dilution onto LB agar plates containing 2% w/v glucose, 50  $\mu$ g/mL ampicillin.

**OPTIONAL** When the cultivating library is XL1-Blue cells, it is recommended to add tetracycline to the culture medium at a concentration of  $50 \ \mu g/mL$ .

- 3. Incubate the plates at 37 °C overnight to allow colony growth.
- 4. Pick individual colonies and inoculate them into 3 mL of LB media supplemented with ampicillin, using Falcon 14 mL round-bottom tubes.
- 5. Incubate these cultures at 37 °C overnight to promote further colony growth.
- 6. Perform a miniprep using the GeneJET Miniprep kit to extract bacterial plasmid.

**OPTIONAL** For higher throughput, consider using a vacuum manifold during the miniprep process.

7. Sequence the extracted DNA using Sanger sequencing with the sequencing primers provided in Appendix A.

# 3.7. Phage Display Affinity Selection (Variable Time Depending on Enrichment, Usually 14–20 Days)

3.7.1. Optimization of Antigen Coating Condition on Different ELISA Plates (1-2 Days)

Optimizing ELISA plate coating concentration is crucial to conserve protein and maximize affinity selection efficiency. Different proteins adhere differently to various types of ELISA plates. For instance, glycoproteins often coat effectively on Nunc MaxiSorp or Immulon 4HBX plates, while lipoproteins may prefer other types [36]. Consult the manufacturer's catalog to choose the optimal ELISA plate for coating your target protein. Additionally, consider the antigen capacity of the plates, but note that this varies with each protein. When selecting antigens for negative selection, consider relevant factors; for instance, use Fc-tagged target proteins, including Fc for negative selection, or BSA if no similar protein is available.

- 1. Select appropriate ELISA plates for optimal coating of the target protein.
- Set various antigen coating concentrations from 1–5 μg/mL based on the manufacturer's guidelines. Target and control antigen proteins should be in similar forms.
- 3. Coat the ELISA plates, 100  $\mu$ L per well, with triplicates of both target and control antigens in PBS (or an appropriate buffer for the antigens) and incubate overnight at 4 °C with gentle shaking.
- 4. Wash the plates three times with PBS-T (0.05% Tween-20 in PBS).
- 5. Block the plates using 200  $\mu$ L of 5% milk in PBS-T, incubating for 1 h at RT with gentle shaking.
- 6. Wash the plates three times with PBS-T.
- Add 100 μL of HRP-conjugated antibody specific to the target or tag of the target protein and control antigens. Follow the manufacturer's guidelines and incubate for 1 h at RT with gentle shaking.

**OPTIONAL** In the absence of an HRP-conjugated antibody specific to the target and control antigens, an indirect ELISA approach can be employed. This involves using a

primary antibody that specifically recognizes the target and control antigens, followed by detection with an HRP-conjugated secondary antibody.

- 8. Wash the plates three times with PBS-T.
- 9. Develop the reaction with 100  $\mu$ L of TMB substrate for 5 to 20 min, then stop the reaction with 50  $\mu$ L of 1 M H<sub>2</sub>SO<sub>4</sub> (or 1 N HCl). Measure the signal intensity at 450 nm.
- 10. Plot a titration curve of OD<sub>450nm</sub> against different antigen coating concentrations to find the optimal concentration for affinity selection.

3.7.2. Induction of Phages Carrying scFv and In Vitro Selection (Multiple Rounds of 3 Days)

This phage display protocol is adapted and modified from Pardon et al. and Lee et al. [35,37], starting with the scFv library stock constructed in Section 3.5.3, "Electroporation into XL1-Blue competent cells."

- 1. Inoculate the library stock into 60 mL of 2xTY-G supplied with 50  $\mu$ g/mL ampicillin (and 50  $\mu$ g/mL tetracycline for XL1-Blue cells). Rinse the tube to collect all bacteria.
- 2. Reamplify the bacteria library culture at 37  $^{\circ}$ C and 200 rpm until OD<sub>600</sub> reaches 0.5.
- 3. Transfer 10 mL of the exponential-phase bacteria to a 50 mL Falcon tube and infect with  $4 \times 10^{10}$  cfu of VCSM13 helper phage. Incubate for 30 min at 37 °C without shaking.

**CRITICAL STEP** Use filtered pipette tips and disinfect the tips before disposal to prevent phage contamination. Used pipette tips should be discarded in a biohazard bag for either dispensing or autoclaving later.

- 4. Centrifuge the cells at  $4000 \times g$  for 10 min at RT. Discard the supernatant, and briefly place the tube on a paper towel.
- 5. Resuspend in 60 mL of 2xTY-I and transfer to a disposable 250 mL baffled flask. Incubate overnight at 37 °C and 200 rpm to induce the phages carrying scFv.
- 6. The next day, start growing naïve TG1 (not carrying the scFv gene) in 10 mL LB media in a 50 mL Falcon tube. Grow until OD<sub>600</sub> reaches 0.5 to 0.6, then keep in ice.

**A CRITICAL STEP** TG1 cells should be first cultured in M9 medium to ensure the presence of F-factor.

- 7. Centrifuge the induced 50 mL culture at  $4000 \times g$  for 30 min at 4 °C.
- 8. Transfer 40 mL of supernatant to a new 50 mL Falcon tube, then add 10 mL of PEG solution. Mix well by inverting the tube 5 times.
- 9. Incubate the mixture in ice for an hour to precipitate the phages.
- 10. Spin down at  $4000 \times g$  for 30 min at 4 °C. Discard the supernatant and place the tube upside down on a paper towel briefly to remove residual supernatant.
- 11. Resuspend the precipitated phage pellet using 1 mL of ice-cold PBS, then transfer the resuspension to a new ice-cold microcentrifuge tube.
- 12. Centrifuge at  $20,000 \times g$  for 1 min at 4 °C to spin down residual bacteria.
- 13. Carefully transfer only the supernatant to a new ice-cold microcentrifuge tube.
- 14. Add 250 µL of ice-cold PEG solution and pipet up and down several times to mix well.
- 15. Incubate on ice for 10 min.
- 16. Spin down the phages at  $14,000 \times g$  for 15 min at 4 °C. Discard the supernatant and resuspend the pellet using 1 mL of ice-cold PBS.
- 17. Centrifuge at 14,000  $\times$  *g* for 1 min at 4 °C to completely remove residual bacteria.
- 18. Transfer the supernatant to a new ice-cold microcentrifuge tube.
- 19. To a round bottom 96-well plate, add 90  $\mu$ L of PBS to a row of 12 wells. Use this PBS to make 12 tenfold serial dilutions of the precipitated phage, beginning with 10  $\mu$ L of precipitated phage. Serial dilution can be alternatively performed using a multichannel pipettor if preferred. In a separate row, add 90  $\mu$ L of naïve TG1 (at OD<sub>600</sub> = 0.5~0.6) and 10  $\mu$ L from each serial dilution to achieve a final volume of 100  $\mu$ L. Incubate the plate at 37 °C for 30 min without shaking.

20. Carefully pipette 5  $\mu$ L drops of the infected TG1 onto LB agar plates containing 2% glucose and ampicillin. As a control, also pipette 5  $\mu$ L of naïve TG1 (uninfected by phages). Allow the drops to air-dry and incubate the plates at 37 °C overnight. Concurrently, inoculate a naïve TG1 into 5 mL of LB media (without antibiotics) in a 50 mL Falcon tube for overnight growth.

**A CRITICAL STEP** When culturing the naïve TG1 cells, also inoculate TG1 into two tubes containing LB with ampicillin or kanamycin to check for potential contamination by either the phage, helper phage, or both. TG1 contaminated with the induced phage will grow in ampicillin-supplemented media, and TG1 contaminated with helper phage will grow in kanamycin-supplemented media.

- 21. For antigen coating on ELISA plates, use the previously optimized concentration (as per Section 3.7.1, "Optimization of antigen coating condition on different ELISA plates." Coat four wells with the control antigen and four non-adjacent wells with the target antigen, using 100  $\mu$ L for each well. Seal the plate with impermeable film and incubate overnight at 4 °C.
- 22. The next day, calculate the colony-forming units (cfu) of the precipitated phage using the following formula:

$$\frac{cfu}{mL} = number of colonies on plate \times dilution factor \times \frac{total \ volume \ (= 1000 \ \mu L)}{plated \ volume \ (= 5 \ \mu L)}$$

**CRITICAL STEP** The cfu/mL should be higher than  $1 \times 10^{12}$  cfu/mL, which would be at least a hundred times greater than the initial bacterial count. If it is lower, repeat the phage production process.

**PAUSE STEP** Store the precipitated phages at 4 °C for weeks. For long-term storage, add glycerol to achieve a final concentration of approximately 40%. The phage in glycerol can then be stored at -80 °C for several months to years.

- 23. Inoculate naïve TG1 into 10 mL of LB media in a 50 mL Falcon tube. Grow at 37 °C and 220 rpm until it reaches  $OD_{600} = 0.5$ . Then, move the Falcon tube into an ice bucket until it is ready for infection with the eluted phages.
- 24. Discard the antigen-coating solution from the ELISA plates and wash three times with PBS-T. Block the plates by adding 200  $\mu$ L of 5% milk in PBS-T to each well, and incubate at RT for 1 h with gentle shaking.
- 25. During the 1 h blocking, pre-block the phages by adding  $5 \times 10^{11}$  cfu of the precipitated phages to 500 µL of 5% milk in PBS-T ( $1 \times 10^{11}$  per well).

**CRITICAL STEP** The actual required phage volume is  $4 \times 10^{11}$  cfu in 400 µL, but we recommend preparing  $5 \times 10^{11}$  cfu in 500 µL to compensate for potential transfer loss. The pre-blocking should be performed for at least 20 min or until the 1 h ELISA plate blocking is complete, with gentle shaking at RT.

26. After blocking the ELISA plates, wash them three times with PBS-T. Transfer 100 μL of the pre-blocked phages to the wells coated with control antigen for negative selection.

**CRITICAL STEP** To prevent the target-antigen-coated wells from drying out, add 100  $\mu$ L of the blocking solution to the wells coated with the target antigen.

- 27. Agitate the plates gently at 700 rpm on a shaking platform for 2 h.
- 28. After incubation, remove the blocking solution from wells that have been coated with the target antigen, then transfer the phages from the wells coated with the control antigen to those coated with the target antigen.

**A CRITICAL STEP** To prevent the controlled-antigen-coated wells from drying out, add 100 µL of the blocking solution to the wells initially coated with the control antigen.

29. Agitate again at 700 rpm for 2 h.

30. After agitation, remove the liquid from all eight wells (four coated with control antigen and four with target antigen). Wash all wells 15 times with 200  $\mu$ L of PBS-T, followed by a single wash with 200  $\mu$ L of PBS. Add 50  $\mu$ L of trypsin (0.25 mg/mL in PBS) to all eight wells and agitate at 700 rpm for 30 min at RT. The addition of trypsin facilitates the dissociation of phages from the target antigen by cleaving the digestion sites between the phage coat protein and the displayed scFv. Phages eluted from both control and target antigens should not be discarded, as they are used in the following steps for enrichment analysis.

**A CRITICAL STEP** Liquid waste containing phage particles must be treated with a disinfectant before being disposed of at a liquid handling station equipped for regular disinfection.

- 31. Prepare two microcentrifuge tubes, each containing 10 μL of AEBSF. Transfer the trypsin-eluted phages from the four wells initially coated with control antigen into the first tube, and those from the target antigen wells into the second tube.
- 32. In a new 50 mL Falcon tube, add 3 mL of exponential-phase naïve TG1 (from step 23). Infect this culture with 50 μL of the eluted phage from the wells initially coated with target antigen. Incubate the mixture at 37 °C without shaking for 30 min. Then, add 7 mL of LB supplemented with ampicillin and glucose to achieve a final concentration of 50 μg/mL ampicillin and 2% glucose in 10 mL. Incubate the 10 mL mixture overnight at 37 °C and 180 rpm.
- 33. Perform tenfold serial dilutions of trypsin with AEBSF mixtures from control and target Ag, referring to Section 3.7.3, "Enrichment analysis." If the enrichment factor is less than 1000 (from the first method) but still increasing (from the second method), proceed to the next round of affinity selection.
- 34. The next day, transfer 500 μL of the overnight-grown 10 mL culture into 50 mL of 2xTY-G (2% glucose, ampicillin) in a baffled disposable flask and incubate at 37 °C, 200 rpm.

**D PAUSE STEP** Aliquot the remaining 10 mL of culture into cryovials, dispensing 1 mL into each vial. Then, add 666  $\mu$ L of 50% glycerol to each vial to make the final 20% glycerol, resulting in a total volume of 1.666 mL per cryovial. These can be stored at -80 °C for future single colony picking or for restarting the affinity selection process in case of error.

#### 3.7.3. Enrichment Analysis (1 Day)

Typically, the affinity selection step in phage display requires 3 to 5 rounds, depending on the behavior of the target antigens. However, excessive rounds of affinity selection increase the chance of clonal focusing; therefore, enrichment analysis is recommended to decide when to conclude the selection process. Most antigens require only up to three rounds. Enrichment analysis, a crucial quality control step, assesses the proportion of antigenspecific populations of scFv after each affinity selection round. This protocol describes two methods for evaluating the enrichment of target antigen-specific scFv-carrying phages based on their titers against control and target antigens. Method 1 compares enrichment factors by comparing phage titers against target Ag and control Ag, and Method 2 calculates the ratio of output phages to input phages after each round. Both methods can be used to quantify affinity selection efficiency and the development of target-specific scFvs.

# Method 1: Enrichment factor calculation

- 1. Prewarm the LB agar plate (2% w/v glucose and 50 µg/mL ampicillin) to 37 °C.
- 2. In a 96-well plate, add 90  $\mu$ L of PBS to two nonadjacent rows of 12 wells each, making 24 wells in total. Use this to prepare 12 tenfold serial dilutions, starting with 10  $\mu$ L of the trypsin + AEBSF mixture for both control and target antigens.
- 3. In two separate rows of 12 wells each, add 90 μL of exponential-phase naïve TG1 (from Section 3.7.2, "Induction of phages carrying scFv and in vitro selection," step 23).

- 4. Add 10  $\mu$ L of each serial dilution to the 90  $\mu$ L of naïve TG1.
- 5. Incubate the plate at 37 °C for 30 min without shaking.
- 6. Using a multichannel pipet, inoculate 5  $\mu$ L from each well of phage-infected TG1 onto multiple LB agar plates (2% glucose and ampicillin). Use three plates each for the control and target antigens to calculate the average values, considering the sticky nature of phages. Include a control spot of 5  $\mu$ L of naïve TG1 not infected with phages as a sentinel against contamination by phage.
- 7. Incubate the plates overnight at 37 °C.
- 8. The next day, count the average number of colonies from control and target antigens. Calculate the enrichment factor as follows:

 $Enrichment \ factor = \frac{colonies \ observed \ from \ target \ Ag}{colonies \ observed \ from \ control \ Ag}$ 

9. An enrichment factor exceeding 1000 is indicative of successful enrichment.

# Method 2: Output/input calculation

- 1. Determine cfu from the LB agar plates with drops of TG1 infected with phages eluted from wells coated with the target antigen. For example, if there are three average colonies at the 2nd dilution point (this is  $10^{-3}$  dilution, making the colony number  $3 \times 10^{3}$ ), cfu is  $\frac{3 \times 10^{3}}{0.005} = 6 \times 10^{5}$ /mL
- 3 × 10<sup>3</sup>), cfu is  $\frac{3 \times 10^3}{0.005} = 6 \times 10^5 / \text{mL}$ 2. Given that 4 × 10<sup>11</sup> cfu in 400 µL was used for input (equivalent to 1 × 10<sup>12</sup>/<sup>mL</sup>), calculate the  $\frac{output}{input}$  ratio as  $\frac{6 \times 10^5}{1 \times 10^{12}}$ .
- 3. This ratio should steadily increase, then plateau or decrease.

**CRITICAL STEP** The optimal round to conclude affinity selection and proceed to monoclonal ELISA is the last round before the plateau or decrease in both methods.

3.7.4. Polyclonal ELISA (2 Days)

- 1. Coat the ELISA plate with control and target antigens at the previously optimized concentration, using 100 μL per well for triplicates. Adjust based on the number of time points being tested in the polyclonal ELISA.
- 2. Incubate the plate overnight at 4 °C.
- 3. Wash the plate three times with PBS-T, then block using 200 µL of 5% milk in PBS-T per well for 1 h at RT with gentle shaking.
- 4. During the blocking period, add  $1 \times 10^{11}$  cfu of phages (based on cfu calculation from Enrichment Analysis Method 2) before affinity selection and after each selection round to 5% milk in PBS-T, preparing 100 µL for each time point.
- 5. After blocking, wash the plates three times with PBS-T and transfer 100  $\mu$ L of phage in milk to each well.
- 6. Incubate for 1 h at RT with gentle shaking.
- Wash the plate three times with PBS-T, add 100 μL of HRP-conjugated anti-M13 coat protein antibody according to the manufacturer's guidelines, and incubate for 1 h at RT with gentle shaking.
- 8. Wash the plate three times with PBS-T, then add 100  $\mu$ L of TMB substrate.
- 9. Develop the color reaction, then add 50  $\mu$ L of 1 M H<sub>2</sub>SO<sub>4</sub> to stop it. Quantify the optical density at 450 nm.
- 10. Draw a titration curve for  $OD_{450}$ . If the  $OD_{450}$  value plateaus or decreases, affinity selection should be halted at that round. The bacterial stock from the round immediately preceding the plateau or decrease can be used for single colony picking.

# 3.8. Screening for Hits by Monoclonal ELISA (4 Days)

1. Review the results of enrichment analyses (output/input calculation, enrichment factor analysis, and polyclonal ELISA). A plateau or decrease in output/input calculations and polyclonal ELISA, along with an enrichment factor over 1000, indicate it is

time to stop affinity selection. For instance, if these criteria are met by the product of Round 3, proceed to pick single colonies from Round 3.

- 2. Thaw a vial from the  $-80 \,^{\circ}\text{C}$  bacteria stock and prepare a serial dilution series from  $10^{-1}$  to  $10^{-8}$ . Inoculate each dilution onto LB agar plates (containing 2% glucose and 50 µg/mL ampicillin) to determine the optimal dilution for picking single colonies. Avoid too many or too few colonies on the plates. Incubate the plates overnight at 37 °C.
- 3. The next day, visually assess the plates to determine which dilution yields a manageable number of colonies for picking.
- 4. Plate multiple (10 to 15) LB agar plates (2% glucose and 50 μg/mL ampicillin). Ensure even spreading and drying of the plates to prevent clumping of colonies.
- 5. Prepare 6 to 12 round-bottom 96-well plates, adding 200 μL of 2xTY-G (2% glucose, 50 μg/mL ampicillin) to each well.
- 6. Using a sterile P200 pipette tip, pick single colonies and inoculate them into the prepared wells.

**OPTIONAL** A1 wells of each plate can be left without any colony being inoculated as a negative control.

7. Secure the round-bottom plates in a plastic box affixed to a shaking platform inside an incubator. Remove the covers from the plates, seal the box with cling film, and incubate the bacteria overnight at 37  $^{\circ}$ C, 250 rpm.

**CRITICAL STEP** To prevent cross-contamination, keep the plates covered while securing them inside the box, and only remove the 96-well plate covers just before sealing the box with cling film.

8. The following day, prepare a new set of 96-well round-bottom plates with 200  $\mu$ L of 2xTY-G. Transfer 5  $\mu$ L from the overnight cultures into the new plates using a multichannel pipette.

**PAUSE STEP** The original 96-well plate with overnight-grown colonies is the 'master plate.' Remove 45  $\mu$ L from the master plate to prevent expansion of the media during freezing, add 100  $\mu$ L of 50% glycerol to achieve a 20% final glycerol concentration, and store at -80 °C for future use.

- 9. After 3 h of growth, add 50  $\mu$ L of 2xTY-G containing 4  $\times$  10<sup>8</sup> cfu of helper phages to each well, gently mixing to ensure even distribution of the phages with the bacteria.
- 10. Incubate at 37 °C for 1 h without shaking to allow for helper phage infection. Then, centrifuge the plates at  $4000 \times g$  for 10 min and discard the supernatant carefully to prevent cross-contamination between wells, by a single vigorous downward motion. Blot the plates on a paper towel and resuspend the bacterial pellets in 200 µL of 2xTY-I, ensuring thorough mixing.
- 11. Return the plates to the incubator and induce phage production overnight at 25  $^{\circ}\mathrm{C}$  and 250 rpm.
- 12. Prior to leaving for the day, coat additional ELISA plates with control and target antigens using the previously optimized conditions. Prepare an equal number of plates for each set of monoclonal phages picked.

**OPTIONAL** To reduce the amount of target antigen used, half-volume microtiter plates can be utilized, which allows for the coating volume to be decreased from 100  $\mu$ L to 50  $\mu$ L.

- 13. The next day, wash the ELISA plates three times with PBS-T and block them with 200  $\mu$ L of 5% milk in PBS-T with gentle shaking for 1 h.
- 14. During the blocking step, centrifuge the phage induction plates at  $4000 \times g$  for 10 min. Transfer 60 µL of the supernatant to new round-bottom 96-well plates containing 180 µL of 5% milk in PBS-T.

**OPTIONAL** A larger volume (up to  $80 \ \mu$ L) can be transferred to increase the volume in the next step, but care should be taken to avoid aspirating bacterial pellets. Also, consider the maximum volume capacity of the round-bottom plates.

- 15. After blocking, add 100 μL of the diluted phage solution to each corresponding well in the ELISA plates. Incubate for 1 h (or 2 h to increase the ELISA signal) at RT with gentle agitation.
- Wash the plates three times with PBS-T, then add 100 μL of HRP-conjugated anti-M13 coat protein antibody, following the manufacturer's dilution guidelines. Incubate for 1 h at RT with gentle shaking.
- 17. Wash the plate five times with PBS-T and add 100  $\mu$ L of TMB substrate to develop the color. Stop the reaction with 50  $\mu$ L of 1 M H<sub>2</sub>SO<sub>4</sub> and quantify the signal of OD<sub>450</sub> using a plate reader.
- 18. Normalize the ELISA results by subtracting the control antigen quantifications from those of the target antigen. Select colonies with normalized values exceeding a predetermined cutoff for further assays specific to the project.

**CRITICAL STEP** High ELISA signals, even after normalization, do not guarantee strong binding in other assays, such as FACS of Western blot using phages carrying scFv. Therefore, it is recommended to use a low cutoff value, such as 'top 33%', for leniency in determining which colonies to further analyze.

19. Concurrently with hit selection, perform sequencing of the selected hits. This can be done by inoculating bacteria from the frozen master plate into 3 to 5 mL of LB media supplemented with ampicillin in a 14 mL tube for miniprep or colony PCR. Standard Sanger sequencing can then be carried out to identify the sequence of the hits, with primer details provided in Appendix A. Identify unique hits to proceed with further specific assays.

# 4. Expected Results

#### 4.1. scFv Library Generation

Our example library was derived from the PBMCs of patients with an active viral infection. Total RNA was extracted and subsequently reverse-transcribed to generate cDNA. This cDNA served as the template for the amplification of the V regions. The amplified products of  $V_H$ -S,  $V_H$ -L,  $V_\kappa$ , and  $V_\lambda$  were approximately 400 bp in size (Figure 3A), which were then precipitated and purified using gel extraction. The extracted  $V_H$  and  $V_L$  constructs were linked via secondary overlapping extension (SOE) PCR, resulting in major bands between 700 and 800 bp, varying based on the linker length (Figure 3B). The scFvs thus generated were electroporated into XL1-Blue competent cells. The test library displayed an estimate size of  $1.8 \times 10^7$  per ligation reaction, with minimal colony formation observed on control #1 (without ligase, Figure 3C) and control #2 plates (without insert, Figure 3D). In contrast, the control ligation involving a linearized vector and a stuffer insert resulted in numerous colonies, indicating successful ligation (Figure 3E). This suggests that the majority of the colonies were products of scFv ligated with the pComb3XSS vector.

The test library, when plated with a dilution series, yielded two and four colonies from duplicate inoculations. The size of the test library was calculated as follows:

*Test library size from single ligation* =

 $(average \ count \ of \ colonies) \times dilution \ factor \times \frac{total \ volume \ of \ elctroporated \ competent \ cells}{inoculation \ volume}$ 

Using our example (Figure 3F), the calculation was:

$$\frac{(2+4)}{2} \times 10^4 \times \frac{3000 \ \mu L}{5 \ \mu L} = 1.8 \times 10^7$$

To ensure a library with higher diversity, we set a goal for the library size to exceed  $1 \times 10^8$  in the large-scale library generation. We conducted 20 ligation reactions and used a

dilution series to determine the final library size (Figure 3G). The size of the final library was calculated as follows:



**Figure 3.** Expected results from scFv library generation. (**A**) Agarose gel electrophoresis of  $V_H$ -S,  $V_H$ -L,  $V_{\kappa}$ , and  $V_{\lambda}$ : This panel shows an example of  $V_H$ -S,  $V_H$ -L,  $V_{\kappa}$ , and  $V_{\lambda}$  analyzed using 1.5% agarose gel electrophoresis. The expected sizes are approximately 400 bp for  $V_H$ -S, 450 bp for  $V_H$ -L, and 350 bp each for  $V_{\kappa}$  and  $V_{\lambda}$ . (**B**) Random scFv fragments: Here, eight randomly selected scFv fragments are analyzed using 1.5% agarose gel electrophoresis. (**C**–**E**) Control ligation results: These panels depict control ligation results for ligation without ligase (**C**), without insert (**D**), and ligation with the vector and stuffer (**E**). (**F**) Test library generation: This panel demonstrates the result of the test library generation, highlighting two colonies in a green box and four colonies in a purple box. (**G**) Large-scale library generation: An example result of large-scale library generation, showing three colonies in a blue box and one colony in a red box.

Upon confirming that the final library's size exceeded our target of  $1 \times 10^8$ , we proceeded to amplify the bacterial library on 245 mm plates. After amplification, the bacteria were scraped from the plates. We then prepared bacterial library stocks, aliquoting them into vials, each containing  $5 \times 10^9$  bacteria.

To assess the diversity of the generated library, we performed Sanger sequencing on randomly selected monoclonal bacterial colonies (as detailed in Table 5). This analysis revealed an even distribution among the  $V_H$ -S- $V_\kappa$ ,  $V_H$ -L- $V_\kappa$ ,  $V_H$ -S- $V_\lambda$ , and  $V_H$ -L- $V_\lambda$  constructs. Such an even distribution is indicative of a well-representative library, encompassing a broad range of the targeted antibody repertoire.

#### 4.2. Phage Display Affinity Selection and Antigen-Specific Antibody Development

Our phage display affinity selection process commenced with a bacterial library containing approximately  $5 \times 10^9$  bacteria. We initiated phage production by infecting the bacterial cells with  $4 \times 10^{10}$  cfu of helper phages. Once the phages were induced, they were precipitated using a PEG solution. Subsequently, the phages carrying scFv repertoires were serially diluted and applied to plates for titration prior to in vitro selection, as demonstrated in Figure 4A.



**Figure 4.** Expected results from phage display screening. (**A**) Input phage titration: This panel shows an example of titrating the input phages. (**B**) Control antigen output phage titration: This panel displays the results of titration for phages bound to the control antigen. (**C**) Target antigen output phage titration: This panel illustrates the results of titration for phages bound to the target antigen. (**D**–**F**) Enrichment analysis results: These panels present the outcomes of the enrichment analysis. Panel D shows the results of enrichment factor calculations; Panel E depicts the output/input phage titre calculations; and Panel F presents data from polyclonal ELISA, all as outlined in Section 3.7.3, "Enrichment analysis," and Section 3.7.4, "Polyclonal ELISA." Blue arrows in these panels indicate the recommended stopping point for affinity selection. R0, R1, R2, R3, and R4 represent phages rescued from the naïve scFv library and after rounds 1, 2, 3, and 4 of affinity selection, respectively.

				1	0 0	0 0	0,			0	5	0 0		
Sample ID	Heavy Chain V Segment Family	Heavy Chain J Segment Family	Light Chain V Segment Family	Light Chain J Segment Family	Sample ID	Heavy Chain V Segment Family	Heavy Chain J Segment Family	Light Chain V Segment Family	Light Chain J Segment Family	Sample ID	Heavy Chain V Segment Family	Heavy Chain J Segment Family	Light Chain V Segment Family	Light Chain J Segment Family
QC1	HUVH3E	Hujh3-S	HUVK6B	hujk1	QC11	HUVH6	Hujh6b-L	HUVK1C	hujk2	QC21	HUVH3C	hujh45	HUVL2C	hujl23
QC2	HUVH3B	Hujh45-S	HUVK3B	hujk2	QC12	HUVH1A	Hujh2-L	HUVL3C	hujl23	QC22	HUVH1A	hujh6b	HUVL2B	hujl23
QC3	HUVH1A	Hujh2-S	HUVL2B	hujl1	QC13	HUVH1C	Hujh6b-L	HUVK3C	hukj2	QC23	HUVH4C	hujh45	HUVK1C	hujk2
QC4	HUVH3C	Hujh1-S	HUVL6	hujl1	QC14	HUVH4A	Hujh3-S	HUVK3A	hujk4	QC24	HUVH1A	hujh45	HUVK1B	hujk3
QC5	HUVH1C	Hujh2-L	HUVK1B	hujk1	QC15	HUVH1C	Hujh1-L	HUVK3A	hujk2	QC25	HUVH1A	Hujh1	HUVL2C	hujl1
QC6	HUVH1A	Hujh3-L	HUVL5A9	hujl1	QC16	HUVH4C	Hujh2-L	HUVK1A	hujk4	QC26	HUVH3E	Hujh3	HUVL8	hujl23
QC7	HUVH1C	Hujh3-S	HUVK2A	hujk1	<b>QC17</b>	HUVH4D	Hujh2-S	HUVL1A	hujl1	QC27	HUVH1C	Hujh3	HUVL6	hujl1
QC8	HUVH1C	Hujh45-L	HUVK3B	hujk3	<b>QC18</b>	HUVH4A	Hujh6a-S	HUVK2B	hujk1	QC28	HUVH4B	hujh6b	HUVL1B	hujl23
QC9	HUVH4D	Hujh45-L	HUVL10	hujl1	QC19	HUVH4B	Hujh6a-S	HUVL2A	hujl7	QC29	HUVH2A	Hujh1	HUVK1A	hujk3
QC10	HUVH1B	Hujh45-S	HUVL10	hujl1	QC20	HUVH1A	hujh2	HUVK6A	hujk4					

**Table 5.** Library quality control (QC) sequence analysis. scFv constructs featuring  $V_H$ - $V_\kappa$  are denoted with an underlined font.  $V_H$ - $V_\lambda$  scFv constructs are indicated in bold font. scFvs incorporating long linkers are highlighted in grey. Those with short linkers are distinguished with a yellow highlight.

The quantification of phage titer was determined using the following formula:

Phage titer as cfu/mL =count of colonies × dillution factor ×  $\frac{total \ volume \ of \ infected \ bacteria}{volume \ of \ inoculated \ bacteria}$ 

In our study (as illustrated in Figure 4A), the calculation was performed as follows:

$$6 \times 10^{11} \times \frac{100 \ \mu L}{5 \ \mu L} = 1.2 \times 10^{13} \frac{cfu}{mL}$$

Initially,  $1 \times 10^{11}$  cfu of induced phages displaying scFv were subjected to control antigen binding (input phage). The unbound phages were then transferred to wells coated with the target antigen. Bound phages for both control and target antigens were eluted separately (output phage, as shown in Figure 4B,C). We quantified the output phage titers from both the control antigen and target antigen across multiple rounds for enrichment analysis, as described in Section 3.7.3, "Enrichment analysis" (Figure 4D,E). Phages eluted from wells initially coated with target antigen were then amplified for subsequent rounds of affinity selection. Notably, both the enrichment factor and the output/input phage titer decreased (Figure 4D,E), while the polyclonal ELISA reached a plateau (Figure 4F). This indicated that maximum enrichment had been achieved, and further rounds of affinity selection were likely to result in the loss of target-specific phages. Therefore, we did not proceed with additional rounds of affinity selection.

We performed serial dilutions of the bacterial stock from Round 3, as the product from the third round of affinity selection was expected to contain the most diverse anti-target antigen scFvs. Single bacterial colonies were then picked into 96-well plates. In our case, we selected a total of 14 plates (1330 single colonies, leaving the A1 well empty as a negative control). Following colony growth and phage induction in each well, the induced phages were used for binding to plates coated with both control antigen (Figure 5A) and target antigen (Figure 5B) in ELISA assays.

The normalized ELISA value was calculated as follows:

#### Normalized ELISA value

#### $= OD_{450}$ from target protein ELISA $- OD_{450}$ from control protein ELISA

We displayed the normalized ELISA values using a color gradient to visually highlight wells where the normalized ELISA exceeded an arbitrary cutoff value (Figure 5C).

So far, this protocol has been successfully employed to generate multiple pathogenspecific libraries and to screen for a variety of antigen-specific binders. Notably, the consistency of the data has been demonstrated across a diverse range of proteins. This includes applications in identifying viral surface glycoproteins for therapeutic antibody development, as well as transmembrane G-protein-coupled receptors (GPCRs) for structural analysis.

In summary, this paper presents a comprehensive and meticulous scFv library generation protocol, along with a detailed description of its subsequent phage display screening procedure. The methods outlined here have been shown to be effective and reproducible, offering significant potential for varied applications in the field of molecular biology and therapeutic development. Δ

		1	2	3	4	5	6	7	8	9	10	11	12
	Α	0.15	0.37	0.38	0.49	0.09	0.1	0.66	0.07	0.12	0.09	0.54	0.09
G	в	0.09	0.06	0.59	0.45	0.81	0.45	0.51	0.46	0.53	0.44	0.08	0.12
	С	0.39	0.12	0.69	0.04	0.53	0.29	1.03	0.04	0.91	0.46	0.79	0.52
	D	0.47	0.9	0.87	0.53	0.87	0.6	0.49	0.94	0.64	0.04	0.64	0.57
	E	0.24	0.4	0.72	0.15	0.18	0.38	0.57	0.04	0.04	0.92	0.05	0.3
2 COODAA AAAAAA	F	0.66	0.71	0.08	0.04	0.06	0.66	0.04	0.41	0.9	0.79	0.21	0.71
	G	0.56	0.49	0.88	0.38	0.17	0.05	0.06	0.68	0.1	0.06	0.4	0.46
	н	0.48	0.5	0.62	0.42	0.99	0.16	0.49	0.06	0.07	0.78	0.1	0.42
No. A CALLAGE CONTRACTOR TO THE POPULATION OF A CALL STORE													

Β.														
	13		1	2	3	4	5	6	7	8	9	10	11	12
		Α	0.1	0.07	0.06	0.09	0.06	0.06	0.06	0.14	0.05	0.05	0.06	0.08
	1 C C C C C C C C C C C C C C C C C C C	в	0.06	0.05	0.05	0.05	0.04	0.05	0.04	0.04	0.05	0.05	0.05	0.07
	CONTRACTOR OF	С	0.06	0.04	0.06	0.04	0.06	0.04	0.04	0.05	0.04	0.04	0.05	0.07
	A CONSIGNOR COLORIDA A	D	0.06	0.04	0.04	0.04	0.04	0.04	0.04	0.04	0.04	0.06	0.04	0.07
		E	0.06	0.05	0.04	0.04	0.04	0.04	0.04	0.04	0.04	0.04	0.04	0.06
	C C C C C C C C C C C C C C C C C C C	F	0.07	0.05	0.25	0.06	0.11	0.09	0.04	0.04	0.1	0.04	0.08	0.09
		G	0.07	0.05	0.04	0.05	0.04	0.04	0.04	0.04	0.04	0.04	0.05	0.09
		Н	0.11	0.07	0.06	0.06	0.06	0.06	0.07	0.06	0.06	0.07	0.08	0.18
	000000000000000													

C.	G													
		1	2	3	4	5	6	7	8	9	10	11	12	
	Α	0.05	0.29	0.32	0.4	0.03	0.04	<u>0.6</u>	-0.06	0.06	0.04	<u>0.48</u>	0.01	
	В	0.04	0.02	<u>0.55</u>	0.4	<u>0.77</u>	0.41	0.47	0.42	<u>0.49</u>	0.4	0.04	0.06	
	С	0.33	0.07	<u>0.63</u>	-0	<u>0.47</u>	0.25	<u>0.99</u>	-0.01	<u>0.87</u>	0.41	<u>0.74</u>	0.44	
	D	0.41	<u>0.86</u>	<u>0.82</u>	<u>0.49</u>	<u>0.83</u>	<u>0.56</u>	0.45	<u>0.9</u>	<u>0.6</u>	-0.02	<u>0.59</u>	<u>0.5</u>	
	Е	0.18	0.36	<u>0.67</u>	0.11	0.14	0.33	<u>0.53</u>	-0	-0	<u>0.88</u>	0	0.24	
	F	<u>0.59</u>	<u>0.67</u>	-0.16	-0.02	-0.05	<u>0.57</u>	-0	0.37	<u>0.8</u>	<u>0.74</u>	0.12	<u>0.63</u>	
	G	<u>0.49</u>	0.44	<u>0.83</u>	0.33	0.13	0.01	0.02	<u>0.63</u>	0.06	0.01	0.35	0.37	
	Н	0.37	0.43	<u>0.56</u>	0.36	<u>0.94</u>	0.1	0.43	0	0	<u>0.71</u>	0.02	0.24	

**Figure 5.** Expected results from monoclonal ELISA. (**A**,**B**) Monoclonal phage ELISA results: (**A**) represents the results using the target antigen as the coating antigen, while (**B**) shows the results with the control antigen used as the coating antigen. (**C**) Normalized monoclonal phage ELISA results: This panel illustrates the outcome of the normalized monoclonal phage ELISA. The normalization was performed by subtracting the OD<sub>450</sub> values obtained from the control antigen ELISA (**B**) from those of the target antigen ELISA (**A**). The results that exceeded our predefined cutoff value (OD<sub>450</sub> = 0.47) are highlighted in red, bold, and underlined font, indicating positive hits.

#### 5. Reagents Setup

A 45% glucose solution: Slowly add glucose to 1 L of deionized water while stirring. Add only 50 g at a time until 450 g is added. Mildly heat the solution up to 50 °C or lower. Filter the solution through a 0.2  $\mu$ m filter.

LB agar plate supplemented with 2% w/v glucose and  $50 \mu g/mL$  ampicillin: Dissolve 20 g of LB broth powder in 900 mL of deionized water. After autoclaving and cooling the mixture to approximately  $60 \degree$ C, add 44.4 mL of filtered 45% glucose solution and ampicillin to the final 50  $\mu$ g/mL. Add autoclaved and deionized water to 1 L. For distribution, pour 20 mL into each 10 cm Petri dish and 150 mL into each 245 mm square plate.

 $H_2SO_4$  ELISA stop solution: Carefully mix 10 mL of 100%  $H_2SO_4$  into 177.7 mL of deionized water, ensuring to add the acid slowly and steadily into the water.

2xTY-Growth medium (2xTY with 2% w/v glucose and ampicillin, 2xTY-G): Dissolve 5 g NaCl, 10 g yeast extract, and 16 g bacto-tryptone in 900 mL of deionized water. After autoclaving, cool the solution to RT, then add 44.4 mL of filtered 45% glucose solution. Add autoclaved and deionized water to 1 L. Just before use, enrich the medium with ampicillin to a final concentration of 50 µg/mL.

2xTY-Induction (2xTY with ampicillin and kanamycin, 2xTY-I) medium: Dissolve 5 g NaCl, 10 g yeast extract, and 16 g bacto-tryptone in 1L of deionized water. After autoclaving, allow the medium to reach RT. Before use, fortify the medium with both ampicillin and kanamycin to reach a final concentration of 50  $\mu$ g/mL each.

PEG solution (20% PEG-6000, 2.5 M NaCl): Dissolve 100 g of PEG-6000 and 73 g of NaCl in 500 mL of deionized water. Filter through a 0.2  $\mu$ m filter for sterility.

PBS-T (0.05% Tween 20): Add 250  $\mu L$  Tween 20 to 500 mL PBS.

AEBSF protease inhibitor: dissolve 4 mg of AEBSF in 1 mL of PBS and store aliquots at  $-20\ ^\circ C$  for later use.

Trypsin stock solution (1 mg/mL): Dissolve 1 g of trypsin powder in 100 mL of PBS. Store this solution in aliquots at -20 °C. When needed, dilute the stock solution with PBS at a ratio of 1:4 to achieve a working trypsin solution of 0.25 mg/mL.

**Author Contributions:** Conceptualization, Z.L. and D.K.; methodology: Z.L., D.K., and S.K.; analysis: Z.L., D.K., and S.K.; writing original draft, Z.L. and D.K.; editing, Z.L., D.K., and J.U.J.; funding acquisition, J.U.J.; project administration, J.U.J.; resources, J.U.J.; supervision, J.U.J.; reviewing, J.U.J. All authors have read and agreed to the published version of the manuscript.

**Funding:** The study was supported by the National Institute of Health (CA251275, AI140718, AI152190, AI171201, AI172252).

**Institutional Review Board Statement:** The study was performed according to Cleveland Clinic Lerner Research Institute Institutional Review Board (IRB) approval #21-581 and Institutional Biosafety Committee (IBC) approval #2030.

Informed Consent Statement: Blood samples were collected from patients upon informed consent.

**Data Availability Statement:** No additional data were created or analyzed in this study. Data sharing is not applicable to this article. The authors can be contacted for any further information regarding the data within the article.

**Conflicts of Interest:** The authors have declared no conflict of interest. The funders had no role in the design of the study; in the collection, analyses, or interpretation of the data; in the writing of the manuscript; or in the decision to publish the results.

# Appendix A Primers for scFv Library Generation

# V<sub>H</sub> F' primers (nineteen in total)

V<sub>H</sub> F1:

5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTGCAGCTGGTGCAGTCTGG-3'

VH F2:

5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTYCAGCTKGTGCAGTCTGG-3'

VH F3:

5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTCCAGCTGGTACAGTCTGG-3'

VH F4:

5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCARATGCAGCTGGTGCAGTCTGG-3'

# VH F5:

5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTSCAGCTGGTGCARTCTGG-3'

VH F6: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGRTCACCTTGAAGGAGTCTGG-3' VH F7: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTCACCTTGAGGGAGTCTGG-3' VH F8: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCSAGGTGCAGCTGGTGGAGTCTGG-3' VH F9: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCGAGGTGCAGCTGTTGGAGTCTGG-3' VH F10: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCGAGGTGCAGCTGGTGGAGWCYGG-3' VH F11: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCGAAGTGCAGCTGGTGGAGTCTGG-3' VH F12· 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTACAGCTGGTGGAGTCTGG-3' VH F13: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGSTGCAGCTGCAGGAGTCGGG-3' VH F14: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTGCAGCTACAGCAGTGGGG-3' VH F15: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGCTGCAGCTGCAGGAGTCCGG-3' VH F16: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTGCAGCTACAACAGTGGGG-3' VH F17: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCGARGTGCAGCTGGTGCAGTCTGG-3' VH F18: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTACAGCTGCAGCAGTCAGG-3' VH F19: 5'-GGTTTCGCTACCGTGGCCCAGGCGGCCATGGCCCAGGTGCAGCTGGTGCAATCTGG-3' VH R' short linker primers (VH-S R', six in total) VH-SR1: 5'-GGAAGATCTAGAGGAACCACC TGAGGAGACGGTGACCAGGGTGCCCTG-3' VH-S R2· 5'-GGAAGATCTAGAGGAACCACC TGAGGAGACAGTGACCAGGGTGCCACG-3'

VH-S R3: 5'-GGAAGATCTAGAGGAACCACC TGAAGAGACGGTGACCATTGTCCCTTG-3'

#### VH-S R4:

5'-GGAAGATCTAGAGGAACCACC TGAGGAGACGGTGACCAGGGTYCCYTG-3'

VH-S R5: 5'-GGAAGATCTAGAGGAACCACC TGAGGAGACGGTGACCGTGGTCCCTTG-3'

VH-S R6: 5'-GGAAGATCTAGAGGAACCACC TGAGGAGACGGTGACCGTGGTCCCTTT-3'

#### VH R' long linker (VH-L R', six in total)

VH-L R'2:

#### VH-L R'3:

#### VH-L R'4:

# VH-L R'5:

#### VH-L R'6:

# Vκ (kappa light chain) F' (twelve in total)

Vκ F'1: 5'-GGTGGTTCCTCTAGATCTTCC GACATCCAGWTGACCCAGTCTCC-3'

Vκ F'2: 5'-GGTGGTTCCTCTAGATCTTCC GCCATCCRGWTGACCCAGTCTCC-3'

Vκ F'3: 5'-GGTGGTTCCTCTAGATCTTCC GTCATCTGGATGACCCAGTCTCC-3'

# Vκ F'4: 5'-GGTGGTTCCTCTAGATCTTCC GATATTGTGATGACCCAGACTCC-3'

Vκ F'5: 5'-GGTGGTTCCTCTAGATCTTCC GATRTTGTGATGACTCAGTCTCC-3'

Vκ F'6: 5'-GGTGGTTCCTCTAGATCTTCC GAAATTGTGTTGACRCAGTCTCC-3' Vκ F'7: 5'-GGTGGTTCCTCTAGATCTTCC GAAATAGTGATGAYGCAGTCTCC-3' Vκ F'8: 5'-GGTGGTTCCTCTAGATCTTCC GAAATTGTAATGACACAGTCTCC-3' Vκ F'9: 5'-GGTGGTTCCTCTAGATCTTCC GACATCGTGATGACCCAGTCTCC-3' Vк F'10: 5'-GGTGGTTCCTCTAGATCTTCC GAAACGACACTCACGCAGTCTCC-3' Vк F′11: 5'-GGTGGTTCCTCTAGATCTTCC GAAATTGTGCTGACTCAGTCTCC-3' Vк F′12∙ 5'-GGTGGTTCCTCTAGATCTTCC GATGTTGTGATGACACAGTCTCC-3' Vκ (kappa light chain) R' (five in total) Vκ R'1: 5'-ATGGTGATGGTGATGGTGCTGGCCGGCCTGGCCACGTTTGATTTCCACCTTGGTCCC-3' Vκ R'2: 5'-ATGGTGATGGTGATGGTGCTGGCCGGCCTGGCCACGTTTGATCTCCAGCTTGGTCCC-3' Vκ R'3: 5'-ATGGTGATGGTGATGGTGCTGGCCGGCCTGGCCACGTTTGATATCCACTTTGGTCCC-3' Vκ R'4: 5'-ATGGTGATGGTGATGGTGCTGGCCGGCCTGGCCACGTTTGATCTCCACCTTGGTCCC-3' Vκ R'5: 5'-ATGGTGATGGTGATGGTGCTGGCCGGCCTGGCCACGTTTAATCTCCAGTCGTGTCCC-3'  $V_{\lambda}$  (lambda light chain) F' (twenty in total) V $\lambda$  F'1: 5'-GGTGGTTCCTCTAGATCTTCC CAGTCTGTGCTGACTCAGCCACC-3'  $V\lambda F'2$ : 5'-GGTGGTTCCTCTAGATCTTCC CAGTCTGTSSTGACGCAGCCGCC-3' Vλ F'3: 5'-GGTGGTTCCTCTAGATCTTCC CAGTCTGTGTTGACGCAGCCGCC-3'  $V\lambda F'4$ : 5'-GGTGGTTCCTCTAGATCTTCC CAGTCTGCCCTGACTCAGCCTCC-3'  $V\lambda F'5$ : 5'-GGTGGTTCCTCTAGATCTTCC CAGTCTGCCCTGACTCAGCCTCG-3' V $\lambda$  F'6:

5'-GGTGGTTCCTCTAGATCTTCC CAGTCTGCCCTGACTCAGCCTGC-3'

Vλ F′7: 5'-GGTGGTTCCTCTAGATCTTCC TCCTATGWGCTGACTCAGCCACC-3'  $V\lambda F'8:$ 5'-GGTGGTTCCTCTAGATCTTCC TCCTATGAGCTGACTCAGCCACT-3' Vλ F'9: 5'-GGTGGTTCCTCTAGATCTTCC TCTTCTGAGCTGACTCAGGACCC-3' Vλ F'10: 5'-GGTGGTTCCTCTAGATCTTCC TCCTATGAGCTGACACAGCYAYC-3' Vλ F'11: 5'-GGTGGTTCCTCTAGATCTTCC TCCTATGAGCTGATGCAGCCAC-3' Vλ F'12· 5'-GGTGGTTCCTCTAGATCTTCC CTGCCTGTGCTGACTCAGCCCCCGT-3' Vλ F'13: 5'-GGTGGTTCCTCTAGATCTTCC CAGCCTGTGCTGACTCAATCATC-3' Vλ F'14: 5'-GGTGGTTCCTCTAGATCTTCC CAGCTTGTGCTGACTCAATCGCC-3' Vλ F'15: 5'-GGTGGTTCCTCTAGATCTTCC CAGCCTGTGCTGACTCAGCCAYC-3' Vλ F'16: 5'-GGTGGTTCCTCTAGATCTTCC CAGGCTGTGCTGACTCAGCCGKC-3' Vλ F'17: 5'-GGTGGTTCCTCTAGATCTTCC AATTTTATGCTGACTCAGCCCCA-3' Vλ F'18: 5'-GGTGGTTCCTCTAGATCTTCC CAGRCTGTGGTGACTCAGGAGCC-3' Vλ F'19: 5'-GGTGGTTCCTCTAGATCTTCC CAGACTGTGGTGACCCAGGAGCC-3' Vλ F'20: 5'-GGTGGTTCCTCTAGATCTTCC CAGGCAGGGCTGACTCAGCCACC-3'  $V_{\lambda}$  (lambda light chain) R' (three in total) V $\lambda$  R'1: 5'-ATGGTGATGGTGCTGGCCGGCCTGGCC ACCTAGGACGGTGACCTTGGTCCCAG-3' Vλ R'2: 5'-ATGGTGATGGTGCTGGCCGGCCTGGCC ACCTAGGACGGTCAGCTTGGTCCCTC-3'

Vλ R'3: 5'-ATGGTGATGGTGCTGGCCGGCCTGGCC ACCGAGGACGGTCAGCTGGGTGCCTC-3'

#### **SOE—Secondary Overlapping Extension (two):**

SOE F': 5'-GAGGAGGAGGAGGAGGAGGGTTTCGCTACCGTGGCCCAGGCGGCC-3'

SOE R':

# Recommended sequencing primer (for pComb3XSS phagemid):

F':

# 5'-ATTTAAAATGAAAAAGACAGCTATCGCG-3'

R′:

# 5'-GCTTAAATTAATTAGCTAGCTTAAGACTCC-3'

#### References

- 1. Hansel, T.T.; Kropshofer, H.; Singer, T.; Mitchell, J.A.; George, A.J.T. The safety and side effects of monoclonal antibodies. *Nat. Rev. Drug Discov.* **2010**, *9*, 325–338. [CrossRef]
- 2. Jin, S.; Sun, Y.; Liang, X.; Gu, X.; Ning, J.; Xu, Y.; Chen, S.; Pan, L. Emerging new therapeutic antibody derivatives for cancer treatment. *Signal Transduct. Target. Ther.* **2022**, *7*, 39. [CrossRef]
- 3. Nahta, R.; Esteva, F.J. Trastuzumab: Triumphs and tribulations. Oncogene 2007, 26, 3637–3643. [CrossRef]
- 4. Zhao, S.; Chadwick, L.; Mysler, E.; Moots, R.J. Review of Biosimilar Trials and Data on Adalimumab in Rheumatoid Arthritis. *Curr. Rheumatol. Rep.* **2018**, *20*, 57. [CrossRef]
- Montgomery, H.; Hobbs, F.D.R.; Padilla, F.; Arbetter, D.; Templeton, A.; Seegobin, S.; Kim, K.; Campos, J.A.S.; Arends, R.H.; Brodek, B.H.; et al. Efficacy and safety of intramuscular administration of tixagevimab–cilgavimab for early outpatient treatment of COVID-19 (TACKLE): A phase 3, randomised, double-blind, placebo-controlled trial. *Lancet Respir. Med.* 2022, 10, 985–996. [CrossRef]
- 6. Wang, S.S.; Yan, Y.S.; Ho, K. US FDA-approved therapeutic antibodies with high-concentration formulation: Summaries and perspectives. *Antib. Ther.* **2021**, *4*, 262–272. [CrossRef]
- 7. Stone, C.A., Jr.; Spiller, B.W.; Smith, S.A. Engineering Therapeutic Monoclonal Antibodies. J Allergy Clin Immunol 2023. [CrossRef]
- 8. Holliger, P.; Hudson, P.J. Engineered antibody fragments and the rise of single domains. *Nat. Biotechnol.* **2005**, *23*, 1126–1136. [CrossRef]
- 9. Bird, R.E.; Hardman, K.D.; Jacobson, J.W.; Johnson, S.; Kaufman, B.M.; Lee, S.-M.; Lee, T.; Pope, S.H.; Riordan, G.S.; Whitlow, M. Single-Chain Antigen-Binding Proteins. *Science* **1988**, *242*, 423–426. [CrossRef]
- 10. Ahmad, Z.A.; Yeap, S.K.; Ali, A.M.; Ho, W.Y.; Alitheen, N.B.M.; Hamid, M. scFv antibody: Principles and clinical application. *J. Immunol. Res.* 2012, 2012, 980250. [CrossRef]
- 11. Baylet, A.; Vyumvuhore, R.; Laclaverie, M.; Marchand, L.; Mainzer, C.; Bordes, S.; Closs-Gonthier, B.; Delpy, L. Transcutaneous penetration of a single-chain variable fragment (scFv) compared to a full-size antibody: Potential tool for atopic dermatitis (AD) treatment. *Allergy Asthma Clin. Immunol.* **2021**, *17*, 73. [CrossRef]
- 12. Duan, Y.; Chen, R.; Huang, Y.; Meng, X.; Chen, J.; Liao, C.; Tang, Y.; Zhou, C.; Gao, X.; Sun, J. Tuning the ignition of CAR: Optimizing the affinity of scFv to improve CAR-T therapy. *Cell. Mol. Life Sci.* **2021**, *79*, 14. [CrossRef]
- KÖHler, G.; Milstein, C. Continuous cultures of fused cells secreting antibody of predefined specificity. *Nature* 1975, 256, 495–497. [CrossRef]
- 14. Schladetsch, M.A.; Wiemer, A.J. Generation of Single-Chain Variable Fragment (scFv) Libraries for Use in Phage Display. *Curr. Protoc.* **2021**, *1*, e182. [CrossRef]
- 15. Safdari, Y.; Farajnia, S.; Asgharzadeh, M.; Khalili, M. Antibody humanization methods—A review and update. *Biotechnol. Genet. Eng. Rev.* **2013**, *29*, 175–186. [CrossRef]
- 16. Ling, W.-L.; Lua, W.-H.; Gan, S.K.-E. Sagacity in antibody humanization for therapeutics, diagnostics and research purposes: Considerations of antibody elements and their roles. *Antib. Ther.* **2020**, *3*, 71–79. [CrossRef]
- 17. Hoogenboom, H.R. Selecting and screening recombinant antibody libraries. Nat. Biotechnol. 2005, 23, 1105–1116. [CrossRef]
- 18. Zhou, H.; Zhang, Y.-L.; Lu, G.; Ji, H.; Rodi, C.P. Recombinant antibody libraries and selection technologies. *New Biotechnol.* **2011**, *28*, 448–452. [CrossRef]
- 19. Pelat, T.; Thullier, P. Non-human primate immune libraries combined with germline humanization. *mAbs* **2009**, *1*, 377–381. [CrossRef] [PubMed]
- Begent, R.H.J.; Verhaar, M.J.; Chester, K.A.; Casey, J.L.; Green, A.J.; Napier, M.P.; Hopestone, L.D.; Cushen, N.; Keep, P.A.; Johnson, C.J.; et al. Clinical evidence of efficient tumor targetting based on single–chain Fv antibody selected from a combinatorial library. *Nat. Med.* 1996, 2, 979–984. [CrossRef] [PubMed]

- 21. Zhao, Q.; Buhr, D.; Gunter, C.; Frenette, J.; Ferguson, M.; Sanford, E.; Holland, E.; Rajagopal, C.; Batonick, M.; Kiss, M.M.; et al. Rational library design by functional CDR resampling. *New Biotechnol.* **2018**, *45*, 89–97. [CrossRef] [PubMed]
- 22. Lerner, R.A. Combinatorial antibody libraries: New advances, new immunological insights. *Nat. Rev. Immunol.* **2016**, *16*, 498–508. [CrossRef] [PubMed]
- 23. Kanekiyo, M.; Bu, W.; Joyce, M.G.; Meng, G.; Whittle, J.R.; Baxa, U.; Yamamoto, T.; Narpala, S.; Todd, J.P.; Rao, S.S.; et al. Rational Design of an Epstein-Barr Virus Vaccine Targeting the Receptor-Binding Site. *Cell* **2015**, *162*, 1090–1100. [CrossRef] [PubMed]
- 24. Smith, G.P. Filamentous Fusion Phage: Novel Expression Vectors That Display Cloned Antigens on the Virion Surface. *Science* **1985**, *228*, 1315–1317. [CrossRef]
- 25. Bazan, J.; Całkosiński, I.; Gamian, A. Phage display—A powerful technique for immunotherapy: 1. Introduction and potential of therapeutic applications. *Hum. Vaccin. Immunother.* **2012**, *8*, 1817–1828. [CrossRef]
- 26. Pande, J.; Szewczyk, M.M.; Grover, A.K. Phage display: Concept, innovations, applications and future. *Biotechnol. Adv.* 2010, *28*, 849–858. [CrossRef]
- Ledsgaard, L.; Kilstrup, M.; Karatt-Vellatt, A.; McCafferty, J.; Laustsen, A.H. Basics of Antibody Phage Display Technology. *Toxins* 2018, 10, 236. [CrossRef]
- 28. Rader, C.; Barbas, C.F. Phage display of combinatorial antibody libraries. Curr. Opin. Biotechnol. 1997, 8, 503–508. [CrossRef]
- 29. Rader, C. Generation of human Fab libraries for phage display. *Methods Mol. Biol.* 2012, 901, 53–79. [CrossRef]
- 30. Barbas, C.F. Phage Display: A Laboratory Manual; Cold Spring Harbor Laboratory Press: Cold Spring Harbor, NY, USA, 2001.
- 31. Janeway, C.A., Jr.; Travers, P.; Walport, M.; Shlomchik, M.J. *Immunobiology: The Immune System in Health and Disease*, 5th ed.; Garland Publishing Inc: New York, NY, USA, 2001.
- 32. Wang, S.; Zheng, C.; Liu, Y.; Zheng, H.; Wang, Z. Construction of multiform scFv antibodies using linker peptide. *J. Genet. Genom.* **2008**, *35*, 313–316. [CrossRef]
- Holliger, P.; Prospero, T.; Winter, G. "Diabodies": Small bivalent and bispecific antibody fragments. *Proc. Natl. Acad. Sci. USA* 1993, 90, 6444–6448. [CrossRef] [PubMed]
- Gu, X.; Jia, X.; Feng, J.; Shen, B.; Huang, Y.; Geng, S.; Sun, Y.; Wang, Y.; Li, Y.; Long, M. Molecular Modeling and Affinity Determination of scFv Antibody: Proper Linker Peptide Enhances Its Activity. *Ann. Biomed. Eng.* 2010, 38, 537–549. [CrossRef] [PubMed]
- Pardon, E.; Laeremans, T.; Triest, S.; Rasmussen, S.G.F.; Wohlkönig, A.; Ruf, A.; Muyldermans, S.; Hol, W.G.J.; Kobilka, B.K.; Steyaert, J. A general protocol for the generation of Nanobodies for structural biology. *Nat. Protoc.* 2014, *9*, 674–693. [CrossRef] [PubMed]
- Thermo Scientific Inc. Thermo Scientific Immunoassay Plate Guide. 2018. Available online: https://assets.thermofisher.com/ TFS-Assets/LCD/Scientific-Resources/Immunoassay\_Plate\_Guide.pdf (accessed on 20 November 2023).
- 37. Lee, C.M.Y.; Iorno, N.; Sierro, F.; Christ, D. Selection of human antibody fragments by phage display. *Nat. Protoc.* 2007, 2, 3001–3008. [CrossRef]

**Disclaimer/Publisher's Note:** The statements, opinions and data contained in all publications are solely those of the individual author(s) and contributor(s) and not of MDPI and/or the editor(s). MDPI and/or the editor(s) disclaim responsibility for any injury to people or property resulting from any ideas, methods, instructions or products referred to in the content.