



Proceeding Paper

# Study of Absorbance and Fluorescence Properties of Laccase and Catechol Solutions in the UV Range <sup>†</sup>

Pauline Conigliaro 1,20, Federica Massaro 3, Marianna Portaccio 4, Maria Lepore 40 and Ines Delfino 1,\*

- Dipartimento di Scienze Ecologiche e Biologiche, Università degli Studi della Tuscia, I-01100 Viterbo, Italy
- Dipartimento di Economia, Ingegneria, Società e Impresa, Università degli Studi della Tuscia, I-01100 Viterbo, Italy
- <sup>3</sup> Dipartimento per l'Innovazione nei Sistemi Biologici, Agroalimentari e Forestali, Università degli Studi della Tuscia, I-01100 Viterbo, Italy
- Dipartimento di Medicina Sperimentale, Università della Campania "L. Vanvitelli", I-80100 Napoli, Italy
- \* Correspondence: delfino@unitus.it
- † Presented at the 9th International Electronic Conference on Sensors and Applications, 1–15 November 2022; Available online: https://ecsa-9.sciforum.net.

Abstract: Laccase is an enzyme belonging to the oxidoreductase class and has copper atoms in the catalytic centre. The catalytic property of this enzyme consists of the enzymatic oxidation of the phenolic compounds, in the corresponding quinones, with the concomitant reduction of molecular oxygen to water. There is growing interest in developing innovative sensing methods for detecting phenolic substrates, such as catechol. Preliminary absorption and fluorescence measurements were carried out in the UV-visible range to evaluate the possibility of using the variations produced in the spectra of laccase and/or catechol to monitor the presence of this substrate. The absorption and fluorescence emission increase upon UV excitation was detected. By monitoring the time course of the fluorescence signal, an evident increase in the signal detected in the UV range is observed until a saturation level is reached. The observed variations in the spectra, in the presence of the catechol, are discussed also in terms of the interactions between the enzyme and the phenolic compound. The results are very promising for the design of new optical detection methods for polyphenol pollutants.

Keywords: laccase; catechol; absorption; fluorescence; sensing



Citation: Conigliaro, P.; Massaro, F.; Portaccio, M.; Lepore, M.; Delfino, I. Study of Absorbance and Fluorescence Properties of Laccase and Catechol Solutions in the UV Range. *Eng. Proc.* **2022**, 27, 32. https://doi.org/10.3390/ ecsa-9-13333

Academic Editor: Jean-marc Laheurte

Published: 1 November 2022

**Publisher's Note:** MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2022 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https://creativecommons.org/licenses/by/4.0/).

## 1. Introduction

Phenolic compounds, such as catechol, bisphenol A, etc., are very dangerous molecules that enter the environment in different ways and are produced in various processes: waste from paper manufacturing, agriculture (pesticides, insecticides, and herbicides), pharmaceuticals, and other processes [1]. Chromatographic techniques are the most used methods for the determination of phenolic compounds, but they require expensive equipment, and specialized operators and do not allow for in situ measurements [2].

There is, therefore, an interest in the development of simple, sensitive, specific, accurate sensing methods for the determination of phenolic compounds. In this frame, optical sensing methods are very promising.

Laccase (*Trametes Versicolor*) is an enzyme belonging to the class of oxidoreductases and has very peculiar optical properties due to the presence of copper atoms in the catalytic centre [2]. They are characterized by low substrate specificity, a characteristic that makes them a key element of numerous biotechnological applications.

The catalytic property of the enzyme consists of the oxidation of phenolic substrates, transforming them into less carcinogenic and mutagenic compounds, the corresponding quinones. Simultaneously, with oxidation, the reduction of  $O_2$  to water takes place. Laccase enzyme has been used for developing a certain number of optical biosensors [1,3], but

Eng. Proc. 2022, 27, 32 2 of 4

further investigation on laccase optical properties can offer new approaches for developing innovative sensing schemes.

As a first step to further evaluate the potential of using the peculiar properties of laccase in realizing biosensors for catechol monitoring, in this work we report on an experimental study of this enzyme's optical absorption and fluorescence properties in the absence and presence of pollutants of environmental interest.

#### 2. Materials and Methods

All chemicals and commercial laccase were purchased from Sigma (Sigma Italia, Milan, Italy). Solutions of free laccase at a concentration of 8  $\mu g/mL$  were obtained using 0.1 M acetate buffer at pH 5.0. The catechol was dissolved in the same buffer at a concentration equal to 0.4 mM.


The absorption spectra were acquired in the 200-900 nm range using a JASCO 630 V spectrophotometer, equipped with a UV Winlab software package. The scan speed was set at 200 nm/min, while the slit was at 5 nm.

Fluorescence emission spectra were acquired in the 350–500 nm range by a PerkinElmer LS55 spectrofluorometer with an excitation wavelength of 320 nm. The slits were set at 5 nm, while the scan speed was at 200 nm/min.

Data analysis was performed with the Origin Pro 7.5 software (OriginLab, Northampton, MA, USA).

#### 3. Results and Discussion

The absorption spectrum of the enzyme laccase Trametes Versicolor in the UV-visible region (Figure 1) shows a peak at 310 nm along with weak features at 330 nm and around 650 nm.




**Figure 1.** Absorption spectrum [Laccase *T. Versicolor*] =  $8 \mu g/mL$  in 0.1 M acetate buffer at pH 5.0. Arrows indicate the presence of the peaks at 310 nm, 330 nm, and 650 nm, respectively.

According to the literature [3,4], the peak at 310 nm can be attributed to the absorption of tryptophan, the weak structure at 330 nm relating to the absorption of pairs of copper atoms of type T3 and relating to the structure at 630 nm with regards to the T1-type copper atom.

The absorption spectrum of the laccase samples was then recorded in the presence of catechol at a fixed concentration (0.4 mM) at different times (t0 = 0 min, 10 min, 25 min, 1 h, 80 min, 100 min, and 2 h after catechol addition). Changes in the absorption spectrum are observed during the time. An increased absorbance is observed in the 300-500 nm range (data not shown). As representative behaviour, the absorbance values detected at 409 nm

Eng. Proc. 2022, 27, 32 3 of 4


(where a signal possibly related to a reaction intermediate could be present) extracted from the recorded spectra are shown in Figure 2 as a function of incubation time. The data show that the absorbance increases with time in the considered time interval.



**Figure 2.** Time course behaviour of absorbance values at a fixed wavelength (409 nm) at different times after catechol addition.

Such modifications could be used to monitor the presence of catechol [5], since they are thought to be related to the laccase–catechol interaction [6].


After measuring the emission spectrum of laccase, we investigated whether the presence of catechol in the solution also produces changes in the fluorescence signal. The fluorescence emission after several incubation periods was detected (Figure 3).



**Figure 3.** Emission spectrum of a solution of laccase and catechol ([laccase] =  $8 \mu g/mL$  and [catechol] = 0.4 mM). Measurements were performed at various times: immediately after catechol addition (t0) (black line), 5 min (red line), 10 min (green line), 50 min (blue line), and 80 min (purple line) after catechol addition.

Eng. Proc. 2022, 27, 32 4 of 4

An increase in the fluorescence emission intensity is observed in the range of 330–420 nm, the most significant effect being detected after 50 min. As an example, we show in Figure 4 the time course behaviour of fluorescence emission intensity at a specific wavelength (360 nm). After 50 min, the fluorescence intensity reaches its saturation level.



**Figure 4.** Time course behaviour of emission values at a fixed wavelength (360 nm). Measurements were performed at various times: immediately after catechol addition, 5 min, 10 min, 50 min and 80 min after catechol addition.

### 4. Conclusions

The results confirm that the use of absorption and fluorescence properties of laccase is promising to build new optical detection methods for the determination of phenol compounds of ecological interest.

**Author Contributions:** Conceptualization, M.L. and I.D.; investigation, P.C., F.M., and M.P.; writing—review and editing, P.C., M.L., and I.D. All authors have read and agreed to the published version of the manuscript.

Funding: This research received no external funding.

Institutional Review Board Statement: Not applicable.

**Informed Consent Statement:** Not applicable.

**Data Availability Statement:** The data presented in this study are available upon request from the corresponding author.

Conflicts of Interest: The authors declare no conflict of interest.

## References

- 1. Delfino, I.; Diano, N.; Lepore, M. Advanced Optical Sensing of Phenolic Compounds for Environmental Applications. *Sensors* **2021**, *21*, 7563. [CrossRef] [PubMed]
- 2. Riva, S. Laccases: Blue enzymes for green chemistry. Trends Biotechnol. 2006, 24, 219–226. [CrossRef] [PubMed]
- 3. Lepore, M.; Portaccio, M. Optical detection of different phenolic compounds by means of a novel biosensor based on sol-gel immobilized laccase. *Biotechnol. Appl. Biochem.* **2017**, *64*, 782–792. [CrossRef]
- 4. Delfino, I.; Viola, D.; Cerullo, G.; Lepore, M. Ultrafast excited-state charge-transfer dynamics in laccase type I copper site. *Biophys. Chem.* **2015**, 200–201, 41–47. [CrossRef]
- 5. Dencheva, N.; Oliveira, S.; Braz, G.; Getya, D.; Malfois, M.; Denchev, Z.; Gitsov, I. Magnetically Responsive PA6 Microparticles with Immobilized Laccase Show High Catalytic Efficiency in the Enzymatic Treatment of Catechol. *Catalysts* **2021**, *11*, 239. [CrossRef]
- 6. Sun, X.; Bai, R.; Zhang, Y.; Wang, Q.; Fan, X.; Yuan, J.; Cui, L.; Wang, P. Laccase-Catalyzed Oxidative Polymerization of Phenolic Compounds. *Appl. Biochem. Biotechnol.* **2013**, *171*, 1673–1680. [CrossRef] [PubMed]