



Review Strategies for Generating Human Pluripotent Stem Cell-Derived-Organoid Culture for Disease Modeling, Drug Screening, and Regenerative Therapy

Zakiya Gania ¹, Syarifah Tiara Noorintan ¹, Ni Putu Diah Pradnya Septiari ¹, Dhea Sandra Fitriany ¹ and Fuad Gandhi Torizal ^{1,2,3},*

- ¹ Department of Biotechnology, Faculty of Science and Technology, Universitas 'Aisyiyah Yogyakarta, Sleman, Yogyakarta 55592, Indonesia
- ² Department of Chemical Systems Engineering, Graduate School of Engineering, The University of Tokyo, Bunkyo City, Tokyo 113-8654, Japan
- ³ Department of Bioengineering, Graduate School of Engineering, The University of Tokyo, Bunkyo City, Tokyo 113-8654, Japan
- * Correspondence: t_gandhi@unisayogya.ac.id

Abstract: Human pluripotent stem cells (hPSCs) have become a powerful tool to generate the various kinds of cell types comprising the human body. Recently, organoid technology has emerged as a platform to generate a physiologically relevant tissue-like structure from PSCs. Compared to an actual human organ, this structure more closely represents a three-dimensional microenvironment than the conventional monolayer culture system for transplantation, disease modeling, and drug development. Despite its advantages, however, the organoid culture system still has various problems related to culture methods, which have become a challenge for attempts to obtain similar physiological properties to their original tissue counterparts. Here, we discuss the current development of organoid culture methods, including the problems that may arise from the currently available culture systems, as well as a possible approach for overcoming their current limitations and improving their optimum utilization for translational application purposes.

Keywords: hPSCs derived-organoids; culture strategy; disease modeling; drug screening; regenerative therapy

1. Introduction

Animal models are a routinely used platform for understanding many biological processes related to human disease. Although they can partially represent human physiological and pathological conditions, these models cannot accurately provide similar features for human tissue physiology and disease mechanisms due to the interspecies genomic divergence which impacts various complex cellular activities [1,2]. In actual translational applications, this discrepancy can reduce accuracy, which potentially contributes to a low success rate in their utilization (for example, clinical trials of newly developed drugs) [3]. On the other hand, commonly used two-dimensional human cell cultures often do not completely represent actual human tissue function and architecture [4].

Organoids are three-dimensional miniature organ-like structures that recapitulate physiologically relevant tissue function through self-organization in vitro. The generation of a human pluripotent stem cells (hPSCs) derived organoid can facilitate the improvement of tissue culture by addressing current limitations in disease modeling or drug screening platforms, as well as providing new insights into their application in regenerative therapy [3].

The development of human-induced pluripotent stem cells (hiPSCs) allows for greater experimental accessibility with fewer ethical constraints than the usage of human embryonic stem cells (hESCs), which require a sacrificial human embryo. Moreover, PSCs can



Citation: Gania, Z.; Noorintan, S.T.; Septiari, N.P.D.P.; Fitriany, D.S.; Torizal, F.G. Strategies for Generating Human Pluripotent Stem Cell-Derived-Organoid Culture for Disease Modeling, Drug Screening, and Regenerative Therapy. *Future Pharmacol.* **2022**, *2*, 360–376. https:// doi.org/10.3390/ futurepharmacol2030025

Academic Editors: Juliana Mara Serpeloni and Colus Ilce Mara

Received: 22 June 2022 Accepted: 30 August 2022 Published: 5 September 2022

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2022 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (https:// creativecommons.org/licenses/by/ 4.0/). be generated either from a healthy or diseased individual. This can be useful for personalized medicine applications, such as disease modeling [5], drug screening [6], or even regenerative therapy [7].

Currently, hPSCs-derived organoid production mostly involves a complex multistep protocol during their differentiation and maturation process, which makes the production of organoids technically challenging. Various methods have been developed to produce organoids from PSCs. However, some important factors related to the organoid culture system need to be considered to achieve the optimum culture conditions for each specific translational purpose. In this review, we highlight the currently available methodology for generating organoids from hPSCs as well as the current challenges involved in and possible alternatives available for improving the generation of hPSCs-derived organoids.

2. The Main Features of hPSCs Derived-Organoids

Organoids have become a promising research tool for elucidating the mechanism of various diseases and for testing the effect of substance exposure, which provides a precise preclinical setting for research [8]. Recently, the concept of organoid transplantation has emerged as a potential method for obtaining insight into alternative partial tissue replacement for regeneration [9]. Nevertheless, possible improvements to the current culture methodologies still need to be explored to achieve better overall reliability for translational applications.

At present, numerous studies of hPSCs-derived organoids have succeeded in representing a broad spectrum of specific organ types for different applications (Table 1). A 3D organoid can enable crosstalk between different cells inside its structure, something which cannot be achieved by a conventional monolayer culture. Therefore, these 3D "miniaturized organs" possess decent potential to become physiologically relevant models and better predictive tools than monolayer cultures, which reflect the cellular interactions between different cells within the organ during their development, healthy or diseased conditions.

Generally, an organoid needs to meet some criteria which confirm the authenticity of the original organ [10]. First, the organoid should be a three-dimensional (3D) structure consisting of multiple cell types that retain the identity of the specific organ. Secondly, this 3D structure should be formed by self-organization according to the similar intrinsic organization principles of the organ. Thirdly, the organoid should recapitulate the key features that represent the functional capability of the actual organ (both in terms of structure and physiological activity) [11,12].

Organoid Type	Culture Method	Targeted Application	Results	Ref	
	Static 96-well V-bottom plates	Disease modeling	The cerebral organoid was successfully delineated defects caused by deletions of the distal tip of chromosome	[13]	
Brain	Static ultra-low attachment 96 well plate and dynamic 24 well plates	Disease modeling	Human cerebral organoids established differential effects in Neurofibromatosis type 1 (NF1) developmental disorder	[14–16	
brain	Static culture low-adhesion 6-cm plates followed by Matrigel embedded	Disease modeling	A patient-specific forebrain organoid was utilized to investigate the pathological changes associated with Miller-Dieker syndrome (MDS), caused by chromosomal deletion	[5]	
	Ultra-low-attachment 96-well plates	Disease modeling	The introduction of hPSCs derived-microglia into hPSCs derived-cerebral organoids reflected an innate immune response system in the central nervous system which showed a phagocytic activity during Zika virus infection.		

Table 1. Selected study of hPSCs derived organoids and their targeted applications.

Organoid Type	Culture Method	Targeted Application	Results	
	Miniaturized spinning bioreactor	Disease modeling and developmental study	Resulted brain organoids were successfully recapitulating key dynamic features of the developing human brain at the molecular, cellular, and structural levels	
	Combined biopolymer scaffold and Matrigel embedded	Developmental study and disease model	The cerebral organoid derived by guided cortical plate formation using a biopolymer scaffold showed a distinctive radial organization of the cerebral cortex and allowed the study of neuronal migration.	[19]
	Matrigel embedded	Transplantation model	Organoid grafts showed a progressive neural maturation and regeneration, followed by neovascularization after they were transplanted into mouse	[20]
	Orbital suspension culture and guided bioprinting	Drug screening and regenerative therapy	Patterned vascularized cortical organoid consisted of neural stem cells, endothelial, and neurons that mimic the vascularized brain showed its potential in drug testing and regeneration model.	[21
	V-bottomed 96-well plate	Drug Screening and disease modeling	Cerebral organoid demonstrated its capability for modeling the teratogenic effects of Zika virus infection as well as testing the therapeutic compounds	[22
	U-bottom ultra-low-attachment 96-well plate	Disease modeling	The presence of vasculature-like structures resulted in an enhanced functional cortical organoids maturation for disease modeling applications	[23
	Ultra-low-attachment 96-well plate	Disease modeling	The 3D system created a reciprocal interaction between the thalamus and cortex by fusing the two distinct region-specific organoids that represent the developing thalamus or cortex.	[24
Esophagus	Matrigel embedded	Disease modeling	Human esophageal organoids were utilized for modeling the human esophageal birth defects, with enhanced survival.	[25
	Matrigel embedded	Disease modeling	The gastric organoid was successfully modeled intestinal metaplasia during Doxorubicin induction	
Gastrointestinal	Static suspension in 24 well plates	Disease modeling	A robust in vitro system was developed to elucidate the mechanisms underlying human stomach development and disease including peptic ulcer disease and gastric cancer	
	Matrigel embedded	Drug Screening and disease modeling	The human fundic-type gastric organoids (hFGOs) platform was reported to be a potential tractable human model system in terms identify and study signaling mechanisms such as normal cellular homeostasis in the fundus and antrum.	
	Matrigel embedded	Disease modeling	Human antral and fundic gastric organoids (hGOs) have been effectively used to study stem cells, recovering gastric cell function and showing a cellular response to injury during H. pylori infection in a patient-specific manner	[29
	Static suspension in a culture insert	Disease modeling	Human intestinal organoids displayed an innate immune response when treated by enteric pathogens	[30
	Matrigel embedded	Disease modeling	Human intestinal organoids (hIOs) were successfully modeled the cystic fibrosis by CFTR alteration as well as the functional recovery after the correction of its mutation	[31

Table 1. Cont.

Organoid Type	Culture Method	Targeted Application	Results		
Thyroid	Hanging drop	Transplantation model	hESC-derived thyroid follicles have successfully produced thyroid hormone in vitro and in vivo after transplantation into thyroid gland ablated mice.		
Thymus	Static transwell-6 well plates	Transplantation model	The thymic organoids were validated for in vitro generation of antigen-specific T cells	[33]	
	Micropatterned	Drug screening	The cardiac organoid successfully represented a decent in vitro pro regenerative drug development and potential reduction of their adverse effects		
	U shape 96 well and Matrigel-embedded	Gene targeting and Drug testing.	Heart-forming organoids were providing simultaneous monitoring of genetic mutation defects, as well as drug screening of teratogenic substances	[34,35	
Heart	2D micropatterned	Developmental study and drug-induced cardiac developmental toxicity.	The development of cardiac organoids presented a platform for comprehensive risk assessment in cardiac developmental toxicity assay, which better predicts thalidomide toxicity in the fetal heart model.	[36]	
	3D micropatterning using PDMS molds and poles	Disease modeling and regeneration	Human cardiac organoids were proposed as the physiologically relevant model of the human heart and completely recover the cardiac function after cryoinjury		
	3D micropattern by agarose hydrogel molds	Drug screening	The cardiac organoids recapitulated 3D tissue-level responses including drug-induced/exacerbated cardiotoxicity and fibrotic effects. This organoid can be applied for studying doxorubicin induced-myocardial infarction	[38]	
Lung	Matrigel embedded	Transplantation model	Human PSCs-derived lung bud tip organoid was able to be fully engrafted into the airways of the mouse model	[39]	
	Suspension culture and Matrigel embedded	Disease modeling	Hepatic organoid consisting of hepatocytes and cholangiocytes was successfully utilized for the disease modeling of liver genetic disease		
Liver	Embedding inside the Matrigel dome	Drug screening and disease modeling	Screening of the drugs that was previously retracted from the market, such as troglitazone, trovafloxacin, and levofloxacin, was successfully employed using liver organoids. Additionally, this organoid was also utilized for modeling hepatic steatosis.	[41]	
	Dynamic suspension culture in Erlenmeyer flask using an orbital shaker	Transplantation, drug screening, and disease modeling	The functional vascularized liver organoid showed liver-like functional features, such as the production of serum proteins and the coagulation factors, as well as supported ureagenesis and bilirubin uptake.		
	Induced pluripotent stem cells	Disease modeling and drug screening	The hepatobiliary organoid was organized as a functional bile canaliculi system, which was disrupted by cholestasis-inducing drugs.		
	Combined static microwell culture, 24-well, 6-well, and 1-well plate.	Transplantation	Vascularized and functional liver organoids generated entirely from iPSCs significantly improved hepatic functionalization and represent functional rescue against acute liver failure via transplantation.		
	Suspension culture in the miniature dialysis device	Drug testing and transplantation	The scalable human liver organoid (hLOs) consisted of functional hepatocytes and cholangiocytes showed a response to drug exposure and showed a bile canaliculi activity	[44]	

Table 1. Cont.

Organoid Type	Culture Method	Targeted Application	Results	Ref	
Pancreas	Microwell chip	Disease modeling	Pancreatic duct-like organoids showed mucins regulations and cystic fibrosis transmembrane conductance.	[45]	
Blood vessel	Coculture techniques in static ultra-low attachment six-well plate	Transplantation model, Modelling diabetic vasculopathy, Drug screening	Transplanted blood vessel organoids which were grown in streptozotocin-treated mice successfully represented a model of diabetic vasculopathy and its drug treatment.	[46]	
Kidney organoids	Ultra-low-attachment 6-well plates followed by spinner flask bioreactor	Disease modeling and drug screening	The organoid successfully modeled the congenital kidney defect that interferes with tubulogenesis and potential drug treatment.	[47,48]	
Assembloid (multiple organoid fusion)	Matrigel embedded	Transplantation model	Transplantation of gut and liver organoids showed an effective therapeutic potential against fulminant liver failure	[49]	
	96-well round bottom low attachment plates	Disease modeling, drug development, and therapeutic transplantation.	Human hepato-biliary-pancreatic organoids represented the therapeutic responses to new potential treatments, modeling biliary atresia and promoting pancreatic fate commitment in Hes1-deficient mice.	[50,51]	

Table 1. Cont.

3. Culture Strategy to Generate hPSCs Derived-Organoids

Currently, several approaches have been developed to generate a miniaturized tissuelike organoid comprising the complex functional cellular components derived from hPSCs. These methods mostly involve a differentiation process via a stepwise treatment using specific combined factors to direct the hPSC cell fates.

This step can be improved by creating or adjusting the invitro culture condition to stimulate subsequent organ-specific lineage differentiation using different available culture platforms. Depending on the particular application and expected organoid type, several factors need to be considered when choosing a culture strategy to obtain the best hPSCs derived-organoid production outcomes (Table 2).

Table 2. A brief comparison of the common culture platform for organoid production.

Consideration	Animal Model	hPSCs	hPSCs Derived-Organoid		
			Matrigel Embedded /Dome	Static Suspension Culture	Dynamic Suspension Culture
Culture maintenance	Not required, but difficult animal handling	Easy	Relatively easy	Relatively Easy	Relatively easy
Represent actual human organ's physiology	Partial, limited by interspecies variability	Poor	good	Good	Good
Tissue complexity	Very good	Poor	good	Good	Good
ECM and cell-cell interaction	Very good, native tissue condition	Poor	good	Good	Good
Ability to represent organogene- sis/developmental biology	No	No	Yes	Yes	Yes
Represent actual human organs physiology	Partial, limited by interspecies variability	Poor	good	Good	Good
Scalability	Not available	Low	Low	High	High
Relative cost production	high	Relatively high	Relatively low	Low	Low
Application for personalized medicine and non-xenogeneic transplantation	Not possible	possible	Not possible	Possible	Possible

Several methodologies were previously developed for generating hPSCs-derived organoids to represent miniaturized specific tissues or organs. The simplest strategy is to directly induce differentiation toward organoid formation from the three-dimensional structure of hPSCs (which are referred to as embryoid bodies (EBs)). This structure can be generated using the dynamic suspension culture, hanging drop techniques, or by using Matrigel (Figure 1). Subsequently, one can employ a stepwise treatment using a cocktail of biomolecule inducers, which mostly consists of lineage-specific signaling factors (e.g., small molecule and/or growth factors) to mimic the cellular signaling during organogenesis in embryo development. This biomolecule formulation differs between the organoid types, depending ultimately on which pathways that needed to direct the hPSCs lineage. For example, hPSCs differentiation into the hepatic lineage requires exposure to a high concentration of Activin A during the endodermal differentiation followed by BMP4 to further direct the endoderm cells into hepatic competent cells. Here, final liver organoid maturation can be induced by HGF and OSM in the differentiation medium [44].




Figure 1. An example of general strategies used to generate hPSCs-derived brain organoids. This organoid can be generated directly from differentiated EBs (upper chart) or by coculture from independently differentiated hPSCs-derived cell types (lower chart).

Direct EBs-based differentiation permits intrinsic self-organization during directed differentiation, which generates three-dimensional structures consisting of multiple cell types comprising a specific organ [52]. Due to their spontaneous morphogenesis during differentiation, this approach can be utilized for developing a model for the developmental study of a certain organ during embryogenesis. However, the emergence of undesirable cell types that did not complement the expected organ is more difficult to control during differentiation. This problem may reduce reliability in other translational applications, such as transplantation or drug testing.

3.2. Coculture of Multiple Differentiated Cells

The problem of cell heterogeneity coming from undesirable differentiated cell types often makes the reproducibility of the resulting organoid quite low. To solve this problem, hPSCs can be differentiated separately into specific cell types, sorted, and then mixed to form an organoid. This technique can be a solution for addressing the uncontrollable heterogeneous cell population problems. The hPSCs-derived cellular components can be achieved by performing an in vitro co-culture of multiple PSCs-derived cells or PSCsderived progenitor cells (Figure 1). For example, the specific progenitor cell type can be individually differentiated from PSCs using a monolayer culture. Then, each cell component is purified and co-cultured together to induce their self-organization into a 3D structure. Since the organoid was assembled from an independent and specific hPSCsderived cell type, the occurrence of other cells which are not normally part of the expected organs can be minimized. Additionally, the cellular diversity inside the organoids can be controlled by mixing each cell type based on the actual cellular composition ratio of the original organ [53]. This approach can be further improved by some engineering techniques, such as 3D cell bioprinting or scaffold based-template to partially direct the arrangement of each cellular component inside the organoid structure [54].

4. General Culture Platform to Generate the hPSCs Derived-Organoids

4.1. Matrigel Embedded Technique

As similar as in vivo organ development, extracellular matrix (ECM) remodeling plays an essential role in organoid morphogenesis. This component is very important (mainly for the branching formation that dictates the functional architecture of certain organs, such as the salivary glands, lungs, or kidneys) [55]. Matrigel is an extracellular matrix extracted from Engelbreth-Holm-Swarm mouse sarcomas which are routinely used in various cell culture applications [56]. During organoid differentiation, the Matrigel dome serves as structural support and provides a three-dimensional niche for the organoid assembly [56–59] (Figure 2a).




Figure 2. The common culture platform for the hPSCs derived-organoid. (**a**) Matrigel embedding technique, (**b**) static suspension culture, and (**c**) dynamic suspension culture.

Despite its benefits in organoid culture, Matrigel also has some disadvantages. The variations in biochemical composition and mechanical variability of this sarcoma-derived-ECM complex may result in single- and batch-to-batch variability during the hPSCs-derived organoids culture [60]. Although this approach can be used for developmental study or disease modeling, this method may not be compatible with therapeutic applications. For example, xenogeneic materials need to be eliminated to achieve safe transplantable organoids. These animal-derived components such as Matrigel may potentially induce antigenicity and transmit pathogens such as viruses into the human recipient [61].

Recently, several efforts have been made to find synthetic alternatives for Matrigel [60,62,63]. These alternative materials may provide a chemically defined and xeno-free niche that can be independently tuned by adjusting various variables (such as their composition) and by using polymerization methods to produce expected mechanical and biological properties [60,62]. Moreover, by using a bioengineering approach, these synthetic matrices can be further utilized to provide a specific boundary guiding tissue morphogenesis [64].

Another challenge of the Matrigel embedded method is related to its application in organoid-based transplantation of large organs. Due to its technical difficulties in scaling up, this technique may be limited to small-scale drug screening and disease modeling.

4.2. Static Suspension Culture

A static suspension culture can be used as an alternative method to produce broad quantities of organoid production, from a small scale (e.g., patterned well plates or hanging drop method), medium-scale (e.g., low attachment culture dish [65,66]), or larger scale (e.g., oxygen-permeable static tissue culture bag [67]) (Figure 2b). This approach can be selected when the usage of ECM embedded is not required and when the interference of hydrodynamic conditions needs to be avoided to generate a specific lineage that is sensitive to shear stress (such as in hepatic differentiation) [68].

Traditionally, the formation of a 3D structure relies on spontaneous aggregation in the static suspension culture. This simple method has been adapted to generate hPSCsderived organoids in medium or large-scale organoid productions (for example, using a low attachment culture dish [65,66] or culture bag [67]). The main drawback of this technique is the failure to control EBs aggregation. Due to the epithelial characteristic of hPSCs, spontaneous aggregation often occurs and results in a random size. These conditions may impact the variability of their 3D-differentiated organoid structure. Therefore, several alternative methods have been developed to improve aggregation control, such as the hanging drop method or patterned culture vessel (Figure 2b). The hanging drop method is a technique used to induce cellular assembly by utilizing gravitational forces. In this method, a cell suspension droplet is hung on the reversed surface, and the EBs are formed from the surface tension properties of the culture medium. This method is relatively simple during cell inoculation, but the medium replacement is very challenging due to the limited droplet volume that can be hung on the culture vessel surface [69]. Patterned well plates (e.g., microwell plates, U or V-bottom-well plates) can be used as a tool to force and control cell aggregation [14,70]. By adjusting the cell number, precise control of the EBs size can be enabled, which significantly reduces the variability of the resulting organoids. This culture system has potential for use in drug screening or disease modeling but is difficult to be applied for the therapeutic purposes of large organs that require a large number of organoids (due to the difficulties involved in their scalability).

4.3. Dynamic Suspension Culture

Dynamic suspension culture is a culture technique that is commonly used to produce a medium to a large number of organoids for various applications (such as transplantation of large organs) [61]. A dynamic culture environment can be achieved by using rotating or stirring mechanisms to promote the formation of a spherical 3D formation and control its initial size. The simple platform that is routinely used for organoid generation is by using the rotational shaker or stirred bioreactor (Figure 2c). These types of suspension cultures enable better control of the size and uniformity of organoids population in the culture system. The size control of initial EBs influences the dynamic profile of their pluripotency and lineage specification during differentiation [71]. The spheroid size and uniformity can be controlled by adjusting the inoculation density and manipulating the medium dynamics (such as rotation/stirring speed, fluidic flow, or medium volume) [61,72].

Hydrodynamic condition is an important biomechanical factor involved in dynamic suspension culture that may affect the differentiation process of hPSCs. Despite its importance in controlling agglomeration and medium mixing, the amount of hydrodynamic shear force needs to be carefully considered. Our previous comparison study showed that diverse hydrodynamic conditions coming from the shear force in various culture vessels which were set up (such as rotational cultures, ring-shaped culture vessels, or spinner flasks) may implicate the differentiation tendencies of hPSCs. The results indicated that higher shear stress may improve the hPSC EBs differentiation tendency toward an ectoderm and mesoderm lineage rather than an endoderm lineage [71]. A study conducted

by Wolfe, et al. revealed that the shear force was potentially inducing early germ specification into ectodermal and mesodermal lineage on stress magnitude, ranging from 1.5 to 15 dynes cm⁻² [73]. This result also corresponds well with another study conducted by Vosough et al. which reports several spontaneously differentiated-mesodermal populations from hepatic differentiation (which is presumably exposed by excessive shear stress in the stirred tank bioreactor) [68].

5. Other Limitations and Possible Improvements of PSCs Derived-Organoids Culture *5.1. Maturation*

Currently, several methods for generating a broad spectrum of hPSCs-derived organoids with specific tissue functionality have been reported, such as nephron filtration [47], beating cardiac organoid [74], insulin-secreting organoid [75], or organoids with drug metabolism capability [41]. However, developing efficient techniques for obtaining a mature and functional organoid remains challenging [48]. PSCs-derived cells are still relatively immature and exhibit the expression profile of fetal or progenitor cells rather than adult tissue cell types [47]. Moreover, these cells are still lacking in both architectural and functional characteristics compared to their original adult tissue counterparts. These maturation problems not only impact their reliability in disease modeling or drug screening but also importantly affect the graft immaturity and safety concerns when they are used in regenerative therapy [76].

A functional organ usually takes several years to be fully mature in vivo [77]. However, current studies in hPSCs-derived organoid generation have mostly been performed only within a couple of weeks or several months. Based on this understanding, extending the culture period using a differentiation induction medium may hypothetically improve their physiological and structural maturation. Giandomenico and his colleagues showed a good example of this approach which was used to achieve further maturation [78]. In the prolonged maturation period, they successfully generated cerebral organoids with the telencephalic identity of later stages of neural development, including axon outgrowth and neuronal maturation systems using micropatterned filaments [78]. Another study also showed a similar result in hPSC-derived-cardiac organoids [79]. The extended culture display (for up to 50 days) improved cardiac specification, survival, and maturation, as well as a good response to cardioactive drugs.

The signaling pathways play a vital role in regulating mature phenotypes during organ development and growth. This mechanism involved several biomolecules, such as cytokines, growth factors, hormones, inhibitory molecules, and other small molecules. By using the screening method, we may identify the specific regulator and optimize the in vitro organoid maturation using the same functional molecule compositions during their growth and development [80,81]. The experimental approach carried out by Huang et al. uncover the potential of this method [82]. A combination of biochemical factors, such as thyroid hormone, dexamethasone, and insulin-like growth factor-1 has been shown to improve the maturation of hiPSC-derived cardiac organoids in a micropatterned culture system [82]. In addition to using an exogenous source, these signaling molecules also can be provided by coculturing the hPSCs-derived cells with other cell types that existed in the original organs. These additional cell types may endogenously secrete molecule components required to promote further maturation of hPSC-derived-organoids. One example showed this by incorporating endothelial and mesenchymal stem cells which secrete a specific factor that plays a role in the functional maturation of hPSCs-derived kidney organoids [83].

Metabolism is a critical process that influences the self-renewal and cell fate specification of hPSCs [84–86]. Based on this knowledge, switching the energy source in the culture medium may direct the cellular metabolic adaptation towards a more mature and functional cell type. For example, the increasing concentration of glucose and an adjustment of palmitate concentration in the culture medium can significantly improve the hPSCs derived-cardiac organoid maturation [87]. The hPSCs-derived-organoids required a complex signal that partially can be mimicked in vitro to become fully mature. However, in vivo niches can provide better complex necessary signals similar to their developmental or maturation period. For instance, Dye et al. [88] and Miller, et al. [39] found that the hPSCs derived-lung organoids that were transplanted in vivo showed improved cellular differentiation of secretory lineages and resulted in airway-like structures that were remarkably similar to the native adult human lung.

5.2. High-Cost Production

Different from disease modeling and drug screening purposes, clinical/translational regenerative therapy applications often require a decent amount of cell numbers. For instance, a billion cells are necessary to regenerate both partial and total organs, such as the heart, pancreas, or liver [89–91]. Therefore, the scalability needs to be performed by a biomanufacturing process that involves a high-cost of production. Among all of the culture components, the high-cost consumption mainly comes from the requirement of expensive recombinant proteins, such as growth factors, that essentially require one to induce the differentiation of PSCs into specific and functional cell types.

Some effort has been made to minimize the requirements of these components, such as the development of growth factor-free differentiation protocols using small molecules [92]. However, differentiated cells often show less functional maturation compared to ones differentiated by the inclusion of growth factors. Since the usage of cytokine or growth factors is still becoming a necessity, possible alternative solutions have been developed to optimize the efficient usage of growth factors during differentiation. For example, researchers have used a specific molecular weight cutoff dialysis membrane barrier to recycle and accumulate high molecular weight endogenous and exogenous growth factors while continuously maintaining the exchange between the small molecular weight nutrition and toxic metabolic byproducts [44,93–96]. This approach may significantly improve the production efficiency per unit of hPSCs-derived organoids for various purposes such as drug screening or transplantation.

5.3. Tissue Complexity

As a model system, hPSCs-derived organoids still show incomplete in vivo features compared to the original human tissue, such as tissue architecture similarity, vascularization, innervation, functional immune regulation, and the representation of multi-organ interactions. Consequently, this miniaturized structure often fails to imitate the complex and integrated interplay between each organ system.

During hPSCs derived-organoids culture, several progenitor cells developed from hPSCs communicate via paracrine factors to arrange the tissue formation in vitro. Their self-organization was not only affected by the microenvironment and cellular interaction inside the structure but also by their external culture environment, which is normally difficult to control, resulting in variability within the same batch or between batches [97]. Several attempts have been made to precisely define the structural tissue-like geometry by niche modification [36], which also potentially reduces batch-to-batch variability during organoid construction. For example, a bioengineering approach called "guided self-organization" which utilized a fiber scaffold to control organoid elongation has been reported [19]. This technique showed improvements in increasing tissue complexity, enhancing differentiation, and creating better reproducibility for their applications such as drug screening and disease modeling.

Vascularization plays a vital role in the exchange of various necessary components for tissue homeostases such as nutrition, oxygen, and signaling molecules in actual human tissue [64]. The absence of this vasculature was becoming one of the major limitations in the current development of organoid technology. During the in vitro culture period, an organoid may excessively grow up to a millimeter in size, causing a diffusion gradient of nutrients, oxygen, signaling molecules, and toxic cellular byproducts inside the structure.

The failure to control this spatiotemporal distribution of these molecules may lead to the development of necrotic tissue inside their three-dimensional structure [70,71]. Although neovascularization can be natively achieved during in vivo tissue engraftment [98,99], the generation of the in vitro vascular structure remained difficult. To address this problem, several techniques can be performed to generate vascularized organoids. A simple method to enable vascularization involves performing a co-culture between the pre-differentiated hPSCs or hPSC derived-progenitor cells with endothelial cells before or during the aggregation period to let them self-organize. The earlier development of this method was introduced by Takebe et al. [100,101] by introducing endothelial cells (HUVECs) and mesenchymal stem cells (MSCs) into hPSCs derived-multipotent hepatic endoderm spheroid to generate a vascularized liver buds organoids. A similar technique was also performed by Pettinato et al. [66] using interlaced human adipose microvascular endothelial cells (HAMEC) within the hiPSCs spheroids followed by hepatic differentiate into endothelial cells cells to form a vascular network [102]

To further improve vascularization inside the organoid, a bioengineering approach can be applied using 3D fabrication techniques. A recent study by Skylar-Scott et al. [21] represents a bioengineering approach used to produce a vascularized organoid by controlling both its composition and structure using tissue bioprinting. In this technique, the hPSCs-derived progenitor cells were individually incorporated within a hydrogel and printed using a layer-by-layer deposition to form the expected structure [54,103,104]. Another strategy to enable vascularization involves utilizing a sacrificial networks method using a degradable- or removable biomaterial. The vascularization was formed by seeding endothelial cells or hPSC-derived-endothelial cells in the space created by the voids from the removal of sacrificial networks (making the specific pattern required inside the organoid) [104]. Although all of these 3D bio-fabrication techniques provide enormous potential to guide the 3D organoid assembly, several challenges still need to be addressed. In addition to the poor survival of hPSCs or hPSC-derived progenitor cells, they are very sensitive to environmental cues due to their embryonic-like nature in response to developmental signals. Moreover, they also tend to spontaneously agglomerate and form random clusters after being printed [54].

The other components that are often lacking in current hPSCs-derived organoids designed to mimic actual organs are the neuronal innervation complex and resident immune system. Neuronal innervation is an essential component that regulates various biological processes. For example, acetylcholine (ACh) is a neurotransmitter that triggers cholinergic signaling and plays a pivotal role in the regulation of the small intestine cells [105]. A study conducted by Workman et al. tried to provide a solution to this problem by introducing hPSCs-derived neural crest cells (NCC) to hPSCs-derived intestinal organoids in a static suspension [106]. Inside the organoid structure, the ECC demonstrated migration and formed a microglia structure. This complex structure showed the functional regulation of intestinal contraction and could be utilized to investigate the mechanisms involved in motility disorders of the human gastrointestinal tract [106].

To demonstrate immune regulation for developmental study and disease modeling applications, several important immune components need to exist inside the organoid or in the culture environment. Several attempts have been carried out to mimic the necessary conditions by adding some cytokines to partially imitate in vivo immune regulation [107]. However, this kind of approach cannot completely reflect the complete regulation by actual immune cells since they also secrete several cytokine complexes [108]. Moreover, immune cells also play a critical role during infection (for example, in phagocytosis) [109]. An experimental approach was conducted by Xu et al. to investigate the possibility of immune cell inclusion by incorporating hPSC-derived microglia into brain organoids in ultra-low-attachment 96-well plates [17]. The results showed that the hPSCs derived-microglia exhibit phagocytic activity and synaptic pruning function inside the cerebral organoid.

Additionally, the secretion of neurotrophic factors and cytokines, such as insulin growth factor 1 (IGF1) and IL-1b, was able to modulate neurogenesis in organoid structures [17].

Most current studies have been concentrated on a single hPSCs-derived organoid type, which means that studies about the complex interactions between different organs have not yet been explored much. A few studies have constructed multi-organ organoids or assembloids via the fusing and interconnecting of multiple organoid types that exhibit a specific function (for example, hepato-biliary-pancreatic organoids [50,51] or retinoganglial organoids [110]). Although organ crosstalk can be simulated, this approach is still limited to the utilization of early progenitor organoid types which do not completely reflect fully mature and functional organs. Another problem that requires further investigation is the development of methods that can be used to precisely control the further maturation of individual organoid after fusion. This is necessary since each organoid type requires specific growth factors and a particular cocktail of small molecules, as well as different culture environment conditions, to maintain their characteristics.

6. Conclusions and Future Perspectives

Organoid culture represents a remarkable technology that may provide an ex vivo platform of miniaturized organs for various applications, such as in vitro models of organogenesis, disease modeling, drug screening, and possibly for regenerative therapy. To realize the best utilization of hiPSCs derived-organoids, the selection of the culture system needs to be carefully considered, depending ultimately on the types of application and the effect of the culture environment during the culture period. Despite their improvement through bioengineering approaches, several challenges still need to be addressed to ensure both their reliability to reconstitute actual complex organ functions and their safety for future regenerative therapy.

To achieve decent organoid production for a specific purpose, the mechanisms related to the culture condition during organoid forming and differentiation from hPSCs into each specific organoid type need to be optimized. By using these optimizations, a standardized methodology can be developed for generating organ-specific-hiPSCs-derived organoids.

A personalized hiPSCs-derived organoid can not only overcome individual variability problems in drug screening and disease modeling but also enable immune-compatible grafts for full- or partial organ transplantation when combined with xeno-free organoid production.

Author Contributions: Conceptualization, F.G.T.; writing—original draft preparation, Z.G., D.S.F., N.P.D.P.S. and S.T.N.; writing—review and editing, F.G.T.; scientific illustration, Z.G. and F.G.T.; supervision, F.G.T. All authors have read and agreed to the published version of the manuscript.

Funding: This research received no external funding.

Institutional Review Board Statement: Not applicable.

Informed Consent Statement: Not applicable.

Data Availability Statement: Not applicable.

Conflicts of Interest: The authors declare no conflict of interest.

References

- Cardoso-Moreira, M.; Sarropoulos, I.; Velten, B.; Mort, M.; Cooper, D.N.; Huber, W.; Kaessmann, H. Developmental Gene Expression Differences between Humans and Mammalian Models. *Cell Rep.* 2020, *33*, 108308. [CrossRef]
- Yue, F.; Cheng, Y.; Breschi, A.; Vierstra, J.; Wu, W.; Ryba, T.; Sandstrom, R.; Ma, Z.; Davis, C.; Pope, B.D.; et al. A comparative encyclopedia of DNA elements in the mouse genome. *Nature* 2014, *515*, 355–364. [CrossRef]
- Schutgens, F.; Clevers, H. Human Organoids: Tools for Understanding Biology and Treating Diseases. Annu. Rev. Pathol. Mech. Dis. 2020, 15, 211–234. [CrossRef]
- Dutta, D.; Heo, I.; Clevers, H. Disease Modeling in Stem Cell-Derived 3D Organoid Systems. *Trends Mol. Med.* 2017, 23, 393–410. [CrossRef]

- Iefremova, V.; Manikakis, G.; Krefft, O.; Jabali, A.; Weynans, K.; Wilkens, R.; Marsoner, F.; Brändl, B.; Müller, F.J.; Koch, P.; et al. An Organoid-Based Model of Cortical Development Identifies Non-Cell-Autonomous Defects in Wnt Signaling Contributing to Miller-Dieker Syndrome. *Cell Rep.* 2017, 19, 50–59. [CrossRef]
- Mills, R.J.; Parker, B.L.; Quaife-Ryan, G.A.; Voges, H.K.; Needham, E.J.; Bornot, A.; Ding, M.; Andersson, H.; Polla, M.; Elliott, D.A.; et al. Drug Screening in Human PSC-Cardiac Organoids Identifies Pro-proliferative Compounds Acting via the Mevalonate Pathway. *Cell Stem Cell* 2019, 24, 895–907.e6. [CrossRef]
- 7. Willenbring, H.; Soto-Gutierrez, A. Transplantable liver organoids made from only three ingredients. *Cell Stem Cell* **2013**, *13*, 139–140. [CrossRef]
- 8. Simian, M.; Bissell, M.J. Organoids: A historical perspective of thinking in three dimensions. J. Cell Biol. 2017, 216, 31–40. [CrossRef]
- 9. Berishvili, E.; Casiraghi, F.; Amarelli, C.; Scholz, H.; Piemonti, L.; Berney, T.; Montserrat, N. Mini-organs forum: How to advance organoid technology to organ transplant community. *Transpl. Int.* **2021**, *34*, 1588–1593. [CrossRef]
- Lancaster, M.A.; Knoblich, J.A. Organogenesis in a dish: Modeling development and disease using organoid technologies. *Science* 2014, 345, 1247125. [CrossRef]
- Sasai, Y. Next-generation regenerative medicine: Organogenesis from stem cells in 3D culture. Cell Stem Cell 2013, 12, 520–530. [CrossRef]
- 12. Sasai, Y. Cytosystems dynamics in self-organization of tissue architecture. Nature 2013, 493, 318–326. [CrossRef]
- Bershteyn, M.; Nowakowski, T.J.; Pollen, A.A.; Di Lullo, E.; Nene, A.; Wynshaw-Boris, A.; Kriegstein, A.R. Human iPSC-Derived Cerebral Organoids Model Cellular Features of Lissencephaly and Reveal Prolonged Mitosis of Outer Radial Glia. *Cell Stem Cell* 2017, 20, 435–449. [CrossRef]
- 14. Anastasaki, C.; Wilson, A.F.; Chen, A.S.; Wegscheid, M.L.; Gutmann, D.H. Generation of human induced pluripotent stem cell-derived cerebral organoids for cellular and molecular characterization. *STAR Protoc.* **2022**, *3*, 101173. [CrossRef]
- Anastasaki, C.; Wegscheid, M.L.; Hartigan, K.; Papke, J.B.; Kopp, N.D.; Chen, J.; Cobb, O.; Dougherty, J.D.; Gutmann, D.H. Human iPSC-Derived Neurons and Cerebral Organoids Establish Differential Effects of Germline NF1 Gene Mutations. *Stem Cell Rep.* 2020, 14, 541–550. [CrossRef]
- Wegscheid, M.L.; Anastasaki, C.; Hartigan, K.A.; Cobb, O.M.; Papke, J.B.; Traber, J.N.; Morris, S.M.; Gutmann, D.H. Patientderived iPSC-cerebral organoid modeling of the 17q11.2 microdeletion syndrome establishes CRLF3 as a critical regulator of neurogenesis. *Cell Rep.* 2021, *36*, 109315. [CrossRef]
- 17. Xu, R.; Boreland, A.J.; Li, X.; Erickson, C.; Jin, M.; Atkins, C.; Pang, Z.P.; Daniels, B.P.; Jiang, P. Developing human pluripotent stem cell-based cerebral organoids with a controllable microglia ratio for modeling brain development and pathology. *Stem Cell Rep.* **2021**, *16*, 1923–1937. [CrossRef]
- Qian, X.; Jacob, F.; Song, M.M.; Nguyen, H.N.; Song, H.; Ming, G.L. Generation of human brain region–specific organoids using a miniaturized spinning bioreactor. *Nat. Protoc.* 2018, 13, 565–580. [CrossRef]
- 19. Lancaster, M.A.; Corsini, N.S.; Wolfinger, S.; Gustafson, E.H.; Phillips, A.W.; Burkard, T.R.; Otani, T.; Livesey, F.J.; Knoblich, J.A. Guided self-organization and cortical plate formation in human brain organoids. *Nat. Biotechnol.* **2017**, *35*, 659–666. [CrossRef]
- 20. Mansour, A.A.; Gonçalves, J.T.; Bloyd, C.W.; Li, H.; Fernandes, S.; Quang, D.; Johnston, S.; Parylak, S.L.; Jin, X.; Gage, F.H. An in vivo model of functional and vascularized human brain organoids. *Nat. Biotechnol.* **2018**, *36*, 432–441. [CrossRef]
- Skylar-Scott, M.A.; Huang, J.Y.; Lu, A.; Ng, A.H.M.; Duenki, T.; Liu, S.; Nam, L.L.; Damaraju, S.; Church, G.M.; Lewis, J.A. Orthogonally induced differentiation of stem cells for the programmatic patterning of vascularized organoids and bioprinted tissues. *Nat. Biomed. Eng.* 2022, *6*, 449–462. [CrossRef]
- Watanabe, M.; Buth, J.E.; Vishlaghi, N.; de la Torre-Ubieta, L.; Taxidis, J.; Khakh, B.S.; Coppola, G.; Pearson, C.A.; Yamauchi, K.; Gong, D.; et al. Self-Organized Cerebral Organoids with Human-Specific Features Predict Effective Drugs to Combat Zika Virus Infection. *Cell Rep.* 2017, 21, 517–532. [CrossRef] [PubMed]
- Cakir, B.; Xiang, Y.; Tanaka, Y.; Kural, M.H.; Parent, M.; Kang, Y.J.; Chapeton, K.; Patterson, B.; Yuan, Y.; He, C.S.; et al. Engineering of human brain organoids with a functional vascular-like system. *Nat. Methods* 2019, *16*, 1169–1175. [CrossRef] [PubMed]
- Xiang, Y.; Tanaka, Y.; Patterson, B.; Kang, Y.J.; Govindaiah, G.; Roselaar, N.; Cakir, B.; Kim, K.Y.; Lombroso, A.P.; Hwang, S.M.; et al. Fusion of Regionally Specified hPSC-Derived Organoids Models Human Brain Development and Interneuron Migration. *Cell Stem Cell* 2017, *21*, 383–398.e7. [CrossRef]
- Trisno, S.L.; Philo, K.E.D.; McCracken, K.W.; Catá, E.M.; Ruiz-Torres, S.; Rankin, S.A.; Han, L.; Nasr, T.; Chaturvedi, P.; Rothenberg, M.E.; et al. Esophageal Organoids from Human Pluripotent Stem Cells Delineate Sox2 Functions during Esophageal Specification. *Cell Stem Cell* 2018, 23, 501–515. [CrossRef]
- 26. Koide, T.; Koyanagi-Aoi, M.; Uehara, K.; Kakeji, Y.; Aoi, T. CDX2-induced intestinal metaplasia in human gastric organoids derived from induced pluripotent stem cells. *Iscience* 2022, 25, 104314. [CrossRef]
- McCracken, K.W.; Catá, E.M.; Crawford, C.M.; Sinagoga, K.L.; Schumacher, M.; Rockich, B.E.; Tsai, Y.H.; Mayhew, C.N.; Spence, J.R.; Zavros, Y.; et al. Modelling human development and disease in pluripotent stem-cell-derived gastric organoids. *Nature* 2014, 516, 400–404. [CrossRef]
- McCracken, K.W.; Aihara, E.; Martin, B.; Crawford, C.M.; Broda, T.; Treguier, J.; Zhang, X.; Shannon, J.M.; Montrose, M.H.; Wells, J.M. Wnt/β-catenin promotes gastric fundus specification in mice and humans. *Nature* 2017, 541, 182–187. [CrossRef]

- Broda, T.R.; McCracken, K.W.; Wells, J.M. Generation of human antral and fundic gastric organoids from pluripotent stem cells. *Nat. Protoc.* 2019, 14, 28–50. [CrossRef]
- Nadkarni, R.R.; Abed, S.; Cox, B.J.; Bhatia, S.; Lau, J.T.; Surette, M.G.; Draper, J.S. Functional Enterospheres Derived In Vitro from Human Pluripotent Stem Cells. Stem Cell Rep. 2017, 9, 897–912. [CrossRef]
- Mithal, A.; Capilla, A.; Heinze, D.; Berical, A.; Villacorta-Martin, C.; Vedaie, M.; Jacob, A.; Abo, K.; Szymaniak, A.; Peasley, M.; et al. Generation of mesenchyme free intestinal organoids from human induced pluripotent stem cells. *Nat. Commun.* 2020, 11, 1–15. [CrossRef] [PubMed]
- Romitti, M.; Fonseca, B.d.F.d.; Doumont, G.; Gillotay, P.; Tourneur, A.; Eski, S.E.; Van Simaeys, G.; Chomette, L.; Lasolle, H.; Monestier, O.; et al. Transplantable human thyroid organoids generated from embryonic stem cells to rescue hypothyroidism. *BioRixv*. 2021. BioRxiv 2021.12.01.470729. Available online: https://www.biorxiv.org/content/10.1101/2021.12.01.470729v1 (accessed on 2 February 2022).
- 33. Montel-Hagen, A.; Seet, C.S.; Li, S.; Chick, B.; Zhu, Y.; Chang, P.; Tsai, S.; Sun, V.; Lopez, S.; Chen, H.C.; et al. Organoid-Induced Differentiation of Conventional T Cells from Human Pluripotent Stem Cells. *Cell Stem Cell* **2019**, *24*, 376–389. [CrossRef]
- Drakhlis, L.; Biswanath, S.; Farr, C.M.; Lupanow, V.; Teske, J.; Ritzenhoff, K.; Franke, A.; Manstein, F.; Bolesani, E.; Kempf, H.; et al. Human heart-forming organoids recapitulate early heart and foregut development. *Nat. Biotechnol.* 2021, 39, 737–746. [CrossRef] [PubMed]
- 35. Drakhlis, L.; Devadas, S.B.; Zweigerdt, R. Generation of heart-forming organoids from human pluripotent stem cells. *Nat. Protoc.* **2021**, *16*, 5652–5672. [CrossRef]
- Hoang, P.; Kowalczewski, A.; Sun, S.; Winston, T.S.; Archilla, A.M.; Lemus, S.M.; Ercan-Sencicek, A.G.; Gupta, A.R.; Liu, W.; Kontaridis, M.I.; et al. Engineering spatial-organized cardiac organoids for developmental toxicity testing. *Stem Cell Rep.* 2021, 16, 1228–1244. [CrossRef] [PubMed]
- 37. Voges, H.K.; Mills, R.J.; Elliott, D.A.; Parton, R.G.; Porrello, E.R.; Hudson, J.E. Development of a human cardiac organoid injury model reveals innate regenerative potential. *Development* 2017, 144, 1118–1127. [CrossRef]
- 38. Richards, D.J.; Li, Y.; Kerr, C.M.; Yao, J.; Beeson, G.C.; Coyle, R.C.; Chen, X.; Jia, J.; Damon, B.; Wilson, R.; et al. Human cardiac organoids for the modelling of myocardial infarction and drug cardiotoxicity. *Nat. Biomed. Eng.* **2020**, *4*, 446–462. [CrossRef]
- Miller, A.J.; Hill, D.R.; Nagy, M.S.; Aoki, Y.; Dye, B.R.; Chin, A.M.; Huang, S.; Zhu, F.; White, E.S.; Lama, V.; et al. In Vitro Induction and In Vivo Engraftment of Lung Bud Tip Progenitor Cells Derived from Human Pluripotent Stem Cells. *Stem Cell Reports* 2018, 10, 101–119. [CrossRef]
- 40. Guan, Y.; Xu, D.; Garfin, P.M.; Ehmer, U.; Hurwitz, M.; Enns, G.; Michie, S.; Wu, M.; Zheng, M.; Nishimura, T.; et al. Human hepatic organoids for the analysis of human genetic diseases. *JCI Insight* **2017**, *2*, e94954. [CrossRef]
- Mun, S.J.; Ryu, J.S.; Lee, M.O.; Son, Y.S.; Oh, S.J.; Cho, H.S.; Son, M.Y.; Kim, D.S.; Kim, S.J.; Yoo, H.J.; et al. Generation of expandable human pluripotent stem cell-derived hepatocyte-like liver organoids. *J. Hepatol.* 2019, *71*, 970–985. [CrossRef]
- Ramli, M.N.B.; Lim, Y.S.; Lim, C.T.; Demircioglu, D.; Tng, W.; Gonzales, K.A.U.; Tan, C.P.; Szczerbinska, I.; Liang, H.; Soe, E.L.; et al. Human Pluripotent Stem Cell-Derived Organoids as Models of Liver Disease. *Gastroenterology* 2020, 159, 1471–1486. [CrossRef] [PubMed]
- Takebe, T.; Sekine, K.; Kimura, M.; Yoshizawa, E.; Ayano, S.; Koido, M.; Funayama, S.; Nakanishi, N.; Hisai, T.; Kobayashi, T.; et al. Massive and Reproducible Production of Liver Buds Entirely from Human Pluripotent Stem Cells. *Cell Rep.* 2017, 21, 2661–2670. [CrossRef] [PubMed]
- 44. Torizal, F.G.; Utami, T.; You, L.Q.; Inamura, K.; Nishikawa, M.; Sakai, Y. Generation of high-density human induced pluripotent stem cell derived-liver organoid by enabling growth factors accumulation in a simple dialysis medium refinement culture platform. *Preprint* **2022**, 1–10. [CrossRef]
- 45. Wiedenmann, S.; Breunig, M.; Merkle, J.; von Toerne, C.; Georgiev, T.; Moussus, M.; Schulte, L.; Seufferlein, T.; Sterr, M.; Lickert, H.; et al. Single-cell-resolved differentiation of human induced pluripotent stem cells into pancreatic duct-like organoids on a microwell chip. *Nat. Biomed. Eng.* 2021, *5*, 897–913. [CrossRef] [PubMed]
- 46. Wimmer, R.A.; Leopoldi, A.; Aichinger, M.; Wick, N.; Hantusch, B.; Novatchkova, M.; Taubenschmid, J.; Hämmerle, M.; Esk, C.; Bagley, J.A.; et al. Human blood vessel organoids as a model of diabetic vasculopathy. *Nature* 2019, 565, 505–510. [CrossRef] [PubMed]
- Przepiorski, A.; Sander, V.; Tran, T.; Hollywood, J.A.; Sorrenson, B.; Shih, J.H.; Wolvetang, E.J.; McMahon, A.P.; Holm, T.M.; Davidson, A.J. A Simple Bioreactor-Based Method to Generate Kidney Organoids from Pluripotent Stem Cells. *Stem Cell Rep.* 2018, 11, 470–484. [CrossRef]
- Sander, V.; Przepiorski, A.; Crunk, A.E.; Hukriede, N.A.; Holm, T.M.; Davidson, A.J. Protocol for Large-Scale Production of Kidney Organoids from Human Pluripotent Stem Cells. STAR Protoc. 2020, 1, 100150. [CrossRef]
- Zhang, R.R.; Koido, M.; Tadokoro, T.; Ouchi, R.; Matsuno, T.; Ueno, Y.; Sekine, K.; Takebe, T.; Taniguchi, H. Human iPSC-Derived Posterior Gut Progenitors Are Expandable and Capable of Forming Gut and Liver Organoids. *Stem Cell Rep.* 2018, 10, 780–793. [CrossRef]
- Koike, H.; Iwasawa, K.; Ouchi, R.; Maezawa, M.; Kimura, M.; Kodaka, A.; Nishii, S.; Thompson, W.L.; Takebe, T. Engineering human hepato-biliary-pancreatic organoids from pluripotent stem cells. *Nat. Protoc.* 2021, 16, 919–936. [CrossRef]

- Koike, H.; Iwasawa, K.; Ouchi, R.; Maezawa, M.; Giesbrecht, K.; Saiki, N.; Ferguson, A.; Kimura, M.; Thompson, W.L.; Wells, J.M.; et al. Modelling human hepato-biliary-pancreatic organogenesis from the foregut–midgut boundary. *Nature* 2019, 574, 112–116. [CrossRef]
- 52. Takebe, T.; Wells, J.M. Organoids by design. Science 2019, 364, 956–959. [CrossRef] [PubMed]
- Liu, C.; Feng, X.; Li, G.; Gokulnath, P.; Xiao, J. Generating 3D human cardiac constructs from pluripotent stem cells. *eBioMedicine* 2022, 76, 103813. [CrossRef] [PubMed]
- 54. Salaris, F.; Rosa, A. Construction of 3D in vitro models by bioprinting human pluripotent stem cells: Challenges and opportunities. *Brain Res.* **2019**, 1723, 146393. [CrossRef] [PubMed]
- 55. Kim, H.Y.; Nelson, C.M. Extracellular matrix and cytoskeletal dynamics during branching morphogenesis. *Organogenesis* **2012**, *8*, 56–64. [CrossRef] [PubMed]
- 56. Hughes, C.S.; Postovit, L.M.; Lajoie, G.A. Matrigel: A complex protein mixture required for optimal growth of cell culture. *Proteomics* **2010**, *10*, 1886–1890. [CrossRef]
- Garreta, E.; Kamm, R.D.; Chuva de Sousa Lopes, S.M.; Lancaster, M.A.; Weiss, R.; Trepat, X.; Hyun, I.; Montserrat, N. Rethinking organoid technology through bioengineering. *Nat. Mater.* 2021, 20, 145–155. [CrossRef]
- 58. Mukhopadhyay, M. Recapitulating early cardiogenesis in vitro. Nat. Methods 2021, 18, 331. [CrossRef]
- Lancaster, M.A.; Renner, M.; Martin, C.A.; Wenzel, D.; Bicknell, L.S.; Hurles, M.E.; Homfray, T.; Penninger, J.M.; Jackson, A.P.; Knoblich, J.A. Cerebral organoids model human brain development and microcephaly. *Nature* 2013, 501, 373–379. [CrossRef]
- 60. Aisenbrey, E.A.; Murphy, W.L. Synthetic alternatives to Matrigel. Nat. Rev. Mater. 2020, 5, 539–551. [CrossRef]
- 61. Torizal, F.G.; Horiguchi, I.; Sakai, Y. Physiological Microenvironmental Conditions in Different Scalable Culture Systems for Pluripotent Stem Cell Expansion and Differentiation. *Open Biomed. Eng. J.* **2017**, *13*, 41–54. [CrossRef]
- 62. Heo, J.H.; Kang, D.; Seo, S.J.; Jin, Y. Engineering the Extracellular Matrix for Organoid Culture. *Int. J. Stem Cells* **2022**, *15*, 60–69. [CrossRef] [PubMed]
- Poudel, H.; Sanford, K.; Szwedo, P.K.; Pathak, R.; Ghosh, A. Synthetic Matrices for Intestinal Organoid Culture: Implications for Better Performance. ACS Omega 2022, 7, 38–47. [CrossRef] [PubMed]
- 64. Schneeberger, K.; Spee, B.; Costa, P.; Sachs, N.; Clevers, H.; Malda, J. Converging biofabrication and organoid technologies: The next frontier in hepatic and intestinal tissue engineering? *Biofabrication* **2017**, *9*, 013001. [CrossRef] [PubMed]
- Pettinato, G.; Ramanathan, R.; Fisher, R.A.; Mangino, M.J. Scalable Differentiation of Human iPSCs in a Multicellular Spheroidbased 3D Culture into Hepatocyte- like Cells through Direct Wnt/β—Catenin Pathway Inhibition. *Sci. Rep.* 2016, *6*, 32888. [CrossRef] [PubMed]
- 66. Pettinato, G.; Lehoux, S.; Ramanathan, R.; Salem, M.M.; He, L.X.; Muse, O.; Flaumenhaft, R.; Thompson, M.T.; Rouse, E.A.; Cummings, R.D.; et al. Generation of fully functional hepatocyte-like organoids from human induced pluripotent stem cells mixed with Endothelial Cells. *Sci. Rep.* **2019**, *9*, 8920. [CrossRef]
- Otsuji, T.G.; Bin, J.; Yoshimura, A.; Tomura, M.; Tateyama, D.; Minami, I.; Yoshikawa, Y.; Aiba, K.; Heuser, J.E.; Nishino, T.; et al. A 3D sphere culture system containing functional polymers for large-scale human pluripotent stem cell production. *Stem Cell Rep.* 2014, 2, 734–745. [CrossRef]
- Vosough, M.; Omidinia, E.; Kadivar, M.; Shokrgozar, M.-A.; Pournasr, B.; Aghdami, N.; Baharvand, H. Generation of functional hepatocyte-like cells from human pluripotent stem cells in a scalable suspension culture. *Stem Cells Dev.* 2013, 22, 2693–2705. [CrossRef]
- 69. Velasco, V.; Shariati, S.A.; Esfandyarpour, R. Microtechnology-based methods for organoid models. *Microsyst. Nanoeng.* 2020, 6, 76. [CrossRef]
- 70. Torizal, F.G.; Kimura, K.; Horiguchi, I.; Sakai, Y. Size-dependent hepatic differentiation of human induced pluripotent stem cells spheroid in suspension culture. *Regen. Ther.* **2019**, *16*, 66–73. [CrossRef]
- Torizal, F.G.; Kim, S.M.; Horiguchi, I.; Inamura, K.; Suzuki, I.; Morimura, T.; Nishikawa, M.; Sakai, Y. Production of homogenous size-controlled human induced pluripotent stem cell aggregates using ring-shaped culture vessel. *J. Tissue Eng. Regen. Med.* 2022, 16, 254–266. [CrossRef]
- 72. Fennema, E.; Rivron, N.; Rouwkema, J.; van Blitterswijk, C.; De Boer, J. Spheroid culture as a tool for creating 3D complex tissues. *Trends Biotechnol.* **2013**, *31*, 108–115. [CrossRef] [PubMed]
- Wolfe, R.P.; Leleux, J.; Nerem, R.M.; Ahsan, T. Effects of shear stress on germ lineage specification of embryonic stem cells. *Integr. Biol.* 2012, 4, 1263–1273. [CrossRef]
- 74. Lee, S.J.; Kim, H.A.; Kim, S.J.; Lee, H.A. Improving Generation of Cardiac Organoids from Human Pluripotent Stem Cells Using the Aurora Kinase Inhibitor ZM447439. *Biomedicines* **2021**, *9*, 1952. [CrossRef] [PubMed]
- 75. Kim, Y.; Kim, H.H.; Ko, U.H.; Oh, Y.; Lim, A.; Sohn, J.W.; Shin, J.H.; Kim, H.H.; Han, Y.M. Islet-like organoids derived from human pluripotent stem cells efficiently function in the glucose responsiveness in vitro and in vivo. *Sci. Rep.* **2016**, *6*, 35145. [CrossRef]
- 76. Nam, S.A.; Seo, E.; Kim, J.W.; Kim, H.W.; Kim, H.L.; Kim, K.; Kim, T.M.; Ju, J.H.; Gomez, I.G.; Uchimura, K.; et al. Graft immaturity and safety concerns in transplanted human kidney organoids. *Exp. Mol. Med.* **2019**, *51*, 1–13. [CrossRef] [PubMed]
- 77. Vreeker, A.; Van Stuijvenberg, L.; Hund, T.J.; Mohler, P.J.; Nikkels, P.G.J.; Van Veen, T.A.B. Assembly of the cardiac intercalated disk during preand postnatal development of the human heart. *PLoS ONE* **2014**, *9*, e94722. [CrossRef]
- 78. Giandomenico, S.L.; Sutcliffe, M.; Lancaster, M.A. Generation and long-term culture of advanced cerebral organoids for studying later stages of neural development. *Nat. Protoc.* **2021**, *16*, 579–602. [CrossRef]

- 79. Ergir, E.; Oliver-De La Cruz, J.; Fernandes, S.; Cassani, M.; Niro, F.; Sousa, D.; Vrbský, J.; Vinarský, V.; Perestrelo, A.R.; Debellis, D.; et al. Generation and Maturation of Human iPSC-derived Cardiac Organoids in Long Term Culture. *BioRxiv* 2022, BioRxiv 2022.03.07.483273. [CrossRef]
- Selvaraj, S.; Mondragon-Gonzalez, R.; Xu, B.; Magli, A.; Kim, H.; Lainé, J.; Kiley, J.; McKee, H.; Rinaldi, F.; Aho, J.; et al. Screening identifies small molecules that enhance the maturation of human pluripotent stem cell-derived myotubes. *Elife* 2019, *8*, e47970. [CrossRef] [PubMed]
- 81. Hergenreder, E.; Zorina, Y.; Zhao, Z.; Munguba, H.; Calder, L.; Baggiolini, A.; Minotti, A.P.; Walsh, R.M.; Levitz, J.; Garippa, R.; et al. Combined small molecule treatment accelerates timing of maturation in human pluripotent stem cell-derived neurons. *BioRxiv* 2022. BioRxiv 2022.06.02.494616. [CrossRef]
- Huang, C.Y.; Peres Moreno Maia-Joca, R.; Ong, C.S.; Wilson, I.; DiSilvestre, D.; Tomaselli, G.F.; Reich, D.H. Enhancement of human iPSC-derived cardiomyocyte maturation by chemical conditioning in a 3D environment. *J. Mol. Cell. Cardiol.* 2020, 138, 1–11. [CrossRef] [PubMed]
- Khoshdel-Rad, N.; Zahmatkesh, E.; Moeinvaziri, F.; Haghparast, N.; Baharvand, H.; Aghdami, N.; Moghadasali, R. Promoting Maturation of Human Pluripotent Stem Cell-Derived Renal Microtissue by Incorporation of Endothelial and Mesenchymal Cells. *Stem Cells Dev.* 2021, 30, 428–440. [CrossRef] [PubMed]
- Tsogtbaatar, E.; Landin, C.; Minter-Dykhouse, K.; Folmes, C.D.L. Energy Metabolism Regulates Stem Cell Pluripotency. *Front. Cell Dev. Biol.* 2020, *8*, 87. [CrossRef] [PubMed]
- Zhang, J.; Nuebel, E.; Daley, G.Q.; Koehler, C.M.; Teitell, M.A. Metabolic regulation in pluripotent stem cells during reprogramming and self-renewal. *Cell Stem Cell* 2012, 11, 589–595. [CrossRef]
- 86. Wu, J.; Ocampo, A.; Belmonte, J.C.I. Cellular Metabolism and Induced Pluripotency. Cell 2016, 166, 1371–1385. [CrossRef]
- Mills, R.J.; Titmarsh, D.M.; Koenig, X.; Parker, B.L.; Ryall, J.G.; Quaife-Ryan, G.A.; Voges, H.K.; Hodson, M.P.; Ferguson, C.; Drowley, L.; et al. Functional screening in human cardiac organoids reveals a metabolic mechanism for cardiomyocyte cell cycle arrest. *Proc. Natl. Acad. Sci. USA* 2017, 114, E8372–E8381. [CrossRef]
- Dye, B.R.; Dedhia, P.H.; Miller, A.J.; Nagy, M.S.; White, E.S.; Shea, L.D.; Spence, J.R. A bioengineered niche promotes in vivo engraftment and maturation of pluripotent stem cell derived human lung organoids. *Elife* 2016, *5*, e19732. [CrossRef]
- 89. Kempf, H.; Andree, B.; Zweigerdt, R. Large-scale production of human pluripotent stem cell derived cardiomyocytes. *Adv. Drug Deliv. Rev.* 2016, *96*, 18–30. [CrossRef]
- 90. Zweigerdt, R. Large Scale Production of Stem Cells and Their Derivatives. *Adv. Biochem. Eng. Biotechnol.* 2009, 123, 127–141. [CrossRef]
- 91. Iansante, V.; Chandrashekran, A.; Dhawan, A. Cell-based liver therapies: Past, present and future. *Philos. Trans. R. Soc. B Biol. Sci.* 2018, 373, 20170229. [CrossRef]
- Siller, R.; Greenhough, S.; Naumovska, E.; Sullivan, G.J. Small-molecule-driven hepatocyte differentiation of human pluripotent stem cells. *Stem Cell Rep.* 2015, *4*, 939–952. [CrossRef] [PubMed]
- Torizal, F.G.; Lau, Q.Y.; Ibuki, M.; Kawai, Y.; Horikawa, M.; Minami, M.; Michiue, T.; Horiguchi, I.; Nishikawa, M.; Sakai, Y. A miniature dialysis-culture device allows high-density human-induced pluripotent stem cells expansion from growth factor accumulation. *Commun. Biol.* 2021, 4, 1316. [CrossRef] [PubMed]
- 94. Torizal, F.G.; Lau, Q.Y.; Ibuki, M.; Kawai, Y.; Horikawa, M.; Minami, M.; Horiguchi, I.; Nishikawa, M.; Sakai, Y. High-density hiPSCs expansion supported by growth factors accumulation in a simple dialysis- culture platform. *Preprint* 2020, 1–17. [CrossRef]
- Choi, H.; Torizal, F.G.; Shinohara, M.; Sakai, Y. Differentiation of Human Induced Pluripotent Stem Cells into Definitive Endoderm Using Simple Dialysis Culture Device. *Methods Mol. Biol.* 2021, 2158, 141–153. [CrossRef]
- 96. Miki, T.; Ring, A.; Gerlach, J.; Ph, D.; Ring, A.; Ph, D.; Gerlach, J. Hepatic differentiation of human embryonic stem cells is promoted by three-dimensional dynamic perfusion culture conditions. *Tissue Eng. Part C. Methods* **2011**, *17*, 557–568. [CrossRef]
- Velasco, S.; Kedaigle, A.J.; Simmons, S.K.; Nash, A.; Rocha, M.; Quadrato, G.; Paulsen, B.; Nguyen, L.; Adiconis, X.; Regev, A.; et al. Individual brain organoids reproducibly form cell diversity of the human cerebral cortex. *Nature* 2019, 570, 523–527. [CrossRef]
- 98. Daviaud, N.; Friedel, R.H.; Zou, H. Vascularization and engraftment of transplanted human cerebral organoids in mouse cortex. *Eneuro* **2018**, *5*, 0219-18. [CrossRef]
- 99. Garreta, E.; Prado, P.; Tarantino, C.; Oria, R.; Fanlo, L.; Martí, E.; Zalvidea, D.; Trepat, X.; Roca-Cusachs, P.; Gavaldà-Navarro, A.; et al. Fine tuning the extracellular environment accelerates the derivation of kidney organoids from human pluripotent stem cells. *Nat. Mater.* **2019**, *18*, 397–405. [CrossRef]
- 100. Takebe, T.; Sekine, K.; Enomura, M.; Koike, H.; Kimura, M.; Ogaeri, T.; Zhang, R.R.; Ueno, Y.; Zheng, Y.W.; Koike, N.; et al. Vascularized and functional human liver from an iPSC-derived organ bud transplant. *Nature* **2013**, *499*, 481–484. [CrossRef]
- Takebe, T.; Zhang, R.-R.; Koike, H.; Kimura, M.; Yoshizawa, E.; Enomura, M.; Koike, N.; Sekine, K.; Taniguchi, H. Generation of a vascularized and functional human liver from an iPSC-derived organ bud transplant. *Nat. Protoc.* 2014, *9*, 396–409. [CrossRef]
- 102. Wörsdörfer, P.; Dalda, N.; Kern, A.; Krüger, S.; Wagner, N.; Kwok, C.K.; Henke, E.; Ergün, S. Generation of complex human organoid models including vascular networks by incorporation of mesodermal progenitor cells. *Sci. Rep.* 2019, *9*, 15663. [CrossRef] [PubMed]
- Zhang, X.; Tang, L.; Yi, Q. Engineering the vasculature of stem-cell-derived liver organoids. *Biomolecules* 2021, 11, 966. [CrossRef]
 [PubMed]
- 104. Grebenyuk, S.; Ranga, A. Engineering organoid vascularization. Front. Bioeng. Biotechnol. 2019, 7, 39. [CrossRef] [PubMed]

- 105. Takahashi, T.; Fujishima, K.; Kengaku, M. Modeling intestinal stem cell function with organoids. Int. J. Mol. Sci. 2021, 22, 10912. [CrossRef] [PubMed]
- 106. Workman, M.J.; Mahe, M.M.; Trisno, S.; Poling, H.M.; Watson, C.L.; Sundaram, N.; Chang, C.F.; Schiesser, J.; Aubert, P.; Stanley, E.G.; et al. Engineered human pluripotent-stem-cell-derived intestinal tissues with a functional enteric nervous system. *Nat. Med.* 2017, 23, 49–59. [CrossRef]
- 107. Jung, K.B.; Lee, H.; Son, Y.S.; Lee, M.O.; Kim, Y.D.; Oh, S.J.; Kwon, O.; Cho, S.; Cho, H.S.; Kim, D.S.; et al. Interleukin-2 induces the in vitro maturation of human pluripotent stem cell-derived intestinal organoids. *Nat. Commun.* **2018**, *9*, 3039. [CrossRef]
- 108. Stanley, A.C.; Lacy, P. Pathways for cytokine secretion. Physiology 2010, 25, 218–229. [CrossRef]
- 109. Rosales, C.; Uribe-Querol, E. Phagocytosis: A Fundamental Process in Immunity. Biomed Res. Int. 2017, 2017, 9042851. [CrossRef]
- Fligor, C.M.; Lavekar, S.S.; Harkin, J.; Shields, P.K.; VanderWall, K.B.; Huang, K.C.; Gomes, C.; Meyer, J.S. Extension of retinofugal projections in an assembled model of human pluripotent stem cell-derived organoids. *Stem Cell Rep.* 2021, 16, 2228–2241. [CrossRef]