

**Supplemental Table S1:** The measurement precision of aCL, aPA, aPC, and aPS (IgG and IgM) ELISAs

	aCL IgG ELISA	aCL IgM ELISA	aPA IgG ELISA	aPA IgM ELISA	aPC IgG ELISA	aPC IgM ELISA	aPS IgG ELISA	aPS IgM ELISA
<i>Within-run CV (%)</i>	8.80	5.40	5.00	7.81	2.39	4.45	4.39	5.92
<i>Between-run CV (%)</i>	17.30	17.20	11.35	12.31	7.20	6.55	15.35	4.09