

TB32047 Clones	Indel Size	PANC1 Clones	Indel Size
1-9	Delete 1bp	1-1	Delete 5 bp
2-9	Insert 1bp/ Delete 92bp	1-10	Insert 265bp
2-10	Delete 2bp/ Delete 26 bp	2-4	Delete 11 bp

TB32047-WT or NC			
NC_000075.7 37155-37174bp			
5'-	GAAGCGGCACGTCTGAGTGG	-3'	
3'-	CTTCGCCGTGCAGACTCACC	-5'	

TB32047-SC1-9			
Mutation: Deletion of 1 bp			
5'-	GAA-CGGCACGTCTGAGTGG	-3'	
3'-	CTT-GCCGTGCAGACTCACC	-5'	

TB32047-WT or NC			
NC_000075.7 24759-24778bp			
5'-	GATATCTGATGTGCTAGAGG	-3'	
3'-	CTATAGACTACACGATCTCC	-5'	

TB32047-SC2-9			
Mutation 1: Insertion of 1bp			
5'-	GATA ^A TCTGATGTGCTAGAGG	-3'	
3'-	CTAT ^T AGACTACACGATCTCC	-5'	

TB32047-SC2-10			
Mutation 2: Deletion of 2bp			
5'-	GATAT--GATGTGCTAGAGG	-3'	
3'-	CTATA--CTACACGATCTCC	-5'	

TB32047-SC2-9			
Mutation 3: Deletion of 92bp			
TB32047-SC2-10			
Mutation 4: Deletion of 26bp			

PANC-1-WT or NC
NC_000006.12 374-393bp

5'-	CGCCCTGCGCTACGATGCCA	-3'
3'-	GCGGGACGCGATGCTACGGT	-5'

PANC-1-SC1-1
Mutation: Deletion of 5 bp

5'-	CGCCCTGCGCTACGATG---	-3'
3'-	GCGGGACGCGATGCTAC---	-5'

PANC-1-SC1-10
Mutation: Insertion of 256 bp

5'-	CGCCCTGCGCTACGATG-CCA	-3'
3'-	GCGGGACGCGATGCTAC-GGT	-5'

PANC-1-WT or NC
NC_000006.12 21295-21314bp

5'-	CCCACCACGACGTCCACAC	-3'
3'-	GGGTGGTGCTGCAGGTGTG	-5'

PANC-1-SC2-4
Mutation: Deletion of 11 bp

5'-	C C C ----- C A C A C	-3'
3'-	G G G ----- G T G T G	-5'
