

Supplementary Materials

Supplementary Table S5:

Different mediums and reagents used in this study

Cell type	Basic medium components	Additional molecules/components	Reason	Reference
RPE	Growth medium: α -MEM 2 mmol/L L-glutamine 1X N1 supplement 1X non-essential amino acids 0.25 mg/mL taurine 0.013 μ g/L Triiodo-thyronine 20 μ g/L Hydrocortisone 1X penicillin-streptomycin 5% fetal bovine serum			[1]
CiPCs (D1-D10)	Photoreceptor induction medium: DMEM/F12 10% knockout serum replacement 2% B27 5 ng/mL Noggin 5 ng/mL IGF-1	0.5 mM VPA 4.8 μ M CHIR99021 2 μ M Repsox 10 μ M Forskolin 10 μ M IWR1 3 nM Sonic hedgehog 100 μ M Taurine 1 μ M Retinoic acid	1.VPA: facilitates photoreceptor differentiation; promotes DNA demethylation; epithelial-mesenchymal transition (EMT) inhibition; regulation of glutamate signaling pathway 2.CHIR99021: promotes neuronal reprogramming 3.Repsox: promotes neuronal reprogramming; assisting mesenchymal-epithelial transition (MET) 4.Forskolin: promotes neuronal reprogramming 5.IWR1: promotes photoreceptor differentiation 6.Sonic hedgehog: facilitates rod photoreceptor differentiation 7.Taurine: promotes retinal cell differentiation 8.Retinoic acid: promotes photoreceptor differentiation	VPA[2-5] CHIR99021[6, 7] Repsox [8, 9] Forskolin [10, 11] IWR1[12] Sonic hedgehog [13] Taurine [14] Retinoic acid [15]

CiPCs (D11~)	Photoreceptor differentiation medium: DMEM/F12 2% B27 1% N2 10 ng/mL Noggin 10 ng/mL IGF1 5 ng/mL bFGF	0.5 mM VPA 4.8 μ M CHIR99021 2 μ M Repsox 10 μ M Forskolin 10 μ M IWR1	9.bFGF: promotes retinal regeneration from RPE	bFGF [16, 17]
HDF	Growth medium: DMEM with glutamine 1X non-essential amino acids 1X penicillin-streptomycin 10% fetal bovine serum			[12]

Reference

1. Maminishkis, A., S. Chen, S. Jalickee, T. Banzon, G. Shi, F. E. Wang, T. Ehalt, J. A. Hammer, and S. S. Miller. "Confluent Monolayers of Cultured Human Fetal Retinal Pigment Epithelium Exhibit Morphology and Physiology of Native Tissue." *Invest Ophthalmol Vis Sci* 47, no. 8 (2006): 3612-24.
2. Mello, M. L. S. "Sodium Valproate-Induced Chromatin Remodeling." *Front Cell Dev Biol* 9 (2021): 645518.
3. Chen, L., A. Alam, A. Pac-Soo, Q. Chen, Y. Shang, H. Zhao, S. Yao, and D. Ma. "Pretreatment with Valproic Acid Alleviates Pulmonary Fibrosis through Epithelial-Mesenchymal Transition Inhibition in Vitro and in Vivo." *Lab Invest* 101, no. 9 (2021): 1166-75.
4. Osakada, F., S. Ooto, T. Akagi, M. Mandai, A. Akaike, and M. Takahashi. "Wnt Signaling Promotes Regeneration in the Retina of Adult Mammals." *J Neurosci* 27, no. 15 (2007): 4210-9.
5. Park, H. J., W. S. Kang, J. W. Paik, and J. W. Kim. "Effect of Valproic Acid through Regulation of Nmda Receptor-Erk Signaling in Sleep Deprivation Rats." *J Mol Neurosci* 47, no. 3 (2012): 554-8.
6. Yang, Y., R. Chen, X. Wu, Y. Zhao, Y. Fan, Z. Xiao, J. Han, L. Sun, X. Wang, and J. Dai. "Rapid and Efficient Conversion of Human Fibroblasts into Functional Neurons by Small Molecules." *Stem Cell Reports* 13, no. 5 (2019): 862-76.
7. Ladewig, J., J. Mertens, J. Kesavan, J. Doerr, D. Poppe, F. Glaue, S. Herms, P. Wernet, G. Kogler, F. J. Muller, P. Koch, and O. Brustle. "Small Molecules Enable Highly Efficient Neuronal Conversion of Human Fibroblasts." *Nat Methods* 9, no. 6 (2012): 575-8.
8. Hu, W., B. Qiu, W. Guan, Q. Wang, M. Wang, W. Li, L. Gao, L. Shen, Y. Huang, G. Xie, H. Zhao, Y. Jin, B. Tang, Y. Yu, J. Zhao, and G. Pei. "Direct Conversion of Normal and Alzheimer's Disease Human Fibroblasts into Neuronal Cells by Small Molecules." *Cell Stem Cell* 17, no. 2 (2015): 204-12.
9. Liu, X., H. Sun, J. Qi, L. Wang, S. He, J. Liu, C. Feng, C. Chen, W. Li, Y. Guo, D. Qin, G. Pan, J. Chen, D. Pei, and H. Zheng. "Sequential Introduction of Reprogramming Factors Reveals a Time-Sensitive Requirement for Individual Factors and a Sequential Emt-Met Mechanism for Optimal Reprogramming." *Nat Cell Biol* 15, no. 7 (2013): 829-38.
10. Li, X., X. Zuo, J. Jing, Y. Ma, J. Wang, D. Liu, J. Zhu, X. Du, L. Xiong, Y. Du, J. Xu, X. Xiao, J. Wang, Z. Chai, Y. Zhao, and H. Deng. "Small-Molecule-Driven Direct Reprogramming of Mouse Fibroblasts into Functional Neurons." *Cell Stem Cell* 17, no. 2 (2015): 195-203.
11. Liu, M. L., T. Zang, Y. Zou, J. C. Chang, J. R. Gibson, K. M. Huber, and C. L. Zhang. "Small Molecules Enable Neurogenin 2 to Efficiently Convert Human Fibroblasts into Cholinergic Neurons." *Nat Commun* 4 (2013): 2183.
12. Mahato, B., K. D. Kaya, Y. Fan, N. Sumien, R. A. Shetty, W. Zhang, D. Davis, T. Mock, S. Batabyal, A. Ni, S. Mohanty, Z. Han, R. Farjo, M. J. Forster, A. Swaroop, and S. H. Chavala. "Pharmacologic Fibroblast Reprogramming into Photoreceptors Restores Vision." *Nature* 581, no. 7806 (2020): 83-88.
13. Levine, E. M., H. Roelink, J. Turner, and T. A. Reh. "Sonic Hedgehog Promotes Rod Photoreceptor Differentiation in Mammalian Retinal Cells in Vitro." *J Neurosci* 17, no. 16 (1997): 6277-88.
14. Altshuler, D., J. J. Lo Turco, J. Rush, and C. Cepko. "Taurine Promotes the Differentiation of a Vertebrate Retinal Cell Type in Vitro." *Development* 119, no. 4 (1993): 1317-28.

15. Kelley, M. W., J. K. Turner, and T. A. Reh. "Retinoic Acid Promotes Differentiation of Photoreceptors in Vitro." *Development* 120, no. 8 (1994): 2091-102.
16. Pittack, C., M. Jones, and T. A. Reh. "Basic Fibroblast Growth Factor Induces Retinal Pigment Epithelium to Generate Neural Retina in Vitro." *Development* 113, no. 2 (1991): 577-88.
17. Luz-Madrigal, A., E. Grajales-Esquivel, A. McCorkle, A. M. DiLorenzo, K. Barbosa-Sabanero, P. A. Tsonis, and K. Del Rio-Tsonis. "Reprogramming of the Chick Retinal Pigmented Epithelium after Retinal Injury." *BMC Biol* 12 (2014): 28.