





Supplementary Figure S1. Confirmation of the purified mAbs.

MAbs (2 μ g) were treated with sodium dodecyl sulfate (SDS) sample buffer containing 2-mercaptoethanol (Nacalai Tesque, Inc.). Proteins were separated on 5%–20% polyacrylamide gel (FUJIFILM Wako). The gel was stained by Bio-Safe CBB G-250 Stain (Bio-Rad Laboratories, Inc.).



Supplementary Figure S2. Body appearance in LN229/PDPN (A), PC-10 (B), and LN319 (C) xenografts-implanted mice treated with control mIgG or PMab-117-mG_{2a} on day 30 (LN229/PDPN) , day 28 (PC-10) or day 22 (LN319) after cell inoculation. Scale bar, 1 cm.