


Supplementary Figure



Supplementary Figure S1. Confirmation of the PCR products corresponding to the data presented in Figure 3. PCR products were amplified from cDNA using Taq polymerase to verify each primer set, separated on a 2% agarose gel, and visualized with SYBR Safe to confirm the expected amplicon sizes. **(a)** Agarose gel image for *C9ORF72* variant 1 (V1). **(b)** Agarose gel image for *C9ORF72* variant 2 (V2). **(c)** Agarose gel image for *C9ORF72* variant 3 (V3). **(d)** Agarose gel image for the *C9ORF72* V1/V3 intron-1 amplicon.