

Supplementary Table S1. Primers used in quantitative real-time PCR and breakpoint analysis

Primer ID	Sequence (5' to 3')
PMP22-QE1-F:	TGGAAGAAGGGGTTACGCTGTT
PMP22-QE1-R:	GGGGGTGGGAAATTAGGCAATT
PMP22-QE2-F:	CTGGTGCTGCTGTTCGTCTC
PMP22-QE2-R:	AACACAGTCCTGAACCAGCA
PMP22-QE3-F:	AGCAATGGATCGTGGGCAA
PMP22-QE3-R:	CAAAACCAGCCTCACCGTTT
PMP22-QE4-F:	GCCACCATGATCCTGTTCGAT
PMP22-QE4-R:	CATTCCGCAGACTTGGATGC
PMP22-QE5-F:	GAGTGGCATCTCAACTCGGAT
PMP22-QE5-R:	CGCGTTTCCGCAAGATCAC
DMD-63-a	CATTCACAGAAATGCAGGCT
DMD-63-b	GTGTGAAACCAAGACAACCC
DMD-53-c	AGTTACCTGGGGAAACTAGC
DMD-53-d	GGCAGCTATTTGAGTGTCTC
DMD-55-e	CCCACTGAGAACTCCACTTGA
DMD-55-f	CTGTGAAGAAAGGCATTGGTAG
DMD-50-g	TTCCATGAAGGCAGTAGGGG
DMD-50-h	AAGCCTTCCATTTTGAGCGG
DMD-3'-i	GGATAAGGAATGGGAGAACTGGT
DMD-3'-j	CAATAGCAGGGAAATGGTGTGC
DMD-55-k	TGCCCAGAGAAATTGTTGAGA
DMD-55-l	TTCCCTTTTCGGGCCAGTAT