



Figure S1. Cytokine release by UniCAR T-cells armed with anti-CD10 TM analyzed using MACSplex technology. UniCAR T-cells were co-cultured for 7 hrs with Nalm-6 Luc cells at 5:1 E:T ratio in the absence or presence of 50 nM of anti-CD10 TM or anti-CD19 TM. The cell-free supernatants were then collected and the cytokines were analyzed by MACSplex cytokine kit. Data is shown as mean \pm SD of three independent donors (* p < 0.0332, ** p < 0.0021, *** p < 0.0002, **** p < 0.0001; comparison to sample w/o TM; One-way ANOVA with Dunnett's multiple comparison test).