

Table S2. The mean and standard deviation of contrast sensitivity under photopic condition, mesopic condition without glare, and mesopic condition with glare at postoperative 3 months.

Photopic condition	3cpd	6cpd	12cpd	18cpd
PanOptix	1.75 ± 0.06	1.95 ± 0.14	1.57 ± 0.15	1.03 ± 0.21
FineVision	1.79 ± 0.08	1.97 ± 0.21	1.67 ± 0.23	1.14 ± 0.18
Symfony	1.69 ± 0.16	1.91 ± 0.16	1.65 ± 0.18	1.10 ± 0.19
MiniWell	1.76 ± 0.10	1.96 ± 0.19	1.73 ± 0.25	1.12 ± 0.20
Mesopic condition without glare	3cpd	6cpd	12cpd	18cpd
PanOptix	1.49 ± 0.18	1.55 ± 0.24	1.08 ± 0.36	0.64 ± 0.32
FineVision	1.60 ± 0.18	1.49 ± 0.18	1.07 ± 0.33	0.54 ± 0.25
Symfony	1.49 ± 0.20	1.41 ± 0.37	0.93 ± 0.29	0.44 ± 0.28
MiniWell	1.51 ± 0.14	1.60 ± 0.21	1.20 ± 0.56	0.70 ± 0.44
Mesopic condition with glare	3cpd	6cpd	12cpd	18cpd
PanOptix	1.69 ± 0.08	1.89 ± 0.20	1.43 ± 0.24	1.00 ± 0.21
FineVision	1.67 ± 0.14	1.78 ± 0.15	1.49 ± 0.17	0.95 ± 0.27
Symfony	1.60 ± 0.20	1.76 ± 0.19	1.43 ± 0.31	0.93 ± 0.27
MiniWell	1.73 ± 0.15	1.88 ± 0.22	1.56 ± 0.31	1.04 ± 0.26

cpd=cycle per degree