

Table S1. The mean and standard deviation of corrected binocular visual acuity measured for defocus curve at postoperative 3 months.

Diopter	+2.0	+1.5	+1.0	+0.5
PanOptix	0.42 ± 0.11	0.27 ± 0.15	0.12 ± 0.13	-0.02 ± 0.13
FineVision	0.15 ± 0.15	0.07 ± 0.11	-0.02 ± 0.09	-0.07 ± 0.11
Symfony	0.31 ± 0.23	0.17 ± 0.20	0.06 ± 0.15	-0.05 ± 0.08
MiniWell	0.13 ± 0.19	0.05 ± 0.16	0.00 ± 0.15	-0.05 ± 0.07
Diopter	0.0	-0.5	-1.0	-1.5
PanOptix	-0.11 ± 0.09	-0.02 ± 0.09	0.00 ± 0.13	-0.06 ± 0.12
FineVision	-0.09 ± 0.10	-0.05 ± 0.10	0.00 ± 0.11	-0.01 ± 0.11
Symfony	-0.10 ± 0.07	-0.11 ± 0.09	-0.06 ± 0.12	-0.01 ± 0.13
MiniWell	-0.09 ± 0.07	-0.09 ± 0.06	-0.06 ± 0.07	-0.07 ± 0.07
Diopter	-2.0	-2.5	-3.0	-3.5
PanOptix	-0.04 ± 0.10	0.02 ± 0.09	0.14 ± 0.13	0.24 ± 0.16
FineVision	-0.05 ± 0.11	0.00 ± 0.12	0.09 ± 0.11	0.24 ± 0.14
Symfony	0.02 ± 0.12	0.11 ± 0.11	0.25 ± 0.14	0.38 ± 0.15
MiniWell	-0.02 ± 0.08	0.05 ± 0.10	0.12 ± 0.11	0.26 ± 0.14