

Supplementary Files

Phosphodiesterase 5 and Arginase Inhibitory Activities of the Extracts from Some Members of Nelumbonaceae and Nymphaeaceae Families

Teerapap Panklai ^{1,2}, Nungruthai Suphrom ³, Prapapan Temkitthawon ¹, Perle Totoson ², Krongkarn Chootip ⁴, Xiao-Liang Yang ⁵, Hui-Ming Ge ⁶, Zhu-Jun Yao ⁵, Nattiya Chaichamnong ⁷, Kornkanok Ingkaninan ^{1,*} and Corine Girard ²

- ¹ Center of Excellence in Cannabis Research, Faculty of Pharmaceutical Sciences and Center of Excellence for Innovation in Chemistry, Naresuan University, Phitsanulok 65000, Thailand; teerapapp62@nu.ac.th (T.P.); prapapantem@gmail.com (P.T.)
- ² Université de Franche-Comté, PEPITE, Besançon 25000, France; perle.totoson@univ-fcomte.fr (P.T.); corine.girard@univ-fcomte.fr (C.G.)
- ³ Department of Chemistry, Faculty of Science and Center of Excellence for Innovation in Chemistry, Naresuan University, Phitsanulok 65000, Thailand; suphrom.n1@gmail.com
- ⁴ Department of Physiology, Faculty of Medical Science and Center of Excellence for Innovation in Chemistry, Naresuan University, Phitsanulok 65000, Thailand; krongkarnc@nu.ac.th
- ⁵ State Key Laboratory of Coordination Chemistry, Jiangsu Key Laboratory of Advance Organic Materials, School of Chemistry and Chemical Engineering, Nanjing University, Nanjing 210023, China; yxlnmr@nju.edu.cn (X.-L.Y.); yaoz@nju.edu.cn (Z.-J.Y.)
- ⁶ State Key Laboratory of Pharmaceutical Biotechnology, Institute of Functional Biomolecules, School of Life Sciences, Nanjing University, Nanjing 210023, China; hmge@nju.edu.cn
- ⁷ Division of Applied Thai Traditional Medicine, Faculty of Public Health, Naresuan University, Phitsanulok 65000, Thailand; orangejussamine@hotmail.com
- * Correspondence: k_inkanan@yahoo.com; Tel.: +66-81-4817350

Supplementary Materials

Table S1: Calibration data, LOD, and LOQ of the flavonoids **1-6** analyzed by HPLC.

Table S2: Intra- and Inter-day precision and accuracy of flavonoids **1-6** analyzed by HPLC.

Table S1. Calibration data, LOD, and LOQ of the flavonoids **1–6** analyzed by HPLC.

Flavonoids	Linearity range ($\mu\text{g}/\text{ml}$)	Regression equation	Correlation coefficient (r^2)	LOD ($\mu\text{g}/\text{ml}$)	LOQ ($\mu\text{g}/\text{ml}$)
1	0.5 - 100	$y = 23100x - 4376.5$	0.9999	0.09	0.5
2	5.0 - 400	$y = 20691x + 69562$	0.9997	0.85	5.0
3	0.5 - 100	$y = 64186x - 48477$	0.9997	0.15	0.5
4	0.5 - 100	$y = 27499x - 5199$	0.9999	0.05	0.5
5	0.5 - 100	$y = 65813x - 10389$	0.9999	0.09	0.5
6	0.5 - 100	$y = 44580x + 4643.1$	0.9999	0.04	0.5

Table S2. Intra- and Inter-day precision and accuracy of flavonoids **1–6** analyzed by HPLC.

Flavonoids	Concentration levels ($\mu\text{g}/\text{ml}$)	Intra-day ($n=3$)		Inter-day ($n=9$)		Accuracy ($n=3$)	
		Measured concentration ($\mu\text{g}/\text{ml}) \pm \text{SD}$	Precision (%RSD)	Measured concentration ($\mu\text{g}/\text{ml}) \pm \text{SD}$	Precision (%RSD)	Recovery (%) Mean \pm SD	Precision (%RSD)
1	3	2.71 \pm 0.003	0.13	2.74 \pm 0.05	1.83	90.04 \pm 1.87	2.08
	25	23.63 \pm 0.56	2.37	25.04 \pm 0.73	3.00	86.69 \pm 1.47	1.69
	65	62.94 \pm 1.01	1.60	64.09 \pm 1.04	1.63	90.05 \pm 2.54	2.82
2	15	13.08 \pm 0.24	1.86	13.01 \pm 0.28	2.12	85.13 \pm 1.53	1.79
	75	75.90 \pm 1.23	1.62	77.15 \pm 1.98	2.57	94.82 \pm 0.73	0.77
	300	309.05 \pm 1.41	0.46	311.04 \pm 8.65	2.78	100.97 \pm 2.44	2.42
3	3	2.94 \pm 0.06	2.06	2.98 \pm 0.09	2.85	106.17 \pm 1.28	1.20
	25	27.01 \pm 0.24	0.89	27.09 \pm 0.61	2.25	98.46 \pm 1.49	1.51
	65	66.27 \pm 1.15	1.74	66.45 \pm 1.64	2.47	109.92 \pm 0.91	0.82
4	3	2.83 \pm 0.07	2.61	2.87 \pm 0.06	2.17	86.21 \pm 2.40	2.78
	25	23.96 \pm 0.56	2.33	24.48 \pm 0.60	2.43	86.14 \pm 2.06	2.40
	65	63.38 \pm 1.26	1.99	64.74 \pm 1.57	2.42	90.23 \pm 0.57	0.63
5	3	2.58 \pm 0.04	1.63	2.65 \pm 0.06	2.18	103.43 \pm 0.72	0.69
	25	25.47 \pm 0.41	1.62	25.50 \pm 0.63	2.46	100.67 \pm 0.65	0.64
	65	64.63 \pm 1.90	2.94	63.16 \pm 1.64	2.49	105.74 \pm 1.13	1.07
6	3	2.69 \pm 0.02	0.90	2.70 \pm 0.03	1.28	89.43 \pm 1.27	1.42
	25	23.76 \pm 0.49	2.06	24.26 \pm 0.55	2.26	96.79 \pm 2.19	2.27
	65	63.06 \pm 1.12	1.77	64.24 \pm 1.56	2.43	92.69 \pm 0.83	0.89