

Supplementary Materials: Single-Domain Antibody Multimers for Detection of Botulinum Neurotoxin Serotypes C, D, and Their Mosaics in Endopep-MS

Michiel M. Harmsen, Jan C. Cornelissen, Fimme J. van der Wal, Jan H. W. Bergervoet and Miriam Koene




Figure S1. Reducing SDS-PAGE analysis of yeast-produced VHH monomers and multimers. Each panel represents an individual gel. Some VHHs were deglycosylated by treatment with endoglycosidase H prior to SDS-PAGE. In the corresponding lanes the position of the endoglycosidase H enzyme (endoH) is indicated (30 kDa; black arrow). N-glycosylated BC24F and N-glycosylated multimers containing BC24 are also indicated (big red arrow). The molecular weight of individual marker proteins is indicated next to each gel.




Figure S3. Analysis of VHH binding to inactivated recombinant BoNT by biolayer interferometry using an Octet Red96 biosensor. Biotinylated VHHs were coupled to SAX sensors. The graphs show the three incubations of VHH loaded sensors, demarcated by arrows, in buffer until $t=0$ (baseline), association phase with BoNT protein and dissociation phase in buffer. (A-F) Affinity determination of individual VHHs binding to different BoNT proteins, as indicated on top. Each panel depicts data for one mono- or multimeric VHH interacting with a specific BoNT protein at different concentrations that were used for global curve fitting to a 1:1 interaction model. Red lines are the fitted curves. BC134BC147F binding to BoNT/DCi (F) resulted in a relatively bad curve fit ($R^2 = 0.96$) while all other affinity measurements resulted in good curve fitting ($R^2 > 0.98$). (G,H) Scouting experiment with 30 nM BoNT/DCi (G) or BoNT/Ci (H) of sensors loaded with 6 biotinylated VHH monomers (1 $\mu\text{g/ml}$) that were also used for construction of multimers. BC126F, BC141F and BC147F do not bind BoNT/DCi.




Figure S4. Effect of 2 M NaCl on inactivated recombinant BoNT dissociation from VHHs analysed by biolayer interferometry using an Octet Red96 biosensor. Biotinylated VHH monomers (A,C,E,G) or multimers (B,D,F,H) were loaded on SAX biosensors and a baseline in PBST was done from -300 to 0 sec. Then 20 nM BoNT/Ci (A,B), BoNT/DCi (C,D) BoNT/Di (E,F) or BoNT/CDi (G,H) analyte was captured. BoNT dissociation was subsequently measured during 4 subsequent incubations of 600 sec in PBST, PBS / 2 M NaCl, HBS-EP and PBST. The sensorgrams show the baseline in PBST for 300 sec obtained after VHH loading, the association of BoNT from 0 to 300 sec, followed by four dissociation phases of 600 sec in four buffers.