

Supplemental Table S9. The list of primers used for Real-time PCR.

Name	Forward sequence: (5' to 3')	Reverse sequence: (5' to 3')
<i>FAXDC2</i>	TGAAGGGAGAAGCTGGAAGG	GGTCCACTCCTTCCCTTCAA
<i>INTS7</i>	CGGATTCTGATGTGCCACTG	AGAGTGCTTTCTGACTGGCT
<i>LIMS2</i>	ATGAACAACAACCTGGCACCC	CGCTGACAGATGTACTTGCC
<i>MBD5</i>	TCCTGCTGTTGTTTCCTTTGC	TCCTTCCAATTCCAGTCCCC
<i>PSEN1</i>	ACGTTGAAATATGGCGCCAA	CAGAATTGAGTGCAGGGCTC
<i>SEPTIN1</i>	CTTGACACAAACGCTGACCA	GCTCGAACTGCTCCTCAATG