

**Supplementary Table S2**

Primers used to perform colony PCR and detection GFP insertion as well as TBRV-pJL89-P1-R1.

<b>GFP-TBRV clone</b>	<b>Name</b>	<b>Sequence (5'-&gt;3')</b>
MP/5/GFP/5CP MP/20/GFP/20/CP	seqMPCPF	CTGGCAACATTCCTGTTGGT
	seqMPCPR	AACACGTGTGTTCGTATTCTTGA
HP/20/GFP/20/MP	R2MJ2F	CTGATGAAGACCTATGTGATGT
	R2MJ2R	CATCCCGCTTCAGCAAGTTTA
pJL89-P1-R2 + 2A CP/2A/GFP	3RACER2PR1	TGACTGTTAGCATAGAGGTTCA
	pJL89seqR	GGGAAATTCGAGCTCTCCCT
TBRV-pJL89-P1-R1	R1KRTBRVF	GGTAAAAGTTCTGGGTGCT
	R1KRTBRVR	GCAAATCCACCTCCTTATCC