


## Supplementary figure 1: Fast QC (quality control) data

Goat ID	Basic statistics		Goat ID	Basic statistics	
<b>X18</b>	Measure	Value	<b>A101</b>	Measure	Value
	Filename	18_1.paired.fq		Filename	A101_1.paired.fq
	File type	Conventional base calls		File type	Conventional base calls
	Encoding	Sanger / Illumina 1.9		Encoding	Sanger / Illumina 1.9
	Total Sequences	18344612		Total Sequences	22740578
	Filtered Sequences	0		Filtered Sequences	0
	Sequence length	50-130		Sequence length	50-105
%GC	52	%GC	55		
<b>A002</b>	Measure	Value	<b>X066</b>	Measure	Value
	Filename	A002_1.paired.fq		Filename	066_1.paired.fq
	File type	Conventional base calls		File type	Conventional base calls
	Encoding	Sanger / Illumina 1.9		Encoding	Sanger / Illumina 1.9
	Total Sequences	20611561		Total Sequences	22153023
	Filtered Sequences	0		Filtered Sequences	0
	Sequence length	50-130		Sequence length	50-105
%GC	51	%GC	52		
<b>H095</b>	Measure	Value	<b>X043</b>	Measure	Value
	Filename	H096_1.paired.fq		Filename	043_1.paired.fq
	File type	Conventional base calls		File type	Conventional base calls
	Encoding	Sanger / Illumina 1.9		Encoding	Sanger / Illumina 1.9
	Total Sequences	19275680		Total Sequences	25440163
	Filtered Sequences	0		Filtered Sequences	0
	Sequence length	50-130		Sequence length	50-105
%GC	51	%GC	50		
<b>X248</b>	Measure	Value	<b>B15</b>	Measure	Value
	Filename	248_1.paired.fq		Filename	B15_1.paired.fq
	File type	Conventional base calls		File type	Conventional base calls
	Encoding	Sanger / Illumina 1.9		Encoding	Sanger / Illumina 1.9
	Total Sequences	22221239		Total Sequences	25789160
	Filtered Sequences	0		Filtered Sequences	0
	Sequence length	50-130		Sequence length	50-105
%GC	51	%GC	50		

Goat ID

Per base sequence quality


X18




Goat ID

Per base sequence quality


A101




A002




X066




H095




X043



X248




B15



Goat ID

Per sequence quality scores


X18




Goat ID

Per sequence quality scores


A101




A002




X066




H095




X043



X248




B15



Goat ID

Per base sequence content


X18




Goat ID

Per base sequence content


A101




A002




X066




H095




X043



X248




B15



Goat ID

Per base GC content


X18




Goat ID

Per base GC content


A101




A002




X066




H095



X043



X248



B15

