

Functionalized collagen/poly(ethylene glycol) diacrylate inter-penetrating network hydrogel enhances beta-pancreatic cell sustenance

Natalia Moreno-Castellanos*¹, Elías Cuartas-Gómez², Oscar Vargas-Ceballos³

¹ Centro de Cromatografía y Espectrometría de Masas, CROM-MASS, Universidad Industrial de Santander, Cra 27 calle 9, Bucaramanga 680002, Colombia

² CICTA Research Group, Department of Basic Sciences, Medicine School, Health Faculty, Universidad

Industrial de Santander, Cra 27 calle 9, Bucaramanga 680002, Colombia

³ GIMAT Research group, Escuela de Ingeniería Metalúrgica y Ciencia de Materiales, Universidad Industrial de Santander, Cra 27 calle 9, Bucaramanga 680002, Colombia

*Correspondence: nrmorcas@uis.edu.co




Figure S1: VEGF cumulative release profile as function of time. The *in vitro* release curve from VEGF (directly loaded) functionalized collagen/PEGDA IPN hydrogels in PBS at 37 °C expressed as a function of initial VEGF loading (n=4).