

Short Note

# 2,3-Bis(4-(3-(trifluoromethyl)-3H-diazirin-3-yl)phenyl)oxirane

Lei Wang, Takuma Yoshida, Masashi Okamoto, Yasuko Sakihama, Yasuyuki Hashidoko and Makoto Hashimoto \*

Division of Applied Bioscience, Graduate School of Agriculture, Hokkaido University; Kita 9, Nishi 9, Kita-ku, Sapporo 060-8589, Japan

\* Author to whom correspondence should be addressed; E-Mail: hasimoto@abs.agr.hokudai.ac.jp.

Received: 20 January 2014 / Accepted: 18 February 2014 / Published: 24 February 2014

**Abstract:** The title compound, which has two photoreactive groups in a molecule, was synthesized by the coupling reaction of 3-(4-(bromomethyl)phenyl)-3-(trifluoromethyl)-3*H*-diazirine in the presence of silver oxide in DMSO.

**Keywords:** diazirine; oxirane; photoaffinity label

Diazirine based photoaffinity labeling is a useful chemical biology method to elucidate interactions between low molecular weight bioactive compounds and biomolecules [1]. Photoreactive compound contained two photophores [2,3] is one of the attractive tools for comprehensive analysis of ligand-biomolecule interactions by photoaffinity labeling.

**Scheme 1.** Synthesis of 2,3-bis(4-(3-(trifluoromethyl)-3*H*-diazirin-3-yl)phenyl)oxirane (2).

$$F_3C$$
 $N$ 
 $N=N$ 
 $CF_3$ 
 $DMSO$ 
 $rt$ 
 $F_3C$ 
 $N=N$ 
 $\mathbf{2}$  (41%)

We have synthesized a new two photophores containing compound, 2,3-bis(4-(3-(trifluoromethyl)-3*H*-diazirin-3-yl)phenyl)oxirane **2**, from 4-[3-(trifluoromethyl)-3*H*-diazirin-3-yl]benzyl bromide **1** with silver oxide in DMSO [4] without decomposition of diazirinyl ring. <sup>1</sup>H- and <sup>13</sup>C-NMR chemical shifts for

Molbank **2014** M816 (Page 2)

oxirane moiety (3.83 and 62.17 ppm) strongly indicated that the purified compound is *trans*- isomer [5,6]. The 2,3-*trans* two-photophores containing oxirane ring is utilized for introduction of **2** to ligand.

## **Experimental**

Synthesis of 2,3-bis(4-(3-(trifluoromethyl)-3H-diazirin-3-yl)phenyl)oxirane (2)

To a solution of 4-[3-(trifluoromethyl)-3*H*-diazirin-3-yl]benzyl bromide (1) [7,8] (0.10 g, 0.36 mmol) in anhydrous DMSO (1 mL), silver oxide (0.09 g, 0.39 mmol) was added. The resulting mixture was stirred overnight at room temperature until TLC showed the reaction to be complete. The reaction was filtered by Celite and treated with H<sub>2</sub>O, which was extracted with diethyl ether. The combined oil layer was dried over MgSO<sub>4</sub> and subjected to rotary evaporation to remove the solvent. The crude product was purified by silica gel column chromatography (CH<sub>2</sub>Cl<sub>2</sub>/hexane, 1:1) to afford **2** (0.06 g, 41%) as a yellow solid. [ $\alpha$ ]<sub>D</sub> –20 (c 0.5, CHCl<sub>3</sub>). UV(CDCl<sub>3</sub>)  $\lambda$ <sub>max</sub> ( $\varepsilon$ ) 250 (2120), 356 (1170). <sup>1</sup>H-NMR (270 MHz, CDCl<sub>3</sub>):  $\delta$  7.37 (4H, d, J = 8.1 Hz,), 7.22 (4H, d, J = 8.1 Hz), 3.83 (2H, s). <sup>13</sup>C-NMR (68 MHz, CDCl<sub>3</sub>):  $\delta$  138.54, 129.48, 126.91, 126.02, 122.14 (q,  ${}^{I}J_{CF}$  = 274.7 Hz), 62.17, 28.24 (q,  ${}^{2}J_{CF}$  = 40.6 Hz). <sup>19</sup>F-NMR (470 MHz, CDCl<sub>3</sub>):  $\delta$  -65.22. HRMS (FD): m/z calcd for C<sub>18</sub>H<sub>10</sub>F<sub>6</sub>N<sub>4</sub>O: 412.0759; found: 412.0744.

## Acknowledgments

This research was partially supported by the Fugaku Foundation for financial support. Part of this work was performed under the Cooperative Research Program of "Network Joint Research Center for Materials and Devices".

#### **Conflicts of Interest**

The authors declare no conflict of interest.

#### **References and Notes**

- 1. Hashimoto, M.; Hatanaka, Y. Recent progress in diazirine-based photoaffinity labeling. *Eur. J. Org. Chem.* **2008**, *2008*, 2513–2523.
- 2. Blencowe, A.; Blencowe, C.; Cosstick, K.; Hayes, W. A carbene insertion approach to functionalized poly(ethylene oxide)-based gels. *React. Funct. Polym.* **2008**, *68*, 868–875.
- 3. Murai, Y.; Wang, L.; Muto, Y.; Sakihama, Y.; Hashidoko, Y.; Hatanaka, Y.; Hashimoto, M. Simple and stereocontrolled preparation of benzoylated phenylalanine using Friedel–Crafts reaction in trifluoromethanesulfonic acid for photoaffinity labeling. *Heterocycles* **2013**, *87*, 2119–2126.
- 4. Wong, F.M.; Chan, Y.M.; Chen, D.X.; Wu, W. Convenient one-pot synthesis of *trans*-1,2-diaryloxiranes from the direct coupling of benzyl halides. *Tetrahedron Lett.* **2010**, *51*, 6649–6650.
- 5. Tao, B.; Lo, M.M.-C.; Fu, G.C. Planar-chiral pyridine *N*-Oxides, a new family of asymmetric catalysts: Exploiting an η<sup>5</sup>-C<sub>5</sub>Ar<sub>5</sub> ligand to achieve high enantioselectivity. *J. Am. Chem. Soc.* **2001**, *123*, 353–354.

Molbank **2014** M816 (Page 3)

6. Okazaki, Y.; Ando, F.; Koketsu, J. Generation of sulfur ylides from sulfonium salts and their reactions. Comparative study of electrochemical reduction with the base method and mechanism elucidation by the MO method. *Bull. Chem. Soc. Jpn.* **2003**, *76*, 2155–2165.

- 7. Brunner, J.; Senn, H.; Richards, F.M. 3-Trifluoromethyl-3-phenyldiazirine. A new carbene generating group for photolabeling reagents. *J. Biol. Chem.* **1980**, *255*, 3313–3318.
- 8. Nakashima, H.; Hashimoto, M.; Sadakane, Y.; Tomohiro, T.; Hatanaka, Y. Simple and Versatile method for tagging phenyldiazirine photophores. *J. Am. Chem. Soc.* **2006**, *128*, 15092–15093.
- © 2014 by the authors; licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution license (http://creativecommons.org/licenses/by/3.0/).