Supplement Information for


An environmentally sensitive silk fibroin/chitosan hydrogel and its drug release behaviors


Zhangpeng Xu¹ · Erni Tang² · Huijing Zhao¹*

1 National Engineering Laboratory for Modern Silk, College of Textile and Clothing Engineering, Soochow University, No. 199 Ren'ai Road, Industrial Park, Suzhou 215123, China

2 China Leather and Footwear Industry Research Institute (Jinjiang) Co.,Ltd, No. 150 Wenhua Road, Hongshan Comprehensive District, Jinjiang 362200, China

Correspondence to zhhj@suda.edu.cn








Fig. SI1. (A) Standard curve of MB in deionized water, (B) Standard curve of MB in PBS (pH=7.4, 37° C), (C) Standard curve of MB in PBS (pH=2.2, 37° C).




Fig. SI2: UV-visible spectroscopy of MB




Fig. SI3: The viscosity of mixed solutions.

The dynamic viscosity of mixed solutions was determined by AR-2000 rheometer (TA Instrument, New Castle, USA), fitted with parallel plate (diameter of 20 mm) kept at a gap distance of 1 mm. The shear velocity ranged from $0\sim100$ rad/s.




Fig. SI4: The thermograms of hydrogels. (A) SF; (B) NSF; (C) CS; (D) SF/CS=10/0; (E) SF/CS=9/1; (F) SF/CS=8/2; (G) SF/CS= 7/3; (H) SF/CS=6/4